3M ESPE

Dental Products

2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Exclusively Distributed by:

3M Unitek
Orthodontic Products

2724 South Peck Road
Monrovia, CA 91016 USA

CE
0123

European Representative:

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Carl-Schurz-Straße 1
D-41453 Neuss – Germany

2015/09

3M and Unitek are trademarks of 3M or
3M Deutschland GmbH. Used under license in Canada.
© 3M 2015. All rights reserved. REF 011-645-2 1509
44-0007-7676-3-A

3M

Unitek™

Temporary Anchorage Device
Dispositif d'ancrage temporaire
Temporäre Verankerungsschraube
Dispositivo di ancoraggio temporaneo
Dispositivo de anclaje temporal
Temporary Anchorage Device (Orthodontisch implantaat
voor tijdelijke verankering)
Miniskruv
Väliaikainen ankkuriruuvi
Προσωρινή οσκευή αγκύρωσης (TAD)
Dispositivo de Ancoragem Temporária
Temporary Anchorage Device
Midlertidig forankringsanordning
Временная опора
Устройство за временно закрепване
Naprava za privremeno sidrište
Ideiglenes rögzítőkészülék
Tymczasowe urządzenie kotwiczące
Dispositiv de ancorare temporară
Pomôcka na dočasné ukotvenie
Začasni pritrdilni pripomoček
Prostředek pro dočasné kotvení
Geçici Ankrāj Cihazı
Ajutine ankurduselement
Pagaidu fiksācijas ierīce
Laikinas tvirtinimo įtaisas
Тимчасова опорна

INDICATIONS:

The Unitek™ Temporary Anchorage Device (TAD) is a 1.8 mm diameter self-tapping titanium screw. Implants are manufactured from Ti 6Al-4V ELI (Grade 23) Titanium alloy, per ASTM F-136. It is indicated to provide a fixed anchorage point for attachment of orthodontic appliances to facilitate the orthodontic movement of teeth. It is used temporarily and is removed after orthodontic treatment has been completed. Implants are intended for single use only. Do not sterilize and re-use. Doing so may result in product failure and potential patient injury.

COMMON LOCATIONS FOR UNITEK TEMPORARY ANCHORAGE DEVICES:

Note: This should only be used as a guide since soft tissue and bone thickness vary from patient to patient.

Length	Unitek™ Temporary Anchorage Device (TAD) Location
6 mm	Facial surface maxillary/mandibular alveolar ridge mesial to 1 st molar, maxillary subANS region, mandibular symphysis
8 mm	Facial surface maxillary/mandibular alveolar ridge distal to 2 nd premolar, parasagittal midpalate
10 mm	Maxillary tuberosity, zygomatic buttress, infrazygomatic crest or posterior lateral palate; mandibular ascending ramus, retromolar region, external oblique ridge

CONTRAINDICATIONS:

The Unitek™ TAD should not be placed if there is insufficient bone width and height to surround the implant. A minimum of 3.3 mm of bone thickness is necessary to ensure optimal performance.

WARNINGS:

The Unitek TAD is a single-use device. Do not re-sterilize or re-use. Doing so may result in product failure and potential patient injury.

Care must be taken to avoid the roots of teeth, major nerves, blood supply and maxillary sinus.

The adequacy of bone should be determined by radiographs, palpations and thorough visual inspection. Determine the location of all avoidable anatomical features prior to initiating any Unitek TAD procedures.

Total forces on any single Unitek TAD should not exceed 300 grams. Applying forces greater than 300 grams may cause the Unitek TAD to pull out.

Due to patient variability, the amount of force will vary depending upon the patient needs. It is not recommended to apply forces to the long axis of the Unitek TAD.

To achieve maximum results the Unitek TAD should be placed where the load is perpendicular to the screw axis (90° angle).

PRECAUTIONS:

It is recommended that Unitek TADs be placed in children over the age of 13. Implants may be used in younger patients in very select cases. Special care must be taken to avoid developing teeth.

NOTES:

Unitek TADs are provided sterile (gamma) and intended for single use only. Contents are considered sterile unless package is opened or damaged or if the expiration date has passed. Do not re-sterilize or use a TAD from a damaged package.

Powder free gloves are recommended when placing a Unitek TAD.

Special consideration should be taken when placing Unitek TADs in patients with any of the following: use of drugs, tobacco or alcohol, oral mucosal pathologies, poor oral hygiene, poor patient compliance, physical handicaps that prevent adequate oral hygiene and or maintenance, and insufficient interradicular/ intraradicular space or para-functional habits.

The titanium alloy used to manufacture MDI implants is not magnetic and thus should be acceptable for MRI procedures. However, it is important that patients advise medical professionals about the presence of dental implants prior to imaging to ensure compatibility with the MRI equipment and to address potential artifacts, depending on the area of imaging.

PATIENT HOME CARE & MAINTENANCE INSTRUCTIONS:

Inform the patient that they need to treat their new implants very much like teeth, which includes the following:

1. Regularly brush the implants to keep them free of debris.
2. Return for regular check-ups.
3. Contact the orthodontist if experiencing soreness or if implant becomes loose.

PLACEMENT PROTOCOL:

1. Remove the vial holding the implant from the outer packaging. Remove the clear vial while holding the white cap in one hand, place either the Unitek™ TAD Straight Driver (REF 504-301) or the Unitek™ TAD Contra Angle Driver (REF 504-315) over the O-Ball and around the square head so that the o-ring in the driver adapter tightly holds the Unitek TAD. The Unitek TAD is unscrewed from the white cap and ready for placement.

2. Gently screw the Unitek TAD into place using one of the Unitek TAD Drivers. **NOTE:** The Unitek TAD should be inserted until the polished collar engages the outer cortex or the square head penetrates the soft tissue by no more than 0.5 mm.

When placing a Unitek™ TAD in dense bone (usually posterior mandible), it may be appropriate to screw the Unitek TAD from 1/2 to 2 complete revolutions until pressure increases considerably, then stop for 10 to 20 seconds, allowing the bone to expand around the Unitek TAD before continuing. This respite should be repeated as often as necessary, and is usually only required for the range between 2.0 to 4.0 mm of the tapered body. After the tapered body is through the cortex and the full 1.8 mm diameter body begins to enter the bone, the bone is no longer required to expand to accommodate the increasing diameter; therefore the respites are usually no longer required.

- Place the Unitek™ TAD O-Cap (REF 504-201) to suppress alveolar/mobile mucosa (optional).
- Load the Unitek TAD by attaching the force module directly to the 0.76 mm holes, the neck or the Unitek TAD Cap. After the Unitek TAD is seated it can be loaded immediately. There is no need to wait to load for either soft tissue or bony healing.
- Remove the Unitek TAD from the bone after it is no longer needed for anchorage or tooth movement. Removal occurs without topical or local anesthetic by unscrewing the Unitek TAD using one of the Unitek TAD Drivers.

ORDERING INFORMATION:

Unitek™ Temporary Anchorage Device (TAD)			
REF	Description	Thread Length	Overall Length
504-006	Unitek TAD 6 mm	6 mm	10 mm
504-008	Unitek TAD 8 mm	8 mm	12 mm
504-010	Unitek TAD 10 mm	10 mm	14 mm

STORAGE AND HANDLING

- Device should be stored at room temperature.
 - All titanium implants and parts must be handled with the appropriate Unitek TAD Driver and/or powder-free gloves to avoid contaminating the product or interfering with its biocompatibility.
 - Refer to individual product labels for special storage and handling conditions.
- For additional information please consult www.3MUnitek.com/TAD

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