



PRODUCT CLINICAL DATA SUMMARY

Product Number 9944

3M Hydrocolloid Adhesive Tape

Effective: December 2002

3M Product Number 9944 has been subjected to the following safety evaluations:

In Vitro Cytotoxicity (Agar Overlay)

Protocol reference: Guess, W. L. et al; "Agar Diffusion Method for Toxicity Screening of Plastics on Cultured Cell Monolayers" J. Pharm. Sci. 54:1545-1547 (1965). (MRFE 05798) **Results:** 1.0/1.0

Acute Primary Skin Irritation on Albino Rabbits

Protocol reference: Draize: Appraisal of the Safety of Chemicals in Food, Drugs and Cosmetics (1965). Published by the Editorial Committee of the Association of Food and Drug Officials of the United States. (MRFE 05798) **Results:** 1.5/8.0.

Acute Intracutaneous Irritation in Albino Rabbits

Protocol reference: Draize: Appraisal of the Safety of Chemicals in Food, Drugs and Cosmetics (1965). Published by the Editorial Committee of the Association of Food and Drug Officials of the United States. (MRFE 05798)

Results: No irritating leachables in test mediums, irritating when injected intracutaneously`

Acute Systemic Toxicity in Albino Mice

Protocol reference: U.S. Pharmacopoeia XXII, 1990, pg. 1499.

<88> Biological Reactivity Tests, In Vivo, Fifth Supplement, USP XXII-NF XVII, 1990. (MRFE 05798) **Results:** No toxic leachables.

Repeated Insult Patch Test (Draize) in Humans

Protocol reference: Draize: Appraisal of the Safety of Chemicals in Food, Drugs and Cosmetics (1965). Published by the Editorial Committee of the Association of Food and Drug Officials of the United States. (MRFE 05794)

Results: No subjects showed evidence of sensitization.

These tests are in accordance with the ISO 10933, Biological Evaluation of Medical Devices Part 1: Evaluation and Testing, as put forth by the FDA. Product Number 9944 has satisfied the requirements for devices in contact with intact skin for short term application (up to 29 days). All laboratory testing was conducted in accordance with the FDA Good Laboratory Practices Regulation of 1978 and in accordance with OECD guidelines.

The use of the term "hypoallergenic" has come to indicate a product that is non-sensitizing to the general public. The hypoallergenic claim for this product is supported by clinical evaluation using the repeated insult patch test in humans, commonly known as the Draize test. This protocol involves repeated application of samples on 200 healthy volunteers for a 2- to 3-week induction period, followed by a 2-week rest period and a challenge application. To be termed hypoallergenic, 3M Medical Specialties products are required to show no evidence of sensitization potential under these test conditions.

It is the responsibility of our customers to determine the final suitability of our products for their application.