
- 
-  **(EN) Dilution Pouch**
 -  **(FR) Poche de dilution**
 -  **(DE) Verdünnungsbeutel**
 -  **(IT) Busta per diluizione**
 -  **(ES) Bolsa de Dilución**
 -  **(NL) Zak voor verdunning**
 -  **(SV) Dilution Pouch**
 -  **(DA) Fortyndingspose**
 -  **(NO) Fortynningspose**
 -  **(FI) Laimennuspussi**
 -  **(PT) Saco de Diluição**
 -  **(EL) Σάκος Αραίωσης**
 -  **(PL) Worek do rozcieńczeń**
 -  **(HU) Hígító tasak**
 -  **(CS) Sáček na ředění**
 -  **(RO) Săculeț pentru diluție**
 -  **(RU) Пакет для разведений**
 -  **(TR) Seyreltme Poşeti**
 -  **(JA) 希釈水 パウチ**
 -  **(ZH) 稀释袋**
 -  **(TH) ไคลูชั่น เพาซ์**
 -  **(KO) 희석 파우치**

Product Instructions

Dilution Pouch

PRODUCT DESCRIPTION AND INTENDED USE

The 3M™ Dilution Pouch is a convenient alternative to the traditional preparation and sterilization of diluent media. The 3M Dilution Pouch prefilled with 99, 225 or 450 mL of specified dilution media consists of a standalone, heavyweight, puncture-resistant pouch which can be used for the dilution of product samples that require stomaching.

Table 1: Product Description and Storage Conditions

Catalog Number	Product Description	Volume	Storage
QDBFD99	3M™ Dilution Pouch with Butterfield's Buffer	99 mL	20-25°C
QDBFD225	3M™ Dilution Pouch with Butterfield's Buffer	225 mL	20-25°C
QDBF225F	3M™ Dilution Pouch with Butterfield's Buffer & Filter	225 mL	20-25°C
QDBFD450	3M™ Dilution Pouch with Butterfield's Buffer	450 mL	20-25°C

3M Food Safety is certified to ISO (International Organization for Standardization) 9001 for design and manufacturing.

SAFETY

The user should read, understand, and follow all safety information in the instructions for the 3M Dilution Pouches. Retain the safety instructions for future reference.

⚠ CAUTION Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury and/or property damage.

⚠ CAUTION

- The 3M Dilution Pouch is intended for use in a laboratory environment by professionals trained in laboratory techniques.
- 3M has not documented the use of this product in industries other than food. For example, 3M has not documented this product for testing water, pharmaceutical, cosmetic, clinical or veterinary samples.
- The 3M Dilution Pouches have not been evaluated with all possible food products, food processes and food processing environments, testing protocols or with all possible strains of bacteria.
- Do not use 3M Dilution Pouches in the diagnosis of conditions in humans or animals.
- 3M has not validated this product using composite samples.
- The user must train its personnel in proper testing techniques, for example: Good Laboratory Practices⁽¹⁾, ISO 7218⁽²⁾, or ISO 17025⁽³⁾.
- For information on documentation of product performance, visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.
- Consult the Safety Data Sheet for additional information.

USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

LIMITATION OF WARRANTIES / LIMITED REMEDY

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1671 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

LIMITATION OF 3M LIABILITY

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

STORAGE

Upon receipt, store unopened 3M Dilution Pouches at the indicated temperatures (See Table 1).

DISPOSAL

After use, 3M Dilution Pouches may contain microorganisms that could be a potential biohazard. Follow current industry standards or local regulations for disposal.

INSTRUCTIONS FOR USE

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Wear appropriate protective apparel and follow standard good laboratory practices GLP¹.

1. Stand the 3M Dilution Pouch on flat surface.
2. Put on gloves.
3. Tear the 3M Dilution Pouch open.
4. Pull red tabs, so the 3M Dilution Pouch stays open.
5. Aseptically add sample to the 3M Dilution Pouch.
6. The 3M Dilution Pouch can be paddle blended utilizing a paddle style blender to homogenize the contents.
7. Fold down the 3M Dilution Pouch to close.
8. Fold ends of blue wires inward.
9. Follow user established procedures.

REFERENCES

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

Refer to the current version of the standard methods listed above.

EXPLANATION OF SYMBOLS

Consult product instructions.

The lot in a box and the hourglass symbols are symbols that represent lot number and expiration date. The hourglass is followed by a year, month, and day which represents the expiration date (year, month, and day: YYYY-MM-DD). The entire line after the hourglass represents the lot number (YYYY-MM AZ).

Store between given temperatures.

3M Food Safety

3M United States

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-563-2921

3M Europe and MEA

3M Deutschland GmbH
Carl-Schurz - Strasse 1
D41453 Neuss/Germany
+49-2131-14-3000

3M Latin America

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-954-340-8263

3M Asia Pacific

No 1, Yishun Avenue 7
Singapore, 768923
65-64508869

3M Japan

3M Health Care Limited
6-7-29, Kita-Shinagawa
Shinagawa-ku, Tokyo
141-8684 Japan
81-570-011-321

3M Australia

Bldg A, 1 Rivett Road
North Ryde, NSW 2113
Australia
61 1300 363 878

3M Health Care

2510 Conway Ave
St. Paul, MN 55144 USA
www.3M.com/foodsafety

© 2014, 3M. All rights reserved.
3M is a trademark of 3M. Used under license in Canada.
34-8715-3330-2

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