


Incognito™ Appliance System and Unitek™ Digital Models



With the 3M™ True Definition Scanner and the Trusted Connection to the Incognito™ Appliance System and Unitek™ Digital Models, you can place customized orders for Incognito™ Appliances or digitally prepared models with added bases and a multitude of treatment planning tools. Access to this Trusted Connection is available to subscribers of 3M's Standard Data Plan.



How does it work?

For Incognito™ Appliances (see below for Unitek™ Digital Models)



1. Scan and Send

Scan the full arch using the 3M™ True Definition Scanner. The large touch screen serves as a virtual digital loupe, enabling you to view, rotate and magnify the 3D model to scrutinize the impression before sending through the Trusted Connection to 3M Unitek™.



2. Create Order

Log in to Unitek™ Treatment Management Portal | TMP to create a new patient then attach the digital impression to the order by selecting the patient's scan to the online prescription form. Complete patient information and add a patient photograph using the file browser. Select the blue arrow under "Create Incognito Appliance System Order."



3. Review and Approve

Review the prescription and cancel, modify or approve.



4. Produce and Place Appliances

When the plan is accepted, the Incognito™ Appliances are manufactured and shipped directly to your office. You can place the appliance, test for fit, and discuss the treatment plan with your patient.

For Unitek™ Digital Models (see above for Incognito™ Appliances)



1. Scan and Send

Scan the full arch using the 3M™ True Definition Scanner. The large touch screen serves as a virtual digital loupe, enabling you to view, rotate and magnify the 3D model to scrutinize the impression before sending through the Trusted Connection to 3M Unitek™.



2. Create Order

Select patient information. Patient information will automatically populate if the patient file was recently created. Add any missing patient information then select the blue arrow under "Create Model Storage Order."



3. Display and Analyze

Once your digital model is uploaded, measurement and analysis tools are available in Unitek™ TMP, enabling you to view your model, perform analysis for treatment planning, or simply store for future reference.



3M ESPE Dental
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3m.com/truedef

Please recycle. Printed in USA © 3M 2015.
All rights reserved. 3M and ESPE are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. 3M, Incognito, Unitek are trademarks of 3M.