- 
-  **(EN) Electronic Pipettor II - 5 mL**
 -  **(FR) Pipette Electronique II - 5 mL**
 -  **(DE) Elektronische Pipette II – 5 ml**
 -  **(IT) Pipettatrice elettronica II - 5 ml**
 -  **(ES) Pipeta electrónica Pipettor II - 5 mL**
 -  **(NL) Elektronisch pipetteerapparaat - 5 mL**
 -  **(PT) Pipeta Eletrônica II - 5 mL**
 -  **(JA) エレクトロニック ピペッター II - 5 mL**

User Manual

Electronic Pipettor II - 5 mL

TABLE OF CONTENTS

1. PRODUCT OVERVIEW.....	1
2. INTENDED USE AND SAFETY INFORMATION.....	1
INTENDED USE.....	1
USER RESPONSIBILITY.....	1
SAFETY INFORMATION.....	1
ENVIRONMENTAL OPERATING CONDITIONS.....	2
3. CONTENTS.....	3
4. FEATURES.....	3
5. INSTRUCTIONS FOR USE.....	3
Pipettor Controls.....	4
6. PIPETTOR PROGRAMMING.....	6
6.1 The Program Menu.....	6
6.2 Mode Selection.....	6
6.3 1:5 Dilution (5 mL Total Volume).....	7
6.4 1:10 Dilution (5 mL Total Volume).....	8
6.5 1:20 Dilution (5 mL Total Volume).....	9
6.6 1:5 Dilution (1 mL Total Volume).....	10
6.7 1:10 Dilution (1 mL Total Volume).....	11
6.8 1 mL Pipetting.....	13
6.9 5 mL Pipetting.....	14
6.10 Open Dilution Programming.....	15
6.11 Speed Selection.....	16
7. TROUBLESHOOTING, CARE AND MAINTENANCE.....	17
7.1 "Error" Display.....	17
7.2. Bubbles on Plate.....	17
7.3. Liquid in Tip Cone/Greasing the Piston.....	17
7.4. Recommended Care and Maintenance.....	18
7.5. Cleaning the Tip Cone.....	19
7.6. Filter Replacement.....	20
7.7. Battery Replacement.....	20
7.8. Calibration Check (Mode d1, 1:5, 5 mL dilution).....	21
8. TECHNICAL DATA.....	22
9. PIPETTOR TIPS.....	22
10. LIMITED WARRANTY.....	23
11. RECYCLING INSTRUCTION (WEEE).....	23

1. PRODUCT OVERVIEW

The 3M™ Electronic Pipettor II - 5 mL is pre-programmed to perform the most common 3M™ Petrifilm™ Plate dilutions and provide high accuracy and precision handling for micro-volumes of liquids. The microprocessor-based system reduces the possibility of human error and instrument contamination by controlling pipetting speed, calibration and volume selection. Its lightweight and ergonomic controls take the effort out of pipetting to help reduce the risk of repetitive strain injuries associated with manual pipetting.

2. INTENDED USE AND SAFETY INFORMATION

Please read, understand, and follow all safety information contained in this Operator's Manual prior to using the 3M Electronic Pipettor II - 5 mL.

Retain these instructions for future reference.

INTENDED USE

The 3M Electronic Pipettor II - 5 mL is intended to be used in food and beverage industry laboratories to deliver liquid food samples onto 3M Petrifilm Plates or other culture media and is intended only to be used by trained personnel. The product has been designed and tested to only be used with the applicable 3M approved wall plug-in transformer. If this product is used in a manner not specified, the protection provided by the 3M Electronic Pipettor II - 5 mL may be impaired.

USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

SAFETY INFORMATION

Explanation of Signal Word Consequences

⚠ WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

⚠ CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE: Indicates a situation which, if not avoided, may result in property damage or damage to the 3M Electronic Pipettor II - 5 mL itself.

Explanation of Safety Related Symbols

	WARNING: Hazardous Voltage
 	WARNING: Read Accompanying Documentation
	WARNING: Biohazard
	The 3M Electronic Pipettor II - 5 mL is subject to the WEEE EU Directive 2002/96/EC

**⚠ WARNING****To reduce the risks associated with exposure to all residual hazards listed in this document:**

- Read, understand and follow all safety information contained in the User's Manual prior to using the pipettor. Retain these instructions for future reference.

To reduce the risk associated with hazardous voltage:

- Do not use transformer if it is damaged. Replace with 3M approved part prior to use.
- Do not self service the pipettor beyond the recommended care & maintenance schedule. There are no user serviceable parts. Use only 3M authorized service personnel.
- Use only the 3M approved transformer specified for the 3M Electronic Pipettor II - 5 mL and certified for the country of use.
- Do not immerse wall plug-in transformer, stand or pipettor in any liquid.

To reduce the risk associated with fire:

- Follow all instructions and recommendations in the ethanol manufacturer's Safety Data Sheet and on the product label.
- Unplug power to the transformer or charging stand when cleaning pipettor or charging stand and allow ethanol to dry completely before restoring power.
- Do not modify the pipettor in any way to possibly short the batteries. Follow 3M instructions for battery replacement. Discharge the battery fully before disposal. Do not dispose of batteries in a fire.

To reduce the risk associated with biohazard exposure:

- Follow all start-of-day, end-of-day and other cleaning and maintenance procedures.

To reduce the risks associated with exposure to biohazard due to a false negative result:

- Ensure correct function of the pipettor prior to use.
- Use only the recommended tips to avoid adversely affecting the accuracy of the instruments.

⚠ CAUTION**To reduce the risk associated with environmental contamination:**

- Dispose of battery and electronic components in accordance with all applicable federal, state and local governmental regulations.
- Do not puncture, incinerate, crush, or disassemble battery cells. May burst or release toxic or hazardous materials. Do not short circuit. May cause burns.
- Follow proper laboratory procedures for the disposal of pipettor filters and tips and all sample and test materials.

To reduce the risk associated with personal injury:

- Never put batteries in mouth. If swallowed, contact your Physician or local Poison Control Center.

To reduce the risk associated with wrist and/or muscle strain:

- Avoid using pipettor continuously without taking occasional breaks.

NOTICE**To reduce the risk associated with cleaning and maintenance of the pipettor:**

- Always follow Good Laboratory Practices.
- Do not autoclave the 3M Electronic Pipettor II - 5 mL.
- Follow proper decontamination/cleaning procedures as indicated in this user manual.
- Never tilt the pipettor more than 15 degrees during the aspiration and dispense procedures. Do not rest the pipettor on its side. Always place the pipettor in the charging stand when not in use.
- Never strike the pipettor tip cone against the tip tray when mounting the tips. This can damage internal components.
- Avoid excessive shock. Do not drop the pipettor, charging stand or transformer.
- Avoid rough handling when using the unit. Moderate keystroke pressure is all that is required for all operating buttons.
- Avoid exposing the unit to extreme changes in temperature, humidity, and dust.
- Use only the recommended grease to avoid causing permanent damage to the piston.
- Use only the recommended tips to avoid adversely damaging the tip cone and ejector mechanism.

ENVIRONMENTAL OPERATING CONDITIONS

Environmental Conditions	Range of Limit
Indoor use only	
Altitude	Maximum 2,500 meters (8,202 Feet) Up to 2,000 meters (6,562 Feet) in Japan and Australia
Storage Temperature	-20 to +40°C (-4 to +104°F)
Mechanical Operating Temperature	+15 to +35°C (+59 to +95°F)
Performance Testing Temperature	+19 to +25°C (+66 to +77°F)
Maximum Humidity	80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C
Pollution Degree	2
Transient Overvoltages	Overvoltage Category II
Mains Supply voltage fluctuations up to ±10% of the nominal voltage.	

3. CONTENTS

- 1 3M Electronic Pipettor II - 5 mL
- 1 3M NiMH battery
- 1 Charging Stand
- 1 Transformer (AC Adaptor)
- 1 Test Performance Certificate
- 1 Important Safety Information
- 1 Quick Reference Card
- 1 Bag replacement Tip Cone Filters (50/bag)

Make sure all items are included and that no damage has occurred during shipment. Also check that the Transformer matches your electrical requirements. (See adapter or Technical Data on Section 5).





4. FEATURES

- Pre-programmed pipetting and dispensing modes to perform common dilutions within the tip
- Open programming mode for user-set pipetting and diluting
- Automatic detection of calibration errors
- Nine speeds for aspiration and dispensing accuracy
- Electronic tip ejection buttons for both right and left handed users
- Soft-touch operating buttons
- LCD screen for viewing speed, mode and volume settings
- Cordless operation
- Recharging stand
- Replaceable NiMH battery
- Ergonomic design

The 3M Electronic Pipettor II - 5 mL meets the specifications outlined in the Standard Methods for the Examination of Dairy Products as well as the FDA 2400 sheet criteria.

5. INSTRUCTIONS FOR USE

ATTENTION: Allow the pipettor to remain in the charging stand for 12 hours to fully charge the battery if the battery is new or has been allowed to fully discharge.

Connecting the Transformer


- Ensure that voltage rating of the transformer correctly corresponds to the voltage from the electrical outlet.
- Do not use electrical outlets that are timed or switched.
- Use only the transformer provided with the unit or an approved replacement from 3M.
- Connect the transformer into the connector jack at the rear of the charging stand's base.
- Plug the transformer into an AC electrical outlet. **When charging, the battery symbol will be displayed with the +/- symbols flashing in the pipettor's display.**
- When the battery is full the battery symbol disappears from the display. If the display shows **RESET** or **PRESS TIP**, press the tip ejection button to reset the pipettor. The pipettor will then be ready for use.

Pipettor Controls

These controls are often referred to this manual.


Keyboard

Select


- Changes aspiration and dispense speeds.

Arrow up/down

- Shifts to the next mode.
- Increases and decreases the open volume or speed.


Enter

- Confirms mode, speed or volume changes.


Operating buttons:

Start Button


- Triggers all piston movements.

Tip ejection buttons

- Press twice to remove tip.

Reset

- Press to reset the pipettor.


Always follow these recommendations for best results:

1. Review Safety Information section, Section 2.
2. Hold the pipettor vertically during use. For greatest accuracy, the pipettor should never be tilted at an angle greater than 15° during aspiration or dispense.
3. During aspiration, consistently place the tip no deeper than 5 mm into the liquid.
4. During dispensing, position the tip at least 5 mm above the surface to prevent bubbles from forming.

Charging the Pipettor Battery

The battery is protected against discharging during shipping with a plastic protector. Remove the plastic protector before use (Fig. 1).

1. Push the buttons simultaneously in order to open the cover.
2. Remove the plastic protector.

- Mount the pipettor into the charging stand. The charging stand operates through the metal contacts in the sides of the 3M Electronic Pipettor II and the head of the charging stand.
- Allow the pipettor to remain in the charging stand for 12 hours to fully charge the battery if the battery is new or has been allowed to fully discharge.
- When the pipettor is not in use it is recommended to keep the pipettor in the charging stand. This ensures that the battery is kept charged and the pipettor is ready for use.

6. PIPETTOR PROGRAMMING

6.1 The Program Menu

1. Preset 3M Petrifilm Plate Dilution Modes

Mode	Dilution Ratio	Diluent Volume	Sample Volume	Total Volume
d1	1:5	4.00 mL	1.00 mL	5.00 mL
d2	1:10	4.50 mL	0.50 mL	5.00 mL
d3	1:20	4.75 mL	0.25 mL	5.00 mL
d4	1:5	0.80 mL	0.20 mL	1.00 mL
d5	1:10	0.90 mL	0.10 mL	1.00 mL

2. Repeat Dispense 1 mL Mode

This mode performs repeat 1 mL dispensing.

Mode	Volume	Repetitions
rd	1.00 mL	1-5

3. Pipette Modes

These modes perform standard 1 mL and 5 mL pipetting with blow-out.

Mode	Volume
P1	1.00 mL
P5	5.00 mL

4. Open Dilution Mode

This mode allows selection of any dilution ratio within the range of the pipettor.

Mode	Diluent Volume Range	Sample Volume Range	Increment	Max. Total Volume
d0	0.10-5.00 mL	0.00-4.90 mL	0.05 mL	5.00 mL

6.2 Mode Selection

Selecting available modes (d1, d2, d3, d4, d5, rd, P1, P5, d0) on the pipettor is done with use of these two keys:


Arrow up/down-key

Enter-key

1. Press Start-button  once to activate the pipettor.

2. Press the Arrow up/down-button  repeatedly to view each of the modes available.


3. Press Enter-key  when the desired mode is displayed. This will enter the mode selection and allows you to use this mode.

 or  or  may not be used in the middle of the aspiration or dispensing cycle.


6.3 1:5 Dilution (5 mL Total Volume)

1. Press  repeatedly until the display reads:


2. Press  to enter selection
3. Attach the pipette tip to the tip cone

The pipettor is now ready to begin the 1:5 dilution (5 mL total volume).
The display indicates that the pipettor is ready to aspirate 4.00 mL.


4. Place the tip in the diluent and press the Start-button. The unit will aspirate the 4.00 mL. 

The following display indicates that the unit is ready to aspirate 1.00 mL.


5. Place the tip in the sample solution and press the Start-button to aspirate 1.00 mL. 

The next display indicates that the unit is ready to dispense the total volume (5 mL) with blow-out.

6. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the full amount from the tip. 


The display now indicates that the unit is ready to repeat the dilution sequence, beginning with the 4.0 mL aspiration.

7. Eject the disposable tip by pressing the tip ejector button twice, and repeat the series beginning with step 3.


6.4 1:10 Dilution (5 mL Total Volume)

1. Press  repeatedly until the display reads:

2. Press  to enter selection.
3. Attach the pipette tip to the tip cone.


The pipettor is now ready to begin the 1:10 dilution (5 mL total volume).
The display indicates that the pipettor is ready to aspirate 4.50 mL.

4. Place the tip in the diluent and press the Start-button. The unit will aspirate the 4.50 mL.


The following display indicates that the unit is ready to aspirate 0.50 mL.

5. Place the tip in the sample solution and press the Start-button to aspirate 0.50 mL.


The next display indicates that the unit is ready to dispense the total volume (5 mL) with blow-out.

6. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the full amount from the tip.


The display now indicates that the unit is ready to repeat the dilution sequence, beginning with the 4.50 mL aspiration.

7. Eject the disposable tip by pressing the tip ejector button twice, and repeat the series beginning with step 3.


6.5 1:20 Dilution (5 mL Total Volume)

1. Press  repeatedly until the display reads:

2. Press  to enter selection.
3. Attach the pipette tip to the tip cone.


The pipettor is now ready to begin the 1:20 dilution (5 mL total volume).
The display indicates that the pipettor is ready to aspirate 4.75 mL.

4. Place the tip in the diluent and press the Start-button. The unit will aspirate the 4.75 mL.


The following display indicates that the unit is ready to aspirate 0.25 mL.

5. Place the tip in the sample solution and press the Start-button to aspirate 0.25 mL.


The next display indicates that the unit is ready to dispense the total volume (5 mL) with blow-out.

6. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the full amount from the tip.

The display now indicates that the unit is ready to repeat the dilution sequence, beginning with the 4.75 mL aspiration.

7. Eject the disposable tip by pressing the tip ejector button twice, and repeat the series beginning with step 3.

6.6 1:5 Dilution (1 mL Total Volume)

1. Press  repeatedly until the display reads:

d 4
1.50 1 mL

2. Press  to enter selection.
3. Attach the pipette tip to the tip cone.

The pipettor is now ready to begin the 1:5 dilution (1 mL total volume).
The display indicates that the pipettor is ready to aspirate 0.80 mL.

d 4
1 0.80 mL
▲
SPEED 5
IN

4. Place the tip in the diluent and press the Start-button. The unit will aspirate the 0.80 mL. 

The following display indicates that the unit is ready to aspirate 0.20 mL.

d 4
2 0.20 mL
▲
SPEED 5
IN

5. Place the tip in the sample solution and press the Start-button to aspirate 0.20 mL. 

The next display indicates that the unit is ready to dispense the total volume (1 mL) with blow-out.

d 4
1.00 mL
▼
SPEED 5
OUT

6. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the full amount from the tip. 


The display now indicates that the unit is ready to repeat the dilution sequence, beginning with the 0.80 mL aspiration.

d 4
1 0.80 mL
▲
SPEED 5
IN

7. Eject the disposable tip by pressing the tip ejector button twice, and repeat the series beginning with step 3.


6.7 1:10 Dilution (1 mL Total Volume)

1. Press  repeatedly until the display reads:


2. Press  to enter selection.
3. Attach the pipette tip to the tip cone.


The pipettor is now ready to begin the 1:10 dilution (1 mL total volume).
The display indicates that the pipettor is ready to aspirate 0.90 mL.


4. Place the tip in the diluent and press the Start-button. The unit will aspirate the 0.90 mL. 

The following display indicates that the unit is ready to aspirate 0.10 mL.


5. Place the tip in the sample solution and press the Start-button to aspirate 0.10 mL. 

The next display indicates that the unit is ready to dispense the total volume (1 mL) with blow-out.

6. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the full amount from the tip. 

The display now indicates that the unit is ready to repeat the dilution sequence, beginning with the 0.90 mL aspiration.

7. Eject the disposable tip by pressing the tip ejector button twice, and repeat the series beginning with step 3.

Repeat Dispensing (1-5 repetitions of 1 mL)

1. Press  repeatedly until the display reads:

rd	
4	-1 mL

2. Press  to enter selection.

3. Press  to select a new number of 1 mL repetitions (1-5). In this case, it will be set to three.

COUNT	
3x	1.00 mL

4. Attach the pipette tip to the tip cone.

The pipettor is now ready to aspirate 3.00 mL (3x1 mL).

rd	
3x	1.00 mL
▲	
SPEED 5	
IN	

5. Place the tip in the sample and press the Start-button. The unit will aspirate the 3.00 mL. 

The following display indicates that the unit is ready to dispense the first of the three 1.00 mL aliquots.

rd	
COUNT	
3x	1.00 mL
▼	
SPEED 5	
OUT	

6. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the first of the three 1.00 mL aliquots. 

The following display indicates that unit is ready to dispense the second of the three 1.00 mL aliquots:

rd	
COUNT	
2x	1.00 mL
▼	
SPEED 5	
OUT	

7. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the second of the three 1.00 mL aliquots. 

The following display indicates that unit is ready to dispense the third of the three 1.00 mL aliquots:

rd	
COUNT	
1x	1.00 mL
▼	
SPEED 5	
OUT	

8. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the third 1.00 mL aliquot. 


The following display indicates that unit is ready to perform a new repeat dispensing series.

rd	
3x	1.00 mL
▲	
SPEED 5	
IN	

9. Eject the disposable tip by pressing the tip ejector button twice, and repeat the series beginning with step 4.


6.8 1 mL Pipetting

1. Press  repeatedly until the display reads:


2. Press  to enter selection.
3. Attach the pipette tip to the tip cone.

The pipettor is now ready to aspirate 1.00 mL.


4. Place the tip in the sample and press the Start-button. The unit will aspirate the 1.00 mL. 

The next display indicates that the unit is ready to dispense the total volume (1 mL) with blow-out.

5. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the full amount from the tip. 


The display now indicates that the unit is ready to repeat the 1 mL pipetting, beginning with the 1.00 mL aspiration.

6. Eject the disposable tip by pressing the tip ejector button twice, and repeat the series beginning with step 3.

6.9 5 mL Pipetting


1. Press  repeatedly until the display reads:

2. Press  to enter selection.


3. Attach the pipette tip to the tip cone.

The pipettor is now ready to aspirate 5.00 mL.


4. Place the tip in the sample and press the Start-button. The unit will aspirate the 5.00 mL. 

The next display indicates that the unit is ready to dispense the total volume (5 mL) with blow-out.

5. Position the pipette tip to dispense onto the 3M Petrifilm Plate and press the Start-button. The unit will dispense the full amount from the tip. 

The display now indicates that the unit is ready to repeat the 5 mL pipetting, beginning with the 5.00 mL aspiration.

6. Eject the disposable tip by pressing the tip ejector button twice, and repeat the series beginning with step 3.

6.10 Open Dilution Programming

1. Press repeatedly until the display reads:

d 0
SPEC mL

2. Press to enter selection.

The display shows the last used setting for the diluent aspiration volume.

1 0.60 mL

3. Press to select a new diluent volume (0.10-5.00 mL). (In this case 4.50 mL).

1 4.50 mL

4. Press to enter selection.

The display shows the last used setting for the sample aspiration volume.

2 0.30 mL

5. Press to select a new sample volume. (In this case 0.50 mL).

2 0.50 mL

6. Press to enter selection.

7. Attach the tip to the tip cone.

The pipettor is now ready to begin the dilution. The display indicates that the pipettor is ready to aspirate the selected diluent volume.

d 0
1 4.50 mL
▲
SPEED 5
IN

8. Place the tip in the diluent and press the Start-button. The unit will aspirate the diluent volume (4.50 mL) selected.


The following display indicates that the unit is ready to aspirate the selected sample volume.

d 0
2 0.50 mL
▲
SPEED 5
IN

9. Place the tip in the sample solution and press the Start-button to aspirate 0.50 mL.


The next display indicates that the unit is ready to dispense the total volume with blow-out.

10. Position the pipette to dispense and press the Start-button. The unit will dispense the full amount from the tip.

The display now indicates that the unit is ready to repeat the dilution sequence, beginning with the diluent volume aspiration.


11. Eject the disposable tip by pressing the tip ejector button twice, and repeat the series beginning with step 7.

6.11 Speed Selection

Speed changes are made using these three keys. The 3M Electronic Pipettor II must be in one of the pre-programmed (d1, d2, d3, d4, d5, rd, P1, P5, d0) modes to modify the speed selection:


Speed “1” is lowest and Speed “9” is fastest.

1. Press  to display the aspiration speed.

2. Press  to select a new aspiration speed.

3. Press  when the desired aspiration speed is displayed.
Dispensing speed is then displayed.

4. Press  to select a new dispensing speed.

5. Press  when the desired aspiration speed is displayed. The display then returns to the diluent aspiration readout.
Speed Selection function may not be used in the middle of any aspiration or dispensing cycle.

7. TROUBLESHOOTING, CARE AND MAINTENANCE

7.1 “Error” Display

The 3M Electronic Pipettor II - 5 mL has a built-in monitoring program to control the pipettor mechanism. A display showing “Error” then “Press Ent.” message indicates that the pipettor has been unable to perform a function properly. Go through the following procedure to clear this display and resume pipetting:

1. Allow the pipettor to charge on the stand or through the transformer for a minimum 30 minutes.
2. Press  to clear the error message. The display will show “reset” and “press tip” messages.
3. Reset the pipettor by pressing either of the tip ejector buttons .

7.2. Bubbles on Plate


The appearance of many small bubbles on the 3M Petrifilm Plate may indicate that the dispense speed is too fast. Lower the dispense speed (see Speed Selection, Section 6.11).

In more viscous samples, a single large bubble may form if the tip is allowed to contact the sample surface during dispensing. To prevent this, position the tip at least 5 mm above the surface during dispensing.


7.3. Liquid in Tip Cone/Greasing the Piston

If liquid is aspirated beyond the tip cone filter, clean the tip cone as follows:

Disassembling:

1. Remove the battery pack (1).
2. Unscrew the tip ejector collar (2) counter-clockwise and remove it.
3. Unscrew the black tip ejector assembly (3) counter-clockwise and remove it.
4. Unscrew the exposed piston (4) counter-clockwise and remove it.
5. Moisten a soft, lint-free wipe with 70% ethanol and gently clean the tip cone. The interior of the tip ejector collar can be cleaned with a cotton swab.
6. Let the parts dry.
7. Apply a thin layer of grease to the surface of the piston (4) on the black seal at the lower part of the piston.

NOTE: Use only the piston grease specially recommended by your 3M Food Safety representative. Other types of grease can permanently damage the piston.

Reassemble the unit and perform a calibration check:

8. Attach the piston (4) by screwing it clock-wise.
9. Attach the black tip ejector assembly (3). Replace the tip cone filter.
10. Attach the tip ejector collar (2) by screwing it clockwise.
11. Replace the battery pack (1). Reset the pipettor and press the start button several times to test the movement. Test the tip ejector operation.

NOTE: This same procedure is used to grease the piston. Monthly greasing of the piston is recommended.

7.4. Recommended Care and Maintenance

The 3M Electronic Pipettor II - 5 mL requires regular cleaning to ensure trouble-free operation. Use a soft cloth lightly moistened with a mild non-abrasive detergent to clean the outer surface of the pipettor. DO NOT AUTOCLAVE.

Recommended Care & Maintenance Schedule


Procedure	Recommended Schedule
Clean Tip Cone	Weekly; more frequently with heavy use. Note: Clean immediately if liquid enters the shaft
Replace Filter	Weekly; more frequently with heavy use. *Note: replace if pipettor is not performing due to a clogged filter, or if sample fluid has contaminated the filter.
Clean & Grease Piston	Monthly; more frequently with heavy use. *Note: if liquid enters the shaft, clean and grease the piston as part of the decontamination procedure.
4-Sample Calibration Check	Monthly; more frequently with heavy use.
10-Sample Calibration Check	Quarterly; more frequently with heavy use.
Replace Battery	**Replace after one year, more frequently with heavy use

****Heavy use is six or more hours of continuous use daily.**


7.5. Cleaning the Tip Cone

Weekly cleaning of the tip cone is recommended. However, with heavy use, clean more frequently. If liquid enters the shaft, clean the tip cone immediately.

Disassembling:

1. Remove the battery pack (1).
2. Unscrew the tip ejector collar (2) counter-clockwise and remove it.
3. Unscrew the black tip ejector assembly (3) counter-clockwise and remove it.
4. Unscrew the exposed piston (4) counter-clockwise and remove it.
5. Moisten a soft, lint-free wipe with 70% ethanol and gently clean the exposed parts. The interior of the tip ejector collar and tip cone can be cleaned with a cotton swab.
6. Let the parts dry.
7. Apply a thin layer of grease to the surface of the piston (4) on the black seal at the lower part of the piston.

Reassembling:

8. Attach the piston (4) by screwing it clockwise.
9. Attach the black tip ejector assembly (3). Replace the tip cone filter.
10. Attach the tip ejector collar (2) by screwing it clockwise.
11. Replace the battery pack (1). Reset the pipettor and press the start button several times to test the movement. Test the tip ejector operation.

7.6. Filter Replacement

Monthly or weekly filter replacement is recommended, depending on use. If the pipettor is not performing due to a clogged filter or if sample fluid has contaminated the filter, the filter should be replaced.

1. Remove the contaminated or clogged filter by pulling it straight out of the end of the tip cone (see Figure 1).
2. Moisten a swab with 70% ethanol and clean the outside surface and the inner filter placement surface of the tip cone (see Figure 2).


Figure 1


Figure 2

3. Place the pipettor in the charging stand for 15 minutes to allow drying of the tip cone.
4. Using dry, clean, gloved hands, place a new filter into the tip cone (see Figure 3).


Figure 3

7.7. Battery Replacement

Replacement of the battery is recommended after 12 months, depending on use. Replace only with the appropriate 3M NiMH battery. If the battery does not hold a sufficient charge for proper operation of the pipettor, follow these steps for replacement of the battery.

1. Hold the pipettor in a horizontal position with the display facing up.
2. Press the two charging contacts simultaneously to release the battery pack.
3. Place the new battery pack into the battery compartment and close the battery cover by pressing it lightly.
4. Dispose of the old battery in accordance with all applicable governmental regulations.

7.8. Calibration Check (Mode d1, 1:5, 5 mL dilution)

Genuine 3M Electronic Pipettor Tips are strongly recommended for verifying pipettor performance. Perform a single-sample measurement following a change of any removable parts. Perform a 4-sample test monthly-or more frequently, depending on use. Perform a 10-sample test quarterly. Record all results and maintain records. Use the procedure indicated by the American Public Health Association (APHA) "Standard Methods for the Examination of Dairy Products".

1. Allow the pipettor, pipettor tips and test liquid to adjust to room temperature (19°-25°C/66°-77°F) for one day before beginning the calibration check.
2. **Note the temperature before starting the check.** You will need this to compare your result to the calibration check chart.
3. Double-rinse the pipettor tip by aspirating water to the 5.0 mL level (4.0 mL+1.0 mL) and then expelling it.
4. Place a 50 mL or smaller vessel containing 20-25 g of distilled water onto an analytical balance. Note its weight. Re-tare balance to 0 (zero).
5. Aspirate 4.0 mL of water from the vessel on the analytical balance. Note its weight. Re-tare balance to 0.
6. Aspirate an additional 1.0 mL from the vessel on the analytical balance. Note its weight. Re-tare balance to 0.
7. Expel the 5.0 mL of water back into the vessel on the analytical balance. Note its weight. Re-tare balance to 0.
8. Repeat steps 4-6, as desired.
9. Calculate the average weight of the 4.0 mL aspirations. Repeat this calculation for the 1.0 mL and 5.0 mL measurements.
10. Refer to Table 1, below. Find the temperature closest to the temperature you recorded during step 1. Confirm that each average weight falls within corresponding weight tolerance range for that temperature.

Table 1

Calibration Check Chart

1:5 Dilution Test Volume>	4.00 mL	1.00 mL	5.00 mL
Volume (mL) Tolerance Range>	3.968-4.032 mL	0.979-1.021 mL	4.945-5.055 mL
Temperature	Weight (g) Tolerance Range*	Weight (g) Tolerance Range*	Weight (g) Tolerance Range*
19°C	3.9585-4.0207 g	0.9767-1.0181 g	4.9332-5.0409 g
20°C	3.9577-4.0199 g	0.9765-1.0179 g	4.9322-5.0399 g
21°C	3.9569-4.0191 g	0.9763-1.0177 g	4.9312-5.0389 g
22°C	3.9557-4.0179 g	0.9760-1.0174 g	4.9297-5.0374 g
23°C	3.9549-4.0171 g	0.9758-1.0172 g	4.9287-5.0364 g
24°C	3.9542-4.0163 g	0.9756-1.0170 g	4.9278-5.0354 g
25°C	3.9530-4.0151 g	0.9753-1.0167 g	4.9263-5.0339 g
26°C	3.9522-4.0139 g	0.9751-1.0164 g	4.9253-5.0324 g
27°C	3.9510-4.0131 g	0.9748-1.0162 g	4.9238-5.0314 g
28°C	3.9498-4.0119 g	0.9745-1.0159 g	4.9224-5.0299 g
29°C	3.9487-4.0107 g	0.9742-1.0156 g	4.9209-5.0283 g
30°C	3.9475-4.0099 g	0.9739-1.0154 g	4.9194-5.0273 g

*Weights adjusted for temperature and pressure effects using the "z-factor".

8. TECHNICAL DATA

Accuracy and Precision

Test Mode: d1 1:5 (5 mL) dilution

Test Action	Precision	Accuracy
Aspirate 4.00 mL	+/- 0.6 %	+/- 0.8%
Aspirate 1.00 mL	+/- 1.6%	+/- 2.1%
Aspirate 5.00 mL	+/- 0.8%	+/- 1.1%

Tolerances

Test Action	Range
Aspirate 4.00 mL	3.968-4.032 mL
Aspirate 1.00 mL	0.979-1.021 mL
Aspirate 5.00 mL	4.945-5.055 mL

Materials

Part	Material
Tip Cone	PEI
Tip Ejector Collar	PVDF
Piston	PEI
Control Handle	ASA/PC

Battery

Rechargeable NiMH battery pack - 4.8 V, 500 mAh
Charging time minimum 12 hours for empty battery

AC-Adaptor

Input voltages and main plug 100-240 VAC 50-60 Hz 135 mA
Output voltage 8 VDC, 740 mA

9. PIPETTOR TIPS

Genuine 3M Electronic Pipettor Tips are strongly recommended for use with 3M Electronic Pipettor II - 5 mL. Other tips may adversely affect the accuracy of the instrument and/or damage the tip cone and ejector mechanism.

Item No.	Description
6489	5 mL, wide pore tip, 500/pkg, racked, sterile
6488	5 mL, wide pore tip, 500/pkg, racked, non-sterile
6487	5 mL, wide pore tip, 250/pkg, bulk, non-sterile

10. LIMITED WARRANTY

LIMITED WARRANTY

3M Food Safety Hardware ("Hardware") includes the 3M Electronic Pipettor II - 5 mL and any associated Software and other components supplied by 3M which are described in the applicable Installation and Use Guides ("User Documentation").

LIMITED WARRANTY AND DISCLAIMER

3M warrants that for one year from the date of shipment that the Hardware will substantially perform in accordance with the User Documentation. THIS WARRANTY IS VOID IF: (A) THE HARDWARE HAS BEEN REPAIRED BY PERSONS NOT AUTHORIZED BY 3M; OR (B) THE HARDWARE HAS BEEN ALTERED, MODIFIED, OR MISUSED; OR (C) THE HARDWARE IS USED WITH PRODUCTS, SUPPLIES, COMPONENTS OR SOFTWARE NOT SUPPLIED BY 3M FOR USE WITH THE HARDWARE; OR (D) THE HARDWARE OR A COMPONENT IS USED FOR OTHER USES (FOR EXAMPLE USE WITH OTHER CIRCUIT BOARDS OR SOFTWARE) OR (E) THE HARDWARE HAS NOT BEEN MAINTAINED OR USED IN ACCORDANCE WITH THE USER DOCUMENTATION. UNLESS PROHIBITED BY LAW, THIS WARRANTY IS MADE IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE, THE IMPLIED WARRANTY OF MERCHANTABILITY, OR ANY IMPLIED WARRANTY ARISING OUT OF A COURSE OF DEALING OR OF PERFORMANCE, CUSTOM OR USAGE OF TRADE. 3M DOES NOT WARRANT THAT THE SOFTWARE WILL FUNCTION ERROR FREE.

If within one year after shipment, the Hardware does not conform to the express warranty set forth above, 3M's sole obligation and User's sole remedy shall be, at 3M's option: 1) to repair or replace the non-conforming component; or, 2) refund the purchase price.

LIMITATION OF 3M LIABILITY

UNLESS PROHIBITED BY LAW, 3M WILL NOT BE LIABLE TO USER OR OTHERS FOR ANY OTHER DIRECT, INDIRECT, CONSEQUENTIAL, INCIDENTAL OR SPECIAL DAMAGES INCLUDING, FOR EXAMPLE, LOST PROFITS, BUSINESS, INVESTMENTS, OR OPPORTUNITIES EVEN IF 3M HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. The parties agree that 3M's total cumulative liability to User for direct damages for all causes shall not exceed One Hundred Dollars, or the price paid for the Hardware, whichever is higher. Some states or countries may have laws which require liability rights different from those stated above. In such state or countries, the minimum required liability terms shall apply. To request service, you must request a Returned Material Authorization (RMA) number from the service center.

Please call 1-800-328-1671

11. RECYCLING INSTRUCTION (WEEE)

In compliance with European Directive, WEEE (2002/96EC) on waste and reduction of hazardous substances in electrical and electronic equipment, this device must not be recycled as unsorted municipal waste. Instead this device must be collected separately in accordance the local recycling regulations. The solid bar used in conjunction with the crossed-out wheeled bin indicates that the product was put on the European market after 13 August 2005.

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3M Deutschland GmbH
3M Health Care Limited