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
Model 505

The Bair Hugger® Total Temperature Management® System was developed by an anesthesiologist to

製品技術 1.  
prevent and/or treat the common but significant problem of hypothermia. Examples of current

applications for Bair Hugger Total Temperature Management Systems are post anesthesia care units

## Bair Hugger® Total Temperature Management System

(PACU), recovery rooms, operating rooms, emergency departments, obstetrical suites, and intensive

care areas. The warming unit draws ambient air through a filter and warms the air to the specified

temperature. It then delivers the warmed air through a hose to the Bair Hugger blanket which is

placed over the patient. When used properly, the Bair Hugger blanket distributes the warm air around

the patient's body, creating a warm environment. This manual includes instructions for operating

and maintaining Bair Hugger warming units, instructions for operating 241® fluid warming sets, and

specifications for both products. For information on how to use the unit with a specific blanket

model or the 241 fluid warming set refer to the instructions included with Bair Hugger blankets and

241 fluid warming sets.

## Bair Hugger® Total Temperature Management® System

The Bair Hugger Total Temperature Management System consists of a disposable blanket and a warming unit. Bair Hugger warming therapy can also include 241® fluid warming therapy, described under *Warming Fluids Using the 241 Fluid Warming Set*.

### Bair Hugger® Blankets

Bair Hugger blankets are made up of long, tubular channels which deliver controlled warmth. The blanket is designed to "hug" the patient. Blankets are designed in various configurations for specific applications. Follow the instructions provided with Bair Hugger blankets for specific information concerning their recommended use.


Figure D. Model 505 warming unit


Figure E. Model 505 warming unit attached to an I.V. pole


Figure F. Model 500/OR and Model 500/OR,E warming unit

### Warming Units

The warming unit uses a high-efficiency motor, a heating element, and a solid-state temperature control to create a continuous flow of warm air to the blanket. Models 500/OR, 500/OR,E and 505 are designed for safe use in all areas, including the operating room.

The folding handle on the Model 500/OR, 500/OR,E warming units can be placed in two positions. The handle is moved forward and down for the folded position; it is pushed up and back for the upright position. In the folded position the warming unit can be rolled out of the way (for example, under the operating table) during warming therapy, and it can be stored more conveniently when warming therapy is completed. The upright position allows the warming unit to be easily transported between clinical areas.

The Model 505 warming unit can be attached to an I.V. pole or to the railing on a bed.

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## Hypothermia: A Common Problem in Patient Care

### Description

Hypothermia occurs when a patient's body temperature drops below 36°C (96.8°F). Sixty to eighty percent of all patients undergoing operative procedures become hypothermic.

### Contributing Factors

Factors contributing to hypothermia include:

- Cold operating room environment
- Anesthetic drug effect
- Administration of cold intravenous fluids
- The opening of the body cavity
- Cold exposure
- Drowning
- Spinal trauma
- Geriatric thermoregulatory processes

### Adverse Consequences

Adverse consequences of hypothermia include:

- Coagulopathy
- Hemodynamic instability
- Immunodepression
- Shivering and patient discomfort
- Altered drug effect
- Postoperative nitrogen wasting
- Cardiac dysfunction

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### Indications for Bair Hugger® Warming Therapy

Examples of indications for the Bair Hugger® Total Temperature Management® system include:

- Body temperature drops below 36°C (96.8°F)
- Patient exhibits shivering
- Patient complains of being uncomfortably cold

Also, to prevent hypothermia, the Bair Hugger Total Temperature Management system should be used whenever conditions exist that could cause patients to become cold.

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## Precautionary Information

### Contraindications


Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

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### Warnings

- Do not use Bair Hugger® warming units with any forced-air blanket or cover other than Bair Hugger blankets. Thermal injury may result.
- Do not warm patients with the warming unit hose alone. Thermal injury may result. Always attach the hose to a Bair Hugger blanket before providing skin surface warming therapy.

Standard length without swivel option  
(replacement only)

Standard length with swivel option  
(replacement only)

Extended length with swivel option  
(note keyway)

Standard length 241 fluid warming hose  
(replacement for 500 series only)

Extended length 241 fluid warming hose  
(note keyway) (standard for 500 series only)

Figure Q. Hose assemblies

#### HOSE ASSEMBLIES

Warming unit hose assemblies are available in five options. The keyway on the extended length hose matches a keyway on warming units specifically calibrated for the extended length hose option, to ensure the delivery of proper temperature to the patient.

Model 505 units are only available with a snap-fit, extended length 241® fluid warming hose assembly, described in *Attaching and Storing the Model 505 warming unit hose*.

## General Maintenance

### Cabinet Cleaning

- 1) Disconnect the warming unit from the power source before cleaning.
- 2) Use a damp soft cloth and a mild detergent to clean the warming unit cabinet. Dry with a separate soft cloth.

#### CAUTION

- Do not use a dripping wet cloth to clean the cabinet. Moisture may seep into the electrical contacts, damaging the components.
- Do not use alcohol or other solvents to clean the cabinet. Solvents may damage the labels and other plastic parts.

## Specifications

### 500/OR

#### Physical Characteristics

##### DIMENSIONS

24 in. high x 16 in. deep x 14 in. wide  
61 cm high x 41 cm deep x 36 cm wide

##### WEIGHT

43 lb; 19.5 kg

### 505

#### Physical Characteristics

##### DIMENSIONS

13 in. high x 10 in. deep x 11 in. wide  
33 cm high x 25 cm deep x 28 cm wide

##### WEIGHT

11.5 lb; 5.2 kg

**500/OR****MOUNTING**

Supplied with swivel casters; floor use only

**RELATIVE NOISE LEVEL**

52 decibels

**HOSE**

Detachable, flexible, washable; compatible with Bair Hugger 241® fluid warming system

**FILTRATION SYSTEM**

0.2µM level

**RECOMMENDED FILTER CHANGE**

Every 6 months or 500 hours of use

**Temperature Characteristics****TEMPERATURE CONTROL**

Electronically controlled using a thermocouple sensor

**HEAT GENERATED**

1800 BTUs (average)

**SYSTEM TIME TO 100°F (37.7°C)**

~35 secs

**OPERATING TEMPERATURES**

*Air temperatures reaching the patient are approximately 2°C lower than the listed temperatures.*

Average temperatures at the end of the hose, assuming the back pressure of an Augustine Medical, Inc. warming blanket, or an Augustine Medical, Inc. temperature test unit:

HIGH:	43° ± 3°C	109.4° ± 5.4°F
MED:	38° ± 3°C	100.4° ± 5.4°F
LOW:	32° ± 3°C	89.6° ± 5.4°F

**Safety System****THERMOSTAT**

Independent bulb and capillary

**OVERCURRENT PROTECTION**

Dual fused input lines

**ALARM SYSTEM**

Over-heat: flashing red light with audible alarm; heater shuts down

**CERTIFICATIONS**

UL 544, CSA C22.2 No. 125, IEC 601-1, IEC 601-1-2, EN55014, AS 3200.1990

**CLASSIFICATION**

Classified under IEC 601-1 Guidelines as Class I, Type BF, Ordinary equipment, Continuous operation

**Electrical Characteristics****BLOWER MOTOR**

Operating speed: 2650 rpm  
Airflow: 28-30 cfm

**POWER CONSUMPTION**

Peak: 1000W  
Average: 500W

✓  
**505****MOUNTING**

IV pole clamp, bed rail hook with safety strap; can be placed on hard surface

**RELATIVE NOISE LEVEL**

53 decibels

**HOSE**

Detachable, flexible, washable; compatible with Bair Hugger 241® fluid warming system

**FILTRATION SYSTEM**

0.2µM level

**RECOMMENDED FILTER CHANGE**

Every 6 months or 500 hours of use

**Temperature Characteristics****TEMPERATURE CONTROL**

Electronically controlled using a thermocouple sensor

**HEAT GENERATED**

1800 BTUs (average)

**SYSTEM TIME TO 100°F (37.7°C)**

~17 secs

**OPERATING TEMPERATURES**

*Air temperatures reaching the patient are approximately 2°C lower than the listed temperatures.*

Average temperatures at the end of the hose, assuming the back pressure of an Augustine Medical, Inc. warming blanket, or an Augustine Medical, Inc. temperature test unit:

HIGH:	43° ± 3°C	109.4° ± 5.4°F
MED:	38° ± 3°C	100.4° ± 5.4°F
LOW:	32° ± 3°C	89.6° ± 5.4°F

**Safety System****THERMOSTAT**

Independent bulb and capillary

**OVERCURRENT PROTECTION**

Dual fused input lines

**ALARM SYSTEM**

Over-heat: flashing red light with audible alarm; heater shuts down

**CERTIFICATIONS**

UL 544, CSA C22.2 No. 125, IEC 601-1, IEC 601-1-2, EN55014, AS 3200.1990

**CLASSIFICATION**

Classified under IEC 601-1 Guidelines as Class I, Type BF, Ordinary equipment, Continuous operation

**Electrical Characteristics****BLOWER MOTOR**

Operating speed: 3150 rpm  
Airflow: 28-30 cfm

**POWER CONSUMPTION**

Peak: 1000W  
Average: 450W

**500/OR**

LEAKAGE CURRENT  
<100 $\mu$ A

HEATING ELEMENT  
850W Resistive

POWER CORD  
15-foot, SJT, 3 cond., 10A  
4.6 meter, HAR, 3 cond., 10A

DEVICE RATINGS  
110-120VAC, 60Hz, 9.5 Amperes, or  
220-240VAC, 50Hz, 4.5 Amperes, or  
100VAC, 50/60 Hz, 10 Amperes

FUSES  
10A, 200mA and 500mA (110 - 120 VAC Units)  
6.3A, 100mA and 500mA (220 - 240 VAC Units)  
10A, 160mA and 400mA (100VAC Units)

DIAGNOSTICS  
Over-heat test can be performed by the bio-  
medical group.

**505**

LEAKAGE CURRENT  
<100 $\mu$ A

HEATING ELEMENT  
850W Resistive

POWER CORD  
15-foot, SJT, 3 cond., 10A  
4.6 meter, HAR, 3 cond., 10A

DEVICE RATINGS  
110-120VAC, 60Hz, 9.5 Amperes, or  
220-240VAC, 50Hz, 4.5 Amperes, or  
100VAC, 50/60 Hz, 10 Amperes

FUSES  
10A, 200mA and 500mA (110 - 120 VAC Units)  
6.3A, 100mA and 500mA (220 - 240 VAC Units)  
10A, 160mA and 400mA (100VAC Units)

DIAGNOSTICS  
Over-heat test can be performed by the bio-  
medical group.

## Definition of Symbols

ON/STANDBY

ON (used on isolation switch)

OFF (used on isolation switch)

ON/OFF push button switch

Temperature Control

Equipotentiality plug (Ground)

Fuse

Warning/Caution (see appropriate documents)

Non Explosion-Proof

Dangerous Voltage

Type BF Equipment (patient applied)

Voltage, Alternating Current (AC)

## Introduction 介紹

### Description of the Total Temperature Management System

The Bair Hugger brand Total Temperature Management system consists of a Model 750 forced-air temperature management unit (with available rolling stand and sheet clip) and Bair Hugger disposable components, including forced-air blankets and the 241\* blood/fluid warming set. You can use the Model 750 temperature management unit in all clinical settings including the operating room to provide patient temperature management.

This manual includes operating instructions and unit specifications for the Model 750 temperature management unit. For information about using Bair Hugger blankets or the 241 blood/fluid warming set with the Model 750 temperature management unit, refer to the "Instructions for Use" included with each of these disposable components.

### Indications for Use 適應症

The Bair Hugger Model 750 Total Temperature Management system is intended to prevent and treat hypothermia and to provide warmth to cold or shivering patients. In addition, this system should be used whenever conditions exist that could cause patients to become cold.

### Contraindication 禁忌症

Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

### Warnings 警告

- Do not warm patients with the temperature management unit's hose alone. Thermal injury may result. Always connect the hose to a Bair Hugger blanket before providing skin-surface warming therapy.
- Use only Bair Hugger disposable components with the Model 750 temperature management unit. This temperature management unit has been designed to operate safely with Bair Hugger disposable components. Use with non-Bair Hugger products may cause thermal injury. To the full extent permitted by law, the manufacturer and/or importer declines all responsibility for thermal injury resulting from the unit being used in conjunction with non-Bair Hugger products.

## Specifications 规格

### Physical Characteristics 物理特性

DIMENSIONS	12" high x 13.5" deep x 10" wide 30 cm high x 34 cm deep x 25 cm wide
WEIGHT	15.5 lb (7kg)
RELATIVE NOISE LEVEL	53 dBA
HOSE	Flexible, compatible with the Bair Hugger brand 241 blood/fluid warming system.
FILTRATION SYSTEM	0.2 $\mu$ m air filter
RECOMMENDED FILTER CHANGE	Every 6 months or as needed.
MOUNTING	Can be clamped to an IV pole, placed on a hard surface, or mounted to the rolling stand accessory.

### Temperature Characteristics 温度特性

RECOMMENDED OPERATING ENVIRONMENT	15°C-25°C
TEMPERATURE CONTROL	Electronically controlled.
HEAT GENERATED	1644 BTU/hr (average), 482 W (average)
OPERATING TEMPERATURES	Average temperatures at the end of the hose: HIGH: 43° ± 1.5°C      109.4° ± 2.7°F MED: 38° ± 1.5°C      100.4° ± 2.7°F LOW: 32° ± 1.5°C      89.6° ± 2.7°F

### Safety System 安全系统

THERMOSTAT	Independent electronic circuit; thermal cutoff shuts the heater OFF at preset high temperature of 53°C ± 3°C at the end of the hose, back-up over-temperature detection at hose inlet.
ALARM SYSTEM	Over-temperature (53° ± 3°C): red <i>Over-temp</i> indicator light flashes, alarm sounds, heater and blower shut down, operating indicator lights turn OFF, control panel becomes unresponsive.  Fault: amber <i>Fault</i> indicator light flashes, alarm sounds, up to three fault codes stored in nonvolatile RAM.
OVERCURRENT PROTECTION	Dual input fused lines.


## Specifications – Continued

### Electrical Characteristics

HEATING ELEMENT	1400 W Resistive
LEAKAGE CURRENT	Meets UL 2601 and EN 60601-1 requirements.
BLOWER MOTOR	Operating speed: approximately 4,000 rpm Airflow: up to 48 cfm or 23 l/s
POWER CONSUMPTION	Peak: 1550 W Average: 800 W
POWER CORD	15 ft., SJT, 3 cond., 13 A 15 ft., SJT, 3 cond., 15 A 4.6 m, HAR, 3 cond., 10 A
DEVICE RATINGS	110-120 VAC, 50/60 Hz, 11.7 A, or 220-240 VAC, 50/60 Hz, 7.2 A, or 100 VAC, 50/60 Hz, 14.5 A
FUSES	12 A (110 - 120 VAC) 8 A (220 - 240 VAC) 15 A (100 VAC)
CERTIFICATIONS	EN 60601-1; EN 60601-1-2; UL 2601-1; CAN/CSA-C22.2, No. 601.1, EN 55011
CLASSIFICATION	Classified under EN 60601-1 Guidelines (and other national versions of the Guidelines) as Class I, Type BF, Ordinary equipment, Continuous operation. Not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide. Classified by Underwriters Laboratories Inc. with respect to electric shock, fire, and mechanical hazards only, in accordance with UL 2601-1 and in accordance with Canadian/CSA C22.2, No. 601.1. Classified under the Medical Device (93/42/EEC) as a Class IIb device.
DIAGNOSTICS	A qualified service technician can perform over-temperature detection system testing, temperature output testing, operating temperature calibration, and fault code troubleshooting.

# Warming units

Bair Hugger® warming units are the foundation of our forced air warming therapy. Each warming unit is individually tested to ensure years of safe, reliable and effective operation.

## Model 505

The standard of excellence in forced air warming units, the Bair Hugger model 505 provides safe, quiet and effective warming to patients around the world.

- Small size gives you more useable work space
- Light weight makes transport and setup a breeze
- Freestanding, or easily attached to IV pole, bedrail or optional rolling stand
- Built-in hour meter makes it easy to monitor usage for preventative maintenance
- Unique snap-fit hose swivels at three points for easy blanket attachment and positioning
- Compatible with 241® blood/fluid warming set — provides two warming therapies with one machine

### *Specifications*

#### **Dimensions:**

13" h x 10" w x 11" d (33 x 25 x 28 cm)

#### **Weight:**

11.5 lb (5.2 kg)

#### **Operating Temperatures:**

High: 43° ± 3°C 109.4° ± 5.4°F

Med: 38° ± 3°C 100.4° ± 5.4°F

Low: 32° ± 3°C 89.6° ± 5.4°F

#### **Leakage Current:**

Meets hospital and regulatory standards for leakage current.

#### **Filter:**

High-efficiency 0.2 µm filter

#### **Device Ratings:**

110-120 VAC, 60 Hz, 9.5 Amperes;

100 VAC, 50/60 Hz, 9.5 Amperes;

220-240 VAC, 50 Hz, 4.5 Amperes


## Bair Hugger® Temperature Management

### Model 750 - Temperature Management Unit

What do you do when you have the standard to which all other forced-air warming units are compared? Make a better one. Bair Hugger introduces the next generation temperature management unit – the Model 750.

- Hose end temperature sensing and microprocessor control for accurate air temperature delivery
- Increased airflow for quicker response
- 3 temperature sensors for increased accuracy and safety
- Built-in hour meter makes it easy to monitor usage for preventative maintenance
- Temperature settings track and report forced-air warming therapy duration
- Over-temperature, calibration and fault code servicing right on the front panel, no need to open the unit
- Compact and lightweight like its predecessor

Bair Hugger temperature management units are the **ONLY** forced-air warming systems that can be used with an integrated, cost-effective method of fluid warming - The [Bair Hugger 241® fluid warming set](#).

#### Specifications

##### Dimensions:

12" h x 10" w x 13.5" d (30 x 25 x 34 cm)

##### Leakage Current:

Meets regulatory standards for leakage current

##### Weight:

15.5 lb (7 kg)

##### Filter:

High-efficiency 0.2 µm filter

##### Certifications

EN 60601-1, UL 2601-1, CSA C22.2 No. 601.1, EN 60601-1-2

##### Operating Temperature:

High: 43° ± 1.5°C (109.4° ± 2.7°F)

Med: 38° ± 1.5°C (100.4° ± 2.7°F)

Low: 32° ± 1.5°C (89.6° ± 2.7°F)

##### Device Rating:

110-120 VAC, 50/60 Hz, 11.7 Amperes

220-240 VAC, 50/60 Hz, 7.2 Amperes

## Model 315

### Multi-Access Blanket

The multi-access blanket offers you access to any part of the patient's body while providing full patient coverage.

- Allows quick access to patient's chest, arms, torso and lower body from six convenient panels
- Non-inflated foot drape minimizes risk of thermal injury in foot and lower leg area

#### Specifications:

**Size:**  
83" x 37" (211 x 94 cm)

**Weight:**  
5.3 oz. (150 g)

Model 315  
Multi-Access

## Model 650

### Sterile Cardiac Blanket

The cardiac blanket is pre-positioned at the waist during prep, unrolled upon closure of the saphenous vein graft incision, and inflated post-bypass

- Warms the lower body to reduce core-to-peripheral heat loss to help you prevent post-bypass after-drop
- Positions easily in the sterile field
- Croin window allows for access to the femoral artery
- Non-inflated foot drape minimizes risk of thermal injury in foot and lower leg area

#### Specifications:

**Size:**  
60" x 36" (152 x 91 cm)

**Weight:**  
8.0 oz. (227 g)

Model 650  
Sterile Cardiac

### Model 522

#### Upper Body Blanket


Designed for use during surgical

procedures on the lower half of the body, the upper body blanket covers the upper torso, arms, neck and head of the patient.

- Commonly used in the supine position
- Can also be used with the patient in the prone and lithotomy position
- Two hose inlet ports provide flexibility in positioning the warming unit on either side of the patient
- Attached clear head drape to keep warm air around the insulated patient's head, and provide for visualization by the clinician (where applicable)

This blanket may be rolled and tucked to accommodate various patient positions.

The straps on each side are provided to keep the blanket in place.

#### Specifications

##### Size:


74" x 24" (188 x 61 cm)

##### Weight:

3.5 oz (99 g) with drape

##### Head Drape:

18" x 18" (45 x 45 cm)

Model 522  
Upper Body

### Model 525


#### Lower Body Blanket

The lower body blanket is designed for use primarily during surgery on the upper half of the body.


- Can be used effectively for patients in the supine, prone or lateral position
- Foot drape minimizes risk of thermal injury in feet and lower leg area

Lower body blanket being used with patient in lateral position

Lower body blanket being used with patient in prone position

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






























645012-E

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## **BAIR HUGGER BLANKETS**

### **PRODUCT DESCRIPTION and INTENDED USE**

The Bair Hugger<sup>®</sup> Blankets are part of the Bair Hugger brand Total Temperature Management System<sup>®</sup>. The Bair Hugger brand Total Temperature Management System is a group of products that are designed to prevent unintended hypothermia. A warming unit forms the base of the system by generating warmed air. The warmed air is transferred from the warming unit through a flexible hose. The hose is connected to an inflatable blanket positioned over, under or around the patient. The warm air is then passed through perforations on the underside of the blanket (that is, the side facing the patient's skin), warming the patient. The blankets are a lightweight, single use and latex free product. All blankets are classified as Class I Normal Flammability as defined by 16 CFR 1610.

Currently, there are 5 categories of blankets:

1. Pre-Op and Outpatient
2. Intraoperative
3. Postoperative
4. Pediatric
5. Specialty

#### **Pre-Op and Outpatient Blankets**

The pre-op and outpatient blankets are typically used for prewarming and to enhance patient comfort during outpatient procedures or before surgery. The blankets conform to the patient in the seated, semi-reclined or supine position. There is ample draping that tucks around the shoulders, torso, and legs for easy positioning and the blanket can follow the patient into recovery for additional thermal comfort. Models in this line are:

1. Model 110 Outpatient with Booties
2. Model 111 Outpatient

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## **Intraoperative**

The intraoperative blankets are typically used for warming patients in the operating room. The blankets are designed to provide even heat distribution in a variety of positions and during various surgical procedures. Some of the features of this line include non-inflated foot drapes (where applicable) to guard against thermal injury and attached clear head drapes (where applicable) to keep warm air around the intubated patient's head. Models in this line are:

1. Model 522 Upper Body
2. Model 523 XL Upper Body
3. Model 525 Lower Body
4. Model 540 Torso
5. Model 542 Dual Port Torso
6. Model 610 Full Body Surgical
7. Model 570 Surgical Access

## **Postoperative**

The postoperative blankets are typically used for warming patients in the Post Anesthesia Care Unit. The blankets are designed to provide even heat distribution in a variety of positions and during various post surgical procedures. All the blankets in this line feature non-inflated foot drapes to guard against thermal injury to the foot and lower leg area. Models in this line are:

1. Model 300 Full Body
2. Model 305 Chest Access
3. Model 310 Pediatric
4. Model 315 Multi-Access

## **Pediatric**

The pediatric blankets are used for warming pediatric patients throughout the perioperative setting. Pediatric blankets provide effective temperature management for the neonate, infant or child and are designed for use in a variety of positions and during various perioperative procedures. Models in this line are:

1. Model 530 Pediatric Long
2. Model 537 Small Lower Body
3. Model 550 Large Pediatric Underbody
4. Model 555 Pediatric Underbody

## **Specialty**


The specialty blankets are typically used in a variety of cardiac procedures and allow easy access for multiple positions. This line of blankets includes sterile blanket models and non-sterile models. Models in this line are:

1. Model 545 Adult Underbody
2. Model 560 Cath Lab
3. Model 630 Sterile Cardiac
4. Model 635 Full Access Underbody
5. Model 645 Sterile Cardiac Access

Handwritten blue ink markings and a circular stamp are visible in the bottom right corner of the page.

# Pre-Op and Outpatient Care

Warming patients for 30 minutes prior to surgery has been shown to add enough heat to maintain normothermia during short (less than two hours) surgical procedures.<sup>4</sup>

## Model 110

### Outpatient Blanket with booties

The Bair Hugger outpatient blanket makes it easy to give your patients the benefit of prewarming and enhance patient comfort during outpatient procedures.

- Conforms to the patient in seated, semi-reclined, or supine position
- Bair Hugger booties to cover the patient's feet are included with each blanket
- Ample draping tucks around shoulders, torso, and legs for easy positioning
- Blanket can follow patient into recovery for additional thermal comfort
- Natural latex-free, lightweight, and disposable
- Non-inflated foot drape minimizes risk of thermal injury to feet and lower leg area


### Specifications

#### Size

84 x 36 in (213 x 91 cm)

#### Weight (with booties)

5.6 oz (157 g)

## Model 111

### Outpatient Blanket

The Bair Hugger outpatient blanket makes it easy to give your patients the benefit of prewarming and enhance patient comfort during outpatient procedures.

- Conforms to the patient in seated, semi-reclined, or supine position
- Ample draping tucks around shoulders, torso, and legs for easy positioning
- Blanket can follow patient into recovery for additional thermal comfort
- Natural latex-free, lightweight, and disposable
- Non-inflated foot drape minimizes risk of thermal injury to feet and lower leg area

### Specifications

#### Size

84 x 36 in (213 x 91 cm)

#### Weight

4.8 oz (134 g)

<sup>4</sup> Steiner DE, Schneider M, Merrifield R, Mandakaw T, Cheng C. Optimal Duration and Temperature of Prewarming. *Anesthesiology*. V. 82, No. 3, Mar 1995

[bairhugger.com](http://bairhugger.com)

### Model 522 Upper Body Blanket

Designed for use during surgical procedures on the lower half of the body. The upper body blanket covers the upper torso, arms, neck, and head of the patient.


- Commonly used in the supine position
- Can also be used with the patient in the prone and lithotomy position
- Two hose inlet ports provide flexibility in positioning the temperature management unit on either side of the patient
- Patented attached clear head drape keeps warm air around the intubated patient's head and provides for visualization by the clinician

#### Specifications

**Size**  
74 x 24 in (188 x 61 cm)

**Weight**  
3.5 oz (99 g)

**Head Drape**  
24 x 24 in (61 x 61 cm)

This blanket may be rolled and tucked to accommodate various patient positions

The strips on each side are provided to keep the blanket in place

### Model 540 Torso Blanket

The torso blanket is designed to affix to the lower abdominal area and to cover the upper half of the body.


- Effectively used for patients in the supine, lateral, or prone position
- Attached clear head drape keeps warm air around the intubated patient's head and provides visualization by the clinician

#### Specifications

**Size**  
42 x 36 in (107 x 91 cm)

**Weight**  
3.2 oz (90 g)

**Head Drape**  
24 x 24 in (61 x 61 cm)

Torso blanket being used on patient in lateral position

### Model 542 Dual Port Torso Blanket


The torso blanket is also available in a two hose port configuration.

#### Specifications

**Size**  
42 x 36 in (107 x 91 cm)

**Weight**  
3.2 oz (90 g)

**Drape**  
24 x 15 in (61 x 41 cm)

## Model 525

### Lower Body Blanket

The lower body blanket is designed for use during surgery on the upper half of the body.

- Effectively used on patients in the supine, lateral, or prone position
- Foot drape minimizes risk of thermal injury in foot and lower leg area

### Specifications

#### Size

60 x 36 in (152 x 91 cm)

#### Weight

4.0 oz (113 g)

Lower body blanket being used on patient in lateral position

Lower body blanket being used on patient in prone position

## Model 610

### Full Body Surgical Blanket

The full body surgical blanket is designed for procedures on the head and neck. This lightweight blanket tapes across the patient's chest, providing a barrier to surgical sites involving the head or neck.

- Foot drape minimizes risk of thermal injury to feet and lower leg area
- Six convenient access panels allow for quick access to patient's chest, arms, and lower body

### Specifications

#### Size

72 x 36 in (183 x 91 cm)

#### Weight

4.8 oz (134 g)

## Model 570

### Surgical Access Blanket

This unique blanket is designed to be used for many surgical procedures, including but not limited to: laminectomies, spinal procedures, abdominal procedures, and transplants.

- Two hose inlet ports provide flexibility in positioning the temperature management unit on either side of the patient
- Attached clear head drape keeps warm air around the intubated patient's head and provides for visualization by the clinician
- Large surgical access window (13.5 x 22 in) (34 x 56 cm)
- Adhesive strips on patient side designed to minimize airflow into sterile field
- Non-inflated foot drape designed to protect feet from thermal injury

### Specifications

#### Size

84 x 36 in (213 x 91 cm)

#### Weight

6.4 oz (181 g)

#### Head Drape

24 x 16 in (61 x 41 cm)

## Model 315

### Multi-Access Blanket

The multi-access blanket offers access to any part of the patient's body while providing full patient coverage.

- Allows quick access to patient's chest, arms, torso, and lower body from six convenient panels
- Non-inflated foot drape minimizes risk of thermal injury to feet and lower leg area


#### Specifications

**Size**

84 x 36 in (213 x 91 cm)

**Weight**

5.3 oz (150 g)

## Model 305

### Chest Access Blanket

The chest access blanket is designed specifically for the thoracic surgery patient.

- Chest drapes provide clear, easy viewing of the chest area
- Non-inflated foot drape minimizes risk of thermal injury to feet and lower leg area

#### Specifications

**Size**


72 x 36 in (183 x 91 cm)

**Weight**

4.7 oz (133 g)

**Drape (2)**

24 x 24 in (61 x 61 cm)

## Model 300

### Full Body Blanket

The full body blanket covers the entire patient.

- Facilitates maximum thermal transfer with a blanket that "hugs" the patient
- Non-inflated foot drape minimizes risk of thermal injury to feet and lower leg area
- Soft and comfortable next to the patient's skin

#### Specifications

**Size**

84 x 36 in (213 x 91 cm)

**Weight**

5.1 oz (145 g)

## Model 310

### Pediatric Blanket

The pediatric blanket is a smaller version of the full body blanket, providing warmth to smaller patients during recovery.

- Non-inflated foot drape minimizes risk of thermal injury to feet and lower leg area


#### Specifications

**Size**

60 x 36 in (152 x 91 cm)

**Weight**

4.8 oz (134 g)

### Model 530

#### Long Pediatric Blanket

This blanket gently surrounds the child with temperature-controlled air. Use of one or two warmed clear plastic drapes keeps warmed air near the skin.

#### Specifications

<b>Size</b>	24 x 8.5 in (188 x 22 cm)
<b>Weight</b>	1.0 oz (28 g)
<b>Drape (2)</b>	24 x 18 in (61 x 46 cm)

### Model 536

#### Short Pediatric Blanket

Similar to model 530, but just 36 inches long and with only one drape.

#### Specifications

<b>Size</b>	37 x 8.5 in (94 x 22 cm)
<b>Weight</b>	0.5 oz (14 g)
<b>Drape</b>	24 x 18 in (61 x 46 cm)

### Model 550

#### Large Pediatric Underbody Blanket

The under-patient design of the large pediatric underbody blanket conveniently warms the large child or small adult patient while allowing full access for surgeons and staff.

- Fluid outlets minimize pooling of fluids on the blanket's surface
- Two drapes (included) help retain the warm air that surrounds the patient
- Hose ports at either end of blanket provide options for hose placement
- Tape and tuck flaps provide stability by securing the blanket to the OR/procedure table

#### Specifications

<b>Size</b>	Un-inflated: 60 x 32 in (152 x 81 cm) Inflated: 51 x 24 in (130 x 61 cm)
<b>Weight</b>	4.8 oz (136 g)
<b>Drape (2)</b>	24 x 24 in (61 x 61 cm)

### Model 555

#### Pediatric Underbody Blanket

The design of the pediatric underbody blanket warms even the tiniest of patients while allowing full access for surgeons and staff.

- Fluid outlets minimize pooling of fluids on the blanket's surface
- Two drapes (included) form a "tent" of warm air that surrounds the intubated patient
- Hose ports at either end of blanket provide options for hose placement
- Tape and tuck flaps provide stability by securing the blanket to the OR/procedure table

#### Specifications

<b>Size</b>	Un-inflated: 36 x 33 in (91 x 84 cm) Inflated: 31 x 29 in (79 x 74 cm)
<b>Weight</b>	3.0 oz (85 g)

#### Drape (2)

24 x 24 in (61 x 61 cm)

### Model 537

#### Small Lower Body Blanket

A smaller version of the model 525 lower body blanket, this blanket is ideal for warming the large child or small adult during upper body surgery. It can be used in the supine, lateral, or prone position.

- Can also be used as a full body blanket on smaller patient
- Non-inflated foot drape minimizes risk of thermal injury to feet and lower leg area

#### Specifications

<b>Size</b>	35 x 24 in (89 x 61 cm)
<b>Weight</b>	1.5 oz (43 g)

## Model 545

### Adult Underbody Blanket

The adult underbody blanket is conveniently positioned on the table before the patient arrives in the cath lab procedure suite. It's ready when you are.

- Fluid outlets minimize pooling of fluids on the blankets surface
- Underbody design warms the patient while allowing full access to the patient
- Unrestricted access to the femoral artery
- Radiolucent materials accommodate imaging requirements
- Non-inflated area at the foot of the blanket protects the feet from thermal injury

### Specifications



#### Size

Un-inflated: 74 x 36 in (188 x 91 cm)

Inflated: 67 x 28 in (170 x 71 cm)

#### Weight

5oz (142 g)

**Model 560** ✓  
Cath Lab Blanket

The cath lab blanket, designed for use in the Cardiac Cath Lab, surrounds the patient with warm air and allows unrestricted access to the femoral artery.

- Radiolucent materials
- Tie strips keep blanket positioned securely around patient
- Foot pad protects feet from thermal injury

**Specifications**

**Size**  
9 x 132 in (23 x 335 cm)

**Weight**  
3.0 oz (85 g)

**Model 646**

Sterile Cardiac Access Blanket with Hose

The cardiac access blanket helps prevent unintended hypothermia and allows unrestricted access to the patient.

- Warms the lower body to reduce core-to-peripheral heat loss
- Provides localized warming to the legs throughout the procedure
- Allows full access to both legs
- Insulated pouch protects feet from thermal injury
- 36 inch sterile hose allows placement of temperature management unit further from sterile field

**Specifications**

**Size**  
9 x 108 in (23 x 274 cm)

**Weight**  
8.0 oz (227 g)

**Hose**  
36 in (91 cm)

**Model 645** ✓  
Sterile Cardiac Access Blanket

The cardiac access blanket helps prevent unintended hypothermia and allows unrestricted access to the patient.

- Warms the lower body to reduce core-to-peripheral heat loss
- Provides localized warming to the legs throughout the procedure
- Allows full access to both legs
- Insulated pouch protects feet from thermal injury

**Specifications**

**Size**  
9 x 108 in (23 x 274 cm)

**Weight**  
4.0 oz (113 g)

# Specialty Blankets - Cardiac Blankets

Help protect your patients from the adverse affects of unintended hypothermia.

## Model 630 ✓

### Sterile Cardiac Blanket

The cardiac blanket is pre-positioned at the waist during prep, unrolled upon closure of the saphenous vein graft incision, and inflated post-bypass.

- Warms the lower body to reduce core-to-peripheral heat loss
- Positions easily in the sterile field
- Groin window allows for access to the femoral artery
- Non-inflated foot drape minimizes risk of thermal injury to feet and lower leg area
- 36 inch sterile hose allows placement of temperature management unit further from sterile field

### Specifications

**Size**  
60 x 36 in (152 x 91 cm)

**Weight**  
8.0 oz (227 g)

**Hose**  
36 in (91 cm)

## Model 635 ✓

### Full Access Underbody Blanket

The full access underbody blanket is positioned on the procedure table before the patient arrives in the room. This revolutionary design delivers the clinical flexibility needed for cardiothoracic and other complex surgical procedures. This advanced blanket design can be used for patients in the supine, lateral or prone positions.

- Pass-through slit design allows for flexible patient positioning and accommodates the use of a drawsheet\*
- Fluid outlets minimize pooling of fluids on the surface of the blanket
- Outer edges of blanket draw up next to the patient to surround them in warmth
- Perforated area at head of blanket can be removed to accommodate a full prone position
- Tape and tuck flaps provide stability by securing the blanket to the OR table
- Hose ports at either end of the blanket provide options for hose placement
- Foot pad protects feet from thermal injury

### Specifications

**Size**

Un-inflated: 84 x 36 in (221 x 91 cm)  
Inflated: 79 x 24 in (201 x 61 cm)

**Weight**

7 oz (198 g)

**Drape**

24 x 24 in (61 x 61 cm)

\* Drawsheet not included

# Bair Hugger® Blanket

for patient temperature management

Preoperative



(01)00608223110003

**Arizant  
Healthcare**  
bright ideas that work

**REF 11000  
with booties**

**REF 11101**

**NOT STERILE**

Model **110/111**  
Outpatient

**Use only with a Bair Hugger temperature management unit.**

**Caution: Federal Law (USA) restricts this device to sale by or on order of a licensed healthcare professional.**

**Do not use if primary packaging has been previously opened or is damaged.**

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99.10.15

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# Bair Hugger® Therapy Temperature Management Unit

Model 775

15

**Quiet. Adjustable. Safe. Effective.**  
It's everything you want from your forced-air warming unit...and more.

# Quiet, adjustable warming options for all Bair Hugger® blanket styles

Arizant Healthcare introduced forced-air warming to the world in 1987 with the launch of Bair Hugger therapy. Since that time, our Bair Hugger warming units and blankets have safely and effectively warmed more than 110 million patients.

Our latest advancement in forced-air warming – **the Bair Hugger Model 775 temperature management unit** – is the future gold standard in patient warming. It offers adjustable airflow, precise temperature delivery and ease of use from a single, quiet warming unit.

Two airflow settings and hose-end temperature sensing ensure precise temperature delivery and compatibility with all Bair Hugger blanket styles – pediatric to adult.

Three temperature sensors and adjustable airflow provide safe, quick temperature response and greater heat transfer.

The Model 775 can track and report forced-air warming therapy duration, over-temperature information, calibration, built-in hour meter and fault codes all via the front panel – no need to open the unit.

Accessories for the Model 775 warming unit allow the unit to be mounted on a hard surface, an IV pole, bed rail or rolling cart.

## Specifications

### Dimensions

13'h x 14'w x 13'd (33h x 35w x 33d cm)

### Weight

16 lb (7.3 kg)

### Operating Temperatures

High:  $43^{\circ} \pm 1.5^{\circ}\text{C}$  ( $109.4^{\circ} \pm 2.7^{\circ}\text{F}$ )

Med:  $38^{\circ} \pm 1.5^{\circ}\text{C}$  ( $100.4^{\circ} \pm 2.7^{\circ}\text{F}$ )

Low:  $32^{\circ} \pm 1.5^{\circ}\text{C}$  ( $89.6^{\circ} \pm 2.7^{\circ}\text{F}$ )

### Fan Speeds and Relative Noise

High Fan Setting: 37 CFM (53 dBA)

Low Fan Setting: 32 CFM (48 dBA)

### Leakage Current

Meets UL 60601-1 and IEC 60601-1 requirements.

EN 60601-1-2; CAN/CSA-C22.2,

No. 601.1, EN 55011; ASTM F2196-02

### Safety Systems:

Audible and visual alarms;

Over-temperature  $\leq 56^{\circ}\text{C}$ ,  $53^{\circ}\text{C} \pm 3^{\circ}\text{C}$  typical

### Filter

High-efficiency 0.2  $\mu\text{m}$  filter

### Device Ratings

110-120 VAC, 50/60 Hz, 11.7 Amperes

220-240 VAC, 50/60 Hz, 7.2 Amperes

100 VAC, 50/60 Hz, 14.5 Amperes

## Ordering Information

Bair Hugger Model 775 Temperature Management Unit

**To order, please contact your Arizant Healthcare Inc. representative or call 1-800-733-7775.**

**Arizant  
Healthcare**

©2009 Arizant Healthcare Inc. All rights reserved.  
10393 West 70th Street, Eden Prairie, MN 55344 USA  
tel 952-947-1200 | 800-733-7775  
fax 952-947-1400 | 800-775-0002  
web arizanthealthcare.com | bairhugger.com

Bair Hugger, Arizant, Arizant Healthcare and the  
Bair Hugger and Arizant logos are trademarks of  
Arizant Healthcare Inc. registered in the U.S.  
Patent & Trademark Office.

303032B 7/09