-  (EN) Enviro Swab Dry
-  (FR) Écouvillon Enaviro Swab
-  (DE) Umfeld-Tupfer Trocken
-  (IT) Tampone ambientale asciutto
-  (ES) Hisopo Enviro Swab Dry
-  (NL) Enviro Droge Swab
-  (SV) Enviro provtagningspinne, torr
-  (DA) Tør Enviro vatpind
-  (NO) Enviro Swab, tørr
-  (FI) Ympäristönäytteenottotikku, kuiva
-  (PT) Enviro Swab Seco
-  (EL) Άνυδρος Στυλεός Enviro
-  (PL) Suchy wacik do pobierania próbek Enviro
-  (HU) Száraz Enviroswab kenetvevő
-  (CS) Tampóny pro stěry z prostředí, suché
-  (RO) Tampon Enviro uscat
-  (RU) Enviro тампон сухой
-  (TR) Enviro Pamuklu Çubuk Kuru
-  (JA) エンバイロスワブ 環境用拭取りスワブ ドライ
-  (ZH) Enviro 干涂抹棒

(English)

Issue Date: 2016-06

Product Instructions

Enviro Swab Dry

PRODUCT DESCRIPTION AND INTENDED USE

A dry swab for environmental testing.

3M assumes that the 3M™ Enviro Swab Dry be utilized by technicians that have been properly trained to follow microbiological techniques.

3M has not documented the 3M Enviro Swab Dry for use in industries other than food and beverage. For example, 3M has not documented the 3M Enviro Swab Dry for testing water, pharmaceuticals or cosmetics. The 3M Enviro Swab Dry has not been evaluated with all possible surfaces or with all possible strains of bacteria.

3M Food Safety is certified to ISO (International Organization for Standardization) 9001 for design and manufacturing.

SAFETY

The user should read, understand, and follow all safety information in the instructions for the 3M Enviro Swab Dry. Retain the safety instructions for future reference.

⚠ WARNING: Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ WARNING

Do not use the 3M Enviro Swab Dry in the diagnosis of conditions in humans or animals.

The user must train its personnel in current proper testing techniques: for example, Good Laboratory Practices or ISO/IEC 17025¹.

To reduce the risks associated with bacterial infection:

- Always follow standard lab safety practices, including wearing appropriate protective apparel while handling reagents and contaminated samples.
- Avoid contact with the contents of the 3M Enviro Swab Dry.
- Dispose of contaminated 3M Enviro Swabs Dry according to applicable Federal, State, and/or local government regulations and/or applicable laboratory procedures.

Consult the Safety Data Sheet for additional information.

For information on documentation of product performance, visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.

USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

LIMITATION OF WARRANTIES / LIMITED REMEDY

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1671 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

LIMITATION OF 3M LIABILITY

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

STORAGE

When not in use, store the 3M Enviro Swab Dry between 2-25°C.

△ DISPOSAL

When testing is complete the 3M Enviro Swab Dry may contain microorganisms that may be a biohazard. Do not use the 3M Enviro Swab Dry past its expiration date. Expiration date and lot number are noted on the outside of each package.






When testing is complete, follow current industry standards for the disposal of contaminated waste. Consult the Safety Data Sheet for additional information and local regulations for disposal.

INSTRUCTIONS FOR USE

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Use a fresh 3M Enviro Swab Dry for each sample area.

1. Label the 3M Enviro Swab Dry with sample identification.
2. Remove the tab from the top of the 3M Enviro Swab Dry (Figure 1).
3. Hold the tube firmly in one hand and release the swab by pushing the top of the tube with your thumb (Figure 2).
4. Carefully remove the swab from its tube and moisten with 0.5 mL of diluent. Swab the area to be sampled (Figure 3).
5. Replace the swab in its tube.
6. Carefully lift the swab from its tube and aseptically pipette sufficient medium into the tube to cover the swab, approximately 20 mL (Figure 4).
7. Push the swab firmly back into the tube to lock the cap. Mix well by vortexing or shaking (Figure 5).
8. Incubate the 3M Enviro Swab Dry according to your method for environmental samples.

For further information refer to the 3M™ Enviro Swab Dry “Interpretation Guide.” If you have questions about specific applications or procedures, please visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.

REFERENCES

1. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

Refer to the current versions of the standard methods listed above.

EXPLANATION OF SYMBOLS

www.3M.com/foodsafety/symbols

3M Food Safety

3M United States

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-563-2921

3M Latin America

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-954-340-8263

3M Europe and MEA

3M Deutschland GmbH
Carl-Shurz - Strasse 1
D41453 Neuss/Germany
+49-2131-14-3000

3M United Kingdom PLC

Morley Street,
Loughborough
Leicestershire
LE11 1EP
United Kingdom
+(44) 1509 611 611

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© 2016, 3M. All rights reserved.
3M is a trademark of 3M. Used under license in Canada.
3M Österreich GmbH
3M Health Care Limited