RelyX™ Fiber Post**RelyX™ Fiber Post 3D**

Glass Fiber Post

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
44097598292

ENGLISH

Product Description
RelyX™ Fiber Post and RelyX™ Fiber Post 3D are radiopaque, translucent, glass fiber reinforced composite posts. RelyX™ Fiber Post 3D posts differ from RelyX™ Fiber Post in that they are outlined with retention rings and are flattened in the coronal area. Both types of posts are available in the four sizes* "0, 1, 2, and 3". Specific drills are used for removal of the root filling and preparation of the root canal; see the table below.

Combination post/drill

Post size/Color code	0/White	1/Yellow	2/Red	3/Blue
Universal drill	Universal	Universal	Universal	Universal

Removal of the root filling

Universal drill	Universal	Universal	Universal
drill	drill	drill	drill

Preparation of the root canal

Drill White	Drill Yellow	Drill Red	Drill Blue
drill	drill	drill	drill

Technical data for the posts

Diameter of apical post end (mm)	0.60	0.70	0.80	0.90
Diameter of coronal post end without retention (mm)	1.10	1.30	1.60	1.90

Taper	2.86° (5%)	3.44° (6%)	4.58° (8%)	5.72° (10%)
Length	20 mm	20 mm	20 mm	20 mm

* Not all sizes are available in all countries.

Use a dual-core composite cement, e.g., RelyX™ U200 Automix, for the cementation of the posts.

The indicated 3M ESPE composite, e.g., Filtek™ Bulk Fill Posterior composite, is used as core-build-up material. For processing details, please see the corresponding Instructions for Use of the cement and the core-build-up material.

* These instructions for Use must be kept for reference for the duration of product use.

Indications

• In case of insufficient residual tooth structure (< 4 mm) the post is needed to support and secure the coronal restoration.

RelyX™ Fiber Post and RelyX™ Fiber Post 3D are not suitable for patients with bruxism or in maxillary anterior with excessive vertical overlap or if the circular coronal tooth structure is less than 1.5 mm.

Precautionary Measures

3M MSDSs can be obtained from www.mmm.com, or contact your local subsidiary.

Precautions During Processing

• Selection of the correct drill and post size as well as their proper use are the responsibility of the user. The universal drill may be used solely to remove a root filling, and the color-coded drills may be used solely to prepare root canals. Do not continue to use drills which are dull or damaged. Replace the drills after 15 applications.

• If there has been a greater degree of destruction of the tooth, a dentin border 2 mm in width must be prepared apically of the core build-up which is covered by the definitive crown.

• Do not use any adhesives or composite cements which are light-cured only to cement the posts.

Preparatory Measures

• Treat the root canal endodontically as usual (root canal filling with gutta-percha).

• Select the post size and the matching drill as indicated by the X-rays and according to the diameter of the prepared root canal.

– To guarantee a precise fit, always use the drills and posts from a single system and do not combine them with products from other manufacturers (see color code).

• Determine the drilling depth with the help of the X-ray and according to the length of the prepared canal. 4 mm of the root filling must remain apically.

– 2/3 of the post should be located in the root canal and 1/3 in the core build-up. This length ratio cannot always be achieved if the canal is crooked.

• Prepare the drills which will be used before every application (see "Cleaning and Disinfection").

• Have a self-adhesive or adhesive resin cement, e.g., RelyX™ U200 Automix with an attached Endo Tip elongation tip, at hand for cementation of the post.

• When other composite cements are used for the adhesive cementation of RelyX™ Fiber Post and RelyX™ Fiber Post 3D, additional preparatory steps or other procedures are required. Please follow the Instructions for Use from the manufacturer in every case.

• We recommend using a cofeeders.

Prepare the Root Canal and the Post for Cementing

• As is standard when preparing the root canal for the insertion of a glass fiber post, begin by removing the gutta-percha filling with a Gates-Glidden, Peeso, or Largo bur.

• Use the universal bur to remove any residual root filling (1000–2000 rpm); see table under Product Description. Avoid any expansion beyond the apex or the perforation of the root.

• Check the fitting of the post in the prepared canal; it should fit the root canal exactly, there must not be any kind of binding, and it must be easily removable.

• Mark the final length of the post or the point where it is to be cut.

• Shorten the post extroradically to the required length using a diamond disk and suction system; wear a mask.

– Do not use scissors, scalpels, or similar tools; they will crush the post, ruining its structure and stability.

– Do not use any post fragments which have been cut off.

• Before inserting the post, clean it with alcohol and dry it with air free of water and oil.

• If you do not use RelyX™ U200 Automix cement, pre-treat the post surface according to the Instructions for Use from the cement manufacturer.

• As the final step before cementing, clean or disinfect the root canal with a 2.5–5.25% sodium hypochlorite solution (NaOCl).

– Do not use hydrogen peroxide (H_2O_2) and other disinfectants and EDTA solutions as the residues can impair the adhesion strength and curing reaction of RelyX™ U200 Automix.

• Rinse immediately with water and dry with paper tips.

Cement the Post and Build Up the Core

• Protect the working area from contamination with water, blood, saliva, and sulci fluid before beginning the cementing and through the setting phase.

• Cement the post as described in the Instructions for Use of the cement, e.g., with RelyX™ U200 Automix and attached Endo Tip elongation tip or with another dual-core self-adhesive or adhesive composite cement.

Pre-treat the remaining post surface and the tooth surface with the bonding which belongs to the core build-up material, e.g., Single Bond Universal Adhesive.

Use the desired core build-up material (e.g., Filtek™ Bulk Fill Posterior composite) according to the Instructions for Use to prepare the core build-up in a direct procedure.

After the treatment

• Carefully check the sulcus of the treated teeth and the surrounding areas; remove any cement residues still remaining.

Revision

• Remove the restoration and the post superstructure.

• Use an X-ray to check the drilling depth and position of the post.

• Drill out the post with a diamond or trepan drill.

– Insert a rubber dam to prevent the patient from inhaling any glass fiber dust. The person performing the treatment should wear a mask and use water spray and a suction unit.

– The transition from the glass fiber post to the root dentin is barely noticeable during the drilling operation; always proceed with extreme caution to avoid perforation of the root. Take X-rays at each step to check progress.

Cleaning the Posts

• The posts are intended for one-time use; their re-use is prohibited.

• Posts may be cleaned solely with alcohol and water; sterilization is not possible.

Preparing the Drills

General

The RelyX™ Fiber Post drill is in non-sterile condition when delivered and must be cleaned, disinfected, and sterilized before first use as well as before every re-use. Effective cleaning and disinfection are absolutely essential requirements for effective sterilization.

If possible, use a mechanical process (cleaning and disinfecting device) for cleaning and disinfecting the drills because the effectiveness and reproducibility are significantly lower for manual processes. Pre-treatment is required regardless of which process is used.

Material Resistance

Make sure that the cleaning and disinfectant agents you have chosen do not contain any of the following materials:

• Organic, mineral, and oxidizing acids (minimum acceptable pH value 5.5)

• Strong bases (maximum acceptable pH value 11); alkaline cleaners of about pH 10–11 (not recommended)

• Organic solvents (e.g., alcohols, ethers, ketones, benzines)

• Oxidation agents (e.g., hydrogen peroxides)

• Halogens (chlorine, iodine, bromide)

• Aromatic/halogenated hydrocarbons

• Oils

Never use metal brushes or steel wool to clean the RelyX™ Fiber Post drills.

The RelyX™ Fiber Post drill may be re-used up to 15 times, provided that it has been handled carefully and is undamaged and clean. Users go beyond this threshold of use damaged and/or contaminated RelyX™ Fiber Post drills solely on their own responsibility.

Pre-Treatment

Immediately after use (within a maximum of 2 hours), remove gross contaminations from the drills.

• Rinse the drills under running water (temperature <35°C/95°F) for at least 1 minute.

• Use a clean, soft, lint-free cloth to manually remove all visible contamination.

• Rinse again for at least 1 min under running water.

Automatic Cleaning and Disinfection

• The chemicals being used must be compatible with the drills (see "Material Resistance").

• Cleaning agents and/or disinfectants must contain alkaline detergents (pH 10–11).

Follow manufacturer's instructions for concentrations, temperatures, development times, and rinsing of the cleaning and disinfectant agents exactly.

• Immerse the RelyX™ Fiber Post drills in the cleaning and disinfecting device (basket for small parts). Make sure that there is no contact between the drills.

• Start the cleaning cycle.

• When the cycle has finished, remove the drills from the cleaning and disinfecting device.

• Check and package the drills as quickly as possible after removing them from the device (see "Packaging"); if needed, after further drying in a clean place.

Manual Cleaning and Disinfection

Alternatively, cleaning and disinfecting can be conducted manually. Information about validated procedures can be obtained from 3M Deutschland GmbH.

Sterilization

Steam sterilization is the only sterilization method which may be used for the RelyX™ Fiber Post. Drill sterilization is prohibited. Follow carefully these instructions for the sterilization.

• Maximum sterilization temperature 134 °C (273 °F), tolerances pursuant to DIN ISO 17665.

• Sterilization time (exposure time at sterilization temperature) at least 5 min or 18 min (prior deactivation) at 132 °C (270 °F)/134 °C (273 °F).

• After sterilization, allow the sterilized instruments to cool down to room temperature.

• Do not use steam sterilization for the sterilization of the RelyX™ Fiber Post.

• Do not use autoclave sterilization for the sterilization of the RelyX™ Fiber Post.

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41453 Neuss - Germany3M, ESPE, Aplicat, Filtek и RelyX are trademarks of 3M or 3M Deutschland GmbH.
© 2015, 3M. All rights reserved.
3M, ESPE, Aplicap, Filtek and RelyX are trademarks of 3M or 3M Deutschland GmbH.
3M, ESPE, Aplicat, Filtek et RelyX are trademarks of 3M or 3M Deutschland GmbH.