ENGLISH			
General Information			
3M™ Filtek™ One Bulk Fill Restorative , is a visible-light activated, restorative composite optimized to create fast and easy restorations. This material provides excellent strength and low wear for durability and improved esthetics. The material can be placed and cured up to 5mm deep, enabled by a stress relieving resin system and optimized optical properties. Filtek One Bulk Fill Restorative serves to enhance the 3M ESPE lineup of restorative materials by improving the esthetic properties of a bulk fill material to allow for broader use in both posterior and anterior restorations. Filtek One Bulk Fill Restorative is in a range of tooth-colored shades. All shades are radiopaque. The fillers are a combination of a non-agglomerated/non-aggregated 20 nm silica filler, a non-agglomerated/non-aggregated 4 to 11 nm zirconia filler, an aggregated zirconia/silica cluster filler (comprised of 20 nm silica and 4 to 11 nm zirconia particles), and a ytterbium trifluoride filler consisting of agglomerate 100 nm particles. The inorganic filler loading is about 76.5% by weight (58.5% by volume). Filtek One Bulk Fill Restorative contains AUDMA, AFM, diurethane-DMA, and 1, 12-dodecane-ADMA. Filtek One Bulk Fill Restorative is applied to the tooth following use of a methacrylate-based dental adhesive, such as manufactured by 3M ESPE, which permanently bonds the restoration to the tooth structure. Filtek One Bulk Fill Restorative is packaged in single-dose capsules and syringes.			
Indications: Filtek One Bulk Fill Restorative is indicated for use in:			
<ul style="list-style-type: none">Direct anterior and posterior restorations (including occlusal surfaces) Base/liner under direct restorations Core build-ups Splinting Indirect restorations including inlays, onlays and veneers Restorations of deciduous teeth Extended fissure sealing in molars and premolars Repair of defects in porcelain restorations, enamel, and temporaries			
Precautionary Information for Patients This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.			
Precautionary Information for Dental Personnel This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed. 3M ESPE SDS information can be obtained from www.3MESPE.com/dental or contact your local subsidiary.			
Instructions for Use Preparation 1. Prophyl: Teeth should be cleaned with pumice and water to remove surface stains.			
2. Shade Selection: Prior to isolation of tooth, select the appropriate shade(s) of Filtek One Bulk Fill Restorative using a standard VITAPAN® classical shade guide.			
3. Isolation: A rubber dam is the preferred method of isolation. Cotton rolls and an evacuator can also be used.			
Directions Direct Restoratives 4. Cavity Preparation: 4.1 Anterior restorations: Use conventional cavity preparations for all Class III, IV and V restorations. 4.2 Posterior restorations: Prepare the cavity. Line and point angles should be rounded. No residual amalgam or other base material should be left in the internal form of the preparation that would interfere with light transmission and therefore, the hardening of the restorative material.			
5. Placement of Matrix: 5.1 Anterior restorations: Mylar strips and crown forms may be used to minimize the amount of material used. 5.2 Posterior restorations: Place a thin dead-soft metal, or a pre-contoured-mylar or a pre-contoured-metal matrix band and insert wedges firmly. Burnish the matrix band to establish proximal contour and contact area. Attach the band to seal the gingival area to avoid overhangs. Note: The matrix may be placed following the enamel etching and adhesive application steps if preferred.			
6. Pulp protection: If a pulp exposure has occurred and the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide on the exposure followed by an application of 3M™ ESPE™ Vitrebond™ or Vitrebond™ Plus Light Cure Glass Ionomer. Vitrebond or Vitrebond Plus liner/bases may also be used to line areas of deep cavity excavation. 7. Adhesive System: To bond Filtek One Bulk Fill Restorative to tooth structure, use of a 3M™ ESPE™ dental adhesive system (for example 3M™ ESPE™ Single Bond Universal Adhesive) is recommended. Refer to adhesive system product instructions for full instructions and precautions for the products. After curing the adhesive, continue to maintain isolation from blood, saliva and other fluids and proceed immediately to placement of Filtek One Bulk Fill Restorative.			
8. Delivery: Dispensing the Composite: 8.1 Syringe: Dispense the necessary amount of restorative material from the syringe onto the mix pad by turning the handle slowly in a clockwise manner. To prevent oozing of restorative when dispensing is completed, turn the handle counterclockwise a half to stop paste flow. Immediately replace syringe cap. If not used immediately, the material should be protected from light.			

8.2 Single-Dose Capsule: Insert capsule into the 3M™ ESPE™ Restorative Dispenser. Use a separate restorative dispenser instructions for full instructions and precautions. Extrude restorative directly into cavity.

9. **Placement:**
 - 9.1 Avoid intense light in the working field. Exposure of paste to intense light may cause premature polymerization.
 - 9.2 Capsule: Start dispensing in the deepest portion of the preparation, holding the tip close to the preparation surface. Withdraw the capsule tip slowly as the cavity is filled, and avoid lifting the tip out of dispensed material while dispensing, to reduce voids. When dispensing has been completed, drag the capsule tip against the cavity wall while withdrawing from the operative field. For proximal areas, hold the tip against the matrix to aid material flow into the proximal box.
 - 9.3 Slightly overfill the cavity to permit extension of composite beyond gingival margins. Contour and shape with appropriate composite instruments.
 10. **Curing:** This product is intended to be cured by exposure to halogen or LED light with a minimum intensity of 550 mW/cm² in the 400-500 nm range. Cure each increment by exposing its entire surface to a high intensity visible light source, such as a 3M ESPE curing light. Hold the light guide tip as close to the restorative as possible during light exposure. Use light cure chart to determine appropriate cure times and conditions for all shades.

Indication	Increment Depth	All halogen lights (with output of 550-1000 mW/cm²)	LED lights (with output 1000-2000 mW/cm²)
Core Build-up and Class II Direct Restorations	5 mm	20 sec occlusal, 20 sec buccal, 20 sec lingual	10 sec occlusal, 10 sec buccal, 10 sec lingual
All indications listed (except Core Build-up and deeper Class II Direct Restorations)	4 mm	40 sec	20 sec
Anterior or shallow Class I Restorations	≤ 3mm	20 sec	10 sec

Note: For Class II restorations, remove the matrix band prior to the buccal and lingual curing steps.

11. **Contouring:** Contour restoration surfaces with fine finishing diamonds, burs or stones. Contour proximal surfaces with 3M™ ESPE™ Sof-Lex™ Finishing Strips.
12. **Adjust Occlusion:** Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.
13. **Finish and Polishing:** Polish with the 3M™ ESPE™ Sof-Lex™ Finishing and Polishing System.

Indirect Procedure for Inlays, Onlays or Veneers

- Dental Operatory Procedure**
1.1 Shade selection: Choose the appropriate shade(s) of Filtek One Bulk Fill Restorative prior to isolation.
1.2 Preparation: Prepare the tooth.
1.3 Impressioning: After preparation is complete, make an impression of the prepared tooth by following the manufacturer's instructions of the impressioning material chosen. An impressioning material, such as manufactured by 3M ESPE, may be used.
2. Laboratory Procedure
2.1 Pour the impression of the preparation with die stone. Place pins at the preparation site at this time if a "triple tray" type of impression was used.
2.2 Separate the cast from the impression after 45 to 60 minutes. Place pins in die and base, then cast as for a typical crown and bridge procedure. Mount or articulate the cast to its counter model on an adequate articulator.
2.3 If a second impression was not sent, pour a second cast using the same impression registration. This is to be used as a working cast.
2.4 Section out the preparation with a laboratory saw and trim away excess or, expose the margins so they can be easily worked. Mark the margins with a red pencil if needed. Add a spacer at this time if one is required.
2.5 Soak the die in water, then with a brush, apply a very thin coat of separating medium to the preparation, let it dry somewhat, and then add another thin layer.
2.6 Add the first increment of composite to the floor of the preparation, stay short of the margins, and follow the cure recommendations described in the Direct Restoration section (Step 10).
2.7 Place and cure additional increments of composite. Allow for the last increment (incisal) to include the contact areas.
2.8 Place the die back into the articulated arch. Add the last increment of composite to the occlusal surface. Overfill very slightly mesially, distally, and occlusally. This will allow for the mesiodistal contacts and the proper occlusal contact when the opposing arch is brought into occlusion with the uncured increment. Light cure for only ten seconds, then remove the die to prevent adhering to adjacent surfaces. Finish the curing process following the cure times in the Direct Restoration section (Step 10).
2.9 With the occlusal contacts already established, begin removing the excess composite from around the points of contact. Develop the inclines and ridges as per remaining occlusal anatomy.
2.10 Care must be taken when removing the prosthesis from the die. Break off small amounts of the die from around the restoration, the die stone should break away cleanly from the cured restoration, until all of the restoration is recovered.

- 2.11 Using the master die, check the restoration for flash, undercuts, and fit. Adjust as necessary, and then polish as noted above in Direct Restorative steps 11-13.
- Dental Operatory Procedure**
 - 3.1 Roughen the interior surfaces of the indirect restoration.
 - 3.2 Clean the prosthesis in a soap solution in an ultrasonic bath and rinse thoroughly.
 - 3.3 Cmentation: Cement the prosthesis using a 3M™ ESPE™ resin cement system.

- #### Storage and Use
- This product is designed to be used at room temperature. If stored in cooler conditions, allow product to reach room temperature prior to use. Shelf life at room temperature is 36 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.
 - Do not expose restorative materials to elevated temperatures, or to intense light.

3. Do not store materials in proximity to eugenol containing products.

Disinfection
The multiple-use syringe dispenser is not intended for direct patient contact. Use new, uncontaminated gloves when handling the syringe dispenser. Directions for cleaning and low level disinfection of the syringe dispenser are provided below:
Step 1 (Cleaning):

Use a CaviWipes™, or equivalent cleaning wipe, and wipe the entire surface of the device thoroughly for at least 30 seconds and until no visible soil remains on the device.

Step 2 (Disinfecting):

Use a new CaviWipes, or equivalent alcohol-quaternary ammonium disinfectant wipe, to disinfect the entire surface of the device by keeping wet for the contact time listed on the disinfectant label.

Disposal
See the Safety Data Sheet (available at www.3MESPE.com/dental or through your local subsidiary for disposal information).

Customer information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a dental professional.

Warranty
3M warrants this product will be free from defects in material and manufacture. 3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the products for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M sole obligation shall be repair or replacement of the 3M product.

Limitation of Liability
Except where prohibited by law, 3M will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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