



3M™ Liqui-Cel™ Membrane Contactors Can Help Meet Pure Steam Requirements in Pharmaceutical Industry

Pure Steam

Pure steam is used in sterilization chambers as a common method to sterilize pharmaceutical products, such as equipment parts, instruments, containers and materials for sterile environments. Although regional Pharmacopeias in Europe, Canada, Japan, India and China have minor differences in limits and measurement techniques, they each have similar specifications on pure steam sterilization and Water for Injection (WFI).

The US Pharmacopeia (USP) defines pure steam as steam that meets the criteria for WFI when condensed. The European and British Pharmacopeias have similar standards with additional specifications for dryness and the noncondensable gas content allowance. The UK National Health document HTM 2010 and the European Standard EN 285 set limits for noncondensable gases at 3.5% v/v and moisture content <5.0%.

The USP does not include a requirement for noncondensable gases but does state, “The level of steam saturation or dryness, and the amount of noncondensable gases are to be determined by the pure steam application.” This places the responsibility of producing the proper quality of steam for a particular application on the user. In many cases pure steam generators adopt the European and British pharmacopeia requirements for noncondensable gases and dryness to assure the steam meets or exceeds the required quality.

Noncondensable Gas

Whenever water and gas come into contact, gas has the tendency to dissolve into the water. The amount of gas that dissolves into the water is governed by Henry’s Law, which dictates that the amount of gas dissolved into a liquid is proportional to the partial pressure of the gas. The partial pressure value of the gas

is dependent on temperature and atmospheric pressure. Henry’s Law shows us that noncondensable gases, such as nitrogen, oxygen and carbon dioxide, will be present in water that comes into contact with the atmosphere.

Any noncondensable gases present in water that is used to produce steam will carryover with the steam when the water is boiled. These noncondensable gases mixed into the steam decreases the steam temperature at a given pressure. This reduces the sterilizer’s ability to properly transfer the required heat to sterilize its contents, so the removal of dissolved gases is a necessary step in the steam generation process.

Conductivity

The conductivity specification for WFI water is 1.1 µS/cm. Dissolved carbon dioxide can contribute to the conductivity of water. Carbon dioxide will freely pass through an RO membrane and dissociate in the permeate water to increase the conductivity.

Meeting the Specification

In order to remove the dissolved gas, many steam generators will preheat the feed water. Heating lowers the solubility of dissolved gases in the water, thereby lowering any noncondensable gas carried

Specification	Pure Steam	HTM 2010/EDN 285
pH	5-7	5-7
TOC	<500 ppb	<500 ppb
Total Bacteria	≤10 CFU/100 ml, pathogen free	≤10 CFU/100 ml, pathogen free
Endotoxin	≤0.5 EU/ml	≤0.5 EU/ml
Conductivity (20°C)	1.1 µS/cm	1.1 µS/cm
Noncondensable Gas	To be determined by application	3.5%
Moisture	To be determined by application	<5.0%

Technical Brief | 3M™ Liqui-Cel™ Membrane Contactors Can Help Meet Pure Steam Requirements in Pharmaceutical Industry

over into the steam. However, this process requires that the hot water be stored to allow time for the gas to be released. This process is energy intensive.

Carbon dioxide removal is often done using a double pass RO with pH adjustment. This requires a higher initial investment and operating cost.




Figure 1. Dissolved gas is removed from water using a 3M™ Liqui-Cel™ Membrane Contactor

3M™ Liqui-Cel™ Membrane Contactors offer a reliable, in-line solution that is far more economical for removing dissolved gas. The feedwater enters the membrane contactor on the outside (shellside) of the microporous hollow fiber. Since the hollow fiber is a hydrophobic membrane, water cannot pass through the small pores in the membrane wall. However, the dissolved gases will pass through the pores to the inside (lumenside) of the hollow fiber and exit the lumenside port.

A vacuum is applied to the lumenside of the hollow fiber. By lowering the gas pressure in contact with water a driving force is generated to remove the dissolved gas from the water. (See Figure 1.)

For more information and system sizing, please contact your 3M representative or visit 3M.com/Liqui-Cel.

Feature	Benefit
In-line	No need for storage tank to allow time for gas to be removed from hot water
Operates under vacuum only	No need to preheat the water - energy savings
Compact	Can fit on existing racks
FDA listed materials	Easily validated



Figure 2. Operating System by Pharmatec GmbH, a Bosch Packaging Technology Company

Technical Information: The technical information, recommendations and other statements contained in this document are based upon tests or experience that 3M believes are reliable, but the accuracy or completeness of such information is not guaranteed.

Product Use: Many factors beyond 3M's control and uniquely within user's knowledge and control can affect the use and performance of a 3M product in a particular application. Given the variety of factors that can affect the use and performance of a 3M product, user is solely responsible for evaluating the 3M product and determining whether it is fit for a particular purpose and suitable for user's method of application.

Warranty, Limited Remedy, and Disclaimer: Unless an additional warranty is specifically stated on the applicable 3M product packaging or product literature, 3M warrants that each 3M product meets the applicable 3M product specification at the time 3M ships the product. 3M MAKES NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR CONDITION OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY IMPLIED WARRANTY OR CONDITION ARISING OUT OF A COURSE OF DEALING, CUSTOM OR USAGE OF TRADE. If the 3M product does not conform to this warranty, then the sole and exclusive remedy is, at 3M's option, replacement of the 3M product or refund of the purchase price.

Limitation of Liability: Except where prohibited by law, 3M will not be liable for any loss or damage arising from the 3M product, whether direct, indirect, special, incidental or consequential, regardless of the legal theory asserted, including warranty, contract, negligence or strict liability.

3M and Liqui-Cel are trademarks of 3M Company. All other trademarks are the property of their respective owners. © 2017 3M Company. All rights reserved.



Separation and Purification Sciences Division
13840 South Lakes Drive
Charlotte, North Carolina
28273 USA
Phone: +1 980 859 5400

3M Deutschland GmbH
Separation and Purification Sciences Division
Öhder Straße 28
42289 Wuppertal Germany
Phone: +49 202 6099 - 0
Fax: +49 202 6099 - 241



ISO 9001

LC-1078
Rev. 01/2017

3M.com/Liqui-Cel