RelyX U200 CLICKER

Cemento de Resina Universal Autoadhesivo

Descripción

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3M Science.
Applied to Life.™

Zirconia Crown Cementation

Simple Steps for 3M™ RelyX™ U200 Self-Adhesive Resin Cement

Step 1 Prepare restoration. Sandblast the restoration with aluminum oxide after trays (Max 2 bar or 20 PSI, particle size $\leq 40\mu\text{m}$).

Step 2 Clean with alcohol and air dry with oil-free air.

Step 3 Remove provisional restoration. Mechanically clean prepared tooth (eg, with pumice paste).

Tip: Make sure any residue (temporary cement, desensitizers, airprimers, disinfectants, etc.) is completely removed. Do not use HCl, EDTA or NaOCl.

Step 4 Rinse and lightly oil. Leave tooth surface moist. Do not over dry the tooth. Over drying can lead to post-operative sensitivity.

Step 5 — For Clicker Delivery:
• Dispense a small amount of cement to ensure even dispensing
• Discard this material.
• Apply two clicks of the cement paste into a mixed.
• Hand mix with a spatula until a homogeneous mix is obtained.
• Wipe the Clicker barrels independently to avoid contamination of the material and premature setting.
• Ensure that the cap is put timely in place (clicking sound).

Tip: For a crown, typically two clicks would be needed.

Step 5 For Automatic Delivery: Dispense a small amount of cement onto the mixed to ensure a perfect mix. Discard this material.

Step 6 For Clicker Delivery: Use the mixing spatula to apply the cement into the crown.

Step 6 For Automatic Delivery: Dispense cement directly into the crown.

Step 7 Firmly seat the crown with finger pressure.

Step 8 Tack cure for 1-2 seconds per surface.

Step 9 Remove excess cement with a scaler while holding the crown in place.

Tip: Do not exceed recommended tool cure time; otherwise clean-up will be difficult. For a combined curing time, use 3M™ Epi-LED™ LED Curing Light or 3M™ Epi-LED™ DeepCure™ LED curing light back curing function.

Step 10 Light cure for 20 seconds per surface or wait 6 minutes from start of mix for dark cure. Finish and polish as needed.

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