3M

Attest™ REF 41482V

- en Super Rapid 5 Steam-Plus Challenge Pack
- fr Paquet-test Super Rapid 5 Steam-Plus
- de Super Rapid 5 Dampfverfahren-Plus-Kontrollpaket
- it Confezione di prova Super Rapid 5 Steam-Plus
- es Paquete de desafío Super Rapid 5 Steam-Plus
- es Paquete de desafío Super Rapid 5 Steam-Plus
- nl PCD met supersnelle aflezing voor stoomsterilisatie plus klasse
- sv Super Rapid 5 Steam-Plus Testpaket
- da Super Rapid 5 testpakke til dampsterilisering
- no Super Rapid 5 Steam-Plus utfordringstest-pakke
- fi Super Rapid 5 Steam-Plus-testipakkaus
- pt Pacote de Desafio a Vapor Superior 5 Super-rápido
- pt Super Rapid 5 Steam Plus Unidade de Teste
- el Πακέτο Ελέγχου Αποστείρωσης με Ατμό Plus 5, Σούπερ Ταχείας Ερμηνείας
- pl Pakiet testowy ultra szybkiego odczytu do monitorowania przebiegu sterylizacji parowej Super Rapid 5 Steam-Plus
- cs Testovací balíček s velice rychlým odečtem pro parní sterilizaci 5 Steam-Plus
- sl Testni paket za nadzorovanje parne sterilizacije z biološkimi indikatorji za super hitro odčitavanje
- ru Комплект сверхбыстрого считывания 5 Steam-Plus для испытания эффективности стерилизации
- hr Superbrzi 5 test paket uspješnosti sterilizacije parom-plus
- tr Süper Rapid 5 Steam-Plus Sorgulama Paketi
- zh 超快 5 蒸汽加强型验证包
- th ชุดตรวจสอบการทำให้เชื้ออ้อยไอน้ำแบบอ่านผลเร็วพิเศษ
- ko Super Rapid 5 Steam-Plus Challenge Pack

3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack 41482V

Product Description

The 3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack 41482V is specifically designed to qualify and monitor dynamic-air-removal steam sterilization processes at 270°F (132°C) and 275°F (135°C) in healthcare facilities. The challenge pack consists of multiple layers of medical index cards, some of which are die-cut to contain the monitoring products. The stacked cards are wrapped with a sterilization wrap. Each challenge pack has a process indicator on the pack label that changes from yellow to brown or darker when exposed to steam. This convenient disposable challenge pack presents a challenge to the sterilization process equivalent to the user-assembled biological indicator (BI) challenge test pack (towel PCD) recommended by the Association for the Advancement of Medical Instrumentation (AAMI). The challenge pack is a single use device.

Each challenge pack contains a 3M™ Attest™ Super Rapid Readout Biological Indicator 1492V (brown cap, hereinafter referred to as a 1492V BI), a 3M™ Attest™ Steam Chemical Integrator, and a record keeping sheet. AAMI recommends that steam sterilization loads containing an implant be monitored with a process challenge device (PCD) containing a biological indicator and an integrating indicator. 3M™ Attest™ Steam Chemical Integrators are Type 5 (Category i5) Integrating Indicators as categorized by ISO 11140-1:2014. 3M™ Attest™ Steam Chemical Integrators are single-use chemical indicators consisting of a paper wick and a steam and temperature sensitive chemical pellet contained in a paper/film/foil laminate. The chemical pellet melts and migrates as a dark color along the paper wick. The migration is visible through a green window marked ACCEPT or a red window marked REJECT; the extent of migration depends on steam, time, and temperature.

The 1492V BI is a self-contained dual readout biological indicator specifically designed for rapid and reliable monitoring of the steam sterilization process when used in conjunction with the 3M™ Attest™ Auto-reader 490 (hereinafter referred to as the 490 Auto-reader), a 3M™ Attest™ Auto-reader 490H having software version 4.0.0 or greater (hereinafter referred to as the 490H Auto-reader), or a 3M™ Attest™ Mini Auto-reader 490M (hereinafter referred to as the 490M Auto-reader). When steam processed, the process indicator on the top of the 1492V BI cap changes color from pink to light brown or darker. 3M™ Attest™ 1492V biological indicator controls are provided with the challenge packs.

The 1492V BI utilizes the α -glucosidase enzyme system, which is generated naturally within growing cells of *Geobacillus stearothermophilus*. The α -glucosidase in its active state is detected by measuring the fluorescence produced by the enzymatic hydrolysis of a non-fluorescent substrate, 4-methylumbelliferyl- α -D-glucoside (MUG). The resultant fluorescent by-product, 4-methylumbelliferone (MU), is detected in the Auto-reader. The presence of fluorescence within the specified incubation time of the 1492V BI in the 490, 490H, or 490M Auto-reader indicates a steam sterilization process failure.

The 1492V BI can also indicate the presence of *G. stearothermophilus* organisms by a visual pH color change reaction. Biochemical activity of the *G. stearothermophilus* organism produces metabolic by-products that cause the media to change color from purple to yellow which also indicates a steam sterilization process failure. Use of this indication method is optional and is typically restricted to special studies.

1492V BIs meet ISO 11138-1:2017 and ISO 11138-3:2017.

Readout Times

The super rapid readout and the optional 48-hour visual pH color change incubation times have been correlated with a 7-day incubation period following the FDA's Reduced Incubation Time protocol. Processed indicators were examined at 48 hours and 7 days for detection of a visual pH color change. The fluorescent result and the 48-hour visual pH color change readings were compared to the 7-day visual pH color change readings to determine the readout time of the indicator. The time to a fluorescent result is determined by the software version programmed on the Auto-reader.

24-minute Fluorescent Result

1492V BIs incubated in a 490 or 490H Auto-reader having software version 4.0.0 or greater or in a 490M Auto-reader have a 24-minute reduced incubation time result that correlates to the 7-day (168 hours) visual readout result \geq 97% of the time.

1-hour Fluorescent Result

1492V BIs incubated in a 490 Auto-reader having a software version less than 4.0.0 have a 1-hour reduced incubation time result that correlates to the 7-day (168 hours) visual readout result \geq 97% of the time.

48-hour Visual pH Color Change Result

1492V BIs incubated in a 490, 490H, or 490M Auto-reader have 48-hour reduced incubation time results that correlate to the 7-day (168 hours) visual readout result \geq 97% of the time.

Due to the high reliability of the fluorescent result there is no advantage to incubating 1492V BIs after the fluorescent result has been determined by the Auto-reader and recorded.

Note: To determine the software version on a 490 or 490H Auto-reader, depress the Information Button once. The Auto-reader will display two lines of text on the LCD display. The top row indicates (left to right):

- Serial number
- Software version (V X.Y.Z)
- Current date and time.

Indications for Use

United States

Use the 3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack 41482V in conjunction with the 3M™ Attest™ Auto-reader 490, the 3M™ Attest™ Auto-reader 490H having software version 4.0.0 or greater, or the 3M™ Attest™ Mini Auto-reader 490M to qualify or monitor dynamic-air-removal steam sterilization cycles of 4 minutes at 270°F (132°C) and 3 minutes at 275°F (135°C).

Outside the United States

Use the 3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack 41482V in conjunction with the 3M™ Attest™ Auto-reader 490, the 3M™ Attest™ Auto-reader 490H having software version 4.0.0 or greater, or the 3M™ Attest™ Mini Auto-reader 490M to qualify or monitor 270°F (132°C) to 275°F (135°C) dynamic-air-removal steam sterilization cycles.

Contraindications

None.

Warnings

There is a glass ampoule inside the plastic vial of the biological indicator. To avoid the risk of serious injury from flying debris due to a ruptured ampoule:

- Allow the biological indicator to cool for the recommended time period before activating. Activating or excessive handling of the BI before cooling may cause the glass ampoule to burst.
- Wear safety glasses when activating the biological indicator.
- Handle the biological indicator by the cap when crushing and flicking.
- Do not use your fingers to crush the glass ampoule.

Precautions

1. To ensure the challenge pack delivers the intended challenge:
 - DO NOT OPEN challenge pack prior to sterilization;
 - DO NOT reuse challenge pack.
2. DO NOT use the challenge pack to monitor sterilization cycles which it is not designed to challenge:
 - a. Gravity-displacement steam sterilization cycles;
 - b. 250°F (121°C) dynamic-air-removal steam sterilization cycles;
 - c. 270°F (132°C) dynamic-air-removal steam sterilization cycles having exposure times <4 minutes or 275°F (135°C) dynamic-air-removal steam sterilization cycles having exposure times <3 minutes;
 - d. Dry heat, chemical vapor, ethylene oxide or other low temperature sterilization processes.
3. To reduce the risk associated with incorrect results:
 - DO NOT place tape or labels on 1492V BI prior to incubation.
 - DO NOT incubate a 1492V BI if, after processing and before BI activation, it is observed to have a broken media ampoule. Retest the sterilizer with a new challenge pack.
 - After processing, allow 1492V BI to cool for 10 minutes before incubation.
 - After 1492V BI activation, ensure media has flowed to the spore growth chamber.
 - DO NOT re-incubate 1492V BIs for which the Auto-reader has already determined a result.

Monitoring Frequency


Follow facility Policies and Procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a best practice and to provide optimal patient safety, 3M recommends that every steam sterilization load be monitored with a biological indicator in an appropriate Process Challenge Device (i.e., BI challenge test pack).

Directions for Use

1. Place a 3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack 41482V flat, with the label side up, in the most challenging area for the sterilant to reach. This is generally on the bottom shelf, over the drain. Do not place objects (e.g. another pack) on top of the challenge pack. This will create too great a challenge for air removal and steam penetration.
2. Process the load according to established procedures.
3. After completion of the cycle, while wearing heat resistant gloves, retrieve the challenge pack.
4. Check to see that the external process indicator on the outside of the challenge pack has changed from yellow to brown or darker. Open the challenge pack and allow the 1492V BI to cool outside the challenge pack for 10 minutes prior to activation.
5. Check the 3M™ Attest™ Steam Chemical Integrator. The dark color should have entered the green ACCEPT window. If the dark color has not entered the green ACCEPT window, this indicates a red REJECT result which means the load was not exposed to sufficient steam sterilization conditions. This load should not be released for use but reprocessed. Record integrator result.
6. Check the process indicator on the top of the 1492V BI cap. A color change from pink to light brown or darker confirms that the biological indicator has been exposed to the steam process. This color change does not indicate that the steam process was sufficient to achieve sterility. If the process indicator is unchanged, check the sterilizer physical monitors.
7. Identify the processed 1492V BI by writing the sterilizer, load number, and processing date on the indicator label. Do not place another label or indicator tape on the biological indicator. While it is best practice to activate and incubate the 1492V BI as soon as it has cooled, the processed BI may be held at room temperature for up to 7 days prior to activation and incubation.
8. If utilized, fill out the required information on the record keeping card. Record the 1492V BI result when available.

9. Discard the challenge pack. Using the challenge pack more than once will invalidate subsequent test results.

10. To activate the biological indicator, while wearing safety glasses, place it in an Attest™ Biological Indicator Activator. Close and squeeze the activator to close the 1492V BI cap and crush the media ampoule. Then remove the BI and flick it (see pictures at right). Visually verify that media has flowed into the growth chamber at the bottom of the vial. If the media hasn't filled the growth chamber, hold the BI by the cap and flick it until media fills the growth chamber. Place the activated 1492V BI:

- in any well of a 490 or 490H Auto-reader having software version 4.0.0 or greater and a black incubation well configuration sticker


or

- in an incubation well of a 490 Auto-reader having software version less than 4.0.0 which is color-coded brown (i.e., configured to incubate 1492V BIs)

or

- in any well of a 490M Auto-reader

and wait for the result. See the Auto-reader Operator's Manual for further information related to its use.

11. Each day that a processed 1492V BI is incubated, activate and incubate at least one non-processed 1492V BI to use as a positive control. Follow the activation instructions provided in Step 10 above. Write a "C" (for "control") and the date on the BI label. The positive control should be from the same lot code as the processed biological indicator. The positive control BI helps confirm:

- correct incubation temperatures are met;
- viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals;
- capability of media to promote rapid growth; and
- proper functioning of the Auto-reader.

12. Incubation and Reading:

Incubate the positive control and steam processed 1492V BIs in a 490 Auto-reader, a 490H Auto-reader having software version 4.0.0 or greater, or a 490M Auto-reader. See the applicable Auto-reader Operator's Manual for the proper use of this equipment. The Auto-reader will indicate a positive result as soon as it is obtained. The final negative 1492V BI reading is made at:

- 24 minutes in 490 or 490H Auto-readers having software version 4.0.0 or greater and in 490M Auto-readers
- 1 hour in 490 Auto-readers having software version less than 4.0.0.

Interpretation of Results

Fluorescent Result

The positive control (unprocessed) 1492V BI must provide a positive fluorescent result (+ on the Auto-reader display). Processed 1492V BI results are not valid until the positive control reads fluorescent positive. If the positive control reads negative (- on the display), check the Auto-reader Operator's Manual Troubleshooting Guide. Retest the Auto-reader with a new positive control.

With processed 1492V BIs, a positive (+ on the display) result indicates a sterilization process failure. A final negative (- on the display) result for the processed 1492V BI after the specified incubation time indicates an acceptable sterilization process.

Optional Visual pH Color Change Result

The 1492V BI is normally discarded after the fluorescent result has been recorded. If, however, special studies are desired, 1492V BIs may be further incubated for a visual pH color change result. After activation and during incubation, the white nonwoven material will absorb the bromocresol purple indicator, the pH-sensitive indicator dye in the growth media, and appear blue. In the case of the positive control BI a yellow color change of the growth media and/or Nonwoven Material will appear within 48 hours. Any observation of a yellow color within the vial indicates a positive result.

In the case of a processed 1492V BI, a media and/or Nonwoven Material color change from purple to yellow indicates a sterilization process failure. A negative pH color change result, i.e., media and Nonwoven Material remain purple/blue, can be assessed at 48 hours.

- Record the processed and control BI lot codes and results. Act immediately on any positive results for processed biological indicators. Determine the cause of the positive biological indicator following facility policies and procedures. Always retest the sterilizer and do not use sterilizer for processing loads until qualification testing yields satisfactory results (typically three consecutive cycles with negative biological indicator results and three consecutive cycles with passing Bowie-Dick test results).

Storage


- Best stored under normal room conditions: 59-86°F (15-30°C), 35-60% relative humidity.
- Store away from direct sunlight. Do not store challenge packs near sterilants or other chemicals.
- After use, the 3M™ Attest™ Steam Chemical Integrator will not change visually within 6 months when stored at above conditions.

Disposal

Dispose of used 1492V BIs according to your health care facility policy. You may wish to steam sterilize any positive biological indicators at 270°F (132°C) for 4 minutes or at 275°F (135°C) for 3 minutes in a dynamic-air-removal steam sterilizer prior to disposal.

For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223, 5.1.1
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Steam indicator		Indicates product is designed for use with steam sterilization processes. Source: ISO 11140-1, 5.6

For more information see, [HCBG/gregulatory.3M.com](https://www.3M.com/HCBG/gregulatory)

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© 2020, 3M. All rights reserved.