3M

Attest™

REF 1492V

- en Super Rapid Readout Biological Indicator
- fr Indicateur biologique à lecture super rapide
- de Bioindikator zur Ultraschnellauslesung
- it Indicatore biologico a lettura super rapida
- es Indicador biológico de lectura superrápida
- es Indicador biológico de lectura superrápida
- nl Biologische indicator met snelle aflezing
- sv Super Rapid Readout Biological Indicator
- da Biologisk indikator med superhurtig udlæsning
- no Superhurtig avlesning biologisk indikator
- fi Super Rapid Readout biologinen pikaindikaattori
- pt Indicador Biológico de Leitura Super-rápida
- pt Indicador Biológico de Leitura Super Rápido
- el Βιολογικός Δείκτης Υπερταχείας Ερμηνείας
- pl Wskaźnik biologiczny o bardzo szybkim odczycie
- cs Biologický indikátor s rychlým odečtem
- sl Biološki indikator za super hitro odčitavanje
- ru Биологический индикатор сверхбыстрого считывания
- hr Super Rapid Readout Biološki indikator
- tr Süper Hızlı Sonuç Veren Biyolojik İndikatör
- zh 超快速判读式生物指示剂
- th ตั๋วบ่งชี้ทางชีวภาพแบบอ่านผลเร็วพิเศษ
- ko Super Rapid Readout Biological Indicator

3M™ Attest™ Super Rapid Readout Biological Indicator 1492V

Product Description

The 3M™ Attest™ Super Rapid Readout Biological Indicator 1492V (brown cap, referred to hereinafter as the 1492V BI) is a self-contained biological indicator specifically designed for rapid and reliable qualification testing and routine monitoring of 270°F (132°C) and 275°F (135°C) dynamic-air-removal (pre-vacuum and steam-flush pressure-pulse (SPPP) steam sterilization processes when used in conjunction with a 3M™ Attest™ Auto-reader 490 (hereinafter referred to as the 490 Auto-reader), a 3M™ Attest™ Auto-reader 490H having software version 4.0.0 or greater (hereinafter referred to as the 490H Auto-reader), or a 3M™ Attest™ Mini Auto-reader 490M (hereinafter referred to as the 490M Auto-reader). The 1492V BI is a single-use device.

A schematic illustrating the design of the 1492V BI is provided in Figure 1. The self-contained design includes a carrier with spores of *Geobacillus stearothermophilus* and a media ampoule containing bacteriological growth medium which meets the requirements for growth promoting ability specified in ISO 11138-1:2017. The spore carrier and media ampoule are contained in a plastic vial topped with a brown cap. A chemical process indicator which changes from pink to light brown or darker upon exposure to steam is located on the top of the cap.

The 1492V BI utilizes the α -glucosidase enzyme system, which is generated naturally within growing cells of *Geobacillus stearothermophilus*. The α -glucosidase in its active state is detected by measuring the fluorescence produced by the enzymatic hydrolysis of a non-fluorescent substrate, 4-methylumbelliferyl- α -D-glucoside (MUG). The resultant fluorescent by-product, 4-methylumbelliferone (MU), is detected in the Auto-reader. The presence of fluorescence within the specified incubation time of the 1492V BI in the 490, 490H, or 490M Auto-reader indicates a steam sterilization process failure.


Figure 1: 3M™ Attest™ Super Rapid Readout Biological Indicator 1492V Design

The 1492V BI can also indicate the presence of *G. stearothermophilus* organisms by a visual pH color change reaction. Biochemical activity of the *G. stearothermophilus* organism produces metabolic by-products that cause the media to change color from purple to yellow which also indicates a steam sterilization process failure. Use of this indication method is optional and is typically restricted to special studies.

1492V BIs meet ISO 11138-1:2017 and ISO 11138-3:2017.

Readout Times

The super rapid readout and the optional 48-hour visual pH color change incubation times have been correlated with a 7-day incubation period following the FDA's Reduced Incubation Time protocol. Processed indicators were examined at 48 hours and 7 days for detection of a visual pH color change. The fluorescent result and the 48-hour visual pH color change readings were compared to the 7-day visual pH color change readings to determine the readout time of the indicator. The time to a fluorescent result is determined by the software version programmed on the Auto-reader.

24-minute Fluorescent Result

1492V BIs incubated in a 490 or 490H Auto-reader having software version 4.0.0 or greater or in a 490M Auto-reader have a 24 minute reduced incubation time result that correlates to the 7 day (168 hours) visual readout result $\geq 97\%$ of the time.

1-hour Fluorescent Result

1492V BIs incubated in a 490 Auto-reader having a software version less than 4.0.0 have a 1 hour reduced incubation time result that correlates to the 7 day (168 hours) visual readout result $\geq 97\%$ of the time.

48-hour Visual pH Color Change Result

1492V BIs incubated in a 490, 490H, or 490M Auto-reader have 48 hour reduced incubation time results that correlate to the 7 day (168 hours) visual readout result $\geq 97\%$ of the time.

Due to the high reliability of the fluorescent result there is no advantage to incubating 1492V BIs after the fluorescent result has been determined by the 490, 490H, or 490M Auto-reader and recorded.

Note: To determine the software version on a 490 or 490H Auto-reader, depress the Information Button once. The Auto-reader will display two lines of text on the LCD display. The top row indicates (left to right):

- Serial number
- Software version (V X.Y.Z)
- Current date and time.

Indications for Use

United States

Use the 3M™ Attest™ Super Rapid Readout Biological Indicator 1492V in conjunction with the 3M™ Attest™ Auto-reader 490, the Attest™ Auto-reader 490H having software version 4.0.0 or greater, or the 3M™ Attest™ Mini Auto-reader 490M to qualify or monitor dynamic-air-removal steam sterilization cycles of:

- 3 minutes at 270°F (132°C)
- 4 minutes at 270°F (132°C)
- 3 minutes at 275°F (135°C).

Outside the United States

Use the 3M™ Attest™ Super Rapid Readout Biological Indicator 1492V in conjunction with the 3M™ Attest™ Auto-reader 490, the Attest™ Auto-reader 490H having software version 4.0.0 or greater, or the 3M™ Attest™ Mini Auto-reader 490M to qualify or monitor 270°F (132°C) to 275°F (135°C) dynamic-air-removal steam sterilization cycles.

Contraindications

None.

Warnings

There is a glass ampoule inside the plastic vial of the biological indicator (BI). To avoid the risk of serious injury from flying debris due to a ruptured BI:

- Allow the BI to cool for the recommended time period before activating. Activating or excessive handling of the BI before cooling may cause the glass ampoule to burst.
- Wear safety glasses and gloves when removing the BI from the sterilizer.
- Wear safety glasses when activating the BI.
- Handle the BI by the cap when crushing or flicking.
- Do not use your fingers to crush the glass ampoule.

Precautions

1. DO NOT use the 1492V BI to monitor sterilization cycles which it is not designed to challenge:
 - a. Gravity-displacement steam sterilization cycles;
 - b. 250°F (121°C) dynamic-air-removal steam sterilization cycles;
 - c. Dry heat, chemical vapor, ethylene oxide or other low temperature sterilization processes.

2. To reduce the risk associated with incorrect results:

- Do not place tape or labels on 1492V BI prior to sterilization or incubation.
- Do not incubate a 1492V BI if, after processing and before BI activation, it is observed to have a broken media ampoule. Retest the sterilizer with a new biological indicator.
- After processing, allow BI to cool for 10 minutes before incubation.
- After BI activation, ensure media has flowed to the spore growth chamber.


Monitoring Frequency

Follow facility Policies and Procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a best practice and to provide optimal patient safety, 3M recommends that every steam sterilization load be monitored with a biological indicator in an appropriate Process Challenge Device (PCD i.e., BI challenge test pack).

Directions for Use

1. Identify the 1492V BI by writing the load number, sterilizer, and processing date on the indicator label. Do not place another label or indicator tape on the vial or on the cap.
2. Place the 1492V BI in a representative tray configuration or Process Challenge Device (PCD) as recommended by professional association guidelines or national standards for healthcare facility practice.
3. Place the PCD in the most challenging area of the sterilizer. This is typically on the bottom shelf, over the drain, however, the sterilizer manufacturer should be consulted to identify the area of the chamber least favorable to sterilization.
4. Process the load according to recommended practices.
5. After completion of the cycle take the PCD out of the sterilizer, and remove the 1492V BI.
6. Allow the 1492V BI to cool for 10 minutes prior to activation.
7. Check the process indicator on the top of the cap of the 1492V BI. A color change from pink to light brown or darker confirms that the 1492V BI has been exposed to the steam process. This color change does not indicate that the steam process was sufficient to achieve sterility. If the process indicator is unchanged, check the sterilizer physical monitors.

8. To activate the biological indicator, place it in an Attest™ Biological Indicator Activator. Close and squeeze the activator to close the 1492V BI cap and crush the media ampoule. Then remove the BI and flick it (see pictures at right). Visually verify that media has flowed into the growth chamber at the bottom of the vial. If the media hasn't filled the growth chamber, hold the BI by the cap and flick it until media fills the growth chamber. Place the activated 1492V BI:

- in any well of a 490 or 490H Auto-reader having software version 4.0.0 or greater


or

- in an incubation well of a 490 Auto-reader having software version less than 4.0.0 which is color-coded brown (i.e., configured to incubate 1492V BIs)

or


- in any well of a 490M Auto-reader

and wait for the result. See the Auto-reader Operator's Manual for further information related to its use.

9. Each day that a processed 1492V BI is incubated, activate and incubate at least one non-processed 1492V BI to use as a positive control. Follow the activation instructions provided in Step 8 above. Write a "C" (for "control") and the date on the BI label. The positive control should be from the same lot code as the processed biological indicator. The positive control BI helps confirm:

- correct incubation temperatures are met,
- viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals,
- capability of media to promote rapid growth, and

- proper functioning of the Auto-reader.

10. Incubation and Reading:

Incubate the positive control and steam processed 1492V BIs in a 490 Auto-reader, a 490H Auto-reader having software version 4.0.0 or greater, or a 490M Auto-reader. See the applicable Auto-reader Operator's Manual for the proper use of this equipment. The Auto-reader will indicate a positive result as soon as it is obtained. The final negative 1492V BI reading is made at:

- 24 minutes in 490 or 490H Auto-readers having software version 4.0.0 or greater and in 490M Auto-readers
- 1 hour in 490 Auto-readers having software version less than 4.0.0.

After the results are displayed and recorded, the 1492V BIs may be discarded.

Interpretation of Results

Fluorescent Results

The positive control (unprocessed) 1492V BI must provide a positive fluorescent result (+ on the Auto-reader display). Processed 1492V BI results are not valid until the positive control reads fluorescent positive. If the positive control reads negative (- on the display), check the Auto-reader Operator's Manual Troubleshooting Guide. Retest the Auto-reader with a new positive control.

With processed 1492V BIs, a positive (+ on the display) result indicates a sterilization process failure. A final negative processed 1492V BI reading (- on the display) after the specified incubation time indicates an acceptable sterilization process.

Act immediately on any positive results for processed BIs. Determine the cause of the positive BI following facility policies and procedures. Always retest the sterilizer and do not use sterilizer for processing loads until three consecutive BI results are negative.

Optional Visual pH Color Change Result

The 1492V BI is normally discarded after the fluorescent result has been recorded. If, however, special studies are desired, 1492V BIs may be further incubated for a visual pH color change result. After activation and during incubation, the white Nonwoven Material will absorb the bromocresol purple indicator, the pH-sensitive indicator dye in the growth media, and appear blue. In the case of the positive control BI a yellow color change of the growth media and/or Nonwoven Material will appear within 48 hours. Any observation of a yellow color within the vial indicates a positive result.

In the case of a processed 1492V BI, a media and/or Nonwoven Material color change from purple to yellow indicates a sterilization process failure. A negative pH color change result, i.e., media and Nonwoven Material remain purple/blue, can be assessed at 48 hours.

Storage


- Best stored in the original box under normal room conditions: 59-86°F (15-30°C), 35-60% relative humidity (RH).
- Do not store 1492V BIs near sterilants or other chemicals.

Disposal

Dispose of used 1492V BIs according to your health care facility policy. You may wish to steam sterilize any positive biological indicators at 132°C (270°F) for 4 minutes or at 275°F (135°C) for 3 minutes in a dynamic-air-removal steam sterilizer prior to disposal.

For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223, 5.1.1
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Steam indicator		Indicates product is designed for use with steam sterilization processes. Source: ISO 11140-1, 5.6
British Standards Institution (BSI) Kitemark™		British Standards Institution (BSI) Kitemark™ License for compliance with BS EN ISO 11138-1:2017 and BS EN ISO 11138-3:2017

For more information see, [HCBGRegulatory.3M.com](https://www.hcbg.com/regulatory)

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3M and Attest are trademarks of 3M.
© 2020, 3M. All rights reserved.