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#### Referenzen:

- 1 Brauer, A., et al. (2004), Conductive Heat Exchange with a Gel-Coated Circulating Water Mattress. *Anaesthesia and Analgesia*, Vol. 6, No. 99, pp. 1742-1746
- 2 Sessler, D.I. (1997), Current Concepts: mild perioperative hypothermia, *New England Journal of Medicine*, No. 336, pp. 1730-1737.
- 3 Karlnoski, R. et. al. Intraoperative Warming with vitalHEAT® during Open Abdominal Surgery. *ASA abstracts*, 2010. A086. vitalHEAT® is a trademark of Dynatherm Medical, Inc..
- 4 Tominaga A, Koitabashi T, Ouchi T, Ban R, Takano E. Efficacy of an underbody forced-air warming blanket for the prevention of intraoperative hypothermia. *Anesthesiology*. Vol 1072007:A91.
- 5 Ouchi, T. et. al. Lithotomy Underbody Air Blanket Can Prevent Intraoperative Redistribution Hypothermia. *ASA abstracts*, 2010. A088.
- 6 Engelen, S, et al. (2010) A Comparison of Under-Body Forced Air and Resistive Heating during Hypothermic, on-pump cardiac surgery. *Anaesthesia*, Vol.66, No.2, pp.104-110.
- 7 Roder, G. et al. (2011), Intra-operative rewarming with Hot Dog® resistive heating and forced-air heating: a trial of lower-body warming, *Anaesthesia*, Vol 66, No. 8, pp. 667-674.
- 8 Sessler DI, Moayeri A. (1990), Skin-surface warming: heat flux and central temperature. *Anesthesiology*, No.73, pp. 218-24.
- 9 Giesbrecht GG, Ducharme MB, McGuire JP. (1994), Comparison of forced-air patient warming systems for perioperative use. *Anesthesiology*, No.80, pp.671-9.
- 10 Hynson JM, Sessler DI. (1992), Intraoperative warming therapies: a comparison of three devices. *Journal of Clinical Anesthesia*, No. 4, pp. 194-9.
11. Kurz A, Kurz M, Poeschl G, Faryniak B, Redl G, Hackl W. (1993), Forced-air warming maintains intraoperative normothermia better than circulating-water mattresses. *Anesthesia Analgesia*, No. 77, pp. 89-95.
12. Borms SF, Englelen SL, Himpe DG, Suy MR, Theunissen WJ. (1994), Bair Hugger forced-air warming maintains normothermia more effectively than thermo-lite insulation. *Journal of Clinical Anesthesia*, No. 6, pp.303-7.
13. Data on file at 3M US, based on internal sales data from Arizant USA, a 3M company 2012.
14. Zink RS, laizzo PA. (1993), Conductive warming therapy does not increase the risk of wound contamination in the operating room. *Anesthesia Analgesia*, No. 76, pp.50-3.
15. Huang JK, Shah EF, Vinodkumar N, Hegarty MA, Greateax RA, (2003), The Bair Hugger patient warming system in prolonged vascular surgery: an infection risk? *Critical Care* No. 7, R13-R16
16. Kurz A., Sessler D.I., Lenhardt R.L. (1996) Perioperative normothermia to reduce the incidence of surgical-wound infection and shorten hospitalization. *New England Journal of Medicine* No. 334, Vol. 19, pp. 1209-1215
17. Sharp RJ, et al. (2002), Do warming blankets increase bacterial counts in the operating field in a laminar-flow theatre? *Journal of Bone Joint Surgery Br*, No. 84, B:486-8.
18. Miyazaki H, Sato M, Okazaki K. (2007), Forced-air warmer did not increase the risk of contamination caused by interference of clean airflow. *Anesthesiology*, No. 107, A1594
19. Sessler DI, Olmsted RN, Kuelpmann R. Forced-air warming does not worsen air quality in laminar flow operating rooms. *Anesthesia & Analgesia*, No. 113, pp. 1416-21
20. Barie PS (2002), *Surgical Site Infections: Epidemiology and Prevention*. *Surgical Infections*. Vol. 3, S-9 – S-21.
21. Tryba M, Leban J, et al. (1996), Does active warming of severely injured trauma patients influence perioperative morbidity? *Anesthesiology*, No. 85, Vol. A283
22. Jeran L. (2001), American Society of PeriAnesthesia Nurses Development Panel. *Clinical Guideline for the Prevention of Unplanned Perioperative Hypothermia*. *Journal of PeriAnesthesia Nursing*. Vol. 16(5), pp. 305-314.
23. Horn, E.F, Bein, B et al. (2012) The effects of short time periods of pre-operative warming in the prevention of peri-operative hypothermia. *Anaesthesia*, Vol. 67, pp.612-617
24. Plowman R, Craves N, Griffin M et al. (2000) *The socio-economic burden of hospital acquired infection*. London: Public Health Service, 1999.
25. National Institute of Clinical Excellence, *Clinical Guideline 65, Inadvertent perioperative hypothermia, The management of inadvertent perioperative hypothermia in adults*, April 2008.
26. Mahoney CB, Odom J. (1999), Maintaining intraoperative normothermia: A meta-analysis of outcomes with costs. *AANA Journal*. No. 67, pp. 155-163.
27. A. Torossian, the TEMMP Study Group, Survey on intraoperative temperature management in Europe, *European Journal of Anesthesiology* 2007; 24:668-675
28. Arizant Healthcare Inc. *Bair Paws Patient Satisfaction Study*. 2006-2007.
29. Siegel JD, Rhinehart E, Jackson M, Chiarello L. Management of multi-drug resistant organisms in healthcare settings, 2006. 2006;1-74.
30. Sharp RJ, et al. Do warming blankets increase bacterial counts in the operating field in a laminar-flow theatre? *J Bone Joint Surg Br* 2002; 84-B:486-8.
31. Miyazaki H, Sato M, Okazaki K. Forced-air warmer did not increase the risk of contamination caused by interference of clean airflow. *Anesthesiology* 2007; 107:A1594
32. Tumia N., Ashcroft GP. Convection warmers – a possible source of contamination in laminar airflow operating theatres? *J Hosp Infect* 2002; 52: 171-4.
33. Memarzadeh F, Zheng Jiang PE. Effect of operation room geometry and ventilation system parameter variations on the protection of the Surgical Site. *IAQ*, 2004: 1-6.

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Bitte recyceln. Gedruckt in Deutschland. © 3M 2015. Alle Rechte vorbehalten.
**3M Deutschland GmbH**