# Filtek™

Supreme XTE
Universal Restorative

Technique Guide

MENU >



## **Technique Guide**

We want you to get the best results using 3M™ Filtek™ Supreme XTE Universal Restorative in your restorative procedures. That's why we've created this interactive Technique Guide. In it, you'll find a list of indications and simple technique guides for the most commonly performed procedures. Each guide has step-by-step illustrations and easy instructions that are important to know in advance of the procedure.

# Menu

Click or tap to jump to the desired technique.

Inc	lications 3
Tech	nique Guides
	Anterior Multishade4
π	Anterior Simple5
1	Posterior Dual Shade6
	Posterior Single7




# **Indications**

- Direct anterior and posterior restorations (including occlusal surfaces)
- Core build-ups
- Splinting
- Indirect restorations including inlays, onlays and veneers





# **Anterior Multishade**\*




## **Preparation:**

- Select tooth shade using the VITAPAN® classical shade guide
- Select three composite shades using the 3M™ Filtek™ Supreme XTE **Universal Restorative** Shade Selector as a guide
- Prepare the tooth
- Bevel enamel margins



## **Bonding procedure:**

 Follow the bonding procedure provided in the Instructions for Use of the adhesive system of choice



Place composite—dentine replacement:

- Place and cure Dentine shade of Filtek Supreme XTE restorative to provide opacity and colour intensity
- Light cure



**Shades** 

Body, Enamel,

Translucent Dentine, A6B

and B5B

Increment

**Depth** 

2 mm

1.5 mm



## Place composite—dentine replacement:

- Place and cure Body shade of Filtek Supreme XTE restorative to the dentine enamel junction
- Light cure



All Halogen and

**LED Lights** 

(with output of 400–1000 mW/cm²)


20 sec.

40 sec.



## Place composite—enamel replacement:

- Fill remaining preparation using Enamel and/or Translucent shades of Filtek Supreme XTE restorative to provide translucency for depth of colour and incisal edge match
- Light cure





3M™ Elipar™ DeepCure-S

**LED Curing Light** 

(with output 1000-2000 mW/cm²)

10 sec.

20 sec.



- Use 3M™ Sof-Lex™ Finishing and Polishing Discs and Strips
- The aesthetic result will be enhanced by creating surface texture similar to the adjacent teeth





# **Anterior Simple**\*





#### **Preparation:**

- Select tooth shade using the VITAPAN® classical shade guide
- Select two composite shades using the 3M™ Filtek™ Supreme XTE **Universal Restorative Shade** Selector as a guide
- Prepare the tooth
- Bevel enamel margins





## **Bonding procedure:**


 Follow the bonding procedure provided in the Instructions for Use of the adhesive system of choice





## Place composite—dentine replacement:

- Place and cure Body (or Dentine) shade of Filtek Supreme XTE restorative in layers to the dentine enamel junction
- Light cure








## Place composite—enamel replacement:

- Fill remaining preparation using the Body (or Enamel) shade of Filtek Supreme XTE restorative
- Light cure





## Finish and polish:


- Use 3M<sup>™</sup> Sof-Lex<sup>™</sup> Finishing and Polishing Discs and Strips
- The aesthetic result will be enhanced by creating surface texture similar to the adjacent teeth





# Posterior Dual Shade<sup>\*</sup>

	Shades	Increment Depth	All Halogen and LED Lights (with output of 400–1000 mW/cm²)	3M <sup>™</sup> Elipar <sup>™</sup> DeepCure-S LED Curing Light (with output 1000-2000 mW/cm²)
	Body, Enamel, Translucent	2 mm	20 sec.	10 sec.
	Dentine, A6B and B5B	1.5 mm	40 sec.	20 sec.





#### Preparation:

- Select shade using the VITAPAN® classical shade guide and the 3M™ Filtek™ Supreme XTE Universal Restorative Shade Selector Wheel
- Prepare the tooth and isolate

## Apply liner/base if desired:

- Mix 3M<sup>™</sup> Vitrebond<sup>™</sup> or 3M<sup>™</sup> Vitrebond<sup>™</sup> Plus
- Apply a thin layer of the liner/base to dentine surfaces
- Light cure





#### **Bonding procedure:**

 Follow the bonding procedure provided in the Instructions for Use of the adhesive system of choice





## Place composite—proximal box:

- To ensure good adaptation, place chosen Body (or Dentine) shade of Filtek Supreme XTE restorative in a 1 mm increment in the proximal box and on the floor of the cavity
- Light cure





## Place composite dentine replacement:

- Continue placing chosen Body (or Dentine) shade of Filtek Supreme XTE restorative in increments to the DEJ
- Use remaining coronal tooth structure as a guide for placement and adaptation
- Light cure





## Place composite enamel replacement:

- Place chosen Enamel (or Translucent) shade of Filtek Supreme XTE restorative in increments
- Use remaining coronal tooth structure as a guide for placement and adaptation
- Overfill slightly
- Light cure





## Finish and polish:

- Finish occlusal surface using an appropriate finishing instrument
- Finish interproximal surfaces with 3M<sup>™</sup>
   Sof-Lex<sup>™</sup> Extra Thin Discs and 3M<sup>™</sup>
   Sof-Lex<sup>™</sup> Strips






#### Occlusion:

- Check lateral and centric occlusion
- Adjust as necessary





# Posterior Single<sup>\*</sup>



#### **Preparation:**

- Select shade using the VITAPAN® classical shade guide and the 3M™ Filtek™ Supreme XTE Universal Restorative Shade Selector Wheel
- Prepare the tooth and isolate


## Apply liner/base if desired:

- Mix 3M<sup>™</sup> Vitrebond<sup>™</sup> or 3M<sup>™</sup> Vitrebond<sup>™</sup> Plus
- Apply a thin layer of the liner/ base to dentine surfaces
- Light cure



## Bonding procedure:

 Follow the bonding procedure provided in the Instructions for Use of the adhesive system of choice




## Place composite—

- To ensure good adaptation, place chosen Body shade of Filtek Supreme XTE restorative in a 1 mm increment in the proximal box and on the floor of the cavity
- Light cure

proximal box:




Increment



## Place composite:

- Continue placing chosen Body shade of Filtek Supreme XTE restorative in increments
- Use remaining coronal tooth structure as a guide for placement and adaptation
- Overfill slightly
- Light cure



All Halogen and

## Finish and polish:

- Finish occlusal surface using an appropriate finishing instrument
- Finish interproximal surfaces with 3M<sup>™</sup> Sof-Lex<sup>™</sup> Extra Thin Discs and 3M<sup>™</sup> Sof-Lex<sup>™</sup> Strips



3M™ Elipar™ DeepCure-S



## Occlusion:

- Check lateral and centric occlusion
- Adjust as necessary

## 3M.com.au/Dental 3M.co.nz/Dental



3M Oral Care

**3M Australia Pty Limited** Building A, 1 Rivett Road North Ryde, NSW 2113 Ph: 1300 363 454 3M.com.au/Dental

3M New Zealand Limited 94 Apollo Drive Rosedale, Auckland 0632 Ph: 0800 80 81 82 3M.co.nz/Dental



f Find us on Facebook

3M, "3M Science. Applied to Life.", Elipar, Filtek, Sof-Lex and Vitrebond are all trademarks of 3M. VITAPAN is a trademark of VITA North America and not of 3M. Published by 3M Oral Care. © 3M 2019. Please always refer to the 3M Product Instructions For Use. All rights reserved. AV011477203