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Referencias

1. Notomi T, et al. Loop-mediated isothermal amplification of DNA. *Nucleic Acids Res.* 2000, 28:E63.
2. Fu S, et al. Applications of loop-mediated isothermal DNA amplification. *Appl. Biochem. Biotechnol.* 2011, 163:845-50
3. Parida M, et al. Loop mediated isothermal amplification (LAMP): a new generation of innovative gene amplification technique; perspectives in clinical diagnosis of infectious diseases. *Rev. Med. Virol.* 2008, 18:407-21.
4. Francois P, et al. Robustness of a loop-mediated isothermal amplification reaction for diagnostic applications. *FEMS Immunol. Med. Microbiol.* 2011, 62:41-48.
5. Kaneko H, et al. Tolerance of loop-mediated isothermal amplification to a culture medium and biological substances. *J. Biochem. Biophys. Methods.* 2007, 70:499-501
6. Kiddle G, et al. GMO detection using a bioluminescent real time reporter (BART) of loop mediated isothermal amplification (LAMP) suitable for field use. *BMC Biotechnology.* 2012, 12:15
7. Plutzer J, and Karanis P. Rapid identification of *Giardia duodenalis* by loop-mediated isothermal amplification (LAMP) from faecal and environmental samples and comparative findings by PCR and real-time PCR methods. *Parasitol. Res.* 2009, 104:1527-1533.
8. Nkouawa A, et al. Evaluation of a loop-mediated isothermal amplification method using fecal specimens for differential detection of *Taenia* species from humans. *J. Clin. Microbiol.* 2010, 48:3350-3352.
9. Opel K, et al. A study of PCR inhibition mechanisms using Real Time PCR. *J. Forensic Sci.* 2010, 55:25-33.
10. Gandelman OA. Novel bioluminescent quantitative detection of nucleic acid amplification in real-time. *PLoS ONE.* 2010, 5: e14155.
11. Rossmannith P, et al. The fluorescence characteristics of enrichment media in the wavelength range of Real-Time PCR thermocycler optical path assignments. *Food Anal. Methods.* 2010, 3:219-224.
12. Rossen L, Nørskov P, Holmstrøm K, Rasmussen OF. 1992. Inhibition of PCR by components of food samples, microbial diagnostic assays and DNA-extraction solutions. *Int J Food Microbiol.* 17:37-45.
13. Hoorfar J, et al. Practical considerations in design of internal amplification controls for diagnostic PCR assays. *J Clin Microbiol.* 2004, 42:1863-1868.
14. ISO 22174:2005. Microbiología de los alimentos para consumo humano y animal. Reacción en cadena de la polimerasa (PCR) para la detección de patógenos en los alimentos. Requisitos generales y definiciones.
15. ISO 16140-2:2016. Microbiología de la cadena alimentaria. Validación de métodos. Parte 2: Protocolo para la validación de métodos alternativos (registrados) frente a los métodos de referencia.

16. Hu L, et al. Evaluation of 3M Molecular Detection System and ANSR Pathogen Detection System for rapid detection of *Salmonella* from egg products. *Poultry Sci.* 2017, 96:1410-1418.
17. Lins P. Detection of *Salmonella* spp. in spices and herbs, *Food Control*, 2017, 83: 61-68.
18. Sarowska J, et al. Detection of *Salmonella* in Foods Using a Reference PN-ISO Method and an Alternative Method Based on Loop-mediated Isothermal Amplification Coupled with Bioluminescence. *Adv Clin Exp Med.* 2016, 25:945-950.
19. Yang Q, et al. Rapid detection of *Salmonella* in food and feed by coupling loop-mediated isothermal amplification with bioluminescent assay in real-time. *BMC Microbiol.* 2016, 16:112.
20. Miks-Krajnik M, et al Loop-mediated isothermal amplification (LAMP) coupled with bioluminescence for the detection of *Listeria monocytogenes* at low levels on food contact surfaces. *Food Control*, 2016, 60:237-240.
21. Abirami N, et al. Evaluation of commercial loop-mediated isothermal amplification based kit and ready-to-use plating system for detection of *Salmonella* in naturally contaminated poultry and their processing environment. *Food Control*, 2016, 70:74-78.
22. Vongkamjan K, et al Various Ready-to-Eat Products from Retail Stores Linked to Occurrence of Diverse *Listeria monocytogenes* and *Listeria* spp. Isolates. *J Food Prot.* 2016, 79:239-45.
23. Lim HS, et al. Evaluation of commercial kit based on loop-mediated isothermal amplification for rapid detection of low levels of uninjured and injured *Salmonella* on duck meat, bean sprouts, and fishballs in Singapore. *J Food Prot.* 2015, 78:1203-1207.
24. Ryan G, et al. Evaluation of Rapid Molecular Detection Assays for *Salmonella* in Challenging Food Matrices at Low Inoculation Levels and Using Difficult-to-Detect Strains. *J Food Prot.* 2015, 78:1632-41.
25. Loff M, et al 3M™ Molecular detection system versus MALDI-TOF mass spectrometry and molecular techniques for the identification of *Escherichia coli* O157:H7, *Salmonella* spp. & *Listeria* spp. *J Microbiol Methods.* 2014, 101:33-43.
26. Fortes ED, et al. Validation of the 3M molecular detection system for the detection of *Listeria* in meat, seafood, dairy, and retail environments. *J Food Prot.* 2013, 76:874-878
27. Bonardi S, et al. Comparison of an isothermal amplification and bioluminescence detection of DNA method and ISO 6579:2002 for the detection of *Salmonella enterica* serovars in retail meat samples. *J Food Prot.* 2013, 76:657-61.
28. Bird P, et al. Evaluation of 3M Molecular Detection Assay (MDA) 2- *E. coli* O157 (including H7) for the Detection of *E. coli* O157 in Select Foods: Collaborative Study, First Action 2017.01. *J AOAC Int.* 2017, In Press.
29. Bird P, et al. Evaluation of 3M Molecular Detection Assay (MDA) 2- *Listeria* for the Detection of *Listeria* Species in Select Foods and Environmental Surfaces: Collaborative Study, First Action 2016.07. *J AOAC Int.* 2017, 100:82-98.
30. Bird P, et al. Evaluation of the 3M™ Molecular Detection Assay (MDA) 2 - *Listeria monocytogenes* for the Detection of *Listeria monocytogenes* in a Variety of Foods and Select Environmental Surfaces: Collaborative Study, First Action 2016.08. *J AOAC Int.* 2017, 100:454-469.
31. Bird P, et al. Evaluation of the 3M™ Molecular Detection Assay (MDA) 2 - *Salmonella* for the Detection of *Salmonella* spp. in Select Foods and Environmental Surfaces: Collaborative Study, First Action 2016.01. *J AOAC Int.* 2016, 99:980-997.
32. Bird P, et al. Evaluation of 3M™ Molecular Detection Assay (MDA) *Listeria* for the Detection of *Listeria* species in Selected Foods and Environmental Surfaces: Collaborative Study, First Action 2014.06. *J AOAC Int.* 2015, 98:993-1002.
33. Bird P, et al. Evaluation of 3M™ Molecular Detection Assay (MDA) *Listeria monocytogenes* for the Detection of *Listeria monocytogenes* in Selected Foods and Environmental Surfaces: Collaborative Study, First Action 2014.07. *J AOAC Int.* 2015, 98:980-992
34. Bird P, et al. Evaluation of Modification of the 3M™ Molecular Detection Assay (MDA) *Salmonella* Method (2013.09) for the Detection of *Salmonella* in Selected Foods: Collaborative Study. *J AOAC Int.* 2014, 97:1329-1342.
35. Bird P, et al. Evaluation of 3M™ Molecular Detection Assay (MDA) *Salmonella* for the detection of *Salmonella* in selected foods: collaborative study. *J AOAC Int.* 2013, 96:1325-1335.

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© 3M 2018. Todos los derechos reservados.