Own the Zone:

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Författare

Al Van Duren, Director of Scientific Affairs, 3M Patient Warming Business

Referenser

- Schroek H, Lyden AK, Benedict WL, Ramachandran SK. Time Trends and Predictors of Abnormal Postoperative Body Temperature in Infants Transported to the Intensive Care Unit. *Anesthesiology Research and Practice*. 2016;7318137.
- Sessler DI. Current concepts: Mild Perioperative Hypothermia. *New Engl J Med*. 1997; 336(24):1730-1737.
- Scott AV, Stonemetz JL, Wasey JO, Johnson DJ, Rivers RJ, Koch CG, et al. (2015). Compliance with Surgical Care Improvement Project for Body Temperature Management (SCIP Inf-10) Is Associated with Improved Clinical Outcomes. *Anesthesiology* 123: 116–125.
- Kurz A, Sessler DI, Lenhardt R. Perioperative normothermia to reduce the incidence of surgical-wound infection and shorten hospitalization. *NEJM*. 9 maj 1996; 334(19): sid. 1209–1216.
- Melling AC, Ali B, Scott EM, Leaper DJ. Effects of preoperative warming on the incidence of wound infection after clean surgery: a randomised controlled trial. *The Lancet*. 15 sep. 2001; 358(9285): sid. 876–880.
- Schmied H, Kurz A, et al. Mild hypothermia increases blood loss and transfusion requirements during total hip arthroplasty. *The Lancet*. 1996;347(8997):289-292.
- Rajagopalan S, et al. The Effects of Mild Perioperative Hypothermia on Blood Loss and Transfusion Requirement. *Anesth*. 2008; 108:71-7.
- Frank SM, Fleisher LA, Breslow MJ, et al. Perioperative maintenance of normothermia reduces the incidence of morbid cardiac events. *JAMA*. 1997;277:1127-1134.
- Bush H Jr., Hydo J, Fischer E, et al. Hypothermia during elective abdominal aortic aneurysm repair: The high price of avoidable morbidity. *J Vasc Surg*. 1995;21(3): 392-402.
- Leslie K, Sessler DI, Bjorksten AR, Moayeri A. Mild hypothermia alters propofol pharmacokinetics and increases the duration of action of atracurium. *Anesth Analg*. 1995;80(5):1007–1014.
- Fritz HG, Bauer R, Walter B, Moertiz KU, Reinhart K. Effects of hypothermia (32°C) on plasma concentration of fentanyl in piglets (abstract). *Anesthesiology*. 1999;91(3A):A444.
- Heier T, Caldwell JE, Sessler DI, Miller RD. Mild intraoperative hypothermia increases duration of action and spontaneous recovery of vecuronium blockade during nitrous oxide-isoflurane anesthesia in humans. *Anesthesiology*. 1991;74(5):815–819.
- Bissonnette B, Sessler DI. Mild hypothermia does not impair postanesthetic recovery in infants and children. *Anesth Analg*. 1993;76(1):168–172.
- Lenhardt R, Marker E, Goli V, et al. Mild Intraoperative Hypothermia Prolongs Postanesthetic Recovery. *Anesth*. 1997; 87(6):1318-1323.
- Sessler, DI. Current Concepts: Mild perioperative hypothermia. *NEJM*, 1997. vol. 336, nr 24, sid. 1730–1737.
- Matsukawa T, Sessler DI, Sessler AM, Schroeder M, Ozaki M, Kurz A, Cheng C. Heat flow and distribution during induction of general anesthesia. *Anesth. Mars* 1995; 82(3): sid. 662–673.
- Sessler DI, Schroeder M, Merrifield B, Matsukawa T, Cheng C. Optimal duration and temperature of prewarming. *Anesth. Mars* 1995 1;82(3): sid. 674–681.
- Sessler DI, Ponte J. Hypothermia during epidural anesthesia results mostly from redistribution of heat within the body, not heat loss to the environment. *Anesth*. 1989;71:A882.
- Arkilic CF, Akça O, Taguchi A, Sessler DI, Kurz A: Temperature monitoring and management during neuraxial anesthesia: An observational study. *Anesth Analg* 2000; 91: sid. 662–666.
- Horn EP, Bein B, Bohm R, Steinfath M, Sahili N, Hocker J. The Effect of Short Time Periods of Pre-Operative Warming in the Prevention of Peri-Operative Hypothermia. *Anaesth*. 2012.67(6).
- Camus Y, Delva E, Sessler DI, Lienhart A. Pre-Induction Skin-Surface Warming Minimizes Intraoperative Core Hypothermia. *J Clinical Anesthesia*. 1995;7:384-388.
- Nelson G, Altman AD, et al. Guidelines for pre- and intra-operative care in gynecologic/oncology surgery: Enhanced Recovery After Surgery (ERAS(R)) Society recommendations - Part I. *Gynecologic Oncology*. 2016;140:313-322.
- American Society of PeriAnesthesia Nurses. Clinical guideline for the prevention of unplanned perioperative hypothermia. *J Perianesth Nurs*. 2001;16:305-314.
- Guideline for prevention of unplanned patient hypothermia. I: Guidelines for Perioperative Practice. Denver, CO: AORN;2017: sid. 567–590.
- Ban KA, Minei JP, Laronga C, Harbrecht BG, Jensen EH, Fry DE, Itani KMF, Dellinger EP, Ko CY, Duane TM. American College of Surgeons and Surgical Infection Society: Surgical Site Infection Guidelines, Uppdatering 2016. *J Am Coll Surg* 2017;224: sid. 59–74.
- Anderson DJ, Podgorny K, et al. Strategies to Prevent Surgical Site Infections in Acute Care Hospitals: Uppdatering 2014. *Infection Control and Hospital Epidemiology*. 2014;35(6)
- Hynson, J.M., et al., The effects of preinduction warming on temperature and blood pressure during propofol nitrous oxide anesthesia. *Anesthesiology*, 1993. 79(2): s. 219–228, diskussion 21A-22A.
- Moayeri, A., et al., Pre-induction skin-surface warming prevents redistribution hypothermia. *Anesthesiology*, 1991. 75 Suppl(3A): sid. A1004.
- Glosten, B., et al., Preanesthetic skin-surface warming reduces redistribution hypothermia caused by epidural block. *Anesth Analg*, 1993. 77(3): sid. 488–493.
- Sessler, DI. Perioperative Heat Balance. *Anesth*. 2000;92:578-596.
- Sessler DI. Perioperative thermoregulation and heat balance. *The Lancet*. 2016;387(10038):2655-2664.
- Eshraghi Y, Nasr V, Parra-Sanchez I, Van Duren A, Botham M, Santoscoy T, Sessler DI. An Evaluation of a Zero-Heat-Flux Cutaneous Thermometer in Cardiac Surgical Patients. *Anesth Analg*. 2014;119(3):543-549.
- Kimberger O. Temperature monitoring in the OR – State of the art och uppdatering från 2012. *Trends Anaesth Crit Care*. 2013;3(1):8-12.
- Torossian A. Survey on intraoperative temperature management in Europe. *European Journal of Anaesthesiology (EJA)*. 2007;24(8):668-675.
- Fossum S, Hays J, Henson MM. A comparison study on the effects of prewarming patients in the outpatient surgery setting. *Journal of PeriAnesthesia Nursing*. 2001;16(3):187-194.
- Wilson L, Kolcaba K. Practical application of comfort theory in the perianesthesia setting. *Journal of PeriAnesthesia Nursing*. 2004;19(3): sid. 164–173.

3M Svenska AB
Medical Solutions Division
Bollstanäsvägen 3
191 89 Sollentuna
08-92 21 00

3M Ireland Limited
The Iveagh Building
The Park, Carrickmines
D18 X015
Irland
+353 (0)1 280 3555

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3M och Bair Hugger är varumärken som tillhör 3M Company. © 3M 2019. Med ensamrätt.