Emphaze™ AEX Hybrid Purifier

**Laboratory Scale Capsules Installation and Operation
Instructions**

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

EN

English Instructions are the Original Instructions. All Other Languages are Translations of the Original Instructions.

Safety Information

Read, understand, and follow all safety information contained in these instructions prior to installation and use of 3M™

Emphaze™ AEX Hybrid Purifier Laboratory Capsules.
Retain these instructions for future reference.

Intended Use:

Single use processing of aqueous based biological pharmaceuticals (drugs) and vaccines strictly following the product operating instructions and cGMP requirements, where applicable. Customers must determine whether the 3M product is suitable for a specific application based on a risk assessment that considers the product leachable characteristics and its impact on drug safety.

Prohibited Use:

For laboratory use only. Not intended for use with materials that will be used on humans or animals. Do not use as a component in a medical device that is regulated by any agency, and/or globally exemplary agencies, including but not limited to: a) FDA, b) European Medical Device Directive (MDD), c) Japan Pharmaceuticals and Medical Devices Agency. Do not use in applications involving permanent implantation into the body, life-sustaining medical applications, applications requiring Food Contact compliance.

Explanation of Signal Word Consequences

 WARNING:	Indicates a hazardous situation which, if not avoided, could result in serious injury or death.
 CAUTION:	Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury and/or property damage.
NOTICE:	Indicates a situation, which, if not avoided, could result in property damage.

Explanation of Safety and Related Symbols

	WARNING: Explosion hazard
	Single use
	Wear PPE
	Storage conditions

 **WARNING**

To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

- Use the installation procedure specified in Operation Instructions.
- The capsule must be subjected to a preconditioning flush prior to use. Use the procedure specified in Operation Instructions.
- At least 54 L/m² preconditioning flush volume.
- Carry out preconditioning flush procedures of 3M™ Zeta Plus™ capsules and 3M™ Emphaze™ AEX Hybrid Purifier capsules separately.
- After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline should be reached (e.g., ± 0.5 pH units from the pH of the buffer used to flush the media).
- Capsules are designed for single use only. Do not reuse.
- Do not store capsules under conditions outside of the recommended storage temperature range.
- Do not use past the expiration date.
- Do not use with organic solvents or flammable liquids.
- Use only with aqueous based process fluids in the pH range of 5 to 9.
- Do not use with or expose this product to hot liquids >40°C (104°F).
- Do not use the product on unlevelled surfaces.
- Do not gamma or ethylene oxide sterilize the capsules.
- Do not steam-in-place sterilize the capsules.
- Only expose the capsule to 1 sterilization or 1 sanitization cycle prior to use.
- Do not autoclave and base sanitize the same capsule prior to use.
- Do not autoclave or base sanitize capsules with the following part numbers BV8 Luer (EMP301AEX020A) and BV8 Sanitary (EMP303AEX020A).

 **WARNING**

To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

- During pre-use autoclave sterilization do not exceed a temperature of 121°C or a duration longer than 30 minutes for numbers BV8R Luer (EMP301AEX020R) and BV8R Sanitary (EMP303AEX020R).
- Do not use with or expose this product to process liquids until capsule temperature goes <40°C (104°F).
- Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1-hour duration for sanitization.
- Do not NaOH sanitize 3M Emphaze AEX Hybrid Purifier and 3M Zeta Plus Encapsulated capsules together in series.
- Use only 70% ethyl and/or isopropyl alcohol to wipe the exterior of the capsules.
- Do not use with or expose this product to pressurized steam or sanitizing chemicals other than NaOH.

 **WARNING**

To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:

- Use the installation procedure specified in Operation Instructions.
- Wear Personal Protective Equipment (PPE) to prevent chemical splash.
- Hardware Qualification testing is recommended prior to use.
- The capsule must be subjected to a preconditioning flush prior to use. Use the procedure specified in Operation Instructions.
- After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline should be reached (e.g., ± 0.5 pH units from the pH of the buffer used to flush the media).
- Capsules are designed for single use only. Do not reuse.
- Do not modify or alter the capsules.
- Do not put this product into service if any damage is occurred to the capsule or any of the connecting parts.
- Do not install where fluid pressure exceeds 2.8 bar (40 psig).
- Do not use in organic solvent based processing.

 **WARNING**




To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:

- Use only with aqueous based process fluids in the pH range of 5 to 9.
- Do not use with or expose this product to hot liquids >40°C (104°F).
- Do not use the product on unlevelled surfaces.
- Do not use the product for more than 8 hours.
- Do not use this product for continuous service with compressed gasses. Use of compressed gas is permissible for post-use integrity testing and blow-down purposes only. Limit post-use integrity test and blow-down gas pressures to less than 1.4 bar (20psig) and temperature less than 25°C for no longer than 30 minutes.
- Do not hang other items on this product.
- Do not gamma or ethylene oxide sterilize the capsules.
- Do not steam-in-place sterilize the capsule.
- Only expose the capsule to 1 sterilization or 1 sanitization cycle.
- Do not autoclave and base sanitize the same capsule prior to use.
- Do not autoclave or base sanitize capsules with the following part numbers BV8 Luer (EMP301AEX020A) and BV8 Sanitary (EMP303AEX020A).
- During pre-use autoclave sterilization do not exceed a temperature of 121°C or a duration longer than 30 minutes for numbers BV8R Luer (EMP301AEX020R) and BV8R Sanitary (EMP303AEX020R).
- Do not use with or expose this product to process liquids until capsule temperature goes <40°C (104°F).
- Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1-hour duration for sanitization.
- Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.
- Use only 70% ethyl and/or isopropyl alcohol to wipe the exterior of the capsules.
- Do not use with or expose this product to pressurized steam or sanitizing chemicals other than NaOH.


⚠ WARNING	
	<p>To reduce the risks associated with chemical exposure due to pressure barrier failure, connector failure, purification media or capsule failure from use of disinfectants within the product:</p> <ul style="list-style-type: none"> • Use up to 1 M Sodium Hydroxide (NaOH) for pre-use and post-use sanitization and up to 5% Bleach (NaClO) for post-use sanitization. • Do not expose the purification capsule to base sanitization chemicals for more than 60 minutes. • Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization. • Wear personal protective equipment (PPE) to prevent chemical splash.

⚠ CAUTION	
	<p>To reduce the risk of burn or exposure injuries associated with autoclave sterilization.</p> <ul style="list-style-type: none"> • Wear personal protective equipment (PPE) to handle hot capsules.


1. Operating Parameters and Product Specifications

Liquid		
Maximum Inlet Pressure	2.8 bar (40 psig)	
Maximum Differential Pressure	2.4 bar (35 psig)	
Maximum Operating Temperature	40°C (104°F)	
Post-Use Compressed Gas for Blow Down Only ¹		
Maximum Inlet Pressure	1.4 bar (20 psig)	
Maximum Gas Temperature	25°C (77°F)	
Required Preconditioning Flush		
Preconditioning Flush Volume	54 L/m ² (15.3 mL flush solution per mL of bed volume)	
Recommended Flow Rate	210 L/m ² /h (1 mL/min per mL of bed volume)	
Product Specifications		
	BV8R Luer	BV8R Sanitary
	EMP301AEX020R	EMP303AEX020R
Dimensions		
Height x Diameter	1 3/4" x 3" (4.5 cm x 7.7 cm)	3 5/16" x 3" (8.8 cm x 7.7 cm)
Bed Volume	8.5 mL	
Filter Surface Area	25 cm ²	
Approximate Weight		
Dry	71 g	77 g
Wet Post Blow Down	80 g	85 g
Material of Construction		
Filter Media	Quaternary Ammonium Functional Polymer Grafted to a Polypropylene Nonwoven	
Membrane	Polyamide	

Product Specifications		
	BV8R Luer	BV8R Sanitary
	EMP301AEX020R	EMP303AEX020R
Material of Construction		
Membrane Support	Polypropylene Nonwoven	
Capsule Shell	Polypropylene	
Edge Seal	Filled Polypropylene	
Luer Cap and Barb Connectors	Polypropylene	
Approximate Volume		
Capsule Fill Volume ²	13 mL	16 mL
Holdup Volume Post Blow Down ³	9 mL	
Pre-use Sterilization and Sanitization		
Autoclave Sterilization (Only 1 autoclave cycle)	121°C for 30 minutes using Pre-vac cycle	
Base Sanitization (Only 1 base sanitization cycle)	1M Sodium Hydroxide (NaOH) for 60 minutes	
Post-use Sterilization and Sanitization		
Autoclave Sterilization	Post use 1 autoclave cycle at 126°C for 30 minutes	
Sanitization	1M Sodium Hydroxide (NaOH) or 5% Bleach (NaClO) for post-use sanitization	
Required Preconditioning Flush Volume⁴		
Flush Volume, No Pre-treatment	130 mL	
Flush Volume, Post-autoclave	130 mL	
Flush Volume, Post Base Sanitization - 50 mM phosphate buffered saline at pH 7.4 ⁵	130 mL	

Product Specifications		
	BV8R Luer	BV8R Sanitary
	EMP301AEX020R	EMP303AEX020R
Miscellaneous		
Connector	Luer	3/4" Sanitary Connectors
Recommended Use Flow Rate	8 mL/min	
Maximum Flow Rate ⁶	24 mL/min	
Storage Conditions ⁷	Controlled indoor temperature: 0 - 30°C (32 - 86°F) in original packaging	
Shelf-Life ⁸	2 years from date of manufacture @ 30°C maximum storage temperature	

⚠ WARNING

¹ Do not use this product for continuous service with compressed gasses. Use of compressed gas is permissible for post-use integrity testing and blow-down purposes only. Limit post-use integrity test and blow-down gas pressures to less than 1.4 bar (20psig) and temperature less than 25°C for no longer than 30 minutes.

² Capsule Fill Volume is defined as the volume of liquid that is required to fill the capsule.

³ Post Blow down Hold up Volume is defined as the volume of the residual liquid after gas (air or nitrogen) blow down.

⁴ The Required Preconditioning Flush is a saline solution, having a conductivity in the same range as that of 25 to 150 mM NaCl (typical 3-16 mS/cm). DO NOT use water to flush the capsule.

⁵ Review Regulatory Support File for recommended alternate buffer selection.

⁶ Do not allow the liquid differential pressure across the capsule to exceed 2.4 bar (35 psid). Maximum flow rate may be limited by system differential pressure.

⁷ Do not allow product to be exposed to freezing or extreme heat conditions.

⁸ USP <88> Class VI compliance and high adsorption capacity for up to 2 years from date of manufacture if stored as indicated. Product labeling includes the expiration date.

Do not steam-in-place sterilize the capsules.

Preconditioning Flush Procedure

The Laboratory Capsule must be flushed prior to use to remove the glycerin, that is present as a media stabilizer. Flush the capsule with 54 L/m² (15.3 mL/bed volume) of a saline solution or buffer having a conductivity in the same range as that of 25 to 150 mM NaCl (3 - 16 mS/cm) at a flow rate no higher than 210 L/m²/h (LMH). DO NOT use water as the flush solution. Flush in the forward flow direction only (see Figure 2). Intermittently bleed air off of the upstream side of the capsule using the Luer Capped Air Vent on the inlet side of the capsule.

 **WARNING**

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

- The capsule must be subjected to a preconditioning flush prior to use.

2. Capsule Use with No Pre-use Sterilization and Sanitization

2.A. Required Preconditioning Flush Procedure

1. Connect the capsule inlet and outlet with the appropriate tubing and process equipment.
2. Connect the capsule inlet to the reservoir containing the pre-conditioning flush solution.
3. Open the capsule vent.
4. Close the outlet valve and open the inlet valve slowly to begin the flow of pre-conditioning flush solution through the capsule inlet.
5. When fluid emerges from the vent, close the vent while simultaneously opening the outlet valve slowly.
6. Pump the appropriate amount of flush solution per capsule size to complete the required preconditioning flush.
7. Close the inlet valve and stop the flow of pre-conditioning flush solution.


IMPORTANT NOTICE: Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.

2.B. Single-stage Purification

1. Connect the inlet and outlet of the capsule to the upstream and downstream process equipment, respectively.
2. Connect the capsule inlet to the reservoir containing the bioprocess fluid.
3. Start the flow of bioprocess fluid and collect purified permeate.

IMPORTANT NOTICE: For optimal performance, run the filtration process at a flux of 160 – 630 LMH (0.75-3.0 Bed Volume per minute).

4. Monitor the efficiency of the bioprocess purification experiment by taking aliquots and analyzing the performance.
5. Stop the flow of bioprocess fluid at the appropriate volume based by on the laboratory study design.

6. The bioprocess fluid remaining on the upstream side of the capsule may then be processed through the system:
 - a. Open the vent to reduce system pressure to 0 bar (0 psig).
 - b. Close the inlet valve and open the outlet to drain.
 - c. Apply 0.3-0.7 bar (5-10 psig) air or nitrogen gas pressure to vent. Using a non-surfing pump or clean air supply to deliver regulated low-pressure air to push the process fluid through the capsule. Do not allow the inlet pressure to pulse or exceed 1.4 bar (20 psig).
 - d. Monitor the capsule outlet. When liquid flow is reduced to a trickle, reduce the system pressure to 0 bar (0 psig).

3. Capsule Use with Pre-use Autoclave Sterilization

3.A. Recommended Pre-use Autoclave Sterilization of Capsules


1. Wrap the capsule connections and vent with sterilization wrap.
2. Place the capsule in the autoclave.
3. Set-up the autoclave cycle.
 - a. Select pre-vac cycle for all capsule sizes.
 - b. Set the sterilization temperature at 121°C.
 - c. Set the autoclave duration for capsules for 30 minutes.
4. Start the autoclave cycle.
5. At the end of autoclave cycle, remove the capsules from the autoclave and let the capsule cool down at ambient temperature for at least 60 minutes.
6. Remove the Sterilization Wrap.

 **CAUTION**

Wear personal protective equipment (PPE) to handle hot capsules.

- Wear PPE to handle hot capsules.

3.B. Required Preconditioning Flush Procedure for the Pre-use Autoclaved Capsules

⚠ WARNING

To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:

Do not use with or expose this product to process liquids until capsule temperature goes 40°C (104°F).

Perform Required Preconditioning Flush procedure as described in Section 2.A.

IMPORTANT NOTICE: Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.

3.C. Single-stage Purification

Perform Single-stage Purification as described in Section 2.B.

4. Capsule Use with Pre-use NaOH Sanitization

4.A. Pre-use NaOH Sanitization of Capsules


Recommended maximum conditions: Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1 hour duration for sanitization.

⚠ WARNING

To reduce the risk of injury from possible end user exposure associated with biological contamination or degradation products released from the purification media, or risk of property damage due to product loss associated with chemical release from purification media, and nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

Do not NaOH sanitize 3M Emphaze AEX Hybrid Purifier and 3M Zeta Plus Encapsulated capsules together in series.

1. Connect the capsule inlet to the reservoir containing the 1 M NaOH solution.
2. Open the capsule vent.

⚠ WARNING

To reduce the risks associated with chemical exposure due to pressure barrier failure, connector failure, purification media or capsule failure from use of disinfectants within the product: Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.

3. When fluid emerges from the vent, close the vent while simultaneously opening the outlet valve slowly.
4. When fluid emerges from the capsule outlet, stop the flow and close the capsule outlet valve.
5. Allow the capsule to incubate with the 1 M NaOH solution. Incubation time should not exceed 1 hour.
6. After the incubation is complete, drain the NaOH solution from the capsule.

Important Notice: Do not blow down the 3M™ Emphaze™ AEX Hybrid Purifier capsules after NaOH sanitization procedure.

4.B. Required Preconditioning Flush Procedure for the Pre-use NaOH Sanitized Capsules

⚠ WARNING

To reduce the risk of Injury from possible end user exposure associated with degradation products released from the purification media; and

To reduce the risk of product loss associated with a) chemical release from purification media and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components.

- The capsule must be subjected to a preconditioning flush prior to use.

Perform Required Preconditioning Flush procedure as described in Section 2.A.

IMPORTANT NOTICE: Do not blow down with a gas the 3M Emphaze AEX Hybrid Purifier capsules after preconditioning flush procedure.

Capsule is ready for immediate use.

IMPORTANT NOTICE: Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.

Perform Single-stage Purification as described in Section 2.B.

5. Disposal

1. After fluid processing is complete, drain the fluid from the capsule until the internal pressure on the capsule is reduced to 0 bar (0 psig).
2. Drain any remaining fluid from the system and remove the external connections from the inlet, outlet, and vent.

CAUTION

To reduce the risk of fluid exposure:

- Evaluate the risk of your process liquids and determine an appropriate protocol to dispose of spent capsules and all other waste.
- Dispose of capsules and all other waste in accordance with federal, state, and local regulations or your local country's laws and regulations.

Intended Use:

Single use processing of aqueous based biological pharmaceuticals (drugs) and vaccines strictly following the product operating instructions and cGMP requirements, where applicable. Customers must determine whether the 3M product is suitable for a specific application based on a risk assessment that considers the product leachable characteristics and its impact on drug safety.

Prohibited Use:

For laboratory use only. Not intended for use with materials that will be used on humans or animals. Do not use as a component in a medical device that is regulated by any agency, and/or globally exemplary agencies, including but not limited to: a) FDA, b) European Medical Device Directive (MDD), c) Japan Pharmaceuticals and Medical Devices Agency. Do not use in applications involving permanent implantation into the body, life-sustaining medical applications, applications requiring Food Contact compliance.

Product Selection and Use:


Many factors beyond 3M's control and uniquely within user's knowledge and control can affect the use and performance of a 3M product in a particular application. As a result, customer is solely responsible for evaluating the product and determining whether it is appropriate and suitable for customer's application, including conducting a workplace hazard assessment and reviewing all applicable regulations and standards (e.g., OSHA, ANSI, etc.). Failure to properly evaluate, select, and use a 3M product and appropriate safety products, or to meet all applicable safety regulations, may result in injury, sickness, death, and/or harm to property.

Warranty, Limited Remedy, and Disclaimer:

Unless a different warranty is specifically stated on the applicable 3M product packaging or product literature (in which case such warranty governs), 3M warrants that each 3M product meets the applicable 3M product specification at the time 3M ships the product. 3M MAKES NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR CONDITION OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR ARISING OUT OF A COURSE OF DEALING, CUSTOM, OR USAGE OF TRADE. If a 3M product does not conform to this warranty, then the sole and exclusive remedy is, at 3M's option, replacement of the 3M product or refund of the purchase price.

Limitation of Liability

Except for the limited remedy stated above, and except to the extent prohibited by law, 3M will not be liable for any loss or damage arising from or related to the 3M product, whether direct, indirect, special, incidental, or consequential (including, but not limited to, lost profits or business opportunity), regardless of the legal or equitable theory asserted, including, but not limited to, warranty, contract, negligence, or strict liability.

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3M, Emphaze and Zeta Plus are trademarks of 3M
All rights reserved./Alle Rechte vorbehalten.
3M Purification Inc.