
**3M™ Emphaze™ AEX Hybrid Purifier Production  
Scale Capsules**

BV800R (EMP710AEX020R)  
BV5600R (EMP770AEX020R)

<b>EN</b>	Safety Information & Installation and Operation Instructions . . . . .	1-16
<b>DE</b>	Sicherheitsinformationen und Installationsund Bedienungshandbuch . . . . .	18-34
<b>ES</b>	Información de seguridad, instrucciones de funcionamiento e instalación . . . . .	36-52
<b>IT</b>	Informazioni sulla sicurezza e istruzioni di installazione e funzionamento. . . . .	54-70
<b>FR</b>	Consignes de sécurité et instructions d'installation et d'utilisation . . . . .	72-88
<b>NL</b>	Veiligheidsinformatie en instructies voor plaatsing en bediening . . . . .	90-106

**EN** English Instructions are the Original Instructions. All Other Languages are Translations of the Original Instructions.

Failure to follow installation, operation, and maintenance instructions may result in system failure and will void warranty.

**Safety Information**

**Read, understand, and follow all safety information contained in these instructions prior to installation and use of the 3M™ Emphaze™ AEX Hybrid Purifier. Retain these instructions for future reference.**


**Intended Use:**

Single use processing of aqueous based biological pharmaceuticals (drugs) and vaccines strictly following the product operating instructions and cGMP requirements, where applicable. Customers must determine whether the 3M product is suitable for a specific application based on a risk assessment that considers the product leachable characteristics and its impact on drug safety.

**Prohibited Use:**

Do not use as a component in a medical device that is regulated by any agency, and/or globally exemplary agencies, including but not limited to: a) FDA, b) European Medical Device Directive (MDD), c) Japan Pharmaceuticals and Medical Devices Agency (PMDA), applications involving permanent implantation into the body, life-sustaining medical applications, applications requiring Food Contact compliance.

<b>Explanation of Signal Word Consequences</b>	
⚠ <b>WARNING:</b>	Indicates a hazardous situation which, if not avoided, could result in serious injury or death.
⚠ <b>CAUTION:</b>	Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury and/or property damage.
<b>NOTICE:</b>	Indicates a situation which, if not avoided, could result in property damage.
<b>Explanation of Safety and Related Symbols</b>	
	<b>WARNING:</b> Explosion hazard
	Single use
	Wear PPE
	Heavy parts
	Storage conditions

**⚠ WARNING**

To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

- Use the installation procedure specified in Operation Instructions.
- The capsule must be subjected to a preconditioning flush prior to use. Use the procedure specified in Operation Instructions.
- At least 54 L/m<sup>2</sup> preconditioning flush volume.
- Carry out preconditioning flush procedures of 3M™ Zeta Plus™ capsules and 3M™ Emphaze™ AEX Hybrid Purifier capsules separately.
- After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline has been reached (e.g., ± 0.5 pH units from the pH of the buffer used to flush the media).
- Capsules are designed for single use only. Do not reuse.
- Do not store capsules under conditions outside of the recommended storage temperature range.
- Do not use past the expiration date.
- Do not use with organic solvents or flammable liquids.
- Use only with aqueous based process fluids in the pH range of 5 to 9.
- Use only alkaline resistant capsule and manifold sets.
- Do not use with or expose this product to hot liquids >40°C (104°F).
- Perform all required maintenance on or before its scheduled date and always use specified replacement parts in the Service & Maintenance Details.
- Do not use the product on unlevelled surfaces.
- Do not gamma or ethylene oxide sterilize the capsules.
- Do not steam-in-place sterilize the capsules.
- Only expose the capsule to 1 sterilization or 1 sanitization cycle prior to use.
- Do not autoclave and base sanitize the same capsule prior to use.
- Do not autoclave or base sanitize polycarbonate based capsules and manifolds with the following part numbers BV800 (EMP710AEX020A) and BV5600 (EMP770AEX020A).
- During pre-use autoclave sterilization do not exceed a temperature of 121°C or a duration longer than 30 minutes BV800R (EMP710AEX020R), and 40 minutes for BV5600R (EMP770AEX020R).
- Do not use with or expose this product to process liquids until capsule temperature goes <40°C (104°F).
- Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1-hour duration for sanitization.
- Do not NaOH sanitize 3M Emphaze AEX Hybrid Purifier and 3M Zeta Plus Encapsulated capsules together in series.
- Use only 70% ethyl and/or isopropyl alcohol to wipe the exterior of the capsules and manifolds.
- Do not use with or expose this product to pressurized steam or sanitizing chemicals other than NaOH.

**⚠ WARNING**

**To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:**

- Use the installation procedure specified in Operation Instructions.
- Follow the installation and operation manual and inspect the o-rings.
- Hardware Qualification testing is recommended prior to use.
- Wear Personal Protective Equipment (PPE) to prevent chemical splash.
- After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline has been reached (e.g.,  $\pm 0.5$  pH units from the pH of the buffer used to flush the media).
- Capsules are designed for single use only. Do not reuse.
- Do not modify or alter the capsules, manifolds and holders.
- Do not put this product into service if any damage is occurred to the capsule or any of the connecting parts.
- Do not install where fluid pressure exceeds 3.4 bar (50 psig).
- Do not use in organic solvent based processing.
- Use only with aqueous based process fluids in the pH range of 5 to 9.
- Use only alkaline resistant capsule and manifold sets.
- Do not use with or expose this product to hot liquids  $>40^{\circ}\text{C}$  ( $104^{\circ}\text{F}$ ).
- Perform all required maintenance on or before its scheduled date and always use specified replacement parts in the Service & Maintenance Details.
- Do not use the product on unlevelled surfaces.
- Do not use the product for more than 8 hours.
- Do not use this product for continuous service with compressed gasses. Use of compressed gas is permissible for post-use integrity testing and blow-down purposes only. Limit post-use integrity test and blow-down gas pressures to less than 1.4 bar (20 psig) and temperature less than  $25^{\circ}\text{C}$  for no longer than 30 minutes.
- Do not hang other items on this product.
- Do not gamma or ethylene oxide sterilize the capsules.
- Do not steam-in-place sterilize the capsules.
- Only expose the capsule to 1 sterilization or 1 sanitization cycle prior to use.
- Do not autoclave and base sanitize the same capsule prior to use.
- Do not autoclave or base sanitize polycarbonate based capsules and manifolds with the following part numbers BV800 (EMP710AEX020A) and BV5600 (EMP770AEX020A).
- During pre-use autoclave sterilization do not exceed a temperature of  $121^{\circ}\text{C}$  or a duration longer than 30 minutes for BV800R (EMP710AEX020R), and 40 minutes for BV5600R (EMP770AEX020R).
- Do not use with or expose this product to process liquids until capsule temperature goes  $<40^{\circ}\text{C}$  ( $104^{\circ}\text{F}$ ).
- Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1-hour duration for sanitization.

**⚠ WARNING**

**To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:**

- Do not expose the purification capsule to base sanitization chemicals more than 60 minutes.
- Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.
- Use only 70% ethyl and/or isopropyl alcohol to wipe the exterior of the capsules and manifolds.
- Do not use with or expose this product to pressurized steam or sanitizing chemicals other than NaOH.

**⚠ WARNING**

**To reduce the risks associated with chemical exposure due to pressure barrier failure, connector failure, purification media or capsule failure from use of disinfectants within the product:**

- Use up to 1 M Sodium Hydroxide (NaOH) for pre-use and post-use sanitization and up to 5% Bleach (NaClO) for post-use sanitization.
- Do not expose the purification capsule to base sanitization chemicals for more than 60 minutes.
- Do not use polycarbonate based capsules and manifolds for base sanitization.
- Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.
- Wear personal protective equipment (PPE) to prevent chemical splash.

**⚠ CAUTION**


**To reduce the risk of burn or exposure injuries associated with autoclave sterilization.**

- Wear personal protective equipment (PPE) to handle hot capsules.

**⚠ CAUTION**

**To reduce the risk of fluid exposure:**

- Evaluate the risk of your process liquids and determine an appropriate protocol to dispose of spent capsules and all other waste.
- Dispose of capsules and all other waste in accordance with federal, state, and local regulations or your local country's laws and regulations.
- The protective caps on the capsule and manifold connections are designed for the protection of capsule parts. The protective caps are not intended to contain process liquids during the disposal of spent capsules or manifolds.


**Table of Contents**

1. Operating Parameters and Product Specifications ..... page 5  
 2. Capsule and Manifold Design ..... page 7  
 3. 3M™ Zeta Plus™ Encapsulated System Holders..... page 7  
 4. Single Stage Purification Instructions for Pilot-Scale (16EZA) and Production-Scale (16EZB) Holder ..... page 7  
 5. Two Stage Purification Instructions for Pilot-Scale (16EZA) and Production-Scale (16EZB) Holder ..... page 10  
 6. Capsule Use with Pre-use Autoclave Sterilization ..... page 13  
 7. Capsule Use with Pre-use NaOH Sanitization..... page 13  
 8. Pilot-Scale (16EZA) and Production-Scale (16EZB) Holder Unloading Instructions ..... page 15  
 Important Notice and Disclaimer ..... page 16

**1. Operating Parameters and Product Specifications**

**Table 1: Operating Parameters**

Liquid		
Maximum Inlet Pressure	3.4 bar (50 psig)	
Maximum Differential Pressure	2.4 bar (35 psid)	
Maximum Operating Temperature	40°C (104°F)	
Post-Use Compressed Gas for Blow Down Only <sup>1</sup>		
Maximum Inlet Pressure <sup>1</sup>	1.4 bar (20 psig)	
Maximum Operating Temperature	25°C (77°F)	
Liquid		
Preconditioning Flush Volume	54 L/m <sup>2</sup> (15.3 mL flush solution per mL of bed volume)	
Recommended Flow Rate	210 L/m <sup>2</sup> /h (1 mL/min per mL of bed volume)	
Product Specifications		
	BV800R	BV5600R
	EMP710AEX020R	EMP770AEX020R
Dimensions		
Height x Diameter	3" x 17 3/4" (7.6 cm x 45 cm)	8 3/4" x 17 3/4" (22.2 cm x 45 cm)
Bed Volume	800 mL	5600 mL
Filter Surface Area	0.23 m <sup>2</sup>	1.6 m <sup>2</sup>
Approximate Weight		
Dry	3.5 kg	9.9 kg
Wet Post Blow Down	4.1 kg	14.6 kg
Material of Construction		
Filter Media	Quaternary Ammonium Functional Polymer Grafted to a Polypropylene Nonwoven	
Membrane	Polyamide	
Membrane Support	Polypropylene Nonwoven	

Capsule Shell and Manifold	Polyphenylene Oxide/Polystyrene	
Separator and Edge Seal	Filled Polypropylene	
<b>Product Specifications</b>		
	<b>BV800R</b>	<b>BV5600R</b>
	EMP710AEX020R	EMP770AEX020R
End Cap	Polypropylene	
O-ring	Silicone	
Handles	N/A	Filled Polyamide
<b>Approximate Volume</b>		
Capsule Fill Volume <sup>2</sup>	3.4 L	10.6 L
Holdup Volume Post Blow Down <sup>3</sup>	0.7 L	4.7 L
<b>Pre-use Sterilization and Sanitization</b>		
Autoclave Sterilization (Only 1 autoclave cycle)	121°C for 30 minutes using Pre-vac cycle	121°C for 40 minutes using Pre-vac cycle
Base Sanitization (Only 1 base sanitization cycle)	1M Sodium Hydroxide (NaOH) for 60 minutes	
<b>Post-use Sterilization and Sanitization</b>		
Autoclave Sterilization	Post use 1 autoclave cycle at 126°C for 30 minutes	
Sanitization	1M Sodium Hydroxide (NaOH) or 5% Bleach (NaClO) for post-use sanitization	
<b>Required Preconditioning Flush Volume<sup>4</sup></b>		
Flush Volume, No Pre-treatment	12 L	85 L
Flush Volume, Post-autoclave	12 L	85 L
Flush Volume, Post Base Sanitization - 50 mM phosphate buffered saline at pH 7.4 <sup>5</sup>	24 L	85 L
<b>Miscellaneous</b>		
Connector	1 ½" Sanitary Connectors	
Recommended Use Flow Rate	0.8 L/min	5.6 L/min
Maximum Flow Rate <sup>6</sup>	2.4 L/min	16.8 L/min
Storage Conditions <sup>7</sup>	Controlled indoor temperature: 0 - 30°C (32 - 86°F) in original packaging	
Shelf-Life <sup>8</sup>	2 years from date of manufacture @ 30°C maximum storage temperature	

**⚠ WARNING**

<sup>1</sup> Do not use this product for continuous service with compressed gasses. The use of compressed gas is permissible for post-use integrity testing and blow-down purposes. Limit post-use integrity test and blow-down gas pressures to less than 1.4 bar (20 psig) and temperature less than 25°C. for no longer than 30 mins.

<sup>2</sup> Capsule Fill Volume is defined as the volume of liquid that is required to fill the capsule.

<sup>3</sup> Post Blow Down Hold up Volume is defined as the volume of the residual liquid after gas (air or nitrogen) blow down.


<sup>4</sup> The Required Preconditioning Flush is a saline solution, having a conductivity in the same range as that of 25 to 150 mM NaCl (typical 3-16 mS/cm). DO NOT use water to flush the capsule.


<sup>5</sup> Review Regulatory Support File for recommended alternate buffer selection.

<sup>6</sup> Do not allow the liquid differential pressure across the capsule to exceed 2.4 bar (35 psi). Maximum flow rate may be limited by system differential pressure.

<sup>7</sup> Do not allow product to be exposed to freezing or extreme heat conditions.

<sup>8</sup> USP <88> Class VI compliance and high adsorption capacity for up to 2 years from date of manufacture if stored as indicated. Product labeling includes the expiration date.

## 2. Capsule and Manifold Design

Two capsule sizes, BV800R and BV5600R, are to be used with the 3M™ Zeta Plus™ Encapsulated System holders. The BV800R capsule has 0.23 m<sup>2</sup> filtration area (800 ml bed volume). The BV5600R capsule has 1.6 m<sup>2</sup> (5600 ml bed volume). In all production capsule configurations, one male and one female adaptor is attached to opposing sides of each BV800R and BV5600R lenticle stacks. The male and female adaptors mate with the adaptors on the manifolds and the neighboring capsules. The o-rings on the adaptors on the outer most lenticles provide a secure seal between the molded mating features of the manifolds and neighboring capsules. The manifolds have imprinted information indicating top/bottom manifolds and inlet/outlet connection ports. The BV5600R capsule shells have a self-guiding locking mechanism for a fast, easy and robust capsule-to-capsule connection. A set of manifolds is required for connecting the capsules to external components of the purification train. The manifold set ensures the entire purification process is fast and clean. An extra pair of manifolds is required to perform two-stage filtration within one holder assembly. Only use Polyphenylene Oxide/Polystyrene based manifold sets.

## 3. 3M™ Zeta Plus™ Encapsulated System Holders


3M offers two holders; Pilot-scale (model 16EZA) and Production-scale (model 16EZB) holders. The pilot-scale holder can accommodate from one to four BV800R capsules or one BV5600R capsule. Either single-stage or two-stage (BV800R only) purification can be performed within the same pilot scale holder. In a two stage purification train, a 3M Zeta Plus filter is used as the first stage. The production scale holder can accommodate up to seven BV5600R capsules. Similarly, either a single-stage or two-stage purification can be performed within the same production scale holder, wherein the 1st stage is 3M Zeta Plus filter(s).

## 4. Single Stage Purification Instructions for Pilot-Scale (16EZA) and Production-Scale (16EZB) Holder

### 4.A. Pilot-Scale Holder (16EZA) Loading Instructions for Single Stage Purification

Loading is performed in front of the holder (two stop pins are located in the back of the holder base).

1. Ensure the holder is on a level surface and lock the four swivel casters.
2. Turn the hand wheel counterclockwise to move the movable end plate upward to desired height.
3. Place a bottom manifold onto the base with the inlet and outlet facing forward, the back side of the manifold touching the two stop pins, and the manifold stops positioned between the support rods. The female connector of the manifold will face upward when the bottom manifold is in correct orientation.
4. Rotate the manifold to desired location within the support rods.
5. Lubricate all o-rings using water or process fluid.
6. Place one BV5600R capsule on top of the manifold with the female connector facing upward.
7. If a BV800R capsule is used, place one BV800R capsule on top of the bottom manifold with the female connector facing upward.

8. The 16EZA holder can accommodate from one to four BV800R capsules or one BV5600R capsule. If more than one BV800R capsule is needed, place the second capsule on top of the first capsule with the female connector facing up.
9. Repeat step 8 until all the capsules are loaded.
10. Place the top manifold on top of the last capsule that is loaded. The male connector of the top manifold will face downward, the upstream and downstream vents will face forward, and the manifold stops will be positioned between the support rods when the top manifold is in the correct orientation.
11. Turn the hand wheel clockwise until you feel the torque limiter slip and moveable end plate is stationary.

#### **4.B. Production-Scale Holder (16EZB) Loading Instructions for Single Stage Purification**

Loading is performed while facing the front of the empty holder with the capsule compression wheel on your right, with the holder in the horizontal position.

1. Ensure holder is on a level surface and engage both floor locks.

**IMPORTANT NOTICE: Holder must be in the horizontal position for loading.**


2. Pull the locking pin out and rotate the locking bar to the unlocked position.
3. Slide the left end plate to the far left.

**IMPORTANT NOTICE: When installing seven (7) capsules, the locking bar must be touching the rubber bumpers**

4. Standing facing the right end of the holder, fully retract the right end plate by rotating the capsule compression wheel counterclockwise.
5. Standing to the front of the holder place the bottom manifold against the right end plate with the female connector facing left, sanitary connectors facing upward and manifold stops positioned between notched rods.
6. Lubricate all o-rings using water or process fluid.

**IMPORTANT NOTICE: Inspect the o-rings and carry out a hardware qualification prior to use.**

7. Place the capsule in the holder so that the capsule stops rest against the front notched rod.
8. Slide the capsule to the right until it touches the bottom manifold.
9. Push the capsule handles away from you until the straight handle stops against the back notched rod.
10. If using more than one capsule, place the next capsule in the holder so that the capsule stops rest on the front notched rod.
11. Slide the capsule to the right until it touches the previously installed capsule.
12. Push the subsequent capsule handles away from you until the straight handle stops against the back notched rod.
13. Repeat steps 10 - 12 until all the required capsule filters are installed.
14. Place the top manifold next to the last capsule installed with the male connector facing right, sanitary connectors facing upward and manifold stops positioned between notched rods.
15. Standing to either the front or left side of the holder, slide the left end plate until the locking bar is positioned over the notch closest to the top manifold and lock.
16. Rotate both manifolds to the desired locations within notched rods.
17. Standing facing the right end of the holder, turn the capsule compression wheel clockwise until it comes to a complete stop (a minimum torque of 100 inch-lbs [11.3 Nm] is required) and there are no visible gaps between all capsules and manifolds.
18. Standing to the back of the holder, turn the rotation wheel clockwise to pivot the holder to the vertical operating position.


**IMPORTANT NOTICE:** Ensure locking pin is in the locked position when holder is in any position other than horizontal.

**IMPORTANT NOTICE:** Perform all required maintenance on or before its scheduled date and always use specified replacement parts in the Service & Maintenance Details in Literature: 34-8722-0459-8

#### 4.C. Single Stage Purification Instructions

The Production Capsule must be flushed prior to use to remove the glycerin that is present as a media stabilizer. Flush the capsule with 54 L/m<sup>2</sup> (15.3 mL/bed volume) of a saline solution or buffer having a conductivity in the same range as that of 25 to 150 mM NaCl (3 - 16 mS/cm) at a flow rate no higher than 210 L/m<sup>2</sup>/h (LMH). DO NOT use water as the flush solution. Flush in the forward flow direction only. Intermittently bleed air off using the Air Vent on the manifold.

**IMPORTANT NOTICE:** Precondition flush the capsule prior to use.

#### 4.C.1 Required

##### Preconditioning Flush Procedure

1. Connect the inlet and outlet located on the bottom manifold to the upstream and downstream process equipment.

**IMPORTANT NOTICE:** Do not rotate the production-scale holder when process equipment is attached.


2. Connect vent valves to the upstream and downstream vent ports located on the top manifold.
3. Connect the capsule inlet to the reservoir containing the pre-conditioning flush solution.
4. Open both vent valves.
5. Close the outlet valve and open the inlet valve slowly to begin the flow of pre-conditioning flush solution through the inlet valve.
6. When fluid emerges from the upstream vent valve, close the upstream vent valve.
7. When fluid emerges from the downstream vent valve, close the valve while simultaneously opening the outlet valve slowly.
8. Pump the appropriate amount of flush solution per capsule size to complete the required preconditioning flush.
9. Close the inlet valve and stop the flow of pre-conditioning flush solution.

**IMPORTANT NOTICE:** Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.

**IMPORTANT NOTICE:** For optimal performance, run filtration process at a flux of 160 – 630 LMH (0.75-3.0 Bed Volume per minute).

#### 4.C.2. Single-stage Purification

1. Connect the inlet and outlet of the capsule to the upstream and downstream process equipment, respectively.
2. Connect the capsule inlet to the reservoir containing the bioprocess fluid.
3. Start the flow of bioprocess fluid and collect purified permeate.

**IMPORTANT NOTICE: For optimal performance, run the filtration process at a flux of 160 - 630 LMH (0.75-3.0 Bed Volume per minute).**

4. Monitor the efficiency of the bioprocess purification experiment by taking aliquots and analyzing the performance.
5. Stop the flow of bioprocess fluid at the appropriate volume based by on the laboratory study design.
6. The bioprocess fluid remaining on the upstream side of the capsule may then be processed through the system:
  - a. Open the vent to reduce system pressure to 0 bar (0 psig).
  - b. Close the inlet valve and open the outlet to drain.
  - c. Apply 0.3-0.7 bar (5-10 psig) air or nitrogen gas pressure to vent. Using a non-surging pump or clean air supply deliver regulated low-pressure air to push the process fluid through the capsule. Do not allow the inlet pressure to pulse or exceed 1.4 bar (20 psig).
  - d. Monitor the capsule outlet. When liquid flow is reduced to a trickle, reduce the system pressure to 0 bar (0 psig).

## 5. Two Stage Purification Instructions for Pilot-Scale (16EZA) and Production-Scale (16EZB) Holder

3M™ Emphaze™ AEX Hybrid Purifier Production Capsules can be used in a series flow path with 3M™ Zeta Plus™ Production Capsules. When using combinations of 3M Zeta Plus and 3M Emphaze AEX Hybrid Purifier Production capsules, the 3M Emphaze capsules should be placed downstream from the 3M Zeta Plus capsules.

**For Series Operation for two production capsules:**


**IMPORTANT NOTICE: Precondition flush each capsule individually prior to use by following the preconditioning flush procedure in Section 4.C.1.**

### 5.A. Pilot-Scale Holder (16EZA) Loading Instructions for Two-stage Purification

The pilot-scale holder is designed to accommodate two-stage filtration. A maximum of one 3M™ Zeta Plus™ (0.23 m<sup>2</sup>) capsule and one BV800R capsule can be used for each stage. A second set of manifolds is required for use between each purification train.

Loading is performed in front of the holder (two stop pins are located in the back of the holder base).

1. Ensure the holder is on a level surface and lock the four swivel casters.
2. Turn the hand wheel counterclockwise to move the movable end plate upward to desired height.
3. Place a bottom manifold onto the base with the inlet and outlet facing forward, the back side of the manifold touching the two stop pins, and the manifold stops positioned between the support rods. The female connector of the manifold will face upward when the bottom manifold is in the correct orientation.
4. Rotate the manifold to desired location within the support rods.
5. Lubricate both o-rings using water or process fluid.
6. Place one BV800R capsule on top of the manifold with the female connector facing upward.
7. Place the top manifold on the capsule. The male connector of the top manifold will face downward, the upstream and downstream vents will face forward, and the manifold stops will be positioned between the support rods when the top manifold is in the correct orientation.
8. Place a bottom manifold on the top of the 2nd Stage assembly with the inlet and outlet facing forward and the manifold stops positioned between the support rods. The female connector of the manifold will face upward when the bottom manifold is in the correct orientation.
9. Repeat steps 4 - 7 with a 3M Zeta Plus (0.23 m<sup>2</sup>) capsule filter. Manually align 1st Stage assembly and 2nd Stage assembly.
10. Turn the hand wheel clockwise until you feel the torque limiter slip and moveable end plate is stationary.

### 5.B. Production-Scale Holder (16EZB) Loading Instructions for Two-stage Purification


The holder is designed to accommodate two-stage filtration. A combination of two grades totaling six 1.6 m<sup>2</sup> capsules can be used within the holder. A second set of manifolds must be used between the 3M Zeta Plus and 3M Emphaze AEX Hybrid Purifier capsules.

Loading is performed while facing the front of the empty holder with the capsule compression wheel on your right, with the holder in the horizontal position.

1. Ensure holder is on a level surface and engage both floor locks.

**IMPORTANT NOTICE: Holder must be in the horizontal position for loading.**

2. Pull the locking pin out and rotate the locking bar to the unlocked position.
3. Slide the left end plate to the far left.
4. Standing to the right side of the holder, fully retract the right end plate by rotating the capsule compression wheel counterclockwise.
5. Standing to the front of the holder, place the bottom manifold against the right end plate with the female connector facing left, sanitary connectors facing upward and manifold stops positioned between notched rods.
6. Lubricate all o-rings using water or process fluid.
7. Place the BV5600R capsule in the holder so that the capsule stops rest against the front notched rod.
8. Slide the capsule to the right until it touches the bottom manifold.
9. Push the capsule handles away from you until the straight handle stops against the back notched rod.
10. If using more than one BV5600R capsule for the 2nd Stage assembly, place the next capsule in the holder so that the capsule stops rest on the front notched rod.
11. Slide the capsule to the right until it touches the previously installed capsule.
12. Push the subsequent capsule handles away from you until the straight handle stops against the back notched rod.
13. Repeat steps 10 - 12 until all the required 2nd Stage capsules are installed.
14. Place the top manifold next to the last capsule installed with the male connector facing right, sanitary connectors facing upward and manifold stops positioned between notched rods.
15. Place the bottom manifold to the left of the 2nd Stage assembly with the female connector facing left, sanitary connectors facing upward and manifold stops positioned between notched rods.
16. Repeat step 7 - 9 with a 3M Zeta Plus 1.6 m<sup>2</sup> capsules.
17. If using more than one capsule for the 1st Stage assembly, place the next capsule in the holder so that the capsule stops rest on the front notched rod.
18. Slide the capsule to the right until it touches the previously installed capsule.
19. Push the subsequent capsule handles away from you until the straight handle stops against the back notched rod.
20. Repeat steps 17 - 19 until all the required 1st Stage capsules are installed.
21. Place the top manifold next to the last capsule installed with the male connector facing right, sanitary connectors facing upward and manifold stops positioned between notched rods.
22. Standing to either the front or left side of the holder, slide the left end plate until the locking bar is positioned over the notch closest to the top manifold and lock.
23. Rotate both manifolds to the desired locations within notched rods.
24. Standing to the right of the holder, turn the capsule compression wheel clockwise until it comes to a complete stop (a minimum torque of 100 inch-lbs [11.3 Nm] is required) and there are no visible gaps between all capsules and manifolds.
25. Standing to the back of the holder, turn the rotation wheel clockwise to pivot the holder to the vertical operating position.

**IMPORTANT NOTICE:** Ensure locking pin is in the locked position when holder is in any position other than horizontal.

#### 5.C. Two-stage Purification

1. Configure the system as shown in Figure 2.
2. Close all valves (V1 - V8).
3. Position V8 to direct flow to the inlet.
4. Open all vent valves (V2, V3, V5 and V6).

**IMPORTANT NOTICE:** Carry out preconditioning flush procedures of 3M™ Zeta Plus™ capsules and 3M™ Emphaze™ AEX Hybrid Purifier capsules separately.

**IMPORTANT NOTICE: 3M Zeta Plus Capsules:** Pre-flushing the system with 54 L/m<sup>2</sup> of buffer or water is required. Maximum pre-flushing flux is 300 Liters/m<sup>2</sup>/hr (LMH). **3M Emphaze AEX Hybrid Purifier Capsules:** The capsules must be flushed prior to use to remove the glycerin that is present as a media stabilizer. Flush the capsule with 15.3 mL/bed volume of buffer or 25-150 mm sodium chloride solution. Do not use water as the flush solution.

5. Open V1 slowly to begin fluid flow.
6. When the fluid exits V3, close V3.
7. When fluid exits V2, close V2 and simultaneously position V4 to direct the fluid to V8.
8. When fluid exits V6, close V6.
9. When fluid exits V5, close V5 and simultaneously open V7.
10. Close V1 and then close V7 to stop flow.
11. Drain fluid from the system if required by turning V1, V4, V7 and V8 to the drain position. Then open V2, V3, V5, and V6.
12. Filter process fluid by repeating steps 2 – 10 with process fluid.

#### WARNING

To reduce the risk of Injury from possible end user exposure associated with degradation products released from the purification media; and


To reduce the risk of product loss associated with a) chemical release from purification media and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components.

- The capsule must be subjected to a preconditioning flush prior to use.

**IMPORTANT NOTICE:** DO NOT blow down 3M Emphaze AEX Hybrid Purifier capsules after pre-conditioning flush procedure.

**IMPORTANT NOTICE:** For optimal performance, run filtration process at a flux of 160-630 LMH (0.75-3.0 Bed Volume per minute).

13. At the end of filtration, the fluid remaining on the upstream side may be filtered through the system:
  - a. Open V3 and reduce the system pressure to zero (0) psi (0 bar).
  - b. Apply a maximum of 20 psi (1.4 bar) air pressure at a temperature less than 25°C to valve V3.
  - c. Monitor the system outlet. When liquid flow is reduced to a trickle from the 1st Stage capsule, reduce the system pressure to zero (0) psi (0 bar) and then close V3 and V8.
  - d. Open valve V6 and apply a maximum of 20 psi (1.4 bar) air pressure at a temperature less than 25°C to blow down the capsules post use.
  - e. Monitor the system outlet. When liquid flow is reduced to a trickle from the 2nd Stage capsule, reduce the system pressure to zero (0) psi (0 bar) and then close V6 and V8.

## 6. Capsule Use with Pre-use Autoclave Sterilization

### 6.A. Recommended Pre-use Autoclave Sterilization of Capsules

1. Remove the sanitary connection protective covers from the capsule and manifolds.
2. Wrap the capsule and manifold connections with sterilization wrap.
3. Place the capsule and manifolds in the autoclave.
4. Set-up the autoclave cycle:
  - a. Select pre-vac cycle for all capsule sizes
  - b. Set the sterilization temperature at 121°C
  - c. Set the autoclave duration:
    - i. BV800R capsules for 30 minutes
    - ii. BV5600R capsules for 40 minutes
5. Start the autoclave cycle.
6. At the end of autoclave cycle, remove the capsule and manifolds from the autoclave and let them cool down at ambient temperature for at least 60 minutes.
7. Remove the Sterilization Wrap.

#### CAUTION

To reduce the risk of burn or exposure injuries associated with autoclave sterilization:

- Wear personal protective equipment (PPE) to handle hot capsules.

### 6.B. Required Preconditioning Flush Procedure for the Pre-use Autoclaved Capsules

#### WARNING

To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components: And

To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:

Do not use with or expose this product to process liquids until capsule temperature goes <40°C (104°F).

Perform Required Preconditioning Flush procedure as described in Section 4.C.1.

**IMPORTANT NOTICE: Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.**

### 6.C. Single-stage Purification

Perform Single-stage Purification as described in Section 4.C.2.

### 6.D. Two-stage Purification

Perform Two-stage Purification as described in Section 5.C.

## 7. Capsule Use with Pre-use NaOH Sanitization

### 7.A. Pre-use NaOH Sanitization of 3M Emphaze AEX Hybrid Purifier Capsules

Recommended maximum conditions: Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1 hour duration for sanitization.

**⚠ WARNING**

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components: Do not NaOH sanitize 3M Emphaze AEX Hybrid Purifier and 3M Zeta Plus Encapsulated capsules together in series.

1. Load the manifolds and capsules to Pilot-Scale (16EZA) or Production-Scale (16EZB) holders as described in Sections 4.A. and 4.B (for Single-stage Purification) and 5.A and 5.B (for Two-stage Purification), respectively.
2. Connect the inlet and outlet located on the bottom manifold to the upstream and downstream process equipment.
3. Connect vent valves to the upstream and downstream vent ports located on the top manifold.
4. Connect the capsule inlet to the reservoir containing the 1M NaOH solution.
5. Open the both vent valves.
6. Close the outlet valve and open the inlet valve slowly to begin the flow of pre-conditioning flush solution through the inlet valve.

**⚠ WARNING**

To reduce the risks associated with chemical exposure due to pressure barrier failure, connector failure, purification media or capsule failure from use of disinfectants within the product: Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.

7. When fluid emerges from the upstream vent valve, close the upstream vent valve.
8. When fluid emerges from the downstream vent valve, close the downstream vent valve while simultaneously opening the outlet valve slowly.
9. When fluid emerges from the capsule outlet, stop the flow and close the capsule.
10. Allow the capsule to incubate with the 1 M NaOH solution. Incubation time should not exceed 1 hour.
11. After the incubation is complete, drain the NaOH solution from the capsule.

**IMPORTANT NOTICE: Do not blow down the 3M™ Emphaze™ AEX Hybrid Purifier capsules after NaOH sanitization procedure.**

**7.B. Required Preconditioning Flush Procedure for the Pre-use NaOH Sanitized Capsules**

**⚠ WARNING**

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline has been reach (e.g.,  $\pm 0.5$  pH units from the pH of the buffer used to flush the media).

Perform Required Preconditioning Flush procedure as described in Section 4.C.1.


**IMPORTANT NOTICE: Do not blow down with a gas the 3M Emphaze AEX Hybrid Purifier capsules after preconditioning flush procedure.**

Capsule is ready for immediate use

**IMPORTANT NOTICE: Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.**

**7.C. Single-stage Purification**

Perform Single-stage Purification as described in Section 4.C.2.

#### 7.D. Two-stage Purification

Perform Two-stage Purification as described in Section 5.C.

### 8. Pilot-Scale (16EZA) and Production-Scale (16EZB) Holder Unloading Instructions

#### 8A. Pilot-Scale Holder (16EZA) Unloading Instructions


1. After filtration is completed, drain the fluid from capsule filters until the internal pressure on the filter capsules is reduced to zero (0) psi (0 bar).
2. Remove the external connections from the inlet, outlet, and both vents.
3. Turn the hand wheel counterclockwise to raise the movable end plate upward to unseal the capsule assembly.
4. To unload the system (including capsules and manifolds), grasp the bottom manifold, slide system forward and lift system off holder from underneath bottom manifold. Repeat with 2nd stage assembly.
5. Evaluate the risk of your process liquids and determine an appropriate protocol to dispose of spent capsules and manifolds. Dispose of capsules and manifolds in accordance with federal, state and local regulations or your local country's laws regulations.

#### 8B. Production-Scale Holder (16EZB) Unloading Instructions

1. After filtration is completed, open the vent valves (V1-V4) to relieve the internal pressure until it is reduced to zero (0) psig (0 bar).
2. Drain any remaining fluid from the system and remove the external connections from the inlets, outlets, and vents.
3. Standing to the back of the holder, turn the rotation wheel counterclockwise to pivot the holder to the horizontal position.
4. Standing to the right of the holder, fully retract the right end plate by rotating the capsule compression wheel counterclockwise.
5. Standing to either the front or left side of the holder, pull the locking pin out and rotate the locking bar to the unlocked position.
6. Slide the left end plate to the far left.
7. To unload the system, work from left to right. Standing to the front of the holder, grasp the handle of the far left capsule and pull towards you while restraining the capsule immediately to the right. Remove the capsule from the system.
8. Repeat step 7 until all capsules are removed.

**IMPORTANT NOTICE: Manifolds may remain attached to capsule(s) during unloading.**

9. Evaluate the risk of your process liquids and determine an appropriate protocol to dispose of spent capsules and manifolds. Dispose of capsule(s) and manifolds in accordance with federal, state and local regulations or your local country's laws and regulations.

**Intended Use:**

Single use processing of aqueous based biological pharmaceuticals (drugs) and vaccines strictly following the product operating instructions and cGMP requirements, where applicable. Customers must determine whether the 3M product is suitable for a specific application based on a risk assessment that considers the product leachable characteristics and its impact on drug safety.

**Prohibited Use:**

Do not use as a component in a medical device that is regulated by any agency, and/or globally exemplary agencies, including but not limited to: a) FDA, b) European Medical Device Directive (MDD), c) Japan Pharmaceuticals and Medical Devices Agency (PMDA), applications involving permanent implantation into the body, life-sustaining medical applications, applications requiring Food Contact compliance.

**Product Selection and Use:**


Many factors beyond 3M's control and uniquely within user's knowledge and control can affect the use and performance of a 3M product in a particular application. As a result, customer is solely responsible for evaluating the product and determining whether it is appropriate and suitable for customer's application, including conducting a workplace hazard assessment and reviewing all applicable regulations and standards (e.g., OSHA, ANSI, etc.). Failure to properly evaluate, select, and use a 3M product and appropriate safety products, or to meet all applicable safety regulations, may result in injury, sickness, death, and/or harm to property.

**Warranty, Limited Remedy, and Disclaimer:**

Unless a different warranty is specifically stated on the applicable 3M product packaging or product literature (in which case such warranty governs), 3M warrants that each 3M product meets the applicable 3M product specification at the time 3M ships the product. 3M MAKES NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR CONDITION OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR ARISING OUT OF A COURSE OF DEALING, CUSTOM, OR USAGE OF TRADE. If a 3M product does not conform to this warranty, then the sole and exclusive remedy is, at 3M's option, replacement of the 3M product or refund of the purchase price.

**Limitation of Liability**

Except for the limited remedy stated above, and except to the extent prohibited by law, 3M will not be liable for any loss or damage arising from or related to the 3M product, whether direct, indirect, special, incidental, or consequential (including, but not limited to, lost profits or business opportunity), regardless of the legal or equitable theory asserted, including, but not limited to, warranty, contract, negligence, or strict liability.

**3M Purification Inc.**

400 Research Parkway  
Meriden, CT 06450  
U.S.A.

(800) 243-6894

(203) 237-5541


Fax (203) 630-4530

[www.3Mpurification.com](http://www.3Mpurification.com)

3M, Emphaze, and Zeta Plus are trademarks of  
3M Company. Please recycle. Printed in U.S.A.

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© 2018 3M Company. All rights reserved.