# Emphaze™ AEX Hybrid Purifier

**Laboratory Capsules Installation and Operation  
Instructions**

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
**EN**

English Instructions Are the Original Instructions. All Other Languages Are Translations of the Original Instructions.

### Safety Information

Read, understand, and follow all safety information contained in

these instructions prior to installation and use of 3M™ Emphaze™ AEX Hybrid Purifier Laboratory Capsules BVO.3R (EMP101AEX020R) and BV1R (EMP201AEX020R). Retain these instructions for future reference.

#### Intended Use:

Single use processing of aqueous based biological pharmaceuticals (drugs) and vaccines strictly following the product operating instructions and cGMP requirements, where applicable.

#### Product Selection and Use:

Many factors beyond 3M's control and uniquely within user's knowledge and control can affect the use and performance of a 3M product in a particular application. As a result, end-user is solely responsible for evaluating the product and determining whether it is appropriate and suitable for end-user's application, including completing a risk assessment that considers the product leachable characteristics and its impact on drug safety, conducting a workplace hazard assessment and reviewing all applicable regulations and standards (e.g., OSHA, ANSI, etc.). Failure to properly evaluate, select, and use a 3M product and appropriate safety products, or to meet all applicable safety regulations, may result in injury, sickness, death, and/or harm to property.

Explanation of Signal Word Consequences	
 <b>WARNING:</b>	Indicates a hazardous situation which, if not avoided, could result in serious injury or death.
 <b>CAUTION:</b>	Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury and/or property damage.
<b>NOTICE:</b>	Indicates a situation which, if not avoided, could result in property damage.
Explanation of Safety and Related Symbols	
	<b>WARNING:</b> Explosion hazard
	Single use
	Wear PPE
	Storage conditions

**Prohibited Use:**

For laboratory use only. Not intended for use with materials that will be used on humans or animals. Do not use as a component in a medical device that is regulated by any agency, and/or globally exemplary agencies, including but not limited to: a) FDA, b) European Medical Device Directive (MDD), c) Japan Pharmaceuticals and Medical Devices Agency. Do not use in applications involving permanent implantation into the body, life-sustaining medical applications, applications requiring Food Contact compliance.

**Table of Contents**

- 1. Operating Parameters and Product Specifications . . . . .page 8
- 2. Single Stage Use Instructions . . . . .page 12
- 3. Capsule Use with Pre-use Autoclave Sterilization. . . . .page 13
- 4. Capsule Use with Pre-use NaOH Sanitization . . . . .page 15
- 5. Disposal . . . . .page 18

 **WARNING**

**To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and**

**To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:**

- Use the installation procedure specified in Operation Instructions.
- The capsule must be subjected to a preconditioning flush prior to use. Use at least 54 L/m<sup>2</sup> preconditioning flush volume.
- Carry out preconditioning flush procedures of 3M™ Zeta Plus™ capsules and 3M™ Emphaze™ AEX Hybrid Purifier capsules separately.
- After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline has been reached (e.g., ± 0.5 pH units from the pH of the buffer used to flush the media).
- Capsules are designed for single use only. Do not reuse.
- Do not store capsules outside of the recommended conditions [0 - 30°C (32 - 86°F)] in original packaging.
- Do not use past the expiration date.
- Do not use with organic solvents or flammable liquids.
- Use only with aqueous based process fluids in the pH range of 5 to 9.
- Do not use with or expose this product to hot liquids >40°C (104°F).
- Do not use the product on unlevelled surfaces.
- Do not gamma or ethylene oxide sterilize the capsules.
- Do not steam-in-place sterilize the capsules.

 **WARNING**

**To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and**

**To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:**

- Only expose the capsule to 1 sterilization or 1 sanitization cycle prior to use.
- Do not autoclave and base sanitize the same capsule prior to use.
- During pre-use autoclave sterilization do not exceed a temperature of 121°C or a duration longer than 30 minutes for numbers BV0.3R (EMP101AEX020R) and BV1R (EMP201AEX020R)
- Do not use with or expose this product to process liquids until capsule temperature goes <40°C (104°F).
- Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1-hour duration for sanitization.
- Use only 70% ethyl and/or isopropyl alcohol to wipe the exterior of the capsules.
- Do not use with or expose this product to pressurized steam or sanitizing chemicals other than NaOH.

 **WARNING**

**To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:**

- Use the installation procedure specified in Operation Instructions.
- Wear Personal Protective Equipment (PPE) to prevent chemical splash.
- The capsule must be subjected to a preconditioning flush prior to use. Use the procedure specified in Operation Instructions.
- After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline has been reached (e.g.,  $\pm 0.5$  pH units from the pH of the buffer used to flush the media).
- Capsules are designed for single use only. Do not reuse.
- Do not modify or alter the capsules.
- Do not put this product into service if any damage is occurred to the capsule or any of the connecting parts.
- Do not install where fluid pressure exceeds 3.4 bar (50 psig).
- Do not use in organic solvent based processing.
- Use only with aqueous based process fluids in the pH range of 5 to 9.
- Do not use with or expose this product to hot liquids  $>40^{\circ}\text{C}$  ( $104^{\circ}\text{F}$ ).
- Do not use the product on unlevelled surfaces.
- Do not use the product for more than 8 hours.

 **WARNING**

**To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:**

- Do not use this product for continuous service with compressed gasses. Use of compressed gas is permissible for post-use integrity testing and blow-down purposes only. Limit post-use integrity test and blow-down gas pressures to less than 1.4 bar (20 psig) and temperature less than 25°C for no longer than 30 minutes.
- Do not hang other items on this product.
- Do not gamma or ethylene oxide sterilize the capsules.
- Do not steam-in-place sterilize the capsule.
- Only expose the capsule to 1 sterilization or 1 sanitization cycle.
- Do not autoclave and base sanitize the same capsule prior to use.
- During pre-use autoclave sterilization do not exceed a temperature of 121°C or a duration longer than 30 minutes for numbers BVO.3R (EMP101AEX020R) and BV1R (EMP201AEX020R).
- Do not use with or expose this product to process liquids until capsule temperature goes <40°C (104°F).
- Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1-hour duration for sanitization.
- Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.
- Do not expose the capsule to base sanitization chemicals more than 60 minutes.
- Use only 70% ethyl and/or isopropyl alcohol to wipe the exterior of the capsules.
- Do not use with or expose this product to pressurized steam or sanitizing chemicals other than NaOH.


**! WARNING**

**To reduce the risks associated with chemical exposure due to pressure barrier failure, connector failure, purification media or capsule failure from use of disinfectants within the product:**

- Use up to 1 M Sodium Hydroxide (NaOH) for pre-use and post-use sanitization and up to 5% Bleach (NaClO) for post-use sanitization.
- Do not expose the purification capsule to base sanitization chemicals for more than 60 minutes.
- Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.
- Wear personal protective equipment (PPE) to prevent chemical splash.

**Figure 1 - Small Laboratory Capsules**

**! CAUTION**

**To reduce the risk of burn or exposure injuries associated with autoclave sterilization.**

- Wear personal protective equipment (PPE) to handle hot capsules.

**⚠ CAUTION**

**To reduce the risk of fluid exposure:**

- Evaluate the risk of your process liquids and determine an appropriate protocol to dispose of spent capsules and all other waste.
- Dispose of capsules and all other waste in accordance with federal, state, and local regulations or your local country's laws and regulations.

## 1. Operating Parameters and Product Specifications

Liquid	
Maximum Inlet Pressure	3.4 bar (50 psig)
Maximum Differential Pressure	2.4 bar (35 psig)
Maximum Operating Temperature	40°C (104°F)
Required Preconditioning Flush	
Preconditioning Flush Volume	15.3 mL flush solution per mL of bed volume
Recommended Flow Rate	1 mL/min per mL of bed volume

Product Specifications		
	BV0.3R	BV1R
	EMP101AEX020R	EMP201AEX020R
Dimensions		
Height x Diameter	4.8 cm x 4.3 cm (1 7/8" x 1 11/16")	5.9 cm x 4.3 cm (2 5/16" x 1 11/16")
Bed Volume	0.4 mL	1.4 mL
Filter Surface Area	1.0 cm <sup>2</sup>	4.0 cm <sup>2</sup>

<b>Product Specifications</b>		
	<b>BV0.3R</b>	<b>BV1R</b>
	<b>EMP101AEX020R</b>	<b>EMP201AEX020R</b>
<b>Approximate Weight</b>		
Dry	9.0 g	14.5 g
Wet Post Blow Down	9.5 g	16.0 g
<b>Material of Construction</b>		
Filter Media	Quaternary Ammonium Functional Polymer Grafted to a Polypropylene Nonwoven	
Membrane	Polyamide	
Membrane Support	Polypropylene Nonwoven	
Capsule Shell	Polypropylene	
Edge Seal	Polypropylene	
Luer Connectors	Polypropylene	
<b>Approximate Volume</b>		
Capsule Fill Volume <sup>2</sup>	2.0 mL	5.5 mL
Holdup Volume Post Blow Down <sup>3</sup>	0.5 mL	1.5 mL
<b>Pre-use Sterilization and Sanitization</b>		
Autoclave Sterilization (Only 1 autoclave cycle)	121°C for 30 minutes using Pre-vac cycle	
Base Sanitization (Only 1 base sanitization cycle)	1M Sodium Hydroxide (NaOH) for 60 minutes	

<b>Product Specifications</b>		
	<b>BVO.3R</b>	<b>BV1R</b>
	<b>EMP101AEX020R</b>	<b>EMP201AEX020R</b>
<b>Post-use Sterilization and Sanitization</b>		
Autoclave Sterilization	Post use 1 autoclave cycle at 121°C for 40 minutes	
Sanitization	1M Sodium Hydroxide (NaOH) or 5% Bleach (NaClO) for post-use sanitization	
<b>Required Preconditioning Flush Volume<sup>4</sup></b>		
Flush Volume, No Pre-treatment	6 mL	22 mL
Flush Volume, Post-autoclave	6 mL	22 mL
Required Preconditioning Flush Volume, Post Base Sanitization - 50 mM phosphate buffered saline at pH 7.4 <sup>5</sup>	6 mL	22 mL
<b>Miscellaneous</b>		
Connector	Luer	
Recommended Use Flow Rate	0.4 mL/min	1.4 mL/min
Maximum Flow Rate <sup>6</sup>	1.2 mL/min	5 mL/min
Storage Conditions <sup>7</sup>	Controlled indoor temperature: 0 - 30°C (32 - 86°F) in original packaging	
Shelf-Life <sup>8</sup>	2 years from date of manufacture @ 30°C maximum storage	

 **WARNING**

**<sup>1</sup>Do not use this product for continuous service with compressed gasses. Use of compressed gas is permissible for post-use integrity testing and blow-down purposes only. Limit post-use integrity test and blow-down gas pressures to less than 1.4 bar (20 psig) and temperature less than 25°C for no longer than 30 minutes.**

- <sup>2</sup> Capsule Fill Volume is defined as the volume of liquid that is required to fill the capsule.
- <sup>3</sup> Post Blow down Hold up Volume is defined as the volume of the residual liquid after gas (air or nitrogen) blow down.
- <sup>4</sup> The Required Preconditioning Flush is a saline solution, having a conductivity in the same range as that of 25 to 150 mM NaCl (typical 3-16 mS/cm). DO NOT use water to flush the capsule.
- <sup>5</sup> Review Regulatory Support File for recommended alternate buffer selection.
- <sup>6</sup> Do not allow the liquid differential pressure across the capsule to exceed 2.4 bar (35 psi). Maximum flow rate may be limited by system differential pressure.
- <sup>7</sup> Do not allow product to be exposed to freezing or extreme heat conditions.
- <sup>8</sup> Product labeling includes the expiration date.

### **Preconditioning Flush Procedure**

**The Laboratory Capsule must be flushed prior to use to remove the glycerin that is present as a media stabilizer.** Flush the capsule with 54 L/m<sup>2</sup> (15.3 mL/bed volume) of a saline solution or buffer having a conductivity in the same range as that of 25 to 150 mM NaCl (3 - 16 mS/cm) at a flow rate no higher than 210 L/m<sup>2</sup>/h (LMH). DO NOT use water as the flush solution. Flush in the forward flow direction only (see Figure 1). Intermittently bleed air off of the upstream side of the capsule using the Luer Capped Air Vent on the inlet side of the capsule.

 **WARNING**

To reduce the risk of Injury from possible end user exposure associated with degradation products released from the purification media; and

To reduce the risk of product loss associated with a) chemical release from purification media and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components.

- The capsule must be subjected to a preconditioning flush prior to use.

## 2. Capsule Use with No Pre-use Sterilization and Sanitization

### 2.A. Required Preconditioning Flush Procedure

1. Connect the capsule inlet and outlet with the appropriate tubing and process equipment.
2. Connect the capsule inlet to the reservoir containing the pre-conditioning flush solution.
3. Open the capsule vent.
4. Close the outlet valve and open the inlet valve slowly to begin the flow of pre-conditioning flush solution through the capsule inlet.
5. When fluid emerges from the vent, close the vent while simultaneously opening the outlet valve slowly.
6. Pump the appropriate amount of flush solution per capsule size to complete the required preconditioning flush.
7. Close the inlet valve and stop the flow of pre-conditioning flush solution.

**IMPORTANT NOTICE: Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.**

### 2.B. Single-stage Purification

1. Connect the inlet and outlet of the capsule to the upstream and downstream process equipment, respectively.

2. Connect the capsule inlet to the reservoir containing the bioprocess fluid.
3. Start the flow of bioprocess fluid and collect purified permeate.

**IMPORTANT NOTICE: For optimal performance, run the filtration process at a flux of 160 – 630 LMH (0.75-3.0 Bed Volume per minute).**

4. Stop the flow of bioprocess fluid at the end of the experiment.
5. The bioprocess fluid remaining on the upstream side of the capsule may then be processed through the system:
  - a. Apply 0.3-0.7 bar (5-10 psig) air or nitrogen gas pressure to vent. Using a non-surgng pump or clean air supply to deliver regulated low-pressure air to push the process fluid through the capsule. Do not allow the inlet pressure to pulse or exceed 1.4 bar (20 psig).
  - b. Monitor the capsule outlet. When liquid flow is reduced to a trickle, reduce the system pressure to 0 bar (0 psig).

### **3. Capsule Use with Pre-use Autoclave Sterilization**

#### **3.A. Recommended Pre-use Autoclave Sterilization of Capsules**

1. Wrap the capsule connections and vent with sterilization wrap.
2. Place the capsule in the autoclave.
3. Set-up the autoclave cycle.
  - a. Select pre-vac cycle for all capsule sizes.
  - b. Set the sterilization temperature at 121°C.
  - c. Set the autoclave duration for capsules for 30 minutes.
4. Start the autoclave cycle.
5. At the end of autoclave cycle, remove the capsules from the autoclave and let the capsule cool down at ambient temperature for at least 60 minutes.
6. Remove the Sterilization Wrap.

**⚠ CAUTION**

**To reduce the risk of burn or exposure injuries associated with autoclave sterilization:**

- Wear personal protective equipment (PPE) to handle hot capsules.

### **3.B. Required Preconditioning Flush Procedure for the Pre-use Autoclaved Capsules**

**⚠ WARNING**

**To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and**

**To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:**

**Do not use with or expose this product to process liquids until capsule temperature goes <40°C (104°F) for no longer than 30 mins.**

#### **3B.1. Required Preconditioning Flush Procedure**

Perform Required Preconditioning Flush procedure as described in Section 2.A.

**IMPORTANT NOTICE: Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.**

#### **3.C. Single-stage Purification**

Perform Single-stage Purification as described in Section 2.B.

## 4. Capsule Use with Pre-use NaOH Sanitization

**Recommended maximum conditions:** Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1 hour duration for sanitization.

### 4.A. Recommended Pre-use NaOH Sanitization

1. Connect the capsule inlet to the reservoir containing the 1 M NaOH solution.
2. Open the capsule vent.

#### WARNING

**To reduce the risks associated with chemical exposure due to pressure barrier failure, connector failure, purification media or capsule failure from use of disinfectants within the product: Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.**

3. When fluid emerges from the vent, close the vent while simultaneously opening the outlet valve slowly.
4. When fluid emerges from the capsule outlet, stop the flow and close the capsule outlet valve.
5. Allow the capsule to incubate with the 1 M NaOH solution. Incubation time should not exceed 1 hour.
6. After the incubation is complete, drain the NaOH solution from the capsule.

**IMPORTANT NOTICE: Do not blow down the 3M™ Emphaze™ AEX Hybrid Purifier capsules after NaOH sanitization procedure.**

### 4.B. Required Preconditioning Flush Procedure for the Pre-use NaOH Sanitized Capsules

**⚠ WARNING**

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range 5-9 has attained and monitor the pH of the capsule effluent to determine when the baseline has been reach (e.g.,  $\pm 0.5$  pH units from the pH of the buffer used to flush the media).

**4.B. Required Preconditioning Flush Procedure for the Pre-use NaOH Sanitized BVO.3R and BV1R Capsules**

Perform Required Preconditioning Flush procedure as described in Section 2.A.

**IMPORTANT NOTICE:** Monitor the pH of the capsule effluent to determine when the baseline has been reach (e.g.,  $\pm 0.5$  pH units from the pH of the buffer used to flush the media).

**IMPORTANT NOTICE:** Do not blow down with a gas the 3M™ Emphaze™ AEX Hybrid Purifier capsules after preconditioning flush procedure.

Capsule is ready for immediate use.

**IMPORTANT NOTICE:** Do not blow down 3M Emphaze AEX Hybrid Purifier capsules after pre-conditioning flush procedure.

**4.C. Single-stage Purification**

Perform Single-stage Purification as described in Section 2.B.

## 5. Disposal

1. After fluid processing is complete, drain the fluid from the capsule until the internal pressure on the capsule is reduced to 0 bar (0 psig).
2. Drain any remaining fluid from the system and remove the external connections from the inlet, outlet, and vent.

### CAUTION

#### To reduce the risk of fluid exposure:

- Evaluate the risk of your process liquids and determine an appropriate protocol to dispose of spent capsules and all other waste.
- Dispose of capsules and all other waste in accordance with federal, state, and local regulations or your local country's laws and regulations.

#### Intended Use:

Single-use processing of aqueous based biological pharmaceuticals (drugs) and vaccines to remove biological contamination strictly following the product operating instructions and cGMP requirements, where applicable.

#### Prohibited Use:

For laboratory use only. Not intended for use with materials that will be used on humans or animals. Do not use as a component in a medical device that is regulated by any agency, and/or globally exemplary agencies, including but not limited to: a) FDA, b) European Medical Device Directive (MDD), c) Japan Pharmaceuticals and Medical Devices Agency. Do not use in applications involving permanent implantation into the body, life-sustaining medical applications, applications requiring Food Contact compliance.

## **Product Selection and Use:**

Many factors beyond 3M's control and uniquely within user's knowledge and control can affect the use and performance of a 3M product in a particular application. As a result, end-user is solely responsible for evaluating the product and determining whether it is appropriate and suitable for end-user's application, including completing a risk assessment that considers the product leachable characteristics and its impact on drug safety, conducting a workplace hazard assessment and reviewing all applicable regulations and standards (e.g., OSHA, ANSI, etc.). Failure to properly evaluate, select, and use a 3M product and appropriate safety products, or to meet all applicable safety regulations, may result in injury, sickness, death, and/or harm to property.

## **Warranty, Limited Remedy, and Disclaimer:**

Unless a different warranty is specifically stated on the applicable 3M product packaging or product literature (in which case such warranty governs), 3M warrants that each 3M product meets the applicable 3M product specification at the time 3M ships the product. 3M MAKES NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR CONDITION OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR ARISING OUT OF A COURSE OF DEALING, CUSTOM, OR USAGE OF TRADE. If a 3M product does not conform to this warranty, then the sole and exclusive remedy is, at 3M's option, replacement of the 3M product or refund of the purchase price.

## **Limitation of Liability**

Except for the limited remedy stated above, and except to the extent prohibited by law, 3M will not be liable for any loss or damage arising from or related to the 3M product, whether direct, indirect, special, incidental, or consequential (including, but not limited to, lost profits or business opportunity), regardless of the legal or equitable theory asserted, including, but not limited to, warranty, contract, negligence, or strict liability.

**3M Purification Inc.**  
**3M Separation and**  
**Purification Sciences Division**  
400 Research Parkway  
Meriden, CT 06450  
U.S.A.  
(800) 243-6894  
(203) 237-5541  
[www.3Mpurification.com](http://www.3Mpurification.com)

3M, Emphaze, and Zeta Plus are trademarks of 3M Company. 3M, Emphaze und Zeta Plus sind Warenzeichen von 3M Company./3M, Emphaze y Zeta Plus son marcas registradas de 3M Company./3M, Emphaze ed Zeta Plus sono marchi di fabbrica di 3M Company./3M, Emphaze et Zeta Plus sont des marques déposées de 3M Company./3M, Emphaze en Zeta Plus zijn handelsmerken van 3M Company.

Please recycle. Printed in U.S.A./Bitte recyceln. Gedruckt in den Vereinigten Staaten./Recycle. Impreso en EE. UU./Si prega di riciclare. Stampato negli U.S.A./Pensez à recycler. Imprimé aux États-Unis./Gelieve deze folder te recyclen. Gedrukt in de VS.  
© 2019 3M Company. All rights reserved./Alle Rechte vorbehalten./Todos los derechos reservados./Tutti i diritti riservati./Tous droits réservés./  
Alle rechten voorbehouden.  
34-8724-3736-2

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**3M Purification Inc.**  
**3M Separation and Purification**  
**Sciences Division 400**  
Research Parkway  
Meriden, CT 06450  
VS  
(800) 243-6894  
(203) 237-5541  
Fax +1 203 630 4530  
[www.3Mpurification.com](http://www.3Mpurification.com)

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<b>Requester:</b> Justin Haas <b>Creator:</b> 25468 <b>File Name:</b> 34872437362.indd <b>Structure #:</b> SS-81078 <b>Date:</b> 07/24/19	<b>Printed Colors – Front:</b> 
	<b>Printed Colors – Back:</b> 
	<b>Match Colors:</b>
<b>Scale:</b>  1 Inch	<b>This artwork has been created as requested by 3M. 3M is responsible for the artwork AS APPROVED and assumes full responsibility for its correctness.</b>