

Saving treatment time in direct posterior procedures.

Tommaso Brunelli, Italy.

3MSM Health Care Academy

Building up a posterior tooth's occlusal anatomy, cusp-by-cusp, can be a laborious and time-consuming task. To simplify this process, it is now possible to replace the numerous steps required with a single procedure using high-viscosity bulk fill materials. This has been successfully demonstrated with 3M Filtek One Bulk Fill Restorative, which offers easy sculptability and a great quality surface. Using this material, the patient's natural occlusal surface structure was restored by following StyleItaliano™'s the essential lines technique. A desirable optical integration was achieved, as well as significant treatment time savings for the practitioner due to the simplicity of the protocol.

Only two specific instruments and one bulk fill composite are required for this technique. Important material properties are easy sculptability and non-stickiness, as they enable the dental practitioner to remove the exact amount of material needed when modelling the surface and integrating the fissures. Due to the simplicity of the protocol and a reduction of the number of treatment steps made possible by an elimination of the need for layering or a cusp-by-cusp modelling, significant treatment time savings are achieved. The following case example is used to demonstrate the simplicity of the procedure. The material of choice was 3M™ Filtek™ One Bulk Fill Restorative, which offers the desired sculpting properties and a great surface quality.



Figure 1: Initial situation with the maxillary first molar affected by caries.




Figure 2: Cavity after caries excavation and preparation with only a small area with exposed dentine.



Figure 5: The essential lines identified by StyleItaliano for maxillary first molars drawn into the composite surface using the Fissura instrument (LM Arte). The best results are obtained if the instrument is not pulled through the material, but repeatedly dipped into it.



Figure 3: Selective enamel etching with 3M™ Scotchbond™ Universal Etchant for 15 seconds. Subsequently, 3M™ Scotchbond™ Universal Adhesive was applied to the tooth surface, rubbed in for 20 seconds and gently air dried to remove the solvents.



Figure 6: Checking the occlusal contacts after light-curing of the restored surface according the protocol recommended by the material manufacturer.



Figure 4: Situation after the application of 3M™ Filtek™ One Bulk Fill Restorative into the cavity in a single layer (the maximum layer thickness is 5 mm). The natural surface anatomy was imitated using the Condensa instrument (LM Arte), which was placed on the surface and rolled along the remaining cusps, thereby following their natural inclination. As it is not yet cured, the composite still appears slightly too translucent.



Figure 7: Treatment result with a natural appearance. The opacity of the bulk fill material has increased after curing and a great optical integration is obtained.



Tommaso Brunelli

Tommaso Brunelli graduated in Dentistry from the University of Brescia in October 2016. Since then, he has been working in different private practices in Bergamo and Brescia, with a focus on restorative dentistry and endodontics. In June 2018, he became the first prize winner of the OneAndOnlyBulkContest powered by StyleItaliano and 3M.

Contact:

👤 Tommaso Brunelli
✉ tbrunelli91@gmail.com

Contact 3M

Australia

3M.com.au/dental

New Zealand

3M.co.nz/dental

Scientific Affairs

Stephen Langdon

Email: sdlangdon@3M.com

3M Health Care Academy



3M Oral Care

3M Australia Pty Limited
Building A, 1 Rivett Road
North Ryde NSW 2113
Ph: 1300 363 454
3M.com.au

3M New Zealand Limited
94 Apollo Drive
Rosedale Auckland 0632
Ph: 0800 80 81 82
3M.co.nz



Find us on Facebook

Disclaimer: 3M brand names mentioned in this clinical case may have been amended from the originally published version to reflect branding variances in Australia and New Zealand. Please always refer to the 3M Product Instructions For Use.

Published by 3M Oral Care. 3M and "3M Science. Applied to Life." are trademarks of 3M. 3M Health Care Academy is a service mark of 3M. Please recycle. © 3M 2019. All rights reserved.