

# Reducing the risk of Surgical Site Infections along the patient pathway




# Did you know?



## SSIs: a costly problem.

Surgical care is an essential component of health care. Yet each year surgical patients are harmed by complications including surgical site infections (SSIs). SSIs place a huge operational and financial burden on health care systems and providers. A reduction in SSIs means improved outcomes and a better experience for patients and their families.



**2–5%**  
of all surgical  
patients develop  
an SSI<sup>1</sup>



Result in  
**2x**  
increased risk  
of mortality<sup>2</sup>



Result in  
**2x**  
increased  
length of stay<sup>3</sup>



**5x**  
more likely to result  
in readmission<sup>1</sup>



Quality of life  
**affected**<sup>4</sup>



Cost at least  
**£3,000**  
per patient<sup>5</sup>

# Did you know?

## SSIs: a complex problem.

The problem of SSIs is a complex and significant one that cannot be solved by a single solution. There are three variables, as demonstrated by the Centres for Disease Control and Prevention (CDC) formula, that contribute to the risk of SSIs.<sup>38</sup>


Although the variables, virulence of bacteria and the patient risk factors are not easily controllable, the risk of SSIs can be minimised by evidence-based best practices related to the dose of bacteria.



**Risk  
of SSI**


=

**Dose  
of bacteria**



**X**

**Virulence  
of bacteria**



**Resistance of the host (patient)**

Bacterial contamination can occur from the:



**Patient**



**Clinical staff**



**Hospital environment**

Controlling these factors is essential in reducing the risk of SSIs.

**60% of SSIs are preventable**<sup>7,40,41</sup>

# Every touch point matters



**When reducing the risk of surgical site infection, every touch point matters.**



## Opportunities to reduce risk

There are many opportunities through the patient's surgical journey to make small changes that can make a big difference to clinical outcomes.



## Evidence and guidelines

All recommendations referenced in this module are supported by clinical evidence and guidelines.



## 3M Solution

We have a wide range of solutions to support clinicians in their mission to protect patients from preventable surgical site infections.




# Every touch point matters



[Learn more ►](#)



Patients should be assessed and provided with relevant information before their surgery.



Opportunity to inform and prepare the patient for successful surgery.

Evidence and guidelines

Guidelines recommend providing the patient with relevant information in advance of their surgery.<sup>8</sup>



All patients undergoing elective surgery should have a preoperative assessment to stratify risk of developing a surgical site infection or wound breakdown, as well as the risk the surgery will pose.<sup>9</sup>

These risks together, influence decision to ensure appropriate therapy or dressing is planned and explained to each patient.



3M Solution

Our clinical experts have prepared a patient information guide to help them take an active role in preparing for their surgery.





Science.  
Applied to Life.™

Take an active role in your surgical recovery.

This useful guide provides some easy to follow steps to improve your surgical recovery and reduce your risk of surgical site infection.



**What is a surgical site infection?**

Germs are all around us. In fact, our bodies are home to an estimated 100 trillion "friendly" germs, which are essential to our survival. Unfortunately, some germs can be harmful and can cause an infection. A surgical site infection occurs when these harmful germs enter the incision (cut) made through your skin for surgery.



**What are the symptoms of a surgical site infection?**

Symptoms of an infection can include redness, swelling and pain around the surgical area, fever and/or discharge (pus) from the surgical wound. This may occur any time from 1-3 days after your operation until the wound has healed (usually 2-3 weeks after your operation).



**How can I prevent getting a surgical site infection?**

Thankfully, surgical site infections are uncommon and most only affect the skin, but some infections can spread to deeper tissue if not treated. Doctors and nurses are trained to ensure that the risk of you getting an infection is low. However, there are things you can do as a patient to help.

Before your surgery

**1 Share your medical history**

Before your surgery, your doctor or nurse will speak with you to check your medical history. This information will help your surgical team to take the care to make sure you have a safe operation and good recovery. You should let your doctor or nurse know if you:

- are diabetic
- smoke
- take steroids
- have been in hospital for a long time before your operation
- have had a blood transfusion before your operation
- have a known history of resistant bacteria
- have been in a hospital abroad
- had an SSI in the past
- have a known allergy for (certain) antibiotics



**2 Stop smoking as soon as possible**

People who smoke are more prone to developing infections and other complications related to surgery. The earlier you stop smoking, the better your chances of a good recovery. Ideally you should stop smoking at least four weeks before your operation. If you are unable to stop smoking immediately, you should look to stop smoking at least 48 hours before your operation. You should also refrain from smoking for as long as possible after surgery, to help your overall recovery. Ask your medical team if you need help to stop smoking, if at least only temporarily.



**3 Make positive lifestyle choices**

Surgery places a lot of stress upon your body. After surgery your immune system will be working overtime to promote the healing of your surgical wound and fight off infection. To increase your immune systems effectiveness, you should consider physically preparing your body by:

- Eating a well-balanced diet by selecting foods that are rich in vitamins, minerals and proteins
- Exercising regularly to build energy and maintain strength
- Avoiding excessive alcohol. Even moderate amounts of alcohol can weaken the immune system and slow down recovery



**4 Don't remove hair**

It may be tempting to shave, wax, epilate or use hair removal creams around your surgical area prior to surgery, but don't! Removing hair can damage the skin prior to surgery, which can weaken the bodies defence against germs. If the incision site is in an area where you regularly remove hair, stop for at least five days prior to surgery. This will provide enough time for tiny cuts or skin injuries to heal.

It is possible that hair may need to be removed before surgery to allow the surgical team to see the surgical area. If this needs to happen a nurse will remove hair using special medical clippers designed for medical use.

This can be accessed [here](#).





An undetected sterilization process failure can put patients at risk. Frequent monitoring can help to reduce risk of SSI.



Opportunity to monitor and check the sterility of instrumentation.

### Evidence and guidelines

Surgical site infections have been attributed to incompletely processed instruments.<sup>10</sup>



Using fully effective sterilization practices in the CSSD can reduce the risk of surgical site infections and their potential impact on hospital and healthcare costs.


Monitoring using sterilization process indicators plays an important role in sterility assurance.

Frequency of monitoring is determined by national guidelines or standards. The World Health Organization recommends testing steam and hydrogen peroxide sterilizers a minimum of daily and testing every EO sterilization cycle.<sup>11</sup>

Due to process variability, international standards e.g. EN ISO 17665 (1-3):2006 require that every cycle be monitored to ensure sterility assurance is guaranteed.<sup>14</sup>

“Patients depend on sterile processing personnel to be monitoring the effectiveness of their sterilization processes and producing a safe product for patient use.” Havrilla (2005)<sup>12</sup>

### 3M Solution



Sterilization process monitoring should include every load and every pack as each cycle should be considered a unique event. 3M offers four opportunities to monitor the process:

- 1 Equipment monitoring** using 3M™ Comply™ Bowie-Dick Test Packs or 3M™ Electronic Test System (ETS).
- 2 Load monitoring** – Steam, vaporized hydrogen peroxide (VH<sub>2</sub>O<sub>2</sub>) and Ethylene Oxide (EO) sterilization cycles can all be monitored using 3M™ Attest™ Super Rapid Readout Biological Indicators.
- 3 Internal pack monitoring** – 3M™ Comply™ SteriGage™ Steam Chemical Indicators placed inside the pack can be inspected when the pack is opened in the OR.
- 4 External pack monitoring** – Packs should have an external process indicator, such as 3M™ Comply™ Indicator Tapes, for visual confirmation that the pack was exposed to the sterilization process.




Preoperative showering, hair removal and nasal decolonization are recommended as part of the patient preparation process.<sup>8,11,15,16</sup>

Opportunity to prepare the patient's skin for surgery.

Evidence and guidelines

Guidelines recommend using a clipper with a single use blade, if hair has to be removed, as close to the time of surgery as possible.<sup>6,8,11,15–19</sup>



Razors should not be used as they cause nicks on the skin which can create an entry point for microorganisms and increase risk of SSI.<sup>23,24</sup>

*S. aureus* is the leading cause of SSI.

**30%** of the population carry *S. aureus* in the nose.

**84%** of the *S. aureus* strains isolated from the nares were identical to those isolated from the surgical site.<sup>20–22</sup>

NICE and CDC recommends nasal decolonization alongside undergoing a chlorhexidine wash or wipe, for identified high risk patients, high risk procedures, and known *S. aureus* carriers.<sup>8,6</sup>

3M Solution




We have a 3M™ Surgical Clipper to meet your needs. Superior design and performance minimizes the potential for nicks or cuts combined with a long lasting battery, short recharge time, battery indicator light, and submersible ergonomically designed handle. In addition the Pivoting clipper has a unique flexible head to conform to patient's skin.

Not all clippers are the same

3M Surgical Clippers minimise nicks and cuts, which helps reduce risk of SSIs. A nick rate of 2–5% was observed in a clinical study.<sup>25</sup>

Per cent (%) of subjects experiencing ≥ 1 nick or cutmean (±SE)



\*Difference is statistically significant (P = 0.004)





Applying active warming to the periphery reduces the effects of redistribution temperature drop and therefore the risk of perioperative hypothermia.


First opportunity to prevent inadvertent perioperative hypothermia.

Evidence and guidelines

Guidelines recommend pre-warming before induction of anaesthesia.<sup>27–34</sup>



After induction of anaesthesia, the thermoregulatory threshold for vasoconstriction is immediately reduced. Vasodilation allows redistribution of the warmer core circulation, to the cooler periphery, which results in a fast decrease in core body temperature, that is difficult to reverse, especially in shorter procedures.<sup>35</sup>



Heat redistribution is responsible for approximately **80%** of core temperature drop, following the induction of anaesthesia.<sup>35</sup>

Pre-warming the patient before induction of anaesthesia increases the temperature of their periphery, reducing the gradient between the core and peripheral temperatures and therefore a significantly smaller drop in core temperature.<sup>26</sup>

3M Solution





3M™ Bair Hugger™ Warming Gown and Blanket Systems can be utilised preoperatively or during the induction of anaesthesia, to narrow the temperature gradient between the core and peripheral temperatures, minimising the impact of heat redistribution.



Did you know?

As little as 10–20 minutes of pre-induction warming with the 3M™ Bair Hugger™ Normothermia System banks heat in the patient's periphery, helping minimise reduction in core temperature caused by redistribution temperature drop.<sup>26</sup>





Continuously and accurately monitor patient core temperature before, during and after surgery to ensure the patient is normothermic.

### Evidence and guidelines

Guidelines recommend that patients should have their core temperature monitored throughout the perioperative period.<sup>27-32</sup>



NICE recommends the measurement of temperature at sites which are a:

- ▶ direct measure of core temperature or
- ▶ direct estimate of core temperature that has been shown in research studies to be accurate to within 0.5°C of direct measurement<sup>27</sup>

Due to inaccuracy, indirect estimates of temperature are not recommended for use in surgical patients.<sup>27</sup>

**Knowing your patient's core temperature can help you respond quickly to help reduce complications of inadvertent perioperative hypothermia, including 3x increased risk of SSI.<sup>35</sup>**



Opportunity to monitor core temperature continuously.

### 3M Solution

**The 3M™ Bair Hugger™ Core Temperature Monitoring System** is a non-invasive continuous core temperature monitoring system, that has been proven to perform as accurately as a pulmonary artery catheter and other commonly used invasive systems.<sup>36</sup>



Just one 3M™ Bair Hugger™ Temperature Monitoring Sensor can accurately monitor the patient's core temperature throughout surgery.

### Did you know?

NICE now recommends zero heat flux technology as an accurate measurement of core temperature and recently published a medical innovation briefing on the Bair Hugger Temperature Monitoring System.<sup>37</sup>







Maintain normothermia during the perioperative period using active forced-air warming and fluid warming to prevent negative consequences of inadvertent perioperative hypothermia.

### Evidence and guidelines

Guidelines recommend actively warming patients using forced air warming to maintain normothermia intraoperatively.<sup>27-34</sup>




Fluids over 500ml should also be warmed in order to prevent the drop in core temperature caused by the infusion of cold fluids.<sup>27-30,33</sup>

### Benefits of maintaining normothermia include:<sup>35</sup>

- ▶ Reduced risk of surgical site infection
- ▶ Increased patient comfort
- ▶ Reduced morbid cardiac events
- ▶ Faster recovery time, shorter hospital stay
- ▶ Reduction in use of blood products
- ▶ Costs savings

Patients should not be discharged from recovery until they are 36°C or above. If required, continue to warm using forced-air warming.<sup>27,28,32,34</sup>



### 3M Solution



The 3M™ Bair Hugger™ Normothermia System can be used to reduce further heat loss, maintain core temperature and treat hypothermia. The interconnecting channel design of Bair Hugger blankets and gowns reduces resistance to airflow; maximising temperature uniformity within the blanket. Heat is transferred to the patient through convection from uniform perforations created in the blanket, conduction and radiation.

The 3M™ Ranger™ Blood and Fluid Warming System uses warming plates that monitor the temperature of fluid through the disposable cassette, ensuring it is warmed and reducing its impact on mean body temperature, irrespective of the flow rate or the ambient temperature of the fluid.

### Did you know?

There are 170+ clinical studies that support the use of the Bair Hugger forced air warming system.<sup>48</sup>

Bair Hugger blankets demonstrate superior heat transfer and even distribution of warm air.<sup>13</sup>



Opportunity to maintain normothermia during surgery.





Prepare healthcare workers, instruments and surgical environment, to exclude microorganisms from the surgical field.

Opportunity to complete the surgical field ready for incision.

Evidence and guidelines

Guidelines recommend hand disinfection should be carried out by the entire OR team using an antiseptic hand scrub.<sup>8,11,15</sup>



Surgical masks should cover the mouth and nose completely, cover all facial hair, and should be replaced for every operation. When a surgical mask is removed, hands should be disinfected.<sup>15</sup>

The operating team should wear sterile gowns in the operating theatre during the operation<sup>8</sup> and sterile drapes should be used to create a physical barrier around the surgical site.<sup>15</sup>



Surgical drapes and gowns should meet or exceed the minimum performance criteria of EN13795.<sup>39</sup>

Wound edge protectors may offer protection against SSI in certain procedures.<sup>11,42</sup>

3M Solution



Ensure a sterile field with a selection of 3M™ Steri-Drape™ 9000 Surgical Drapes and 3M™ Steri-Drape™ Wound Edge Protectors.

Choosing the correct personal protective equipment (PPE) is vital to protect both clinicians and patients, we have a range of surgical masks and respirators for all situations.

OR staff can complete final checks on instrument sterility using a 3M™ Comply™ SteriGage™ Chemical Indicator that has been placed inside the pack by the sterilization department.

Products you can trust

All 3M surgical drapes and gowns either meet or exceed the minimum performance criteria of EN13795.





Disinfect skin and create a chemical and physical barrier to eliminate microorganisms that could lead to SSI.



Opportunity to reduce skin recolonization and eliminate microorganisms.

### Evidence and guidelines




Most SSIs are caused by the contamination of a surgical incision with microbes from the patient's own body during surgery. Infection caused by microorganisms from an outside source following surgery is less common.<sup>8</sup>

Whilst skin preps reduce microbes on the skin surface, bacteria in the deeper skin layers will remain. Over time these microbes can recolonise the skin surface.<sup>43</sup>

Disinfect skin using an alcohol-based antiseptic solution.<sup>1,3,4,5,6,7,9</sup>

Guidelines recommend if an incise drape is required, it should be iodophor impregnated.<sup>8,15</sup>

Do not use non-iodophor-impregnated incise drapes routinely for surgery as they may increase the risk of surgical site infection.<sup>8,15</sup>



### 3M Solution



Create a chemical and physical barrier around the incision site with a 3M™ Ioban™ 2 Antimicrobial Incise Drape to reduce skin recolonization and therefore the risk of SSI.<sup>44,45</sup>

Evidence demonstrates that iodine released from 3M™ Ioban™ 2 Antimicrobial Incise Drapes is able to penetrate these deeper skin layers at a concentration required for microbial death.<sup>43</sup>

### Not all incise drapes are created equal



3M™ Ioban™ 2 Incise Drapes are classified as Class III medical, because the iodine in the incise drape is a drug which works in the deeper layers of a patient's skin to reduce the risk of SSI.<sup>43-45</sup>

Class II iodophor incise drapes do not offer the same level of assurance. Their drug component has not been assessed for safety or efficacy and does not act on the human body.





Close, protect and cover the incision; choosing the most appropriate post-operative dressing or therapy.

Monitor incision after surgery and during recovery.



Opportunity for successful wound healing and monitoring.

### Evidence and guidelines



Poor incision healing can lead to surgical site complications (SSCs) such as surgical site infection (SSI), dehiscence, haematoma and seroma. Each of these complications can lead to delayed healing, increased length of stay, readmissions, re-operation and increased morbidity and mortality.<sup>37</sup>

Surgical incisions should be covered with an appropriate interactive dressing at the end of the operation.<sup>8</sup> Dressing should not be removed until at least 48 hours post-op, unless there is a reason to inspect the wound.<sup>6,8,15,16</sup>

NICE (2019) indicate that the first dressing can be removed at 48 hours post operatively. However, if the incision can be viewed the dressing can be left until stitch/staple removal or the first follow up appointment depending on the wear time of the dressing.<sup>8</sup>

WHO (2016) suggests the use of prophylactic negative pressure wound therapy in adult patients on primarily closed surgical incisions in high-risk wounds for the purpose of the prevention of SSI, while taking resources into account.<sup>11</sup>

### 3M Solution

#### 3M™ Tegaderm™ Absorbent Clear Acrylic Dressings offer

a completely transparent barrier – allowing active monitoring for signs of infection and early intervention without disturbing the incision. Wear time of 7–28 days gives patients extended protection beyond their initial hospital stay,<sup>46</sup> and waterproof barrier supports early showering without removing the dressing, allowing patients to return more quickly to daily activities.<sup>47</sup>




**PREVENA™ Incision Management System** is the first powered negative pressure therapy designed specifically for closed surgical incisions. Applying continuous negative pressure therapy for up to 7 days, PREVENA helps to promote healing by holding the incision edges together, drawing fluid and exudate out of the wound, reducing oedema and stimulating perfusion. It is recommended that the therapy should be considered as an alternative to standard of care in people at risk of developing surgical site complications.







Patients should be discharged with the relevant information to care for their surgical wound.



Opportunity to empower the patient to improve their quality of life.

Evidence and guidelines



NICE (2019) guidelines state that patients and carers should be given information and advice about how to recognise a surgical site infection and who to contact if they are concerned.<sup>8</sup>

Patients should be prepared about the care of the incision site post operatively.

Well prepared patients can actively contribute to their own after care and report early signs of complications to a healthcare professional. This can also offer assurance of what they should expect to see during the normal healing of a surgical incision.



3M Solution



Our patient information document provides a guide for what to expect when wearing a 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing. Patients can monitor their own incision following surgeon’s instructions and contact a healthcare professional with any concerns.



Science.  
Applied to Life.™



### Dressing Guide

Your surgeon placed a 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing over your incision.

Follow your surgeon's instructions for incision monitoring. Use this document as a guide for what to expect when wearing your dressing. Contact your surgeon with questions about your incision or dressing.

#### Frequently asked questions

##### What might I see when I wear the dressing?

Your dressing has a clear absorbent pad with a waterproof cover dressing. This allows you to see your incision. As you wear the dressing, the pad may get stained with blood and fluid.

**This is okay, and does not mean you need to remove your dressing.** In fact, keeping the dressing on provides a healing environment.<sup>1</sup>

##### Can I shower or bathe with the dressing?

**Yes. You may shower** with this dressing if it is firmly in place and the edges have not lifted or rolled. Do not soak in the bathtub or whirlpool while wearing the dressing.

Follow dressing-removal recommendations from your surgeon/provider:

#### What might I expect while wearing the dressing?

The dressing is designed to stay in place longer than you may have experienced with other dressings. The information below can help you understand what to expect while wearing the dressing.

If you are concerned about your incision or dressing, call your surgeon's office.



**Bubbles**

Bubbles are expected and represent the pad swelling to absorb and capture fluid.

**Keep the dressing on.**



**Fluid in pad or cracked appearance**

Fluid within the pad appears or there is a cracked appearance, but the edges of the dressing remain sealed to the skin. This is expected.

**Keep the dressing on.**



**Rolled edges**

Rolled edges appear, but the dressing and pad are still in place. This is okay and does not necessarily require reinforcement.

**Keep the dressing on.**



**Fluid present in the dressing border**

Fluid or staining is seen outside the pad area, but it is not actively leaking. This is expected for some types of incisions.

**Keep the dressing on.**



Medical Solutions Division  
3M Health Care  
2050 Conway Ave.  
St. Paul, MN 55144 USA  
Phone 1-800-228-3957  
Web 3M.com/Medical

3M data on file, EM-05-001534  
3M and Tegaderm are marks and/or registered marks of 3M. Unauthorized use prohibited. Please recycle. Printed in USA.  
© 3M 2020. All rights reserved. GL\_IN\_70-2011-7937-7

# Summary



## Opportunities to reduce risk

There are several opportunities along the patient pathway to make small changes that can make a big difference to patient outcomes.



## Evidence and guidelines

Evidence and guidelines support proactive measures to reduce the risk of preventable surgical site infections.



## 3M Solution

3M can support clinicians with evidence-based solutions and practical resources to deliver positive patient outcomes.

# References

- 1 Smyth ET *et al.* (2008) Four Country Healthcare Associated Infection Prevalence Survey 2006: Overview of the results. *Journal of Hospital Infection*; 69:230–48.
- 2 Kirkland *et al.* (1999) The impact of surgical-site infections in the 1990s: attributable mortality, excess length of hospitalization, and extra costs. *Infect Control Hosp Epidemiol*; 20(11): 725–730.
- 3 Coello R, *et al.* (2005) Adverse impact of surgical site infections in English hospitals *J. Hosp. Infect* 60: 93–103.
- 4 Whitehouse *et al.* (2002) The impact of surgical-site infections following orthopedic surgery at a community hospital and a university hospital: adverse quality of life, excess length of stay, and extra cost. *Infect Control*
- 5 Plowman R, Graves N, Griffin M *et al.* (1999) The socio-economic burden of hospital acquired infection. London: Public Health. Laboratory Service.
- 6 CDC (2017) Centers for Disease Control and Prevention, Guideline for the Prevention of Surgical Site Infection. Published August 2017.
- 7 Meeks DW, Lally KP, Carrick MM *et al.* Compliance with guidelines to prevent surgical site infections: As simple as 1-2-3? *Am J Surg* 2011; 201(1):76–83.
- 8 NICE (2019) Surgical site infections: prevention and treatment, Clinical guideline [NG125]. Published April 2019.
- 9 AHSN Network (2019) National wound care strategy programme 2018/19 <https://www.ahsnnetwork.com/about-academic-health-science-networks/national-programmes-priorities/national-wound-care-strategy-programme>: accessed 30 March 2020.
- 10 Dancer, S.J., Stewart, M., Coulombe, C., *et al.*, Surgical site infection linked to contaminated surgical instruments, *Journal of Hospital Infection*, 2012;81(4):231–238.
- 11 WHO (2018) Global guidelines for the prevention of surgical site infection, second edition. Geneva: World Health Organization. Published 2018.
- 12 Havrilla, G., Hicks, R., Larson, D., *et al.*, Troubleshooting steam sterilization process failures – A series of unfortunate events, *Managing Infection Control*, October 2005;82–95.
- 13 Brauer A, *et al.* What Determines the Efficacy of Forced-Air Warming? A Manikin Evaluation with Upper Body Blankets. *Anesth Analg* 2009;108:192–8.
- 14 EN ISO 17665 (1–3):2006 specifies requirements for the development, validation and routine control of a moist heat sterilization process for medical devices. Moist heat sterilization processes covered by ISO 17665-1:2006 include but are not limited to: saturated steam venting systems; saturated steam active air removal systems; air steam mixtures; water spray; water immersion. <https://www.iso.org/standard/43187.html>
- 15 RKI (2018) Prevention of postoperative wound infections: Commission recommendation for hospital hygiene and infection prevention (KRINKO) at Robert Koch Institute, Published April 2018.
- 16 WIP (2011) Prevention of post-operative wound infections, Work-group Infection Prevention. Published May 2006, Last updated May 2011.
- 17 AOS (2016) American College of Surgeons and Surgical Infection Society: Surgical Site Infection Guidelines. Updated 2016.
- 18 SF2H (2013) Gestion préopératoire du risque infection. Published October 2013.
- 19 NASIC (2018) National Association of Specialists on Infection Control (Russian Federation). Clinical Recommendations. Surgical Site Infections Prophylaxis. Published 2018.
- 20 Sievert DM, *et al.* Antimicrobial-Resistant Pathogens Associated with Healthcare-Associated Infections: Summary of Data Reported to the National Healthcare Safety Network at the Centers for Disease Control and Prevention, 2009–2010. *Infect Control Hosp Epidemiol* 2013; 34(1):1–14.
- 21 Kuehnert MJ, *et al.* Prevalence of Staphylococcus aureus nasal colonization in the United States, 2001–2002. *J of Inf Diseases* 2006;193:172–179. National Center for Infectious Diseases and National Center for Health Statistics, Centers for Disease Control and Prevention, Atlanta, Georgia.
- 22 Perl TM., *et al.* Intranasal mupirocin to prevent postoperative Staphylococcus aureus infections. *NEJM*. 2002;346(24):1871–77.
- 23 Alexander JW Fischer JE, *et al.* The influenced hair removal methods on wound infections archives of Surgery 1983; 118357–352.
- 24 Van Pelt-Koops, *et al.* Preventing surgical site infections: an update. Hospital Healthcare Europe, 2012.
- 25 3M data on file 05-013632.
- 26 Horn EP, Bein B, Bohm R, Steinfath M, Sahili N, Hocker J. The Effect of Short Time Periods of Pre-Operative Warming in the Prevention of Peri-operative Hypothermia. *Anaesth.* 2012. 67(6).
- 27 NICE (2016) Hypothermia prevention and management in adults having surgery, Clinical guideline [CG65] Published date: April 2008, Last updated: December 2016.
- 28 Di Marco P, Canneti A (2017)SIAARTI Clinical Best Practice: Perioperative Normothermia, Published March 2017.
- 29 Torossian A, Bräuer A, *et al.*, (2014) S3 German and Austrian Guideline: Preventing Inadvertent Perioperative Hypothermia, Published May 2014.
- 30 SFAR (2018) Prevention of inadvertent perioperative hypothermia in adults, Societe Francaise d’Anesthesie et de Reanimation, published September 2018.
- 31 SEDAR (2018) A guide to clinical practice of unintentional perioperative hypothermia. Rev Esp Anesthesiol Reanim. published 2018.
- 32 TARD The Turkish Anaesthesiology and Reanimation Society Guidelines for the prevention of inadvertent perioperative hypothermia.
- 33 CFKR Center For Kliniske Retningslinjer Danish PW.
- 34 PTPALIO Polish Society of Anaesthesia and Intensive Care Nurses – PW.
- 35 Sessler, DI. Perioperative Heat Balance. *Anesth.* 2000; 92: 578–596.
- 36 Eshraghi Y *et al.* An Evaluation of a Zero-Heat-Flux Cutaneous Thermometer in Cardiac Surgical Patients. *Anesth Analg.* 2014 Sep;119(3):543–9.
- 37 World Union of Wound Healing Societies [WUWHS], Closed surgical incision management: understanding the role of NPWT, 2016.
- 38 Mangram, A.J. *et al.* (2019) Guideline for Prevention of Surgical Site Infection. Centers for Disease Control and Prevention (CDC). Published April 1999.
- 39 EN 13795-1 Surgical clothing and drapes – Requirements and test methods – Part 1: Surgical drapes and gowns.
- 40 Umscheid CA, Mitchell MD, Doshi JA *et al.* Estimating the proportion of healthcare-associated infections that are reasonably preventable and the related mortality and costs. *Infect Control Hosp Epidemiology* 2011; 32(2):101–114.
- 41 Anderson, Deverick J., *et al.* Strategies to Prevent Surgical Site Infections in Acute Care Hospitals: 2014 Update. *Infection Control and Hospital Epidemiology.* 2014; 35(6): 605–627.
- 42 Mihaljevic, A.L. *et al.* (2014) Multi-center Double-Blinded Randomized Controlled Trial of Standard Abdominal Wound Edge Protection With Surgical Dressings Versus Coverage With a Sterile Circular Polyethylene Drape for Prevention of Surgical Site Infections A CHIR-Net Trial (BaFO; NCT01181206) *Annals of Surgery.* Volume 260, Number 5, November 2014.
- 43 Elliott *et al.* Antimicrobial activity and skin permeation of iodine present in an iodine-impregnated surgical incise drape. *J. Antimicrobial Chemotherapy.* 2015.
- 44 Yoshimura *et al.* Plastic iodophor drape during liver surgery operative use of the iodophor impregnated adhesive drape to prevent wound infection during high risk surgery. *World J. Surgery.* 2003; 27:685–688.
- 45 Bejko *et al.* Comparison of efficacy and cost of iodine impregnated drape vs. standard drape in cardiac surgery:Study in 5100 patients. *J Cardiovasc Trans. Res.* 2015;8:431–437.
- 46 3M data on file. EM-05-014692. EM-05-014725. May be left in place for up to 28 days depending on the condition on the wound/incision and the surrounding skin, or as indicated by clinical practice. Can be worn until it leaks, loses adhesion, needs to be removed for wound/incision inspection or other clinical needs.
- 47 3M data on file. EM-05-014684.
- 48 3M™ Bair Hugger™ System Research Compendium. Available at: <https://engage.3m.com/normothermiacompendiumuk>





# Disclaimer

Please note, the design or ‘flow’ of the pathway is not intended to be linear. The steps are not necessarily carried out in this order and many will take place at the same time. The steps are separated to demonstrate to evidence and guidelines associated with each step of the journey.

All recommendations are taken from clinical studies and/or published guidelines and are not the opinion of 3M.

**3M United Kingdom PLC**  
Charnwood Campus  
10 Bakewell Road  
Loughborough  
LE11 5RB  
01509 611 611  
[www.3M.co.uk/medical](http://www.3M.co.uk/medical)

**3M Ireland Limited**  
The Iveagh Building  
Carrickmines Park  
Carrickmines  
Dublin 18  
1 800 320 500  
[www.3M.ie](http://www.3M.ie)

3M, Bair Hugger, Ranger, Comply, Attest, SteriGage, Tegaderm, Ioban and Steri-drape are trademarks of 3M Company. PREVENA is a trademark of KCI, a 3M Company.  
© 3M 2020. All rights reserved. OMG49713.

