Description 说明

The following discussion is intended to help you differentiate standard versions and surgical versions of N95 particulate filtering facepiece respirators.

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NIOSH-Certified N95 Respirators NIOSH 认证的 N95 呼吸器

Particulate respirators are designed to help reduce the wearer's exposure to airborne particulate hazards. In the U.S., respirators are tested and certified by the U.S. National Institute of Occupational Safety and Health (NIOSH). NIOSH tests and certifies respirators based on their physical and performance characteristics, including filtration efficiency. For example, N95-rated filtering facepiece respirators have a filtration efficiency of at least 95% against non-oily particles when tested using the NIOSH criteria. The particles used to test the filtration are in a size range that is considered the most penetrating. Therefore, the test methods ensure that the filter media can filter particles of all sizes with at least 95% efficiency.

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Surgical masks, in contrast, are designed to be worn by healthcare professionals during surgery and nursing, to help prevent contamination of the surgical field and/or the patient by capturing liquid droplets that are expelled by the wearer. Surgical masks are cleared for use as medical devices by the U.S. Food and Drug Administration (FDA), or equivalent agencies outside the U.S. That clearance is based on data and proposed claims provided by the manufacturer to the FDA for review, which the FDA evaluates and then "clears" for those products that meet the requirements. Because surgical masks are meant for use during surgeries, a key performance requirement is fluid resistance – the ability of masks to resist penetration by high-pressure streams of liquid, such as those that might result from a human artery being punctured during surgery¹.

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Surgical N95 Respirators 医用 N95 防护口罩

Surgical respirators are both certified by NIOSH as an N95 respirator and also cleared by the FDA as a surgical mask. These products are frequently referred to as medical respirators, healthcare respirators, or surgical N95s.

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Comparing Standard N95s to Surgical N95s 对比普通 N95 和医用 N95 防护口罩

Putting this all together will help you differentiate between a standard NIOSH-certified N95 respirator and a surgical N95 respirator. While similar in appearance, the key difference is the fluid resistance and the resulting FDA clearance of surgical N95s. But when is that fluid resistance necessary?

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Many tasks performed by healthcare workers – such as patient intake and non-emergency patient evaluation – pose little to no risk of generating high-pressure streams of liquid and are not conducted in a sterile field. For workers performing such tasks, a primary potential hazard to consider is airborne viruses and bacteria, such as those generated by coughs and sneezes, which are effectively filtered by an N95 respirator.

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
1. ASTM F1862 is a standard test method for resistance of medical facemasks to penetration by synthetic blood. This test is required because during certain medical procedures, a blood vessel may occasionally be punctured, resulting in a high-velocity stream of blood impacting a protective medical facemask. The test procedure specifies that a mask or respirator is conditioned in a high-humidity environment to simulate human use and is placed on a test holder. Synthetic blood (2cc) is shot horizontally at the mask at a distance of 30 cm (12 inches). Surgical masks and respirators are tested on a pass/fail basis at three velocities corresponding to the range of human blood pressure (80, 120, and 160 mmHg). The inside of the mask is then inspected to see if any synthetic blood has penetrated to the inside of the facemask. Fluid resistance according to this test method is when the device passes at any level.

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Therefore, if a healthcare facility is prioritizing respirator use – due to, for example, limited supply during a health emergency – they may want to consider prioritizing use of surgical N95 respirators for those healthcare workers requiring respiratory protection while performing surgery or other tasks that may expose them to high pressure streams of bodily fluid or conducting work in a sterile field. For other workers who will not be performing such surgical procedures or do not need to maintain a sterile field, a standard non-surgical N95 (or equivalent) respirator can be worn to help reduce those workers' exposure to patient-generated airborne viruses and bacteria.

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The following chart demonstrates some key similarities and differences between three respirator

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a. **ASTM F1862** is a standard test method for resistance of medical facemasks to penetration by synthetic blood. This test is required because during certain medical procedures, a blood vessel may occasionally be punctured, resulting in a high-velocity stream of blood impacting a protective medical facemask. The test procedure specifies that a mask or respirator is conditioned in a high-humidity environment to simulate human use and is placed on a test holder. Synthetic blood (2cc) is shot horizontally at the mask at a distance of 30 cm (12 inches). Surgical masks and respirators are tested on a pass/fail basis at three velocities corresponding to the range of human blood pressure (80, 120, and 160 mmHg). The inside of the mask is then inspected to see if any synthetic blood has penetrated to the inside of the facemask. Fluid resistance according to this test method is when the device passes at any level.

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For a list of 3M medical facemasks and surgical respirators, see this [3M Masks and Respirators Brochure](#).

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Personal Safety Division
3M Center, Building 235-2W-70
St. Paul, MN 55144-1000

3M PSD products are
occupational use only.

In United States of America
Technical Service: 1-800-243-4630
Customer Service: 1-800-328-1667
3M.com/workersafety
In Canada
Technical Service: 1-800-267-4414
Customer Service: 1-800-364-3577
3M.ca/Safety

© 3M 2020. All rights reserved.
3M is a trademark of 3M Company and its affiliates.
Used under license in Canada. All other trademarks
are property of their respective owners.
Please recycle.

3M