

STEP-BY-STEP



3M[™] RelyX[™] Universal



System Components

3M[™] Scotchbond[™] Universal Plus Adhesive



- Self-etch, selective-etch and total-etch adhesive
- MDP primer for zirconia and metal
- Silane primer for glass ceramics, composites, fiber reinforced posts

3M[™] RelyX[™] Universal Resin Cement

- Dual-cure two paste resin cement
- Adhesive and self-adhesive indications



 $3M^{\mbox{\tiny M}}$ RelyX $^{\mbox{\tiny M}}$ Universal Micro Mixing Tips with Elongation Tips

|--|









Self-sealing automix syringe



To ensure proper function, the valve must be kept in correct closed starting position.





Full rotation until stop.



Apply moderate centric force without bending the plunger.



Do not pull back the plunger of the syringe. This might clog the syringe permanently.



Discard a small amount of cement.







Store syringe without mixing tip.



Don't insert any instrument to rotate or disassemble the valve. Don't rotate the valve during cleaning.

Restoration Pre-Treatment Overview

Fiber-reinforced Composite Post

Priming with adhesive not required for 3M[™] RelyX[™] Fiber Posts

Zirconia, Alumina, Metal

Priming with adhesive optional, recommended for veneer, table top, adhesive bridge

Composite/Hybrids

Alternative primer: Apply a Silane,

e.g. 3M[™] RelyX[™] Ceramic Primer

Glass-Ceramic

20 sec. Alcoho Al₂O₃ 30-50µm Alcoho 20 sec. 5 Al₂O₃ 30-50µm 20 sec. 5 sec. HE 20 sec.

See section



Follow the instructions for use of the restoration material for specific instructions.

Tooth Pre-Treatment Options

3M[™] RelyX[™] Universal Resin Cement



Restoration **Pre-Treatment Recommendation**



Sec

Priming optional, recommended for veneer, table top, adhesive bridge.







Follow the instructions for use of the restoration material for specific instructions.

Self-Adhesive Post Cementation



 $3M^{\mathbb{M}}$ RelyX $^{\mathbb{M}}$ Fiber Posts only. For pretreatment of other posts see section 1.









To unlock syringe fully rotate mixer until stop.



Do not use high pressure. Check mixer position (image 5) if extrusion force is high.

-25



No lentulo!













Self-Adhesive Cementation

Pre-treat restoration according to section 1









To unlock syringe fully rotate mixer until stop.



Do not use high pressure. Check mixer position (image 5) if extrusion force is high.



8







Tack cure approx. 2-3 sec per surface (>1200 mW/cm²).





Light cure 10 sec per surface (>1200 mW/cm²) or wait 6 minutes.





(Selective-Etch) Adhesive Cementation

Pre-treat restoration according to **section 1**



Selective-etch of uncut enamel. Optional for cut enamel.











To unlock syringe fully rotate mixer until stop.

8



Do not use high pressure. Check mixer position (image 7) if extrusion force is high.

9







Tack cure approx. 2-3 sec per surface (>1200 mW/cm²).



Light cure 10 sec per surface (>1200 mW/cm²) or wait 6 minutes.





13

Total-Etch Adhesive Cementation

Pre-treat restoration according to **section 1**



Total-etch.











To unlock syringe fully rotate mixer until stop.

7

8



Do not use high pressure. Check mixer position (image 7) if extrusion force is high.

9







Tack cure approx. 2-3 sec per surface (>1200 mW/cm²).



Light cure 10 sec per surface (>1200 mW/cm²) or wait 6 minutes.





Please refer to product instructions for use for complete information.



3M Deutschland GmbH Health Care Business Carl-Schurz-Str. 1 41453 Neuss - Germany 3M, RelyX and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 2021, 3M. All rights reserved.

3M, RelyX et Scotchbond sont des marques de commerce de 3M ou 3M Deutschland GmbH. Utilisées sous licence au Canada. © 2021, 3M. Tous droits réservés.

