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Science. Applied to Life.™

Dressing Guide

Your surgeon placed a 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing over your incision.

Follow your surgeon's instructions for incision monitoring. Use this document as a guide for what to expect when wearing your dressing. Contact your surgeon with questions about your incision or dressing.

Frequently asked questions

What might I see when I wear the dressing?

Your dressing has a clear absorbent pad with a waterproof cover dressing. This allows you to see your incision. As you wear the dressing, the pad may get stained with blood and fluid.

This is okay, and does not mean you need to remove your dressing. In fact, keeping the dressing on provides a healing environment.¹

Can I shower or bathe with the dressing?

Yes. You may shower with this dressing if it is firmly in place and the edges have not lifted or rolled. Do not soak in the bathtub or whirlpool while wearing the dressing.

Follow dressing-removal recommendations from your surgeon/provider:

What might I expect while wearing the dressing?

The dressing is designed to stay in place longer than you may have experienced with other dressings. The information below can help you understand what to expect while wearing the dressing.

If you are concerned about your incision or dressing, call your surgeon's office.

Bubbles

Bubbles are expected and represent the pad swelling to absorb and capture fluid.

Keep the dressing on.

Fluid in pad or cracked appearance

Fluid within the pad appears or there is a cracked appearance, but the edges of the dressing remain sealed to the skin. This is expected.

Keep the dressing on.

Rolled edges

Rolled edges appear, but the dressing and pad are still in place. This is okay and does not necessarily require reinforcement.

Keep the dressing on.

Fluid present in the dressing border

Fluid or staining is seen outside the pad area, but it is not actively leaking. This is expected for some types of incisions.

Keep the dressing on.

Medical Solutions Division
3M Health Care
2050 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-228-3957
Web 3M.com/Medical

3M data on file, EM-05-001534
3M and Tegaderm are marks and/or registered marks of 3M. Unauthorized use prohibited. Please recycle. Printed in USA.
© 3M 2020. All rights reserved. GL_IN_70-2011-7937-7

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Viitteet

- 1 Smyth ET *et al.* (2008) Four Country Healthcare Associated Infection Prevalence Survey 2006: Overview of the results. *Journal of Hospital Infection*; 69:230–48.
- 2 Kirkland *et al.* (1999) The impact of surgical-site infections in the 1990s: attributable mortality, excess length of hospitalization, and extra costs. *Infect Control Hosp Epidemiol*; 20(11): 725–730.
- 3 Coello R, *et al.* (2005) Adverse impact of surgical site infections in English hospitals *J. Hosp. Infect* 60: 93–103.
- 4 Whitehouse *et al.* (2002) The impact of surgical-site infections following orthopedic surgery at a community hospital and a university hospital: adverse quality of life, excess length of stay, and extra cost. *Infect Control*
- 5 Plowman R, Graves N, Griffin M *et al.* (1999) The socio-economic burden of hospital acquired infection. London: Public Health. Laboratory Service.
- 6 CDC (2017) Centers for Disease Control and Prevention, Guideline for the Prevention of Surgical Site Infection. Julkaistu elokuussa 2017.
- 7 Meeks DW, Lally KP, Carrick MM *et al.* Compliance with guidelines to prevent surgical site infections: As simple as 1-2-3? *Am J Surg* 2011; 201(1):76–83.
- 8 NICE (2019) Surgical site infections: prevention and treatment, Clinical guideline [NG125]. Julkaistu huhtikuussa 2019.
- 9 AHSN Network (2019) National wound care strategy programme 2018/19 <https://www.ahsnnetwork.com/about-academic-health-science-networks/national-programmes-priorities/national-wound-care-strategy-programme>: noudettu 30.3.2020.
- 10 Dancer, S.J., Stewart, M., Coulombe, C., *et al.*, Surgical site infection linked to contaminated surgical instruments, *Journal of Hospital Infection*, 2012;81(4):231–238.
- 11 WHO (2018) Global guidelines for the prevention of surgical site infection, second edition. Geneve: Maailman terveysjärjestö. Julkaistu 2018.
- 12 Havrilla, G., Hicks, R., Larson, D., *et al.*, Troubleshooting steam sterilization process failures – A series of unfortunate events, *Managing Infection Control*, October 2005;82–95.
- 13 Brauer A, *et al.* What Determines the Efficacy of Forced-Air Warming? A Manikin Evaluation with Upper Body Blankets. *Anesth Analg* 2009;108:192–8.
- 14 EN ISO 17665 (1–3):2006 specifies requirements for the development, validation and routine control of a moist heat sterilization process for medical devices. Moist heat sterilization processes covered by ISO 17665-1:2006 include but are not limited to: saturated steam venting systems; saturated steam active air removal systems; air steam mixtures; water spray; water immersion. <https://www.iso.org/standard/43187.html>
- 15 RKI (2018) Prevention of postoperative wound infections: Commission recommendation for hospital hygiene and infection prevention (KRINKO) at Robert Koch Institute, julkaistu huhtikuussa 2018.
- 16 WIP (2011) Prevention of post-operative wound infections, Work-group Infection Prevention. Julkaistu toukokuussa 2006, viimeksi päivitetty toukokuussa 2011.
- 17 AOS (2016) American College of Surgeons and Surgical Infection Society: Surgical Site Infection Guidelines. Päivitetty 2016.
- 18 SF2H (2013) Gestion préopératoire du risque infection. Julkaistu lokakuussa 2013.
- 19 NASIC (2018) National Association of Specialists on Infection Control (Russian Federation). Clinical Recommendations. Surgical Site Infections Prophylaxis. Julkaistu 2018.
- 20 Sievert DM, *et al.* Antimicrobial-Resistant Pathogens Associated with Healthcare-Associated Infections: Summary of Data Reported to the National Healthcare Safety Network at the Centers for Disease Control and Prevention, 2009–2010. *Infect Control Hosp Epidemiol* 2013; 34(1):1–14.
- 21 Kuehnert MJ, *et al.* Prevalence of Staphylococcus aureus nasal colonization in the United States, 2001–2002. *J of Inf Diseases* 2006;193:172–179. National Center for Infectious Diseases and National Center for Health Statistics, Centers for Disease Control and Prevention, Atlanta, Georgia.
- 22 Perl TM., *et al.* Intranasal mupirocin to prevent postoperative Staphylococcus aureus infections. *NEJM*. 2002;346(24):1871–77.
- 23 Alexander JW Fischer JE, *et al.* The influenced hair removal methods on wound infections archives of Surgery 1983; 118357–352.
- 24 Van Pelt-Koops, *et al.* Preventing surgical site infections: an update. Hospital Healthcare Europe, 2012.
- 25 3M:n asiakirja-aineisto 05-013632.
- 26 Horn EP, Bein B, Bohm R, Steinfath M, Sahili N, Hocker J. The Effect of Short Time Periods of Pre-Operative Warming in the Prevention of Peri-operative Hypothermia. *Anaesth*. 2012. 67(6).
- 27 NICE (2016) Hypothermia prevention and management in adults having surgery, Clinical guideline [CG65] Julkaisujankohhta: huhtikuu 2008, viimeksi päivitetty: joulukuu 2016.
- 28 Di Marco P, Canneti A (2017) SIAARTI Clinical Best Practice: Perioperative Normothermia, julkaistu maaliskuussa 2017.
- 29 Torossian A, Bräuer A, *et al.*, (2014) S3 German and Austrian Guideline: Preventing Inadvertent Perioperative Hypothermia, julkaistu toukokuussa 2014.
- 30 SFAR (2018) Prevention of inadvertent perioperative hypothermia in adults, Societe Francaise d’Anesthesie et de Reanimation, julkaistu syyskuussa 2018.
- 31 SEDAR (2018) A guide to clinical practice of unintentional perioperative hypothermia. Rev Esp Anesthesiol Reanim. Julkaistu 2018.
- 32 TARD The Turkish Anaesthesiology and Reanimation Society Guidelines for the prevention of inadvertent perioperative hypothermia.
- 33 CFKR Center For Kliniske Retningslinjer Danish PW.
- 34 PTPALIO Polish Society of Anaesthesia and Intensive Care Nurses – PW.
- 35 Sessler, DI. Perioperative Heat Balance. *Anesth*. 2000; 92: 578–596.
- 36 Eshraghi Y *et al.* An Evaluation of a Zero-Heat-Flux Cutaneous Thermometer in Cardiac Surgical Patients. *Anesth Analg*. 2014 Sep;119(3):543–9.
- 37 World Union of Wound Healing Societies [WUWHS], Closed surgical incision management: understanding the role of NPWT, 2016.
- 38 Mangram, A.J. *et al.* (2019) Guideline for Prevention of Surgical Site Infection. Centers for Disease Control and Prevention (CDC). Julkaistu huhtikuussa 1999.
- 39 EN 13795-1 Surgical clothing and drapes – Requirements and test methods – Part 1: Surgical drapes and gowns.
- 40 Umscheid CA, Mitchell MD, Doshi JA *et al.* Estimating the proportion of healthcare-associated infections that are reasonably preventable and the related mortality and costs. *Infect Control Hosp Epidemiology* 2011; 32(2):101–114.
- 41 Anderson, Deverick J., *et al.* Strategies to Prevent Surgical Site Infections in Acute Care Hospitals: 2014 Update. *Infection Control and Hospital Epidemiology*. 2014; 35(6): 605–627.
- 42 Mihaljevic, A.L. *et al.* (2014) Multi-center Double-Blinded Randomized Controlled Trial of Standard Abdominal Wound Edge Protection With Surgical Dressings Versus Coverage With a Sterile Circular Polyethylene Drape for Prevention of Surgical Site Infections A CHIR-Net Trial (BaFO; NCT01181206) *Annals of Surgery*. Volume 260, Number 5, November 2014.
- 43 Elliott *et al.* Antimicrobial activity and skin permeation of iodine present in an iodine-impregnated surgical incise drape. *J. Antimicrobial Chemotherapy*. 2015.
- 44 Yoshimura *et al.* Plastic iodophor drape during liver surgery operative use of the iodophor impregnated adhesive drape to prevent wound infection during high risk surgery. *World J. Surgery*. 2003; 27:685–688.
- 45 Bejko *et al.* Comparison of efficacy and cost of iodine impregnated drape vs. standard drape in cardiac surgery:Study in 5100 patients. *J Cardiovasc Trans. Res*. 2015;8:431–437.
- 46 3M:n asiakirja-aineisto. EM-05-014692. EM-05-014725. Voidaan jättää paikalleen 28 päiväksi haavan ja ympäröivän ihon kunnosta riippuen; noudatetaan kliinistä käytäntöä. Sidosta voi käyttää, kunnes se alkaa vuotaa, liima menettää tehonsa tai side on irrotettava haavan tarkastamiseksi tai muun hoitotarpeen vuoksi.
- 47 3M:n asiakirja-aineisto. EM-05-014684.
- 48 3M™ Bair Hugger™ System Research Compendium. Saatavilla osoitteessa: <https://engage.3m.com/normothermiacompendiumuk>

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Suomen 3M Oy