


A man with dark hair is lying in a hospital bed, eyes closed, resting his head on a white pillow. A white, circular, non-invasive forehead sensor with blue patterns is attached to his forehead. A blue cable extends from the sensor across the pillow. The background is a blurred hospital room.

**3M** Science.  
Applied to Life.™

**Consistent. Accurate.  
Non-invasive.**

**The future of temperature monitoring.**



# Core body temperature is a critical vital sign that should be monitored.

Proactively monitoring temperature with a consistent, accurate, and non-invasive system is the only true way to own the normothermic temperature zone (36.0°C – 37.5°C)<sup>1</sup> and protect patients from unintended perioperative hypothermia.



## The importance of measuring and monitoring core temperature.

Inadvertent perioperative hypothermia is defined as a core body temperature of less than 36.0°C, which can:



Increase the rate of SSIs<sup>2</sup>



Increase surgical blood loss<sup>3</sup>



Lead to increased mortality<sup>4</sup>



Extend recovery time<sup>5</sup>




Cause patient discomfort<sup>6</sup>

Inadequate monitoring of core temperature can increase the risk of death associated with malignant hyperthermia (MH).<sup>7,8</sup>




When temperature was not monitored, 30% of patients with an MH episode died.<sup>7,8</sup>



If core temperature monitoring was used, the risk of death was reduced to 2% (P = 0.0012).<sup>7,8</sup>

## Temperature monitoring trade-offs.

Most current technologies are unable to non-invasively and accurately measure core body temperature. The 3M™ Bair Hugger™ Temperature Monitoring System can do both, allowing you to improve active warming practices. You can only manage what you can accurately measure.



3M™ Bair Hugger™ Temperature Monitoring System

### Global guidelines

AORN<sup>9</sup> recognises zero-heat-flux as a core temperature modality

ASPAN<sup>10</sup> recommends the use of one modality throughout the perioperative journey

NICE<sup>11</sup> recommends continuous temperature monitoring

DGAJ<sup>12</sup> recommends continuous core temperature monitoring with a consistent method


Teunissen LPJ. Measurement and manipulation of body temperature in rest and exercise. Amsterdam, Vrije Universiteit - PhD Thesis; 2012.

## The future of non-invasive temperature monitoring is here.

The Bair Hugger temperature monitoring system is the first accurate, non-invasive solution that allows you to easily monitor core body temperature consistently throughout the perioperative journey.

Clinically equivalent to invasive methods<sup>13</sup>


Displays current and 2-hour temperature trends



Eliminates the potential of trauma from invasive modalities<sup>14,15</sup>


One modality eliminates variability

## Technology that's more than skin deep.



Above image shows the inverse direction of how the product appears when in-use.


Zero-heat-flux technology produces near-perfect insulation, so heat cannot escape from the skin surface, creating an isothermal pathway.



Directly beneath the sensor, an isothermal pathway forms, allowing the deep tissue temperature to rise to the skin surface.

## Monitor with confidence.

The Bair Hugger temperature monitoring system makes it possible to get accurate readings throughout the entire perioperative journey with a single sensor.



The sensor remains on the patient throughout the perioperative journey but may be disconnected from and reconnected to the control unit as needed.

## One sensor, one proven solution.

Now that there is one accurate, consistent temperature monitoring solution, you have the power to Own the Zone.

Learn more about the Bair Hugger temperature monitoring system at [3M.co.uk/medical](http://3M.co.uk/medical) or contact your 3M representative.



Fight SSIs from every angle.

## References

- 1 Schroeck H, Lyden AK, Benedict WL, Ramachandran SK. Time Trends and Predictors of Abnormal Postoperative Body Temperature in Infants Transported to the Intensive Care Unit. *Anesthesiology Research and Practice*. 2016;7318137.
- 2 Kurz A, Sessler DI, Lenhardt R. Perioperative normothermia to reduce the incidence of surgical-wound infection and shorten hospitalization. *NEJM*. 1996 May 9;334(19):1209-16. Melling AC, Ali B, Scott EM, Leaper DJ. Effects of preoperative warming on the incidence of wound infection after clean surgery: a randomised controlled trial. *The Lancet*. 2001 Sep 15;358(9285):876-80. Scott AV, Stonemetz JL, Wasey JO, Johnson DJ, Rivers RJ, Koch CG, et al. (2015) Compliance with Surgical Care Improvement Project for Body Temperature Management (SCIP Inf-10) Is Associated with Improved Clinical Outcomes. *Anesthesiology* 123: 116–125.
- 3 Schmied H, Reiter A, Kurz A, Sessler DI, Kozek S. Mild hypothermia increases blood loss and transfusion requirements during total hip arthroplasty. *The Lancet*. 1996 Feb 3;347(8997):289-92.
- 4 Bush H Jr., Hydo J, Fischer E, et al. Hypothermia during elective abdominal aortic aneurysm repair: The high price of avoidable morbidity. *J Vasc Surg*. 1995;21(3): 392-402.
- 5 Lenhardt R, Marker E, Goll V, et al. Mild intraoperative hypothermia prolongs postanesthetic recovery. *Anesth*. 1997;87:1318-23.
- 6 Sessler, DI. Current Concepts: Mild perioperative hypothermia. *NEJM*, 1997. Vol. 336, No. 24, pp. 1730-1737.
- 7 Larach MG, Brandom BW, Allen GC, Gronert GA, Lehman EB. Malignant Hyperthermia Deaths Related to Inadequate Temperature Monitoring, 2007–2012: A Report from The North American Malignant Hyperthermia Registry of the Malignant Hyperthermia Association of the United States. *Anesthesia & Analgesia*. 2014;119(6):1359-1366.
- 8 Shafer SL, Dexter F, Brull SJ. Deadly Heat: Economics of Continuous Temperature Monitoring During General Anesthesia. *Anesthesia & Analgesia*. 2014;119(6):1235-1237.
- 9 Association of periOperative Registered Nurses.
- 10 American Society of PeriAnesthesia Nurses.
- 11 National Institute for Health and Care Excellence.
- 12 German Society of Anaesthesiology and Intensive Care Medicine.
- 13 <0.28°C error between 3M™ Bair Hugger™ Temperature Monitor and a nasopharyngeal probe. <0.28°C error between Bair Hugger temperature monitor and a pulmonary artery catheter.
- 14 Insler SR, Sessler DI. Perioperative thermoregulation and temperature monitoring. 2006. *Anesthesiology Clinics*, (24): 823-837.
- 15 Cereda M, & Maccioli GA. Intraoperative temperature monitoring. 2004. *International Anesthesiology Clinics*, (42)2: 41-54.

### 3M United Kingdom PLC

Charnwood Campus  
10 Bakewell Road  
Loughborough  
LE11 5RB

+44 (0)1509 611 611

### 3M Ireland Limited

The Iveagh Building  
The Park, Carrickmines  
D18 X015  
Ireland

+353 (0)1 280 3555

[www.3M.co.uk/medical](http://www.3M.co.uk/medical)

J431686

3M and Bair Hugger are trademarks of the 3M company.  
© 3M 2018. All rights reserved.

