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ICU Patient Rooms	ICU Common Areas
Bedside supply cart handle	Cardiac monitor control panel
IV pump monitor	Telephone
Bed rail	Workstation keyboard
Ventilator control panel	Medicine dispensing station
Cardiac monitor control panel	Crash Carts

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- [Clean trace Environment surface sample guide](#)
- [Clean trace Environment implementation guide](#)
- [Clean trace Surgical instrument implementation guide](#)
- [Clean trace surgical instrument sample guide](#)
- [Improving the quality of Flexible GI Endoscope reprocessing](#)

References:

- 1 World Health Organization. Home care for patients with suspected novel coronavirus (nCoV) infection presenting with mild symptoms and management of contacts (Interim guidance, 17 March 2020 Publication). [https://www.who.int/publications-detail/home-care-for-patients-with-suspected-novel-coronavirus-\(ncov\)-infection-presenting-with-mild-symptoms-and-management-of-contacts](https://www.who.int/publications-detail/home-care-for-patients-with-suspected-novel-coronavirus-(ncov)-infection-presenting-with-mild-symptoms-and-management-of-contacts)
- 2 Centers for Disease Control and Prevention. Interim guidance for persons who may have 2019 Novel Coronavirus (2019-nCoV) to prevent spread in homes and residential communities. https://www.cdc.gov/coronavirus/2019-ncov/hcp/guidance-prevent-spread.html#First_heading
- 3 Centers for Disease Control and Prevention. Interim Infection Prevention and Control Recommendations for Patients with Confirmed 2019 Novel Coronavirus (2019-nCoV) or Patients Under Investigation for 2019-nCoV in Healthcare Settings. February 3, 2020. <https://www.cdc.gov/coronavirus/2019-nCoV/hcp/infection-control.html>
- 4 World Health Organization. Infection prevention and control during health care when novel coronavirus (nCoV) infection is suspected. January 25, 2020. [https://www.who.int/publications-detail/infection-prevention-and-control-during-health-care-when-novel-coronavirus-\(ncov\)-infection-is-suspected-20200125](https://www.who.int/publications-detail/infection-prevention-and-control-during-health-care-when-novel-coronavirus-(ncov)-infection-is-suspected-20200125)
- 5 Ong SWX, Tan YK, Chia PY, Lee TH, Ng OT, Wong MSY, Marimuthu K.; Air, Surface Environmental, and Personal Protective Equipment Contamination by Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) From a Symptomatic Patient; JAMA. 2020 Mar 4. [Epub ahead of print]
- 6 van Doremalen N, Bushmaker T, Morris DH, Holbrook MG, Gamble A, Williamson BN, Tamin A, Harcourt JL, Thornburg N, Gerber SI, Lloyd-Smith JO, de Wit E, Munster VJ; Aerosol and Surface Stability of SARS-CoV-2 as Compared with SARS-CoV-1.; N Engl J Med. 2020 Mar 17. [Epub ahead of print]
- 7 Kampf G, Todt D, Pfaender S, Steinmann E.; Persistence of coronaviruses on inanimate surfaces and their inactivation with biocidal agents.; J Hosp Infect. 2020 Mar;104(3):246-251.
- 8 Cleaning Hospital Room Surfaces to Prevent Health Care–Associated Infections: A Technical Brief; Jennifer H. Han, MD, MSCE *; Nancy Sullivan, BA *; Brian F. Leas, MS, MA *; David A. Pegues, MD; Janice L. Kaczmarek, MS; Craig A. Umscheid, MD, MSCE; Ann Intern Med. 2015;163(8):598-607.
- 9 Snyder GM, Holyoak AD, Leary KE, Sullivan BF, Davis RB, Wright SB. Effectiveness of visual inspection compared with non-microbiologic methods to determine the thoroughness of post-discharge cleaning; Antimicrob Resist Infect Control 2013 2(1):26.
- 10 Mulvey, D. et al., Finding a benchmark for monitoring hospital cleanliness; Journal of Hospital Infection, 2011 77(1):25-30.
- 11 Luick L, Thompson PA, Loock MH, Vetter SL, Cook J, Guerrero DM. Diagnostic assessment of different environmental cleaning monitoring methods; Am J Infect Control. 2013 41(8):751-2.
- 12 Smith PW, Gibbs S, Sayles H, Hewlett A, Rupp ME, Iwen PC. Observations on hospital room contamination testing; Healthc Infect 2013 18:10-3.
- 13 Malik RE, Cooper RA, Griffith CJ. Use of audit tools to evaluate the efficacy of cleaning systems in hospitals; Am J Infect Control 2003 31:181-7.
- 14 Branch-Elliman W, Robillard E, McCarthy G Jr, Gupta K. Direct feedback with the ATP luminometer as a process improvement tool for terminal cleaning of patient rooms; Am J Infect Control 2014 42:195-7.

15 Smith PW, Beam E, Sayles H, Rupp ME, Cavalieri RJ, Gibbs S, et al. Impact of adenosine triphosphate detection and feedback on hospital room cleaning; *Infect Control Hosp Epidemiol* 2014 35:564-9.

16 Sitzlar, B., Deshpande, A., Fertelli, D., Kundrapu, S., Sethi, A., & Donskey, C. Sitzlar, B., Deshpande, A., Fertelli, D., Kundrapu, S., Sethi, A., Donskey, C. An environmental disinfection odyssey: Evaluation of sequential interventions to improve disinfection of clostridium difficile isolation rooms. *The role of the environment in Infection Prevention* May 2013. *Infection Control and Hospital Epidemiology* 34(4):459-465.

17 Katherine Hardy, Gill Abbott, Sarah Bashford, Helen Bucior, Jane Codd, Madelaine Holland, Mandy Reynolds, Avril Simms, Diane Thomlinson Can measuring environmental cleanliness using ATP aid in the monitoring of wards with periods of increased incidence of Clostridium difficile?; *J Infect Prev.* 2014 Jan; 15(1): 31–35

18 Deshpande A, Donskey CJ; Practical Approaches for Assessment of Daily and Post-discharge Room Disinfection in Healthcare Facilities; *Curr Infect Dis Rep.* 2017 Sep;19(9):32.

19 Deshpande A, Sitzlar B, Fertelli D, Kundrapu S, Sunkesula VC, Ray AJ, Donskey CJ; Utility of an adenosine triphosphate bioluminescence assay to evaluate disinfection of Clostridium difficile isolation rooms; Infect Control Hosp Epidemiol. 2013 Aug;34(8):865-7.

20 Boszczowski, Icaro; Alfa, Michele; Hopman, Joost; Donskey, Curtis J. et al.; Comparison of ATP, Colony Count and Surface Marker as a Measure of Environmental Cleaning Compliance for Intensive Care Discharge Rooms; American Journal of Infection Control, Volume 44, Issue 6, S77

21 Deshpande A, Dunn AN, Fox J, Cadnum JL, Mana TSC, Jencson A, Fraser TG, Donskey CJ, Gordon SM; Monitoring the effectiveness of daily cleaning practices in an intensive care unit (ICU) setting using an adenosine triphosphate (ATP) bioluminescence assay. Am J Infect Control. 2019 Dec 26. [Epub ahead of print]