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1. Sessler DI. Perioperative thermoregulation and heat balance. The Lancet. 2016;387(10038):2655-2664.
2. AORN. Updated Clinical Information to Support Nurses & Healthcare Teams During COVID-19 Pandemic. <https://www.aorn.org/guidelines/aorn-support>. Published 2020. Accessed 5/19/2020.
3. ANSI/ASHRAE/ASHE. Standard 170-2013: Ventilation of health care facilities. In: ASHRAE; 2013. Accessed 5/19/2020.
4. Clark RP, de Calcina-Goff ML. Some aspects of the airborne transmission of infection. Journal of the Royal Society Interface. 2009;6(Suppl 6): S767-S782.
5. Tang JW, Li Y, Eames I, Chan PKS, Ridgway GL. Factors involved in the aerosol transmission of infection and control of ventilation in healthcare premises. Journal of Hospital Infection. 2006;64(2):100-114.
6. ANSI/ASHRAE Standard 52.2-2012: Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. Published 2017. Accessed 5/19/2020.
7. Chow TT, Yang XY. Ventilation performance in operating theatres against airborne infection: review of research activities and practical guidance. Journal of Hospital Infection. 2004;56(2):85-92.
8. ASHRAE. HVAC design manual for hospitals and clinics. Atlanta: ASHRAE; 2003.
9. CDC. Guidelines for environmental infection control in health-care facilities <https://www.cdc.gov/infectioncontrol/guidelines/environmental/appendix/air.html>. Published 2003. Accessed 5/19/2020.
10. American Society of Anesthesiologists. COVID-19 FAQs. <https://www.asahq.org/about-asa/governance-and-committees/asa-committees/committee-on-occupational-health/coronavirus/clinical-faqs>. Published 2020. Accessed 5/19/2020.
11. CDC. Infection Control Guidance for Healthcare Professionals about Coronavirus (COVID-19). <https://www.cdc.gov/coronavirus/2019-nCoV/hcp/infection-control.html>. Published 2020. Accessed 5/19/2020.
12. CDC. Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings (2007). <https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html>. Published 2020. Accessed 5/19/2020.
13. CDC. Interim Infection Prevention and Control Recommendations for Patients with Suspected or Confirmed Coronavirus Disease 2019 (COVID-19) in Healthcare Settings. <https://www.cdc.gov/coronavirus/2019-ncov/hcp/infection-control-recommendations.html>. Published 2020. Accessed 5/19/2020.

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Please Recycle. Printed in Japan.
© 3M 2020. All Rights Reserved.
HPM-967-A(0920)IT

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