Acelity™

ENGLISH

SNAP PLUS™ 125mmHg THERAPY CARTRIDGE

Patient / Caregiver Guide

Rx Only

An Acelity Company

CE 0086

KCI USA, Inc.
12930 IH 10 West
San Antonio, TX 78249 USA
1-800-275-4524
www.acelity.com

KCI Manufacturing
Unlimited Company
IDA Business & Technology Park
Dublin Road, Athlone,
Co. Westmeath, Ireland
www.kci-medical.com

Refer to the SNAP PLUS™ Therapy System Instructions for Use provided with each SNAP™ Dressing Kit for complete safety information and instructions for use.

This is not a complete set of instructions.

SAY HELLO TO THE SNAP PLUS™ THERAPY CARTRIDGE

We have designed the SNAP PLUS™ Therapy Cartridge to be easy to use, lightweight and completely silent. The device delivers negative pressure wound therapy, which may promote wound healing through the removal of excess exudate, infectious material and tissue debris.

Visually inspect system at least every 8 hours.

THE SNAP PLUS™ THERAPY CARTRIDGE IS ACTIVE IF:

- **Green** capacity indicator is visible and stationary
- Dressing has a sucked down appearance
- Dressing is firm to the touch

WHEN SHOULD I CHANGE MY DEVICE?


The SNAP PLUS™ Therapy Cartridge should be changed when the **red** pressure discharge indicator is visible at the top of the SNAP PLUS™ Therapy Cartridge and further treatment with the SNAP™ Therapy System is desired or indicated by your clinician.

INSTRUCTIONS FOR CHANGING THE SNAP PLUS™ THERAPY CARTRIDGE

If the **red** pressure discharge indicator is visible at the top of the cartridge, the SNAP PLUS™ Therapy Cartridge will need to be replaced with a new cartridge following the steps below if further treatment with the SNAP™ Therapy System is desired or indicated by your clinician:

1. Remove the cartridge from the patient.

2. Squeeze the release buttons on the tube fitting and pull the tube fitting out of the cartridge.
3. Insert the tube fitting into the new cartridge.
4. Push down the activation key to release and remove the key from the new cartridge.
5. Secure cartridge to patient.

NOTE: When the SNAP™ Therapy Strap is placed around an extremity, take care to ensure that the strap is not placed too tightly as this may cause discomfort or potentially decrease blood flow to the extremity.

GENERAL PRECAUTIONS

- Contact your clinician if the dressing leaks. If you cannot see your clinician within 24 hours, remove the cartridge, dressing and other dressing materials and cover the wound with sterile moist gauze or other dressing as instructed by your clinician.
- Immediately contact your clinician if you develop signs of infection such as increased swelling, fever, pain or redness at the wound area. Remove the cartridge, dressing, and other dressing materials and cover the wound with sterile moist gauze or other dressing as instructed by your clinician.
- Patients may shower or bathe with the SNAP™ Therapy System. If bathing, do not submerge the SNAP™ Therapy System under water.
- BIOLOCK™ Technology will gel exudate in the SNAP™ Therapy Cartridge.

EXPLANATION OF SYMBOLS USED

Consult Instructions
For Use

Manufacturer

Single Use Only

Use By

Do not use if package
is damaged or open

Keep Dry

Non-Sterile

Lot number

Date of
Manufacture

MR Unsafe

Catalog Number

Rx Only

CAUTION: Federal (US) law restricts this device to sale/
rental by or on the order of a physician

Authorized Representative in the
European Community

All trademarks designated herein are proprietary to KCI Licensing, Inc., its affiliates and / or licensors.

©2018 KCI Licensing, Inc. All rights reserved.

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Alle andre varemerker som nevnes her, tilhører KCI Licensing, Inc., deres datterselskaper og/eller lisensgivere.
©2018 KCI Licensing, Inc. Med enerett.