




Prevena™
Incision Therapy

Evidence Brochure

Expanded FDA Indication
3M™ Prevena™ 125 Therapy
Unit and 3M™ Prevena™
Plus 125 Therapy Unit



Incisions can be complicated

Certain surgical procedures and patient conditions can make healing difficult

Surgical procedures that most commonly lead to complications include sternotomies, C-sections, open hysterectomies, hip and knee arthroplasties, open reduction fractures, lower extremity bypasses, femoropopliteal bypasses, renal transplants, and breast reconstruction.¹

Risk Factors that may compromise healing²⁻⁴

- Obesity
- Nicotine use
- Diabetes—poor control
- Radiation therapy
- Age >65
- Wound infection
- Pulmonary disease
- Peripheral vascular disease
- Hemodynamic instability
- Ostomies
- Hypoalbuminemia
- Systemic infection
- Uremia
- Hyperalimentation
- Ascites
- Malignancy
- Hypertension
- Length and depth of incision
- Anemia
- Jaundice
- Type of injury
- Steroid use
- Malnutrition

Incisions can be costly

There are

8 Million People

at risk for healthcare-associated infections annually.⁵

Post-surgical complications lead to significant costs

Surgical site infections (SSIs) are

21.8%

of all healthcare-associated infections⁶

Of the top 5 Healthcare Acquired Infections (HAIs), SSI is

33.7%

of the **\$9.8 Billion** cost to the US healthcare system*⁵

SSIs increase average length of hospital stay by

9.58 days

at an additional cost of **\$38,656**⁷

Other common complications include **dehiscence**, **hematoma** and **seroma** formation²⁻⁴

Consequences extend beyond discharge

Patients with an SSI are **6 times** more likely to have a

30-day readmission

than patients without an SSI⁷

Patients with SSIs have an ICU length of stay that is

2.2 times greater

than patients without SSIs⁷

Postoperative dehiscence increases average length of hospital stays by

9.42 days

and average costs by **\$40,323**^{2,8-10}

The Centers for Medicare & Medicaid Services emphasize the need to decrease costs and improve care by identifying hospital-acquired conditions that will not be reimbursed, including 3 SSIs¹¹:



1. Mediastinitis following coronary artery bypass graft (CABG)
2. SSIs following certain orthopedic procedures
3. SSIs following bariatric surgery for obesity

*Top five HAIs are central line-associated bloodstream infections (CLABSI), ventilator-associated pneumonia (VAP), Clostridium difficile infection (C diff), Surgical Site Infections (SSI) and catheter-associated urinary tract infection (CAUTI).

How 3M™ Prevena™ Therapy can help

The FDA granted the following Indications for Use:*


3M™ Prevena™ 125 Therapy Unit and 3M™ Prevena™ Plus 125 Therapy Unit manage the environment of closed surgical incisions and remove fluid away from the surgical incision via the application of -125mmHg continuous negative pressure. When used with legally marketed compatible dressings, Prevena 125 Therapy Unit and Prevena Plus 125 Therapy Unit are intended to aid in reducing the incidence of seroma and, in patients at high risk for post-operative infections, aid in reducing the incidence of superficial surgical site infection in Class I and Class II wounds.

Contraindication:

Sensitivity to Silver.

Optimum Use:†

For maximum benefit, the 3M™ Prevena™ Incision Management System should be applied immediately post surgery to surgically closed incisions. It is to be continuously applied for a minimum of 2 days up to a maximum of 7 days. It can transition home with the patient.




Clinical evidence supporting the new indication

A systematic literature review and associated meta-analysis were used to support the safety and effectiveness of 3M™ Prevena™ Therapy over closed incisions in reducing the incidence of surgical site infections (SSIs) and seromas versus conventional wound dressings.

- Out of 426 studies in the initial search, ultimately, sixteen (16) prospective studies were included in this meta-analysis for SSI characterization
- A total of up to 6,187 evaluable patients were included in this meta-analysis for SSI with 1,264 in the Prevena Therapy (treatment) group and 4,923 in the conventional wound dressing (control) group
- 9 randomized controlled trials (RCTs) were included in a subgroup analysis for SSI in high risk patients

Prevena Therapy demonstrated the greatest benefit in reducing SSIs in high risk patients

Forest Plot of Meta-Analysis on Surgical Site Infection




***The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at 3M/Medical.com**

[†]Refer to the 3M™ Prevena™ Incision Management System Clinician Guide for additional information relating to Optimum Use, Indications and Contraindications, Warnings and Precautions, and Important Safety Information.

3M™ Prevena™ Therapy manages and protects surgical incisions utilizing unique 3M™ Prevena™ Dressings by:



Delivering continuous
-125mmHg up
to 7 days



Helping to hold
incision edges
together



Protecting the incision
from external
infectious
sources



Removing fluids
and infectious
materials



Reducing
edema





3M™ Prevena™ Therapy utilizes reticulated open cell foam technology and -125mmHg pressure

Passive Therapy



3M™ Prevena™ Therapy



Under -125mmHg of negative pressure, the Reticulated Open Cell Foam dressing collapses to its geometric center. This brings the incision edges together, reduces lateral tension, and also allows for improved fluid management.¹²⁻¹⁴


- Contours in 3M™ Prevena™ Dressings allow for even distribution of negative pressure
- Adhesive film creates a barrier to external contaminants
- Designed to conform to articulating joints to allow movement
- Skin interface layer contains 0.019% ionic silver, which reduces bacterial colonization in the fabric
- Multiple sizes and configurations
- 3M™ Prevena™ 125 Therapy Unit and Prevena Dressings are shower friendly*

*See Prevena Therapy Patient and Clinician Guides for additional details


Use of closed incisional negative pressure wound therapy after revision total hip and knee arthroplasty in patients at high risk for infection: A prospective, randomized clinical trial¹⁵

- The purpose of the Newman study was to compare the use of 3M™ Prevena™ Therapy to a sterile antimicrobial dressing (Aquacel® Ag) in revision arthroplasty patients who were at high risk to develop wound complications.
- 160 patients undergoing elective revision arthroplasty were prospectively randomized to receive either 3M™ Prevena Therapy or Aquacel® Ag in a single institution.
- Patients were included if they had at least 1 risk factor for developing wound complication.
- Study endpoints included wound complications (such as SSI, drainage, and cellulitis) readmission, and reoperation rates were collected at 2, 4, and 12 weeks postoperatively.
- **The postoperative wound complication rate was significantly higher in the Aquacel® Ag compared to the Prevena Therapy group (19 [23.8%] vs 8 [10.1%], $p = 0.022$).**
- There was no significant difference between the Aquacel® Ag and Prevena Therapy cohorts in terms of readmissions (19 [23.8%] vs 16 [20.3%], $p = 0.595$).
- **Reoperation rate was higher in Aquacel® Ag patients compared to Prevena Therapy patients (10 [12.5%] vs 2 [2.5%], $p = 0.017$).**
- After adjusting for the history of a prior periprosthetic joint infection and inflammatory arthritis, **the Prevena Therapy cohort had a significantly decreased wound complication rate (odds ratio 0.28, 95% confidence interval 0.11-0.68).**

Wound Complications (wks. 2, 4, and 12)¹⁵



Reoperation Rate¹⁵



Note: Although the authors reported use of Prevena Therapy for a mean of 3.6 days (ranging from 2 to 15 days), this mean time of application is outside the recommendations for Optimum Use as stated in the 3M™ Prevena™ Incision Management System Clinician Guide Instructions for Use: "The Prevena Incision Management System is to be continuously applied for a minimum of 2 days up to a maximum of 7 days." Use for greater than 7 days is not recommended or promoted by 3M.

A hypothetical cost model applied to the clinical results of the Newman study shows potential cost savings per patient of \$1,959 with the use of Prevena Therapy.

Economic Model

Revision Hip (THA) and Knee (TKA) Surgery Hypothetical Economic Model	3M™ Prevena™ Therapy (n = 79)	Aquacel® Ag (n = 80)
Number of Reoperations at 2, 4, and 12 weeks (a)	2	10
Average Estimated Cost of Reoperation* (b)	\$24,200	\$24,200
Total Reoperation Cost (a*b)	\$48,400	\$242,000
Per Patient Cost of Reoperation (a*b)/n	\$613	\$3,025
Per Patient Cost of Therapy [§]	\$495	\$42
Total Cost per Patient	\$1,108	\$3,067

*Kurtz, Steven M. et al. Economic Burden of Periprosthetic Joint Infection in the United States. *J Arthroplasty*. 2012 Sep;27(8 Suppl):61-5.e1.


§3M estimate based on price of 3M™ Prevena™ Peel and Place System Kit and Aquacel® Ag; individual prices may vary.

The hypothetical economic model uses select study data to provide an illustration of estimates of costs for use of Prevena Therapy or Aquacel® Ag. This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.


Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy¹⁶

- A prospective, single center clinical trial evaluated the use of 3M™ Prevena™ Therapy compared to standard post-operative dressings (Control) for the prevention of wound infection within 90 days after median sternotomy procedures in 150 consecutive obese (BMI ≥ 30) patients.
- **Patients treated with Prevena Therapy developed fewer wound infections (3/75 [4%] vs 12/75 [16%], $p=0.0266$) than patients treated with standard post-operative dressings.**
- Wound infections with Gram-positive skin flora were found in only 1 patient in the Prevena Therapy group compared with 10 patients in the control group ($p=0.0090$).
- A hypothetical cost model applied to the clinical results of this study shows a **potential cost savings per patient of \$4,024 with the use of Prevena Therapy.***

Infection Rate¹⁶



Total Cost Per Patient*



Economic Model

Post-sternotomy Hypothetical Economic Model	3M™ Prevena™ Therapy (n = 75)	Control (n=75)
Number of Infections (a)	3	12
Percent of Infections	4.0%	16.0%
Cost per Infection ⁶⁶ (b)	\$37,513	\$37,513
Cost of Infection Per Patient (a*b)/n)	\$1,501	\$6,002
Cost of Therapy Per Patient [†]	\$495	\$18
Total Cost Per Patient	\$1,996	\$6,020

*The hypothetical economic model uses select study data to provide an illustration of estimates of costs for use of Prevena Therapy or standard post-operative dressings (Control). This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results.


[†]3M estimate based on price of 3M™ Prevena™ Peel and Place System Kit and Control therapy (gauze) changed once a day at \$18 a week.

The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.


Using closed incision negative pressure therapy (ciNPT) in high-risk general surgery patients following laparotomy: A retrospective study¹⁷

- The aim of this study was to compare wound complication (skin wound dehiscence or deep incisional infection) occurrence in patients with multiple comorbidities who received 3M™ Prevena™ Therapy or conventional care (Control) within 30 days following laparotomy.
- Inclusion criteria were obesity (BMI ≥ 35 kg/m²), or two or more of the following risk factors: malignancy, history of smoking, immunosuppression, malnutrition, emergency surgery, diffuse atherosclerotic disease.
- Compared to the control group, **patients treated with Prevena Therapy had significantly fewer wound complications (2/69 [2.9%] vs 23/112 [20.5%],** respectively; $p < 0.0009$).
- The relative risk of a wound complication in the Prevena Therapy group was 0.12 (95% CI 0.03-0.51) compared with the control group of 0.14 (95% CI 0.03-0.58); suggesting that infection is less likely to occur in Prevena Therapy treated incisions.
- A hypothetical cost model applied to the clinical results of this study shows **potential cost savings per patient of \$2,777 with the use of Prevena Therapy.***

Total Wound Complications¹⁷



Total Per Patient Cost*



Economic Model

Laparotomy Hypothetical Economic Model	3M™ Prevena™ Therapy (n = 69)	Control (n=112)
Deep Incisional Infection	1	23
Skin Wound Dehiscence	1	0
Total Wound Complications	2	23
Total Infection Cost (Incremental cost of infection = \$17,995 per patient) ⁶⁶	\$17,995	\$413,885
Total Dehiscence Cost (Incremental cost of dehiscence = \$12,407 per patient) ⁶⁵	\$12,407	\$0
Total Complication Cost (a)	\$30,402	\$413,885
Complication Cost Per Patient (a/n)	\$441	\$3,695
Cost of Therapy Per Patient [†]	\$495	\$18
Total Cost Per Patient	\$936	\$3,713

*The hypothetical economic model uses select study data to provide an illustration of estimates of costs for use of Prevena Therapy or standard post-operative dressings (Control). This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results.

[†]3M estimate based on price of 3M™ Prevena™ Peel and Place System Kit and Control therapy (gauze) changed once a day at \$18 a week.

The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.


Reduction of wound groin complications in vascular surgery patients using closed incision negative pressure therapy: A prospective, randomized single-institution study¹⁸

- The aim of this prospective, randomized, single-institution study was to investigate the effectiveness of 3M™ Prevena™ Therapy compared to conventional adhesive dressing (Control) on groin incisions after vascular surgery.
- The Prevena Therapy group had 43 patients and 58 groin incisions and the control group consisted of 57 patients and 71 groin incisions.
- Wound evaluation based on the Szilagyi classification (Grade I, II, and III) took place postoperatively on days 5–7 and 30.
- In this study, patients with cutaneous wound dehiscence, skin necrosis and single local infection signs were classified as grade I. Wound dehiscence in the subcutaneous layer, hematoma, lymphatic fistula, lymphocele, seroma, single local infection signs and systemic infection parameters were classified as grade II.


All classical local infection signs (pain, swelling, redness and hyperaemia, warmth, dysfunction), systemic infection parameters and arterial graft infections were classified as grade III.

- **Prevena Therapy significantly reduced the incidence of local infection** compared to the conventional dressing (1/43 [2.3%] vs 10/57 [17.5%], respectively; $p=0.022$).
- Compared to the control group, the Prevena Therapy group showed a significant reduction in wound complications after both evaluation periods (5/58 [8.62%] vs 30/71 [42.3%], $p<0.0005$).
- Prevena Therapy showed a significant reduction in revision surgeries (1/58 [1.7%] vs 10/71 [14.1%], respectively; $p=0.022$) until 30 days postoperatively compared to the control group.
- A hypothetical cost model applied to the clinical results of this study shows a **potential cost savings per patient of \$2,694 with the use of Prevena Therapy.***


Local Infection Rate¹⁸



Wound Complications At Both Evaluation Periods¹⁸



Total Per Patient Cost*



Economic Model

Vascular Groin Hypothetical Economic Model	3M™ Prevena™ Therapy	Control
Number of Patients (n)	43	57
Number of Local Infections (a)	1	10
Percent of Local Infections	2.3%	17.5%
Cost Per Local Infection ⁶⁶ (b)	\$20,842	\$20,842
Cost of Local Infection Per Patient (a*b)/n)	\$485	\$3,656
Cost of Therapy Per Patient [†]	\$495	\$18
Total Cost Per Patient	\$980	\$3,674

*The hypothetical economic model uses select study data to provide an illustration of estimates of costs for use of Prevena Therapy or standard post-operative dressings (Control). This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results.

[†]3M estimate based on price of 3M™ Prevena™ Peel and Place System Kit and Control therapy (gauze) changed once a day at \$18 a week.

The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

The negative pressure incision management system with the most published clinical evidence*

3M™ Prevena™ Therapy helps manage and protect closed surgical incisions utilizing a unique 3M™ Prevena™ Dressing by:

- Protecting the incision from external infectious sources
- Delivering continuous negative pressure (-125 mmHg) for up to 7 days
- Helping to hold incision edges together
- Removing fluids and infectious materials
- Reducing edema

*Among negative pressure based incision management systems

Ordering information

Item #	Product Name	Qty	Item #	Product Name	Qty
PRE1001US	3M™ Prevena™ Peel and Place Incision Management System – 20cm	1	PRE1121US	3M™ Prevena™ Duo Incision Management System with Peel and Place Dressings - 13cm/13cm	1
PRE1055US	3M™ Prevena™ Peel and Place Dressing – 20cm	5	PRE3321US	3M™ Prevena™ Duo Incision Management System with Peel and Place Dressings - 13cm/20cm	1
PRE1101US	3M™ Prevena™ Peel and Place Incision Management System – 13cm	1	PRE3021US	3M™ Prevena™ Duo Incision Management System with Peel and Place Dressings - 20cm/20cm	1
PRE1155US	3M™ Prevena™ Peel and Place Dressing – 13cm	5	PRE4000US	3M™ Prevena™ Plus 125 Therapy Unit	1
PRE3201US	3M™ Prevena™ Plus Peel and Place Incision Management System – 35cm	1	PRE4010	3M™ Prevena™ Plus 125 Therapy Unit (14 Day)	1
PRE3255US	3M™ Prevena™ Plus Peel and Place Management System 35cm Dressing	5	PRE1095	3M™ Prevena™ 45ml Canister	5
PRE4001US	3M™ Prevena™ Plus Customizable Incision Management System	1	PRE4095	3M™ Prevena™ Plus 150ml Canister	5
PRE4055US	3M™ Prevena™ Plus Customizable Dressing	5	PRE9090	3M™ Prevena™ Therapy V.A.C.® Connector	10

To learn more about the benefits of 3M™ Prevena™ Therapy, contact your 3M representative.

References

1. Shrestha BM, Nathan VC, Delbridge MS, et al. Vacuum-assisted closure (VAC) therapy in the management of wound infection following renal transplantation. *Kathmandu Univ Med J.* 2007;5:4-7. 2. Riou JP, Cohen JR, Johnson H Jr. Factors influencing wound dehiscence. *Am J Surg.* 1992;163:324-330. 3. Wilson JA, Clark JJ. Obesity: impediment to postsurgical wound healing. *Adv Skin Wound Care.* 2004;17:426-435. 4. Abbas SM, Hill AG. Smoking is a major risk factor for wound dehiscence after midline abdominal incision; case-control study. *ANZ J Surg.* 2009;79:247-250. 5. Zimlichman E, Henderson D, Tamir, et al. Health care-associated infections a meta-analysis of costs and financial impact on the U.S. health care system. *JAMA Intern ed.* 2013;173(22):20-46. 6. Magill SS, Edwards JR, Bamberg W, et al. Multistate point-prevalence survey of health care-associated infections. *N Engl J Med.* 2014;370:1198-208. 7. Shepard J, Ward W, Millstone A, et al. Financial impact of surgical site infections on hospitals. The hospital management perspective. *JAMA Surg.* 2013;148(10):907-914. doi:10.1001/jamasurg.2013.2246 Published online August 21, 2013. 8. Olsen K. Prevention of surgical site infections: improving compliance with the surgical care improvement project measures. <http://www.medscape.com/viewarticle/705366>. Accessed September 20, 2010. 9. Klevens RM, Edwards JR, Richards CL Jr, et al. Estimating health care-associated infections and deaths in U.S. hospitals, 2002. *Public Health Rep.* 2007;122:160-166. 10. National Nosocomial Infections Surveillance report, data summary from October 1986-April 1996, issued May 1996. A report from the national nosocomial infections surveillance system. *Am J Infect Control.* 1996;24:380-388. 11. US Department of Health and Human Services. Hospital-acquired conditions and present on admission indicator reporting provision. ICN 901046. <https://www.cms.gov/Outreach-and-Education/Medicare-Learning-Network-MLN/MLNProducts/downloads/wPOAFactSheet>. Published September 2014. Accessed June 10, 2015. 12. Wilkes RP, Kilpadi DV, Zhao Y, et al. Closed incision management with negative pressure wound therapy (CIM): biomechanics. *Surgical Innovation.* 2012;19(1):67-75. 13. Kilpadi DV, Cunningham MR. Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system. *Wound Repair and Regeneration.* 2011;19:588-596. 14. Glaser DA, Farnsworth CL, Varley ES, et al. Negative pressure therapy for closed spine incisions: a pilot study. *Wounds.* 2012;24(11):308-316. 15. Newman JM, Siqueira MBP, Klika AK, Molloy RM, Barsoum WK, Higuera CA. Use of closed incisional negative pressure wound therapy after revision total hip and knee arthroplasty in patients at high risk for infection: a prospective, randomized clinical trial. *Journal of Arthroplasty.* 2018;34:554-559. 16. Grauhan O, Navasardyan A, Hofmann M, Muller P, Stein J, Hetzer R. Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy. *J Thorac Cardiovasc Surg.* 17. Zaidi A, El-Masry S. Closed-incision negative-pressure therapy in high-risk general surgery patients following laparotomy: a retrospective study. *Colorectal Dis.* 2017 Mar;19(3):283-287. 18. Pleger SP, Nink N, Elzien M, Kunold A, Koshty A, Böning A. Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study. *Int Wound J.* 2018 Feb;15(1):75-83. doi: 10.1111/iwj.12836. Epub 2017 Oct 25.



3M Company
2510 Conway Ave
St. Paul, MN 55144 USA

Phone 1-800-275-4524 (NPWT products)
1-800-228-3957
Web 3M.com/medical

The effectiveness of 3M™ Prevena™ Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at 3M.com.

Note: Specific indications, limitations, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.

© 2021 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.
70-2013-1105-0 PRA-PM-US-01710 (06/21)