





This guide may be used in deciding which wound dressing to use in conjunction with V.A.C. VERAFLOR™ Therapy

	V.A.C. VERAFLOR™ Dressings	V.A.C. VERAFLOR™ Large Dressing	V.A.C. VERAFLOR CLEANSE™ Dressing	V.A.C. VERAFLOR CLEANSE CHOICE™ Dressing
Wound Types				
Wound Characteristics	Open wounds, including wounds with shallow undermining or tunnel areas where the distal aspect is visible	Large open wounds, including wounds with shallow undermining or tunnel areas where the distal aspect is visible	Cavity wounds or wounds with complex geometries	Wounds with thick fibrinous exudate, slough, infectious material and other wound bioburden
Key Goal(s) of Therapy	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR™ Therapy, to help facilitate the removal of infectious material and other wound bioburden.</li> <li>Generation of robust granulation tissue</li> </ul>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR™ Therapy, to help facilitate the removal of infectious material and other wound bioburden.</li> <li>Generation of robust granulation tissue in large wounds</li> </ul>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR™ Therapy, to initiate therapy and to help facilitate the removal of infectious material and other wound bioburden.</li> <li>Easy application into tunneling and undermining</li> </ul>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR™ Therapy, to initiate therapy and to help facilitate the removal of infectious material such as thick fibrinous exudate, slough, and other wound bioburden.</li> <li>Provides a wound cleansing option for clinicians when surgical debridement must be delayed or is not possible or appropriate.<sup>1</sup></li> </ul>

	V.A.C. VERAFLOR™ Dressings	V.A.C. VERAFLOR™ Large Dressing	V.A.C. VERAFLOR CLEANSE™ Dressing	V.A.C. VERAFLOR CLEANSE CHOICE™ Dressing
SKU	Small - ULTVFL05SM / Medium - ULTVFL05MD	ULTVFL05LG	ULTVCL05MD	Medium - ULTVCC05MD / Large - ULTVCC05LG
Shape	Spiral-cut foam	Block foam pre-slit into two layers	Tubular shape	Block foam pre-slit into three layers
Application Characteristics	<p><b>Easy application:</b></p> <ul style="list-style-type: none"> <li>• Size without scissors</li> <li>• Precut area for pad application when used for bridging</li> <li>• Single pad application</li> </ul>	<p><b>Easy application:</b></p> <ul style="list-style-type: none"> <li>• Ideal for large surfaces areas with shallow depths</li> <li>• Provided with V.A.C. VERAT.R.A.C. DUO™ Tube set for extended surface area coverage</li> </ul>	<p><b>Application flexibility:</b></p> <ul style="list-style-type: none"> <li>• Ideal for addressing wounds with complex geometries (eg, tunnels, undermining and circumferential applications)</li> <li>• Single pad application</li> </ul>	<p><b>Easy application with flexibility:</b></p> <ul style="list-style-type: none"> <li>• Thin layers for improved conformability</li> <li>• Multiple layers provide application options for wounds with varying depths</li> <li>• Ideal for dirty wounds needing active therapy</li> <li>• Single or duo pad application</li> </ul>
Dressing Characteristics	<p><b>Small: 7.7 x 11.2 x 1.75cm per piece (2 pieces spiral)</b></p> <ul style="list-style-type: none"> <li>• Each piece 24cm uncoiled</li> <li>• Each piece 18cm long bridge</li> </ul> <p>- 2x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Pad - 1x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>Medium: 14.7 x 17.4 x 1.75cm per piece (2 pieces spiral)</b></p> <ul style="list-style-type: none"> <li>• Each piece 80cm uncoiled</li> <li>• Each piece 74cm long bridge</li> </ul> <p>- 3x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Pad - 2x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>FILL ASSIST SET POINT RANGE</b> Small: 12-80mL (1 piece) 26-160mL (2 pieces) Medium: 38-240mL (1 piece) 80-450mL (2 pieces)</p>	<p>25.0 x 15.0 x 1.6 cm per piece (oval - 2 pieces)</p> <p>- 5x V.A.C.® Advanced Drape - 1x V.A.C. V.A.C. VERAT.R.A.C. Duo™ Pad - 4x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>FILL ASSIST SET POINT RANGE</b> 55-350mL (1 piece) 110-500mL (2 pieces)</p>	<p>1 each coiled rod (24" rod; 1.25" diameter)</p> <p>- 3x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Pad - 2x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>FILL ASSIST SET POINT RANGE</b> 65-130mL (1/2 rod) 130-250mL (full rod)</p>	<p><b>Medium: 18.0cm x 12.5cm x Varying Thicknesses (0.8cm – 3.2cm)</b></p> <p>- 4x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Pad - 4x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>Large: 25.6cm x 15.0cm x Varying Thicknesses (0.8cm – 3.2cm)</b></p> <p>- 6x V.A.C.® Advanced Drape - 1x V.A.C. VERAT.R.A.C.™ Duo Pad - 5x 3M™ Cavilon™ No Sting Barrier Film</p> <p><b>FILL VOLUME START POINTS</b></p> <p><b>Medium:</b> 85mL (1.6cm Cover Layer); 42mL (0.8cm Cover Layer); 24mL (0.8cm Wound Contact Layer)</p> <p><b>Large:</b> 150 mL (1.6cm Cover Layer); 75mL (0.8cm Cover Layer); 42mL (0.8cm Wound Contact Layer)</p>

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a physician and product instructions for use prior to application. This material is intended for healthcare professionals.

© 2021 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.

PRA-PM-IN-00186 (01/21)