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# Closed incision negative pressure therapy in oncological breast surgery: comparison with standard care dressings.<sup>12</sup>

Ferrando PM, Ala A, Bussone R, Bergamasco L. *Plast Reconstr Surg Glob Open*. Juni 2018;6(6):e1732.

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

# Economic analysis based on the use of closed-incision negative-pressure therapy after postoperative breast reconstruction.<sup>7</sup>

Gabriel A, Maxwell GP. *Plast Reconstr Surg.* 2019; 143: 36S-40S.

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

# Closed-incision negative-pressure therapy decreases complications in ventral hernia repair with concurrent panniculectomy.<sup>13</sup>

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
# Closed incision negative pressure therapy over incisions following hernia repair and panniculectomy.

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## Niveau van klinisch bewijs.

| Geciteerde bron                                                                                                                                                                                                                                                                                                                                                                          | Type wond/operatie                                      | Niveau van klinisch bewijs* |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|-----------------------------|
| Muller-Sloof E, de Laat HEW, Hummelink SLM, Peters JWB, Ulrich DJO. The effect of postoperative closed incision negative pressure therapy on the incidence of donor site wound dehiscence in breast reconstruction patients: DEhiscence PREvention Study (DEPRES), pilot randomized controlled trial. <i>Journal of Tissue Viability</i> . 2018;27(4):262-266.                           | Borstreconstructie                                      | 1b ●                        |
| Ferrando PM, Ala A, Bussone R, Bergamasco L, Actis Perrinetti F, Malan F. Closed Incision Negative Pressure Therapy in Oncological Breast Surgery: Comparison with Standard Care Dressings. <i>Plastic and Reconstructive Surgery - Global Open</i> . Juni 2018;6(6):e1732.                                                                                                              | Borstreconstructie                                      | 2 ●                         |
| Papp AA. Incisional negative pressure therapy reduces complications and costs in pressure ulcer reconstruction. <i>International Wound Journal</i> . December 2018. doi:10.1111/iwj.13045.                                                                                                                                                                                               | decubitus-vorming door ruggenmergletsel                 | 2 ●                         |
| Renno I, Boos AM, Horch RE, Ludolph I. Changes of perfusion patterns of surgical wounds under application of closed incision negative pressure wound therapy in postbariatric patients. <i>Clinical Hemorheology and Microcirculation</i> . Januari 2019. doi:10.3233/CH-180450.                                                                                                         | Abdominoplastiek                                        | 2 ●                         |
| Swanson EW, Cheng HT, Susarla SM, Lough DM, Kumar AR. Does negative pressure wound therapy applied to closed incisions following ventral hernia repair prevent wound complications and hernia recurrence? A systematic review and meta-analysis. <i>Plastic Surgery</i> . Zomer 2016;24(2):113-8.                                                                                        | Ventrale hernia-correctie                               | 2 ●                         |
| Chowdhry SA, Wilhelmi BJ. Comparing Negative Pressure Wound Therapy with Instillation and Conventional Dressings for Sternal Wound Reconstructions. <i>Plastic and Reconstructive Surgery - Global Open</i> . 2019;7(1). doi:10.1097/gox.0000000000002087.                                                                                                                               | Spierflap-reconstructie van sternale wondcomplicaties   | 3 ●                         |
| Conde-Green A, Chung TL, Holton LH 3rd, Hui-Chou HG, Zhu Y, Wang H, Zahiri H, Singh DP. Incisional negative-pressure wound therapy versus conventional dressings following abdominal wall reconstruction: a comparative study. <i>Annals of Plastic Surgery</i> . Oktober 2013;71(4):394-7.                                                                                              | Abdominale herniacorrectie                              | 3 ●                         |
| Jorgensen MG, Toyserkani NM, Thomsen JB, Sorensen JA. Prophylactic incisional negative pressure wound therapy shows promising results in prevention of wound complications following inguinal lymph node dissection for Melanoma: A retrospective case-control series. <i>J Plast Reconstr Aesthet Surg</i> . 2 mrt 2019.                                                                | Dissectie lymfeknoten in de lies                        | 3 ●                         |
| Jorgensen MG, Toyserkani NM, Thomsen JB, Sorensen JA. Prophylactic incisional negative pressure wound therapy shows promising results in prevention of wound complications following inguinal lymph node dissection for Melanoma: A retrospective case-control series. <i>Journal of Plastic, Reconstructive &amp; Aesthetic Surgery</i> . 2019;000:1-6. doi:10.1016/j.bjps.2019.02.013. | Dissectie lymfeknoten in de lies                        | 3 ●                         |
| Lo Torto F, Monfrecola A, Kaciulyte J, Ciudad P, Casella D, Ribuffo D, Carlesimo B. Preliminary result with incisional negative pressure wound therapy and pectoralis major muscle flap for median sternotomy wound infection in a high-risk patient population. <i>Int Wound J</i> . Dec 2017;14(6):1335-1339.                                                                          | Grote-borstspierflap voor wondinfecties bij sternotomie | 3 ●                         |
| Gabriel A, Sigalove S, Storm-Dickerson T, Rice J, Maxwell P, Griffin L. The Impact of Closed Incision Negative Pressure Therapy on Postoperative Breast Reconstruction Outcomes. <i>Plastic and Reconstructive Surgery</i> . Global Open. Aug 2018;6(8):e1880.                                                                                                                           | Borstreconstructie                                      | 3 ●                         |
| Diaconu SC, McNichols CHL, Ngaage LM, Liang Y, Ikkeloe E, Bai J, Grant MP, Nam AJ, Rasko YM. Closed-incision negative-pressure therapy decreases complications in ventral hernia repair with concurrent panniculectomy. <i>Hernia</i> . 17 dec 2018. [Epub voorafgaand aan gedrukte publicatie]                                                                                          | Ventrale herniacorrectie                                | 3 ●                         |
| Abatangelo S, Saporiti E, Giatsidis G. Closed Incision Negative-Pressure Therapy (ciNPT) Reduces Minor Local Complications in Post-bariatric Abdominoplasty Body Contouring: a Retrospective Case. <i>Obese Surg</i> . Juli 2018;28(7):2096-2104.                                                                                                                                        | Abdominoplastiek                                        | 3 ●                         |

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## Literatuur

- 1 Starnoni M, Pinelli M, De Santis G. Surgical Wound Infections in Plastic Surgery: Simplified, Practical, and Standardized Selection of High-risk Patients. *Plast Reconstr Surg Glob Open*. Apr 2019; 7(4): e2202.
- 2 Alderman AK, Wilkins EG, Kim HM, Lowery JC. Complications in postmastectomy breast reconstruction: two-year results of the michigan breast reconstruction outcome study. *Plast Reconstr Surg*. 2002;109:2265-2274. Geciteerd door: Gabriel A1, Maxwell GP. Economic analysis based on the use of closed-incision negative-pressure therapy after postoperative breast reconstruction. *Plast Reconstr Surg*. Jan 2019;143(1S Management of Surgical Incisions Utilizing Closed-Incision Negative-Pressure Therapy):36S-40S.
- 3 Phillips BT, Bishawi M, Dagum AB, Khan SU, Bui DT. A systematic review of antibiotic use and infection in breast reconstruction: what is the evidence? *Plast Reconstr Surg*. 2013;131:1-13. Geciteerd door: Gabriel A1, Maxwell GP. Economic analysis based on the use of closed-incision negative-pressure therapy after postoperative breast reconstruction. *Plast Reconstr Surg*. Jan 2019;143(1S Management of Surgical Incisions Utilizing Closed-Incision Negative-Pressure Therapy):36S-40S.
- 4 Sullivan SR, Fletcher DR, Isom CD, Isik FF. True incidence of all complications following immediate and delayed breast reconstruction. *Plast Reconstr Surg*. 2008;122:19-28.
- 5 Olsen M, Nickel KB, et al. Incidence of Surgical Site Infection Following Mastectomy With and Without Immediate Reconstruction Using Private Insurer Claims Data. *Inf Control & Hosp Epidemiol*. 2015;36(8):907-914.
- 6 Bennett KG, Qi J, Kim HM, et al. Comparison of 2-year complication rates among common techniques for postmastectomy breast reconstruction. *JAMA Surg*. 2018;153(10):901-908. doi:10.1001/jamasurg.2018.1687.
- 7 Gabriel A, Maxwell GP. Economic analysis based on the use of closed-incision negative-pressure therapy after postoperative breast reconstruction. *Plast Reconstr Surg*. Jan 2019;143(1S Management of Surgical Incisions Utilizing Closed-Incision Negative-Pressure Therapy):36S-40S. doi: 10.1097/PRS.0000000000000531.
- 8 Wilkes RP, Kilpadi DV, Zhao Y, et al. Closed incision management with negative pressure wound therapy (CIM): biomechanics. *Surgical Innovation*. 2012;19(1):67-75.
- 9 Kilpadi DV, Cunningham MR. Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system. *Wound Repair and Regeneration*. 2011;19:588-596.
- 10 Glaser DA, Farnsworth CL, Varley ES, et al. Negative pressure therapy for closed spine incisions: a pilot study. *Wounds*. 2012;24(11):308-316.
- 11 Federal Drug Administration. De Novo Classification Request for PREVENA 125 and PREVENA PLUS 125 Therapy Units. De Novo Summary (DEN180013), 2019. [https://www.accessdata.fda.gov/cdrh\\_docs/reviews/DEN180013.pdf](https://www.accessdata.fda.gov/cdrh_docs/reviews/DEN180013.pdf)
- 12 Ferrando PM, Ala A, Bussone R, Bergamasco L. Closed incision negative pressure therapy in oncological breast surgery: comparison with standard care dressings. *Plast Reconstr Surg Glob Open*. Juni 2018;6(6):e1732. doi:10.1097/GOX00000000000001732.
- 13 Diaconu SC, McNichols CHL, Ngaage LM, et al. Closed-incision negative-pressure therapy decreases complications in ventral hernia repair with concurrent panniculectomy. *Hernia*. Online gepubliceerd op 17 december 2018. <https://doi.org/10.1007/s10029-018-1865-2>.
- 14 Willy, C., Agarwal, A., Andersen, C.A., Santis, G.D., Gabriel, A., Grauhan, O., Guerra, O.M., Lipsky, B.A., Malas, M.B., Mathiesen, L.L., Singh, D.P. and Reddy, V.S. (2017), Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. *Int Wound J*, 14: 385-398. doi:10.1111/iwj.12612.

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© 2020 3M. Alle rechten voorbehouden. 3M en de andere weergegeven merken zijn merken of geregistreerde merken. Ongeoorloofd gebruik verboden. PRA-PM-NL-00061 (07/20). OMG102183.