



Filtek™

Supreme
Flowable Restorative



**Here's a tip to virtually
eliminate bubbles**

Key features

- Virtually no bubbles or material “run-on”
- Easy-to-bend cannula improves access
- Ergonomic design is easy to hold and inject
- Good adaptation, polish retention and wear resistance
- Versatile for a wide range of techniques including injection molding



3M's popular flowable is made with TRUE nanotechnology. This patented technology makes the nanoparticles more uniform in size, resulting in excellent wear resistance and polished retention.



A shade-matched palette


3M™ Filtek™ Supreme Flowable Restorative shades are formulated to be an identical match to 3M™ Filtek™ Z350 XT Universal Restorative (A2, A3, OA3, A4).

A2

**3M™ Filtek™ Supreme
Flowable Restorative**

A2B

**3M™ Filtek™ Z350 XT
Universal Restorative**



Warm with confidence

Why warm 3M Composites?

- Warmed composites are more fluid, for easier handling
- Warming may improve adaptation of 3M composite to tooth structure¹
- Warming lowers extrusion force²

158°F*


Get performance with peace of mind.

Physical Properties Unchanged:³

- Diametral tensile strength
- Flexural modulus
- Flexural strength
- Color stability
- Depth of cure

Capsules


3M™ Filtek™ Z350 XT Universal Restorative,
3M™ Filtek™ One Bulk Fill Restorative



Heating Lowers Capsule Extrusion Force by 75-80%

Flowable Syringe

3M™ Filtek™ Supreme Flowable Restorative
3M™ Filtek™ Bulk Fill Flowable Restorative



Heating Lowers Flowable Syringe Extrusion Force by 46-52%

*Recommended heating temperature

1 Based on a 3M sponsored in vitro study.

11 dentists placed 88 Class II MOD restorations. Teeth were microscopically examined for flaws, defects and voids. Comparisons made between techniques and operators.

2 3M Internal Data.

3 When warmed according to our IFU (up to 70°C or 158°F for up to 1 hour).

Tips for success

You can safely warm Filtek Restoratives in as little as 5 minutes! Follow these general guidelines, and refer to the instructions for use for more information.

DO

DO heat flowable syringes up to 70°C (158°F) for up to 1 hour, up to 25 times.



DO consult the Instructions for Use for more information.

DO only use commercially available composite warmer devices.

DON'T

Do **NOT** use heat from a light bulb, microwave, mug warmer or other unapproved heating devices. Consult the Instructions for Use for more information.





Did you know

A bottle of capsule refill* generates

70% MORE plastic waste

for the dentist than a syringe refill**

*Capsule refill of 3M™ Filtek™ Supreme Ultra Flowable (U.S. only) includes: 20 capsules with caps
**Refill of 3M™ Filtek™ Supreme Flowable Restorative includes: 2 syringes, 20 tips, 2 caps

Same great material, new name and syringe

3M™ Filtek™ Supreme Flowable Restorative

Syringe refills include: 2 – 20g syringes;
20 dispensing tips

SHADE	ITEM NO.
A1	7032A1
A2	7032A2
A3	7032A3
A3.5	7032A3.5
A4	7032A4
B1	7032B1
B2	7032B2

ITEM NO.	PRODUCT INFO.
3700T	3M™ Filtek™ Supreme Flowable Restorative Dispensing Tips, 20g, Gray (20)



New syringe for bulk fill flowable

3M™ Filtek™ Bulk Fill Flowable Restorative

Syringe refills include: 2 – 2g syringes; 20 dispensing tips

- Injectable, one-step placement up to 4 mm
- Excellent adaptation and low shrinkage stress
- Can be safely warmed up to 70°C for one hour up to 25 times in commercially available composite warmer devices

SHADE	ITEM NO.
A1	4862A1
A2	4862A2
A3	4862A3
Universal	4862U

ITEM NO.	PRODUCT INFO.
3740T	3M™ Filtek™ Bulk Fill Flowable Restorative Dispensing Tips, 19-gauge, Green (20)



A billion smiles around the world.

As the global leader in restorative and esthetic technology, 3M™ Filtek™ Dental Restoratives has been put to the test in over 1 billion restorations worldwide.

We pioneered the nano-composite category and continue to raise the performance bar to meet restorative challenges around the world.

Simplify your procedures, your inventory and your practice with a single brand – 3M Filtek Dental Restoratives.



3M

3M Malaysia Sdn Bhd
Level 8, Block F, Oasis Square
No.2, Jalan PJU 1A/7A, Ara Damansara,
47301 Petaling Jaya, Selangor, Malaysia

Tel +603 7884 2888
Fax +603 7884 2902
Web www.3M.com.my

This collateral is intended for Healthcare Professionals only, not for general public. Registered under Act 737. MDA registration numbers: GB27042394117, GB35649394017, GB77092433717

3M and Filtek are trademarks of 3M Company. © 2019. All rights reserved.