



Filtek™

Supreme  
Flowable Restorative

Giving you the upper hand  
**with our new syringe design**



# Now our most popular flowable is easier to inject

Designed for improved comfort and injection molding, our 3M™ Filtek™ Supreme Flowable Restorative is now available for expanded indications. You can use this popular, clinically proven flowable for all classes of direct restorations. Plus, take advantage of our updated syringe, with its ergonomic design and finer, bubble-reducing tip.

## Expanded indications

Versatile for an expanded range of indications and techniques, Class I - V.

## Reduce bubbles and waste

Virtually no bubbles or material “run-on” during dispensing gives you better control.

## Easy to hold and inject

Triangle finger plate and plunger make it quick and easy to personalize your grip.

## No more guessing

The blue syringe barrel indicates the remaining material volume. When your syringe is empty, only the white plunger shows.

## Bendable cannula

Easy-to-bend without kinking for better access.

## Good adaptation


Good adaptation, polish retention and wear resistance.

## Versatile

A wide range of indications and techniques.



3M's popular flowable is made with TRUE nanotechnology. This patented technology makes the nanoparticles more uniform in size, resulting in excellent wear resistance and polished retention.






# Warm with confidence

## Why warm 3M Composites?

- Warmed composites are more fluid, for easier handling
- Warming may improve adaptation of 3M composite to tooth structure<sup>1</sup>
- Warming lowers extrusion force<sup>2</sup>

158°F\*



## Get performance with peace of mind

Physical Properties Unchanged:<sup>3</sup>

- Diametral tensile strength
- Flexural modulus
- Flexural strength
- Color stability
- Depth of cure



## Tips for success

You can safely warm Filtek Restoratives in as little as 5 minutes! Follow these general guidelines, and refer to the instructions for use for more information.

### DO

**DO** heat flowable syringes up to 70°C (158°F) for up to 1 hour, up to 25 times.

**DO** only use commercially available composite warmer devices.

**DO** consult the Instructions for Use for more information.

### DO NOT

**Do NOT** use heat from a light bulb, microwave, mug warmer or other unapproved heating devices. Consult the Instructions for Use for more information.

\*Recommended heating temperature

<sup>1</sup> Based on a 3M sponsored in vitro study. 11 dentists placed 88 Class II MOD restorations. Teeth were microscopically exemplified for flaws, defects and voids. Comparisons made between techniques and operators.

<sup>2</sup> 3M Internal Data

<sup>3</sup> When warmed according to our IFU (up to 70°C or 158°F for up to 1 hour).



# Clinical cases

3M™ Filtek™ Supreme Flowable Restorative offers a smooth transition from material to tooth for esthetic restorations.



Before




Flowable placement




After

Photos courtesy of Dr. Gunnar Reich, Munich, Germany

Filtek Supreme Flowable Restorative adapts easily and seals the margin.



Before



Flowable placement



After

Photos courtesy of Dr. Giuseppe Chiodera, Brescia, Italy

**Expanded Indications**  
(Incremental up to 2 mm)

Class I    Class II    Class III    Class IV

Class V    Liner    Repair of Indirect (small defects)

Repair<sup>1</sup>    Pit & Fissure    Base

Refer to Instructions for Use for full list of indications.  
1. Repair of resin and acrylic temporary materials

## Ordering information



**3M™ Filtek™ Supreme Flowable Restorative Refills**  
Syringe refills include: 2 – 2g syringes;  
20 dispensing tips.

| SHADE | ITEM NO. |
|-------|----------|
| A1    | 7032A1   |
| A2    | 7032A2   |
| A3    | 7032A3   |
| A3.5  | 7032A3.5 |
| A4    | 7032A4   |
| B1    | 7032B1   |
| B2    | 7032B2   |

| ITEM NO. | PRODUCT INFO.  |
|----------|--|
| 3700T    | 3M™ Filtek™ Supreme Flowable Restorative Dispensing Tips, 20g, Gray (20) |



**3M Malaysia Sdn Bhd** (7251 V)  
Level 8, Block F, Oasis Square  
No. 2, Jalan PJU 1A/7A, Ara Damansara  
47301 Petaling Jaya, Selangor, Malaysia  
Tel: +603 7884 2888  
Fax: +603 7884 2902  
www.3M.com.my

3M and Filtek are trademarks of 3M.

This collateral is intended for Healthcare Professionals only, not for general public.  
Registered under Act 737.  
MDA registration number:  
GB27042394117