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Conventional cementation possible?	Pre-treatment of material for adhesive bonding	Posterior crown	Fabrication effort / cost	Anterior crown	Fabrication effort / cost
------------------------------------	------------------------------------------------	-----------------	---------------------------	----------------	---------------------------

no	acid-etch with 5 % hydrofluoric acid	+/-	*****	+	*****
yes	acid-etch with 5 % hydrofluoric acid	+	****	+	****

yes	sandblasting with alumina ¹ (grain size max. 30-50 µm, max. 2 bar)	+	***	+/-	***
yes	sandblasting with alumina ¹ (grain size max. 30-50 µm, max. 2 bar)	+	***	hand-veneered	*****
yes	sandblasting with alumina ¹ (grain size max. 30-50 µm, max. 2 bar)	hand-veneered	*****	hand-veneered	*****

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


3. Treatment result with a veneer made of feldspathic ceramic and a handveneered crown with coping made of 3MTM LavaTM Frame Zirconia.

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3M™ Lava™ Esthetic
Fluorescent Full-Contour
Zirconia* (5Y-TZP)
(5Y-TZP)

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Beier US, Kapferer I, Burtscher D, Dumfahrt H. Clinical Performance of Porcelain Laminate Veneers for Up to 20 Years. *Int J Prosthodont* 2012; 25: 79-85.

Burke FJ, Murray MC, Shortall AC. Trends in indirect dentistry: 6. Provisional restorations, more than just a temporary. *Dent Update.* 2005 Oct; 32(8): 443-4, 447-8, 450-2.

Burns DR, Beck DA, Nelson SK; Committee on Research in Fixed Prosthodontics of the Academy of Fixed Prosthodontics. A review of selected dental literature on contemporary provisional fixed prosthodontic treatment: report of the Committee on Research in Fixed Prosthodontics of the Academy of Fixed Prosthodontics. *J Prosthet Dent.* 2003 Nov; 90(5): 474-97.

Castelnuovo J, Tjan AH, Phillips K, Nicholls JI, Kois JC. Fracture load and mode of failure of ceramic veneers with different preparations. *J Prosthet Dent.* 2000; 83(2): 171-80.

Christensen GJ. The state of fixed prosthodontic impressions: room for improvement. *J Am Dent Assoc.* 2005 Mar; 136(3): 343-6.

Cobb CM, Daubert DM, Davis K, Deming J, Flemmig TF, Pattison A, Roulet JF, Stambaugh RV. Consensus Conference Findings on Supragingival and Subgingival Air Polishing. *Compend Contin Educ Dent.* 2017 Feb; 38(2): e1-e4.

Das AK, Muddugangadhar BC, Amarnath GS, Garg A, Kumar U, Rao TR. Comparative Evaluation of Push Out Bond Strength of a Fiber Post System using Four Different Resin Cements: An In-Vitro Study. *J Int Oral Health.* 2015; 7(Suppl 1): 62-7.

Dias WR, Pereira PN, Swift EJ Jr. Effect of bur type on microtensile bond strengths of self-etching systems to human dentin. *J Adhes Dent.* 2004 Autumn; 6(3): 195-203.

Dogan S, Raigrodski AJ, Zhang H, Mancl LA. Prospective cohort clinical study assessing the 5-year survival and success of anterior maxillary zirconia-based crowns with customized zirconia copings. *J Prosthet Dent.* 2017 Feb; 117(2): 226-232. doi: 10.1016 / j. prosdent.2016.07.019. Epub 2016 Oct 17.

Donovan TE, Chee WW. A review of contemporary impression materials and techniques. *Dent Clin North Am.* 2004 Apr; 48(2): vi-vii, 445-70.

Ebrahimi Chaharom ME, Ajami AA, Bahari M, Rezazadeh H. Effect of smear layer thickness and pH of self-adhesive resin cements on the shear bond strength to dentin. *Indian J Dent Res.* 2017 Nov-Dec; 28(6): 681-686. doi: 10.4103 / ijdr.IJDR_12_16.

Edelhoff D, Liebermann A, Beuer F, Stimmelmayer M, Güth JF. Minimally invasive treatment options in fixed prosthodontics. *Quintessence Int.* 2016 Mar; 47(3): 207-16. doi: 10.3290 / j. qi.a35115.

Figueiredo FE, Martins-Filho PR, Faria-E-Silva AL. Do metal post-retained restorations result in more root fractures than fiber post-retained restorations? A systematic review and meta-analysis. *J Endod*. 2015 Mar; 41(3): 309-16. doi: 10.1016/j.joen.2014.10.006. Epub 2014 Nov 11.

Fradeani M, Redemagni M, Corrado M. Porcelain Laminate Veneers: 6-to-12 Year Clinical Evaluation. A Retrospective Study. *Int J Perio Rest Dent* 2005; 25: 9-17.

Galal RM, Omar N, Nabil H, Aly Y. Efficacy of Fiber Post Bonding To Root Dentin after Different Obturation Techniques and Cementation Timings: In Vitro Study. *Open Access Maced J Med Sci*. 2018 Sep 21; 6(9): 1707-1711. doi: 10.3889/oamjms.2018.343. eCollection 2018 Sep 25.

Galvão Ribeiro BR, Galvão Rabelo Caldas MR, Almeida AA Jr, Fonseca RG, Adabo GL. Effect of surface treatments on repair with composite resin of a partially monoclinic phase transformed yttrium-stabilized tetragonal zirconia. *J Prosthet Dent*. 2018 Feb; 119(2): 286-291. doi: 10.1016/j.prosdent.2017.02.014. Epub 2017 May 20.

Gratton DG, Aquilino SA. Interim restorations. *Dent Clin North Am*. 2004 Apr; 48(2): vii, 487-97.

Gresnigt M, Magne M, Magne P. Porcelain veneer post-bonding crack repair by resin infiltration. *Int J Esthet Dent* 2017; 12(2): 156-170.

Huang C, Somar M, Li K, Mohadeb JV. To cord or not to cord? That is still a question. *Evidence-Based Dentistry* (2017) 18: 21-2.

Huang C, Somar M, Li K, Mohadeb JVN. Efficiency of Cordless Versus Cord Techniques of Gingival retraction: A Systematic Review. *J Prosthodont* 2017 Apr; 26(3): 177-85.

Jorgensen KD. The relationship between retention and convergence angle in cemented veneer crowns. *Acta Odontol Scand* 1956;13:35-40.

Keul C, Köhler P, Hampe R, Roos M, Stawarczyk B. Glass Fiber Post / Composite Core Systems Bonded to Human Dentin: Analysis of Tensile Load vs Calculated Tensile Strength of Various Systems Using Pull-out Tests. *J Adhes Dent*. 2016;18(3):247-56. doi: 10.3290/j.jad.a36136.

Kurtzman GM, Strassler HE. Provisional fixed restorations. *Dental Economics*. 2006; 3 (Suppl): 1-12.

Langeland K, Langeland LK. Pulp reactions to crown preparation, impression, temporary crown fixation, and permanent cementation. *J Prosthet Dent*. 1965 Jan-Feb;15: 129-43.

Laxe L, Marchiori RH, De Goes MF and Sabrosa CE. Bond Strength of Different Cements to a Resin-nano-ceramic CAD-CAM Material. *J Dent Res* 93 (Spec Iss A): 1134, 2014.

Laxe L, Salina L, Sartori BT, Silva PA, Possidonio L, Marchiori RH, Sabrosa CE. Light energy transmission through various shades of a CAD-CAM Material. *J Dent Res* 94 (Spec Iss A): 2220, 2015.

Lee H, So JS, Hochstedler JL, Ercoli C. The accuracy of implant impressions: a systematic review. *J Prosthet Dent*. 2008 Oct;100(4): 285-91. doi: 10.1016/S0022-3913(08)60208-5.

Lin J, Matinlinna JP, Shinya A, Botelho MG, Zheng Z. Effect of fiber post length and abutment height on fracture resistance of endodontically treated premolars prepared for zirconia crowns. *Odontology*. 2018 Apr; 106(2): 215-222. doi: 10.1007/s10266-017-0320-7. Epub 2017 Dec 14.

Ma S, Fenton A. Screw versus cement-retained implant prostheses: a systematic review of prosthodontic maintenance and complications. *Int J Prosthodont*. 2015 Mar-Apr; 28(2): 127-45.

Magne P; Belser U. Bonded porcelain restorations in the anterior dentition: a biomimetic approach. Quintessence Pub. Co., 2002.

Mamoun J. Post and core build-ups in crown and bridge abutments: Bio-mechanical advantages and disadvantages. *J Adv Prosthodont*. 2017 Jun; 9(3): 232-237. doi: 10.4047/jap.2017.9.3.232. Epub 2017 Jun 19.

Michalakakis KX, Bakopoulou A, Hirayama H, Garefis DP, Garefis PD. Pre- and post-set hydrophilicity of elastomeric impression materials. *J Prosthodont*. 2007 Jul-Aug; 16(4): 238-48. Epub 2007 Jun 9.

Miragaya L, Vasconcelos L, Sabrosa CE. Hydrophilicity of unset impression materials. *J Dent Res* 90 (Spec Iss A): 3280, 2011. Miragaya LM, Maia LC, Sabrosa CE, Goes MF, Silva EM. Evaluation of self-adhesive cement bond strength to yttria-stabilized zirconia ceramic (Y-TZP) using four surface treatments. *J Adhes Dent*, 2011; 13(5): 473-80.

Novaes SA, Laxe LAC, Marchiori RH, Sartori BT and Sabrosa CE. Light energy transmission through various thicknesses of a CAD-CAM ceramic Material. *J Dent Res* 93 (Spec Iss B): 913, 2014.

Oliveira SS, Pugach MK, Hilton JF, Watanabe LG, Marshall SJ, Marshall GW Jr. The influence of the dentin smear layer on adhesion: a self-etching primer vs. a total-etch system. *Dent Mater*. 2003 Dec;19(8):758-67.

Ostlund LE. Cavity design and mathematics: Their effect on gaps at the margins of cast restorations. *Oper Dent* 1985; 10: 122-37.

Pascoe DF. Analysis of the geometry of finishing lines for full crown restorations. *J Prosthet Dent* 1978;40:157-62.

Peumans M, De Munck J, Fieuws S, Lambrechts P, Vanherle G, Van Meerbeek B. A Prospective Ten-year Clinical Trial of Porcelain Veneers. *J Adhes Dent*, 2004; 6: 65-75.

Piovesan EM, Demarco FF, Cenci MS, Pereira-Cenci T. Survival rates of endodontically treated teeth restored with fiber-reinforced custom posts and cores: a 97-month study. *Int J Prosthodont*. 2007 Nov-Dec; 20(6): 633-9.

Podhorsky A, Rehmann P, Wöstmann B. Tooth preparation for full-coverage restorations-a literature review. *Clin Oral Investig*. 2015 Jun; 19(5): 959-68.

Poggio CE, Ercoli C, Rispoli L, Maiorana C, Esposito M. Metal-free materials for fixed prosthodontic restorations. *Cochrane Database Syst Rev*. 2017 Dec 20;12:CD009606. doi: 10.1002 / 14651858.CD009606.pub2. Review.

Powers JM, Sakaguchi RL. Impression materials. *Craig's restorative dental materials*. Elsevier Mosby, 2018

Radz GM. Minimum thickness anterior porcelain restorations. *Dent Clin N Am*, 2011;55:353-370.

Rosner D. Function, placement and reproduction of bevels for gold castings. *J Prosthet Dent* 1963; 13: 1160-6.


Sabrosa CE, Morgano SM. Effect of different rotary instruments on the surface of prepared dentin and enamel. ACP Meeting 1997.


Sabrosa CE, Sartori BT, Andrade P, Salina L, Possidonio L, Machado KC. Long-term stability of bisacrylic-composite crowns fabricated chairside after 36 months. *J Dent Res* 93(Spec Iss B): 934, 2014.

Sabrosa CE, Sartori BT, Silva PA, Possidonio L, Rocha Jr MA, Ferreira K, Felix C, Deacon C. Light transmission through anterior teeth in vivo. *J Dent Res* 94 (Spec Iss A): 2219, 2015.

Sabrosa CE, Miragaya L, Nascimento R, Andreiuolo R, Vasconcelos L, Alves L, DeGoes MF. Flowability on different working time of polyether and VPS. *J Dent Res* 89 (Spec Iss A): 238, 2010.

- Sailer I, Balmer M, Hüsler J, Hämmerle CHF, Känel S, Thoma DS. 10-year randomized trial (RCT) of zirconia-ceramic and metal-ceramic fixed dental prostheses. *J Dent*. 2018 Sep; 76: 32-39. doi: 10.1016 / j.jdent.2018.05.015. Epub 2018 May 25.
- Sartori BT, Andrade P, Marchiori RH, Felix C and Sabrosa CE. Irradiance quantification from three LED LCUs at various distances. *J Dent Res* 93 (Spec Iss A): 318, 2014.
- Shiratori FK, Valle AL, Pegoraro TA, Carvalho RM, Pereira JR. Influence of technique and manipulation on self-adhesive resin cements used to cement intraradicular posts. *J Prosthet Dent*. 2013 Jul; 110(1): 56-60. doi: 10.1016 / S0022-3913(13)60341-8.
- Skurrow HM, Nevins M. The rationale of the preperiodontal provisional biologic trial restoration. *Int J Periodontics Restorative Dent*. 1988; 8(1): 8-29.
- Stawarczyk B, Liebermann A, Eichberger M, Güth JF. Evaluation of mechanical and optical behavior of current esthetic dental restorative CAD/CAM composites. *J Mech Behav Biomed Mater*. 2015 Mar; 55:1-11. doi: 10.1016 / j.jmbbm.2015.10.004. Epub 2015 Oct 19.
- Stewardson DA. Trends in indirect dentistry: 5. Impression materials and techniques. *Dent Update*. 2005 Sep; 32(7): 374-6, 379-80, 382-4 passim.
- Tabassum S, Adnan S, Khan FR. Gingival Retraction Methods: A Systematic Review. *J Prosthodont*. 2017 Dec; 26(8): 637-643. doi: 10.1111 / jopr.12522. Epub 2016 Jul 28.
- Tamura Y, Takamizawa T, Shimamura Y, Akiba S, Yabuki C, Imai A, Tsujimoto A, Kurokawa H, Miyazaki M. Influence of air-powder polishing on bond strength and surface-free energy of universal adhesive systems. *Dent Mater J*. 2017 Nov 29; 36(6): 762-769. doi: 10.4012/dmj.2016-185. Epub 2017 Jul 12.
- Teichmann M, Wienert AL, Rückbeil M, Weber V, Wolfart S, Edelhoff D. Ten-year survival and chipping rates and clinical quality grading of zirconia-based fixed dental prostheses. *Clin Oral Investig*. 2018 Nov; 22(8): 2905-2915. doi: 10.1007 / s00784-018-2378-1. Epub 2018 Mar 8.
- Tiu J, Al-Amleh B, Waddell JN, Duncan WJ. Clinical tooth preparations and associated measuring methods: A systematic review. *J Prosthet Dent* 2015 Mar; 113(3): 175-84
- Tjan AH, Castelnovo J, Shiotsu G. Marginal fidelity of crowns fabricated from six proprietary provisional materials. *J Prosthet Dent*. 1997 May; 77(5): 482-5.
- Ubalini ALM, Benetti AR, Sato F, Pascotto RC, Medina Neto A, Baesso ML, Peutzfeldt A. Challenges in luting fibre posts: Adhesion to the post and to the dentine. *Dent Mater*. 2018 Jul; 34(7): 1054-1062. doi: 10.1016 / j.dental.2018.04.001. Epub 2018 May 1.
- Vasconcelos L, Miragaya LM, Maia LC, Al-Harbi FA, Sabrosa CE. Flexural strength of resins used to fabricate provisional restorations. *J Dent Res* 90 (Spec Iss A): 2006, 2011.
- 3MTM ESPETM EspertiseTM Publication. 3rd Edition (2008). Impressioning Compendium - A Guideline for Excellent Impressions in Theory and Practice.

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