3M™ Comply™ 1251 EO Chemical Indicator Strip

Product Description

The 3M™ Comply™ 1251 EO Chemical Indicator Strip is a paper strip 1,5 cm wide by 20 cm long (0.625 inches wide by 8 inches long) printed with a chemical indicator ink that turns from red to green when exposed to the critical parameters of the ethylene oxide (EO) sterilization process (i.e., gas concentration, RH, time and temperature). It is designed to indicate whether EO sterilization conditions have been met at the point of placement of the strip, which should be placed in the area determined to be the least accessible to EO penetration. Each strip is perforated at the center if a shorter internal chemical indicator is preferred. An internal chemical indicator should be used inside each package to be sterilized.

Indications for Use

The Comply™ 1251 EO Chemical Indicator Strip is designed to monitor all warm and cool cycles of 100% EO and EO/HCFC (hydrochlorofluorocarbon) gas mixture sterilizer cycles.

Contraindications

None.

Precautions

Do not use Comply™ 1251 EO Chemical Indicator Strip to monitor steam, dry heat, EO/carbon dioxide mixtures or other low temperature sterilization processes.

Directions for Use

- Carefully open package. Remove only the number of 3M™ Comply™ EO Chemical Indicators 1251 needed. Reseal the package.
- Place a Comply™ 1251 EO Chemical Indicator Strip inside each pack, peel pouch, container system, or tray to be EO sterilized in the area determined to be the least accessible to EO penetration.
- After processing, the red chemical indicator line will turn green when exposed to ethylene oxide.

Storage

Best stored in a dry (<50% relative humidity) condition at room temperature (15-30°C [59-86°F]) and protect from light. Do not store near strong alkaline or acidic products such as cleaning/disinfection agents.

Explanation of Symbols

 Caution, see instructions for use

 Do not reuse

 Use by date

 Batch code

 Manufacturer

 Date of manufacture

 Product is designed for use with ethylene oxide sterilization cycles.

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3M and Comply are trademarks of 3M.
© 2018, 3M. All rights reserved.