Prevena^{MC}
Thérapie de l'incision

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Citation	Type de plaie/ d'intervention chirurgicale	Niveau d'évidence clinique*
Higuera-Rueda, C.; Emara, A.K.; Nieves-Malloure, Y.; Klika, A.K.; Cooper, H.J.; Cross, M.B.; Guild, G.N.; Nam, D.; Nett, M.; Scuderi, G.R.; Cushner, F.D.; Piuze, N.S.; Silverman, R.P. « The Effectiveness of Closed Incision Negative Pressure Therapy versus Silver-Impregnated Dressings in Mitigating Surgical Site Complications in High-Risk Patients after Revision Knee Arthroplasty: The PROMISES Randomized Controlled Trial. » <i>J. Arthroplasty</i> (2021), doi : https://doi.org/10.1016/j.arth.2021.02.076	Révision d'arthroplastie totale du genou	1 ●
Newman, J.M.; Siqueira, M.B.P.; Klika, A.K.; Molloy, R.M.; Barsoum, W.K.; Higuera, C.A. Use of Closed Incisional Negative Pressure Wound Therapy After Revision Total Hip and Knee Arthroplasty in Patients at High Risk for Infection: A Prospective, Randomized Clinical Trial. <i>Journal of Arthroplasty</i> . 2018 Nov 17. [Publication électronique avant impression]	Arthroplastie totale du genou et de la hanche	1b ●
Crist, B.D.; Oladeji, L.O.; Khazzam, M.; Della Rocca, G.J.; Murtha, Y.M.; Stannard, J.P. Role of acute negative pressure wound therapy over primarily closed surgical incisions in acetabular fracture ORIF: A prospective randomized trial. <i>Injury</i> . 2017 April 27.pii: S0020-1383(17)30283-8.	Fractures acétabulaires	1b ●
Pauser, J.; Nordmeyer, M.; Biber, R.; Jantsch, J.; Kopschina, C.; Bail, H.J.; Brem, M.H. Incisional negative pressure wound therapy after hemiarthroplasty for femoral neck fractures - reduction of wound complications. <i>International Wound Journal</i> . 2014; 13(5): 663-667.	Hémiarthroplastie pour les fractures du col du fémur	1b ●
Manoharan, V.; Grant, A.; Harris, A.; Hazratwala, K.; Wilkinson, M.; McEwen, P. Closed Incision Negative Pressure Wound Therapy vs Conventional Dry Dressings After Primary Knee Arthroplasty: A Randomized Controlled Study. <i>J Arthroplasty</i> . 2016 April 28. pii: S0883-5403 (16)30083-3.	Arthroplastie du genou	1b ●
Howell, R.D.; Hadley, S.; Strauss, E.; Pelham, F.R. Blister formation with negative pressure dressings after total knee replacement. <i>Current Orthopaedic Practice</i> . 2011 March; 22(2): 176-179.	Arthroplastie du genou	1b ●
Stannard, J.P.; Robinson, J.T.; Anderson, E.R.; McGwin, G. Jr.; Volgas, D.A.; Alonso, J.E. Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma. <i>Journal of Trauma</i> . 2006 June; 60(6): 1301-6.	Fracture des extrémités inférieures	1b ●
Stannard, J.P.; Volgas, D.A.; McGwin, G.; Stewart, R.L.; Obremesky, W.; Moore, T.; Anglen, J.O. Incisional negative pressure wound therapy after high-risk lower extremity fractures. <i>Journal of Orthopedic Trauma</i> . 2012 January; 26(1): 37-42.	Fracture des extrémités inférieures	1b ●
Stannard, J.P.; Volgas, D.A.; Stewart, R.; McGwin, G. Jr.; Alonso, J.E. Negative pressure wound therapy after severe open fractures: a prospective randomized study. <i>Journal of Orthopedic Trauma</i> . 2009 September; 23(8): 552-7.	Fracture des extrémités inférieures	1b ●
Pachowsky, M.; Gusinde, J.; Klein, A.; Lehl, S.; Schulz-Drost, S.; Schlechtweg, P.; Pauser, J.; Gelse, K.; Brem, M.H. Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty. <i>International Orthopaedics</i> . 2012 April; 36(4): 719-22.	Arthroplastie totale de la hanche	1b ●
Redfern, R.E.; Cameron-Ruetz, C.; O'Drobinak, S.; Chen, J.; Beer, K.J. Closed incision negative pressure therapy effects on postoperative infection and surgical site complication after total hip and knee arthroplasty. <i>J. Arthroplasty</i> 2017 November; 32(11): 3333-3339.	Arthroplastie du genou et de la hanche	2 ●
Reddix, R.N. Jr.; Leng, X.I.; Woodall, J.; Jackson, B.; Dedmond, B.; Webb, L.X. The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery. <i>Journal of Surgical Orthopaedic Advances</i> . 2010 June; 19(2): 91-7.	Arthroplastie de la hanche	3 ●
Cooper, H.J.; Roc, G.C.; Bas, M.A.; Berliner, Z.P.; Hepinstall, M.S.; Rodriguez, J.A.; Weiner, L.S. Closed incision negative pressure therapy decreases complications after periprosthetic fracture surgery around the hip and knee. <i>Injury</i> . 2018 February; 49(2): 386-391. doi: 10.1016/j.injury.2017.11.010. Publication électronique Nov 2017 14.	Chirurgie de fracture périprothétique	3 ●
Cooper, H.J.; Bas, M.A. « Closed-Incision Negative-Pressure Therapy Versus Antimicrobial Dressings After Revision Hip and Knee Surgery: A Comparative Study. <i>J. Arthroplasty</i> . 2016 May; 31(5): 1047-52.	Révision du genou et de la hanche	3 ●
Anatone, A.J.; Shah, R.P.; Jennings, E.L.; Geller, J.A.; Cooper, J. A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty. <i>Arthroplasty Today</i> . 2018; 4(4): 493-498. doi:10.1016/j.artd.2018.09.004.	Arthroplastie du genou et de la hanche	3 ●
Curley, A.J.; Terhune, E.B.; Velott, A.T.; Argintar, E.H. Outcomes of Prophylactic Negative Pressure Wound Therapy in Knee Arthroplasty. <i>Orthopedics</i> . 2018; 41(6): e837-e840. doi:10.3928/01477447-20181010-02.	Arthroplastie du genou	3 ●
Reddix, R.N.; Tyler, H.K.; Kulp, B.; Webb, L.X. Incisional vacuum-assisted wound closure in morbidly obese patients undergoing acetabular fracture surgery. <i>The American Journal of Orthopedics</i> . Sep 2009; 38(9) : 32-5.	Fractures acétabulaires	4 ●
Hansen, E.; Durinka, J.B.; Costanzo, J.A.; Austin, M.S.; Deirmengian, G.K. Negative pressure wound therapy is associated with resolution of incisional drainage in most wounds after hip arthroplasty. <i>Clinical Orthopaedics and Related Research</i> . Oct 2013; 471(10) : 3230-6.	Arthroplastie de la hanche	4 ●

Citation	Type de plaie/ d'intervention chirurgicale	Niveau d'évidence clinique*
Stannard, J.P.; Atkins, B.Z.; O-Malley, D.; Singh, H.; Bernstein, B.; Fahey, M.; Masden, D.; Attinger, C.E. Use of negative pressure therapy on closed surgical incisions: A case series. <i>Ostomy Wound Management</i> . 2009 August; 55(8): 58-66.	Fractures des extrémités inférieures	4 ●
Gomoll, A.H.; Lin, A.; Harris, M.B. Incisional vacuum-assisted closure therapy. <i>Journal of Orthopaedic Trauma</i> . 2006 November-December; 20(10): 705-9.	Traumatisme orthopédique	4 ●
Stannard, J.P.; Gabriel, A.; Lehner, B. Use of Negative Pressure Wound Therapy Over Clean, Closed Surgical Incisions. <i>International Wound Journal</i> . 2012; 9:32-39.	Traumatisme orthopédique	4 ●
Brem, M.H.; Bail, H.J.; Biber, R. Value of Incisional Negative Pressure Wound Therapy in Orthopedic Surgery. <i>International Wound Journal</i> . 2014 June; 11(Suppl 1): 3-5.	Mixte	5 ●
Berkowitz, M.J. Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures. <i>Techniques in Foot and Ankle Surgery</i> . 2013 December; 12(4): 172-174.	Fractures du calcaneum	5 ●
Suleiman, L.I.; Mesko, D.R.; Nam, D. Intraoperative Considerations for Treatment/Prevention of Prosthetic Joint Infection. <i>Current Reviews in Musculoskeletal Medicine</i> . 2018:1-8.	Arthroplastie du genou et de la hanche	5 ●
Chotanaphuti, T.; Courtney, P.M.; Fram, B.; Kleef, N.J.; Kim, T.K.; Kuo, F.C.; Lustig, S.; Moojen, D.J.; Nijhof, M.; Oliashirazi, A.; Poolman, R.; Purtill, J.J.; Rapisarda, A.; Rivero-Boschert, S.; Veltman, E.S. Hip and Knee Section, Treatment, Algorithm: Proceedings of International Consensus on Orthopedic Infections. <i>The Journal of Arthroplasty</i> . 34(2S):S393-S397. doi: 10.1016/j.arth.2018.09.024.	Arthroplastie du genou et de la hanche	5 ●
DeCarbo, W.T.; Hyer, C.F. Negative-Pressure Wound Therapy Applied to High-Risk Surgical Incisions. <i>Journal of Foot and Ankle Surgery</i> . 2010 May; 49(3): 299-300.	Traumatisme orthopédique	5 ●
Nam, D.; Sershon, R.A.; Levine, B.R.; Della Valle, C.J. The Use of Closed Incision Negative-Pressure Wound Therapy in Orthopaedic Surgery. <i>J. Am. Acad. Orthop. Surg</i> . 2018:1-8. doi: 10.5435/JAAOS-D-17-00054.	Chirurgie orthopédique	5 ●
Al-Hourabi, R.K.; Aalirezaie, A.; Adib, F.; Anoushiravani, A.; Bhashyam, A.; Binlaksar, R.; Blevins, K.; Bonanzinga, T.; Chih-Kuo, F.; Cordova, M.; Deirmengian, G.K.; Fillingham, Y.; Frenkel, T.; Gomez, J.; Gundtoft, P.; Harris, M.A.; Harris, M.; Heller, S.; Jennings, J.A.; Jimenez-Garrido, C.; Karam, J.A.; Khlopas, A.; Klement, M.R.; Komnos, G.; Krebs, V.; Lachiewacz, P.; Miller, A.O.; Mont, M.A.; Montanez, E.; Romero, C.A.; Schwarzkopf, R.; Shaffer, A.; Sharkey, P.F.; Smith, B.M.; Sodhi, N.; Thienpont, E.; Villanueva, A.O.; Yazdi, H. General Assembly, Prevention, Wound Management: Proceedings of International Consensus on Orthopedic Infections. <i>The Journal of Arthroplasty</i> . 2019;34(2):S157-S168. doi:10.1016/j.arth.2018.09.066.	Infections orthopédiques	5 ●
Agarwal, A. Management of Closed Incisions Using Negative-Pressure Wound Therapy in Orthopedic Surgery. <i>Plastic and reconstructive surgery</i> . 2019;143(1 Management of Surgical Incisions Utilizing Closed Incision Negative Pressure Therapy): 21S-26S.	Chirurgie de traumatisme orthopédique	5 ●

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Références :

- ¹ De Lissovoy, G.; Fraeman, K.; Hutchins, V.; Murphy, D.; Song, D.; Vaughn, B.B. Surgical site infection: incidence and impact on hospital utilization and treatment costs. *Am. J. Infect. Control.* 2009 June; 37(5): 387-97.
- ² Merkow, R.; et al. Underlying reasons associate with hospital readmission following surgery in the US. *JAMA.* 2015; 313(5): 483-95.
- ³ Kurtz, S.M.; Lau, E.; Watson, H.; Schmier, J.K.; Parvizi, J. Economic burden of periprosthetic joint infection in the United States. *J. Arthroplasty.* 2012 September; 27(8 Suppl.): 61-5.e1. doi: 10.1016/j.arth.2012.02.022. Publication électronique 2 mai 2012.
- ⁴ Thakore, R.V.; et al. « Surgical site infection in orthopedic trauma: A case-control study evaluating risk factors and cost. » *Journal of Clinical Orthopaedics and Trauma.* 2015; (6): 220-226.
- ⁵ Wilkes, R.P.; Kilpadi, D.V.; Zhao, Y.; et al. « Closed Incision Management With Negative Pressure Wound Therapy (CIM): Biomechanics. *Surgical Innovation.* 2012; 19(1): 67-75.
- ⁶ Kilpadi, D.V.; Cunningham, M.R. Evaluation of Closed Incision Management with Negative Pressure Wound Therapy (CIM): Hematoma/Seroma and Involvement of the Lymphatic System. » *Wound Repair and Regeneration.* 2011; 19: 588-596.
- ⁷ Glase, D.A.; Farnsworth, C.L.; Varley, E.S.; et al. Negative pressure therapy for closed spine incisions: A pilot study. *Wounds.* 2012; 24(11): 308-316.
- ⁸ Anatone, A.J.; Shah, R.P.; Jennings, E.L.; Geller, J.A.; Cooper, J. A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty. *Arthroplasty Today.* 2018; 4(4): 493-498. doi:10.1016/j. artd.2018.09.00.
- ⁹ Willy, C.; Agarwal, A.; Andersen, C.A.; Santis, G.; Gabriel, A.; Grauhan, O.; Guerra, O.M.; Lipsky, B.A.; Malas, M.B.; Mathiesen, L.L.; Singh, D.P. Reddy VS. Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. *Int. Wound J.* 2017 April; 14(2): 385-398. doi: 10.1111/iwj.12612. Publication électronique 2016 May 12. PMID: 27170231.
- ¹⁰ Higuera-Rueda, C.; Emara, A.K.; Nieves-Malloure, Y.; Klika, A.K.; Cooper, H.J.; Cross, M.B.; Guild, G.N.; Nam, D.; Nett, M.; Scuderi, G.R.; Cushner, F.D.; Piuze, N.S.; Silverman, R.P. The Effectiveness of Closed Incision Negative Pressure Therapy versus Silver-Impregnated Dressings in Mitigating Surgical Site Complications in High-Risk Patients after Revision Knee Arthroplasty: The PROMISES Randomized Controlled Trial. *J. Arthroplasty* (2021), doi : <https://doi.org/10.1016/j.arth.2021.02.076>
- ¹¹ Sullivan, D.; Chung, K.C.; Eaves, F.F.; Rohrich, R.J. The Level of Evidence Pyramid: Indicating Levels of Evidence in Plastic and Reconstructive Surgery Articles. *Plast. Reconstr. Surg.* 2011; 128(1): 311-314.

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© 2021, 3M. Tous droits réservés. 3M et les autres marques indiquées sont des marques et/ou des marques déposées. Toute utilisation non autorisée est interdite. Utilisées sous licence au Canada. PRA-PM-CA-00223 (05/21) 2105-20672 E
KCI USA, Inc., une société 3M
KCI Medical Canada Inc.