3M

PREVENA™

Incision Management System

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
引用	切開・術式	臨床エビデンス・レベル*
Newman JM, Siqueira MBP, Klika AK, Molloy RM, Barsoum WK, Higuera CA. Use of Closed Incisional Negative Pressure Wound Therapy After Revision Total Hip and Knee Arthroplasty in Patients at High Risk for Infection: A Prospective, Randomized Clinical Trial. <i>Journal of Arthroplasty</i> . 2018 Nov 17. [Epub Ahead of Print]†	人工股関節・人工膝関節置換術	1b
Crist BD, Oladeji LO, Khazzam M, Della Rocca GJ, Murtha YM, Stannard JP. Role of acute negative pressure wound therapy over primarily closed surgical incisions in acetabular fracture ORIF: A prospective randomized trial. <i>Injury</i> . 2017 Apr 27;pii: S0020-1383(17)30283-8.	臼蓋骨折	1b
Pauser J, Nordmeyer M, Biber R, Jantsch J, Kopschina C, Bail HJ, Brem MH. Incisional negative pressure wound therapy after hemiarthroplasty for femoral neck fractures - reduction of wound complications. <i>International Wound Journal</i> . 2014;13(5):663-667.	大腿骨頸部骨折に対する大腿骨頭置換術	1b
Manoharan V, Grant A, Harris A, Hazratwala K, Wilkinson M, McEwen P. Closed Incision Negative Pressure Wound Therapy vs Conventional Dry Dressings After Primary Knee Arthroplasty: A Randomized Controlled Study. <i>J Arthroplasty</i> . 2016 Apr 28. pii: S0883-5403(16)30083-3.	人工膝関節置換術	1b
Howell RD, Hadley S, Strauss E, Pelham FR. Blister formation with negative pressure dressings after total knee replacement. <i>Current Orthopaedic Practice</i> . 2011 Mar;22(2):176-179.	人工膝関節置換術	1b
Stannard JP, Robinson JT, Anderson ER, McGwin G Jr, Volgas DA, Alonso JE. Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma. <i>Journal of Trauma</i> . 2006 Jun;60(6):1301-6.	下肢の骨折	1b
Stannard JP, Volgas DA, McGwin G, Stewart RL, Obremsky W, Moore T, Anglen JO. Incisional negative pressure wound therapy after high-risk lower extremity fractures. <i>Journal of Orthopedic Trauma</i> . 2012 Jan;26(1):37-42.	下肢の骨折	1b
Stannard JP, Volgas DA, Stewart R, McGwin G Jr, Alonso JE. Negative pressure wound therapy after severe open fractures: a prospective randomized study. <i>Journal of Orthopedic Trauma</i> . 2009 Sep;23(8):552-7.	下肢の骨折	1b
Pachowsky M, Gusinde J, Klein A, Lehl S, Schulz-Drost S, Schlechtweg P, Pauser J, Gelse K, Brem MH. Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty. <i>International Orthopaedics</i> . 2012 Apr;36(4):719-22.	人工股関節置換術	1b
Redfern RE, Cameron-Ruetz C, O'Drobinak S, Chen J, Beer KJ. Closed incision negative pressure therapy effect on postoperative infection and surgical site complication after total hip and knee arthroplasty. <i>J Arthroplasty</i> 2017 Nov;32(11):3333-3339.†	人工股関節・人工膝関節置換術	2
Reddy RN Jr, Leng XI, Woodall J, Jackson B, Dedmond B, Webb LX. The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery. <i>Journal of Surgical Orthopaedic Advances</i> . 2010 Jun;19(2):91-7.	人工股関節置換術	3
Cooper HJ, Roc GC, Bas MA, Berliner ZP, Hepinstall MS, Rodriguez JA, Weiner LS. Closed incision negative pressure therapy decreases complications after periprosthetic fracture surgery around the hip and knee. <i>Injury</i> . 2018 Feb;49(2):386-391. doi: 10.1016/j.injury.2017.11.010. Epub 2017 Nov 14.	人工関節周囲骨折に対する手術	3
Cooper HJ, Bas MA. Closed-Incision Negative-Pressure Therapy Versus Antimicrobial Dressings After Revision Hip and Knee Surgery: A Comparative Study. <i>J Arthroplasty</i> . 2016 May;31(5):1047-52.	人工股関節・人工膝関節再置換術	3
Anatone AJ, Shah RP, Jennings EL, Geller JA, Cooper J. A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty. <i>Arthroplasty Today</i> . 2018;4(4):493-498. doi:10.1016/j.artd.2018.09.004.	人工股関節・人工膝関節置換術	3
Curley AJ, Terhune EB, Velott AT, Argintar EH. Outcomes of Prophylactic Negative Pressure Wound Therapy in Knee Arthroplasty. <i>Orthopedics</i> . 2018;41(6):e837-e840. doi:10.3928/01477447-20181010-02.	人工膝関節置換術	3

引用	症例・術式	臨床エビデンス・レベル*
Reddix RN, Tyler HK, Kulp B, Webb LX. Incisional vacuum-assisted wound closure in morbidly obese patients undergoing acetabular fracture surgery. <i>The American Journal of Orthopedics</i> . 2009 Sep;38(9):32-5.	臼蓋骨折	4
Hansen E, Durinka JB, Costanzo JA, Austin MS, Deirmengian GK. Negative pressure wound therapy is associated with resolution of incisional drainage in most wounds after hip arthroplasty. <i>Clinical Orthopaedics and Related Research</i> . 2013 Oct;471(10):3230-6.	人工股関節置換術	4
Stannard JP, Atkins BZ, O-Malley D, Singh H, Bernstein B, Fahey M, Masden D, Attinger CE. Use of negative pressure therapy on closed surgical incisions: A case series. <i>Ostomy Wound Management</i> . 2009 Aug;55(8):58-66.	下肢の骨折	4
Gomoll AH, Lin A, Harris MB. Incisional vacuum-assisted closure therapy. <i>Journal of Orthopaedic Trauma</i> . 2006 Nov-Dec;20(10):705-9.	整形外傷	4
Stannard JP, Gabriel A, Lehner B. Use of Negative Pressure Wound Therapy Over Clean, Closed Surgical Incisions. <i>International Wound Journal</i> . 2012;9:32-39.	整形外傷	4
Berkowitz MJ. Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures. <i>Techniques in Foot and Ankle Surgery</i> . 2013 Dec;12(4):172-174.	踵骨骨折	5
Brem MH, Bail HJ, Biber R. Value of Incisional Negative Pressure Wound Therapy in Orthopedic Surgery. <i>International Wound Journal</i> . 2014 Jun;11(Suppl 1):3-5.	整形全般	5
Suleiman LI, Mesko DR, Nam D. Intraoperative Considerations for Treatment/Prevention of Prosthetic Joint Infection. <i>Current Reviews in Musculoskeletal Medicine</i> . 2018:1-8.	人工股関節・人工膝関節置換術	5
Chotanaphuti T, Courtney PM, Fram B, Kleef N.J., Kim TK, Kuo FC, Lustig S, Moojen DJ, Nijhof M, Oliashirazi A, Poolman R, Purtil JJ, Rapisarda A, Rivero-Boschert S, Veltman ES. Hip and Knee Section, Treatment, Algorithm: Proceedings of International Consensus on Orthopedic Infections. <i>The Journal of Arthroplasty</i> . 34(2S):S393-S397. doi: 10.1016/j.arth.2018.09.024.	人工股関節・人工膝関節置換術	5
DeCarbo WT, Hyer CF. Negative-Pressure Wound Therapy Applied to High-Risk Surgical Incisions. <i>Journal of Foot and Ankle Surgery</i> . 2010 May;49(3):299-300.	整形外傷	5
Nam D, Sershon RA, Levine BR, Della Valle CJ. The Use of Closed Incision Negative-Pressure Wound Therapy in Orthopaedic Surgery. <i>J Am Acad Orthop Surg</i> . 2018:1-8. doi: 10.5435/JAAOS-D-17-00054.	整形外科手術	5
Al-Houraihi RK, Aalirezaie A, Adib F, Anoushiravani A, Bhashyam A, Binlaksar R, Blevins K, Bonanzinga T, Chih-Kuo F, Cordova M, Deirmengian GK, Fillingham Y, Frenkel T, Gomez J, Gundtoft P, Harris MA, Harris M, Heller S, Jennings JA, Jimenez-Garrido C, Karam JA, Khlopas A, Klement MR, Komnos G, Krebs V, Lachiewicz P, Miller AO, Mont MA, Montanez E, Romero CA, Schwarzkopf R, Shaffer A, Sharkey PF, Smith BM, Sodhi N, Thienpont E, Villanueva AO, Yazdi H. General Assembly, Prevention, Wound Management: Proceedings of International Consensus on Orthopedic Infections. <i>The Journal of Arthroplasty</i> . 2019;34(2):S157-S168. doi:10.1016/j.arth.2018.09.066.	整形外科領域の感染症	5
Agarwal A. Management of Closed Incisions Using Negative-Pressure Wound Therapy in Orthopedic Surgery. <i>Plastic and reconstructive surgery</i> . 2019;143(1 Management of Surgical Incisions Utilizing Closed Incision Negative Pressure Therapy):21S-26S.	整形外傷	5

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- Public Health England. Surveillance of surgical site infections in NHS hospitals in England April 2018 to March 2019. Published December 2019.
- Jenks, P.J. Clinical and economic burden of surgical site infection (SSI) and predicted financial consequences of elimination of SSI from an English hospital, Volume 86 (2014), Issue 1, pg 24–33.
- Merkow R, . Underlying reasons associate with hospital readmission following surgery in the US. . 2015;313(5):483–95.
- M. Nobile, P. Navone, A. Orzella, et al. Developing a model for analysis the extra costs associated with surgical site infections (SSIs): an orthopaedic and traumatological study run by the Gaetano Pini Orthopaedic Institute, 4 (2015), p. P68.
- Wilkes RP, Kilpadi DV, Zhao Y, . Closed Incision Management With Negative Pressure Wound Therapy (CIM): Biomechanics. 2012;19(1):67–75.
- Kilpadi DV, Cunningham MR. Evaluation of Closed Incision Management with Negative Pressure Wound Therapy (CIM): Hematoma/Seroma and Involvement of the Lymphatic System. . 2011;19:588–596.
- Glaser DA, Farnsworth CL, Varley ES. Negative pressure therapy for closed spine incisions: A pilot study. 2012;24(11):308–316.
- Federal Drug Administration. De Novo Classification Request for PREVENA 125 and PREVENA PLUS 125 Therapy Units. De Novo Summary (DEN180013), 2019. https://www.accessdata.fda.gov/cdrh_docs/reviews/DEN180013.pdf
- Newman JM, Siqueira MBP, Klika AK, Molloy RM, Barsoum WK, Higuera CA. Use of Closed Incisional Negative Pressure Wound Therapy After Revision Total Hip and Knee Arthroplasty in Patients at High Risk for Infection: A Prospective, Randomized Clinical Trial. *J Arthroplasty*. 2019;34(3):554–559.
- Cooper HJ, Bas MA. Closed-Incision Negative-Pressure Therapy Versus Antimicrobial Dressings After Revision Hip and Knee Surgery: A Comparative Study. *J Arthroplasty*. 2016;31(5):1047–1052.
- Anatone AJ, Shah RP, Jennings EL, Geller JA, Cooper J. A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty. *Arthroplasty Today*. 2018;4(4):493-498. doi:10.1016/j.artd.2018.09.004.
- de Lissoyov G, Fraeman K, Hutchins V, Murphy D, Song D, Vaughn BB. Surgical site infection: incidence and impact on hospital utilization and treatment costs. *Am J Infect Control*. 2009 Jun;37(5):387-97.
- Kurtz SM, Lau E, Watson H, Schmier JK, Parvizi J. Economic burden of periprosthetic joint infection in the United States. *J Arthroplasty*. 2012 Sep;27(8 Suppl):61-5.e1. doi: 10.1016/j.arth.2012.02.022. Epub 2012 May 2.
- Thakore RV, et al. Surgical site infection in orthopedic trauma: A case-control study evaluating risk factors and cost. *Journal of Clinical Orthopaedics and Trauma*. 2015;(6):220-226.
- Reddy HV, Ujwala J, Swetha M, Ramya SB. Seroma: an interesting case report. *Int J Reprod Contracept Obstet Gynecol*. 2014;3(1):254-257.
- Carlson MA. Acute wound failure. *Surg Clin Nor Am*. 1997;77(3):607–36.
- Son D, Harijan A. Overview of surgical scar prevention and management. *J Korean Med Sci*. 2014;29:751-57.
- Mangram AJ, Horan TC, Pearson ML, et al. Guidelines for prevention of surgical site infection, 1999. *Inf Control Hosp Epidemiol*. 1999;20(4):247-78.
- Culver DH, Horan TC, Gaynes RP, et al. Surgical wound infection rates by wound class, operative procedure, and patient risk index. *Am J Med*. 1991;91 (suppl 3B):152-157.
- Jones ME, Hardy CJ, Ridgway JM. Head and neck keloid management: a retrospective early review on a new approach using surgical excision, platelet rich plasma and in-office superficial photon X-ray radiation. *Eduorium J Otolaryngol*. 2015;2:14-19.
- Willy, C., Agarwal, A., Andersen, C.A., Santis, G.D., Gabriel, A., Grauhan, O., Guerra, O.M., Lipsky, B.A., Malas, M.B., Mathiesen, L.L., Singh, D.P. and Reddy, V.S. (2017), Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. *Int Wound J*, 14: 385–398. doi:10.1111/iwj.12612.

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Please Recycle. Printed in Japan. ©3M 2021, All Right Reserved. KCI-114-A(1021)IT