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Cita	Tipo de herida o cirugía	Nivel de evidencia clínica*
Engelhardt M, Rashad NA, Willy C, Müller C, Bauer C, Debus S, Beck T. Closed-incision negative pressure therapy to reduce groin wound infections in vascular surgery: a randomized controlled trial. <i>Int Wound J.</i> 2018;doi: 10.1111/iwj.12848	Heridas inguinales	1b ●
Gombert A, Babilon M, Barbati ME, Keskei A, et al. Closed incision negative pressure therapy reduces surgical site infections in vascular surgery: a prospective randomised trial (AIMS Trial). <i>Eur J Vasc Endovasc Surg.</i> 2018;doi:10.1016/j.ejvs.2018.05.018	Heridas inguinales	1b ●
Kwon J, Staley C, McCullough M, Goss, S, et al. A randomized clinical trial evaluating negative pressure therapy to decrease vascular groin incision complications. <i>J Vasc Surg.</i> 2018;1-9	Heridas inguinales	1b ●
Lee AJ, Sheppard CE, Kent WD, Mewhort H, Sikdar KC, Fedak PW. Safety and efficacy of prophylactic negative pressure wound therapy following open saphenous vein harvest in cardiac surgery: a feasibility study. <i>Interact Cardiovasc Thorac Surg.</i> 2017;24(3):324-328.	Extracción abierta de la vena safena en cirugía cardíaca	1b ●
Lee K, Murphy PB, Ingves MV, Duncan A, DeRose G, Dubois L, Forbes TL, Power A. Randomized clinical trial of negative pressure wound therapy for high-risk groin wounds in lower extremity revascularization. <i>J Vasc Surg.</i> 2017;66(6):1814-19.	Heridas inguinales	1b ●
Pleger SP, Nink N, Elzien M, Kunold A, Koshty A, Böning A. Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study. <i>Int Wound J.</i> 2018;15(1):75-83. doi: 10.1111/iwj.12836.	Incisiones inguinales	1b ●
Colli A. First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients. <i>J Cardiothoracic Surg.</i> 2011;6(1):160.	Esternotomía	2 ●
Weir G. The use of a surgical incision management system on vascular surgery incisions: a pilot study. <i>Int Wound J.</i> 2014;11 Suppl 1:10-2.	Bypass vascular	2 ●
Grauhan O, et al. Effect of surgical incision management on wound infections in a poststernotomy patient population. <i>Int Wound J.</i> 2014;11 Suppl 1:6-9.	Esternotomía	2b ●
Grauhan O, et al. Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy. <i>J Thorac Cardiovasc Surg.</i> 2013;145(5):1387-92.	Esternotomía	2b ●
Matatov T, et al. Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients. <i>J Vasc Surg.</i> 2013;57(3):791-5.	Bypass vascular	3 ●
Santarpino G, Gazdag L, Sirch J, Vogt F, Ledwon M, Fischlein T, Pfeiffer S. A Retrospective Study to Evaluate Use of Negative Pressure Wound Therapy in Patients Undergoing Bilateral Internal Thoracic Artery Grafting. <i>Ostomy Wound Manage.</i> 2015;61(12):26-30.	Injerto de arteria torácica	3 ●
Simon K, et al. [Use of Negative Pressure Wound Therapy on Surgical Incisions (Prevena™) after Surgery of Pectus Deformities Reduces Wound Complications.]. <i>Zentralblatt für Chirurgie.</i> 2014. [en alemán]	Esternotomía	3 ●
Atkins BZ, Tetterton JK, Petersen RP, et al. Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications? <i>Surg Innov.</i> 2009;16(2):140-6.	Esternotomía	4 ●
Atkins BZ, et al. Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy. <i>Int Wound J.</i> 2011;8(1):56-62.	Esternotomía	4 ●
Reddy VS. Use of Closed Incision Management with Negative Pressure Therapy for Complex Cardiac Patients. <i>Cureus.</i> 2016;8(2):e506.	Esternotomía	4 ●
Simon K, et al. [Use of Negative Pressure Wound Therapy on Surgical Incisions (Prevena™) after Surgery of Pectus Deformities Reduces Wound Complications.]. <i>Zentralblatt für Chirurgie.</i> 2014. [en alemán]	Esternotomía	3 ●
Chopra K, Tadisina KK, Singh DP. The 'French Fry' VAC Technique: Hybridization of Traditional Open Wound NPWT with Closed Incision NPWT. <i>Int Wound J.</i> 2016;13(2):216-9.	Linfedema masivo localizado	5 ●

Cita	Tipo de herida o cirugía	Nivel de evidencia clínica*
Dohmen PM, et al. Can post-sternotomy mediastinitis be prevented by a closed incision management system? <i>GMS Hyg Infect Control</i> . 2014;9(3):Doc19.	Esternotomía	5 ●
Dohmen PM, et al. Use of incisional negative pressure wound therapy on closed median sternal incisions after cardiothoracic surgery: clinical evidence and consensus recommendations. <i>Med Sci Monit</i> . 2014;20:1814-25.	Esternotomía	5 ●
Haghshenas Kashani A, Varcoe RL. A New Negative Pressure Dressing (Prevena) to Prevent Wound Complications Following Lower Limb Distal Arterial Bypass. <i>Br J Diabetes Vasc Dis</i> . 2011;11(1):21-4.	Bypass vascular	5 ●
Wu RT, Sumpio BJ, Miller S, Sumpio BE. Use of closed-incision negative pressure therapy: cardiothoracic and vascular surgery. <i>Plast Reconstr Surg</i> . 2019;143:31S.	Esternotomía	5 ●

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Referencias:

1. Zimlichman E, Henderson D, Tamir O, et al. Health care-associated infections. A meta-analysis of costs and financial impact on the US health care system. *JAMA Intern Med.* 2013;173(22):2039-2046. doi:10.1001/jamainternmed.2013.9763.
2. Magill SS, Edwards JR, Bamberg W, et al. Multistate point-prevalence survey of health care-associated infections. *N Engl J Med.* 2014;370(13):1198-1208.
3. Zhan C, Miller MR. Excess Length of Stay, Charges, and Mortality Attributable to Medical Injuries During Hospitalization. *JAMA.* 2003;290:1868-1874.
4. Urban JA. Cost analysis of surgical site infections. *Surg Infect.* 2006;17(Suppl 1):S19-S22.
5. Fowler VG Jr, O'Brien SM, Muhlbaier LH, et al. Clinical predictors of major infections after cardiac surgery. *Circulation.* 2005;112(9 Suppl.):1358-65.
6. de Lissoyoy G, Fraeman K, Hutchins V, et al. Surgical site infection: Incidence and impact on hospital utilization and treatment costs. *Am J Infect Control.* 2009;37(5):387-397
7. Giles KA, Hamdan AD, Pomposelli FB, Wyers MC, Siracuse JJ, Schermerhorn ML. Body mass index: surgical site infections and mortality after lower extremity bypass from the national surgical quality improvement program 2005-2007. *Ann Vasc Surg.* 2007;24:48-56.
8. Davenport DL, Zwischenberger BA, Xenos ES. Analysis of 30-day readmission after aortoiliac and infrainguinal revascularization using the American College of Surgeons national surgical quality improvement program data set. *J Vasc Surg.* 2011;60:1266-74.
9. Kuy S, Dua A, Desai S, Dua A, Patel B, Tondravi N, et al. Surgical site infections after lower extremity revascularization procedures involving groin incisions. *Ann Vasc Surg.* 2014;28:53-8.
10. Ozaki CK, Hamdan AD, Barshes NR, Wyers M, Hevelone ND, Belkin M, et al. Prospective, randomized, multi-institutional clinical trial of a silver alginate dressing to reduce lower extremity vascular surgery wound complications. *J Vasc Surg.* 2015;61:419-27.e1.
11. Wiseman JT, Fernandes-Taylor S, Barnes ML, Saunders RS, Saha S, Havlena J, et al. Predictors of surgical site infection after hospital discharge in patients undergoing major vascular surgery. *J Vasc Surg.* 2012;62:1023-31.e5.
12. Nguyen L, Leya GA, Hevelone ND, et al. Prospective cost analysis and implications of wound complications in lower extremity vascular surgery procedures. *J Vasc Surg.* 2014;60(3):813.
13. Wilkes RP, Kilpadi DV, Zhao Y, et al. Closed incision management with negative pressure wound therapy (CIM): Biomechanics. *Surg Innov.* 2012;19(1):67-75.
14. Kilpadi DV, Cunningham MR. Evaluation of Closed Incision Management with Negative Pressure Wound Therapy (CIM): Hematoma/Seroma and Involvement of the Lymphatic System. *Wound Repair Regen.* 2011;19:588-596.
15. Glaser DA, Farnsworth CL, Varley ES, et al. Negative pressure therapy for closed spine incisions: A pilot study. *Wounds.* 2012;24(11):308-316.
16. Antoniou GA, Onwuka CC, Antoniou SA, Russell D. Meta-analysis and trial sequential analysis of prophylactic negative pressure therapy for groin wounds in vascular surgery. *J Vasc Surg.* 2019; doi 10.016/jvs.2019.01.083.
17. Svensson-Björk R, Zarrouk M, Ascitutto G, Hasselmann J, Acosta S. Meta-analysis of negative pressure wound therapy of closed groin incisions in arterial surgery. *BJS.* 2018; doi:10.1002/bjs.11100.
18. Grauhan O, Navasardyan A, Hofmann M, et al. Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy. *J Thorac Cardiovasc Surg.* 2013;145:1387-1392.
19. Pleger SP, Nink N, Elzien M, Kunold A, Koshty A, Böning A. Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study. *Int Wound J.* 2018;15(1):75-83.
20. Reddy HV, Ujwala J, Swetha M, Ramya SB. Seroma: an interesting case report. *Int J Reprod Contracept Obstet Gynecol.* 2014;3(1):254-257.
21. Carlson MA. Acute wound failure. *Surg Clin Nor Am.* 1997;77(3):607-36.
22. Son D, Harijan A. Overview of surgical scar prevention and management. *J Korean Med Sci.* 2014;29:751-57.
23. Mangram AJ, Horan TC, Pearson ML, et al. Guidelines for prevention of surgical site infection, 1999. *Inf Control Hosp Epidemiol.* 1999;20(4):247-78.
24. Culver DH, Horan TC, Gaynes RP, et al. Surgical wound infection rates by wound class, operative procedure, and patient risk index. *Am J Med.* 1991;91 (suppl 3B):152-157.
25. Jones ME, Hardy CJ, Ridgway JM. Head and neck keloid management: a retrospective early review on a new approach using surgical excision, platelet rich plasma and in-office superficial photon X-ray radiation. *Edorium J Otolaryngol.* 2015;2:14-19.
26. Public Health England. Surveillance of surgical site infections in NHS hospitals in England, 2014-15. Diciembre de 2015. Disponible en: www.gov.uk/phe
27. Willy C, Agarwal A, Andersen CA, et al. Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. *Int Wound J.* 2017;14:385-398. doi:10.1111/iwj.12612.
28. Stannard J, Atkins B, O'Malley D, et al. Use of negative pressure therapy on closed surgical incisions: a case series. *Ostomy Wound Manage.* 2009;55:58-66.
29. Documento de consenso de la Unión Mundial de Sociedades de Cuidado de Heridas (WUWHs). Closed surgical incision management: understanding the role of NPWT. *Wounds International.* 2016.

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3M y las otras marcas que aparecen son marcas o marcas registradas. © 3M 2021. Todos los derechos reservados.
KCI, USA, Inc.