
Ranger™

en	Blood/Fluid Warming Unit, Model 245 Operator's Manual	01
fr	Manuel de l'utilisateur de l'unité de réchauffement pour sang/liquide, modèle 245	10
de	Erwärmungseinheit für Blut/Flüssigkeit, Modell 245, Bedienungsanleitung	20
it	Manuale d'uso dell'Unità riscaldante per sangue/fluidi, modello 245	30
es	Manual del operario de la unidad de calentamiento de sangre/fluidos (Modelo 245)	40
nl	Bedieningshandleiding verwarmingsunit bloed/vloeistof, model 245	50
sv	Blod-/vätskevärmeeinheit, modell 245 Bruksanvisning	60
da	Varmeeinheit til blod/væske, model 245, betjeningsmanual	70
no	Varmeeinheit for blod/væske, modell 245 Brukerhåndbok	80
fi	Veren/nesteen lämmitysyksikkö malli 245 Käyttöopas	90
pt	Manual do Operador da Unidade de Aquecimento de Fluidos/Sangue, Modelo 245	100
el	Μονάδα Θέρμανσης Αίματος/Υγρών, Μοντέλο 245 Εγχειρίδιο Χειριστή	110
pl	Urządzenie do podgrzewania, model 245 — instrukcja obsługi	121
hu	Vér-/folyadékmelegítő egység, 245-ös modell, Kezelői kézikönyv	131
cs	Návod k obsluze ohřívací jednotky pro krev/tekutiny, model 245	141
sk	Používateľská príručka jednotky na ohrev krvi/tekutín, model 245	151
sl	Enota za ogrevanje krvi/tekočin, model 245 Uporabniški priročnik	161
et	Vere/vedeliku soojendussüsteemi mudeli 245 kasutusjuhend	171
lv	Asins/šķīdruma sildītājs, modelis 245 – lietotāja rokasgrāmata	181
lt	Kraujo / skysčio šildymo įrenginio (245 modelio) naudotojo vadovas	191
ro	Unitate de încălzire a sângelui/fluidului, Model 245 Manualul operatorului	201
ru	Инструкция по эксплуатации устройства для подогрева крови/растворов, модель 245	211
hr	Grijalica krvi/tekućina, model 245, operaterov priručnik	222
bg	Ръководство за оператора на устройство за затопляне на кръв/течност, модел 245	232
tr	Kan/Sıvı Isıtma Ünitesi, Model 245 Kullanım Kılavuzu	244
zh	245 型输血/输液加温器操作手册	254
ar	دليل المشغل لوحدة احتراز السوائل/الدم، الطراز 245	264

Section 1: Technical Service and Order Placement	2
Technical Service and Order Placement	2
USA	2
Outside of the USA.	2
In-warranty repair and exchange	2
When you call for technical support	2
Section 2: Introduction.	2
Product Description	2
Indications for use.	2
Patient Population and Settings.	2
Explanation of Signal Word Consequences	2
WARNING:	2
CAUTION:	3
NOTICE:	3
Overview and Operation	3
Ranger blood/fluid warming set example	3
Model 245 Product safety features	4
Section 3: Instructions for Use	4
Preparation and setup of the Ranger blood/fluid warming unit	4
Removing the warming set from the Ranger blood/fluid warming unit	4
Transferring the warming set from one Ranger warming unit to another	4
Section 4: Troubleshooting	5
Section 5: General Maintenance and Storage	5
Cleaning Instructions	5
To clean the exterior of the warming unit:	5
To clean the heating plates:	6
To clean resistant, dried-on fluids:	6
Storage	6
Servicing	6
Symbol Glossary	6
Section 6: Specifications.	7
Physical Characteristics	9
Classifications	9
Leakage Current	9
Environmental Conditions	9

Section 1: Technical Service and Order Placement

USA: TEL: 1-800-228-3957 (USA Only)

Outside of the USA: Contact your local 3M representative.

3M assumes no responsibility for the reliability, performance, or safety of the temperature warming unit if the following events occur:

- Modifications or repairs are not performed by a qualified, medical equipment service technician who is familiar with good practice for medical device repair.
- The unit is used in a manner other than that described in the Operator's or Preventive Maintenance Manual.
- The unit is installed in an environment that does not provide grounded electrical outlets.
- The warming unit is not maintained in accordance with the procedures described in the Preventive Maintenance Manual.

In-warranty repair and exchange

To return a device to 3M for service, first obtain a Return Authorization (RA) number from a customer service representative. Please use the (RA) number on all correspondence when returning a device for service. A shipping carton will be delivered to you at no charge, if needed. Call your local supplier or sales representative to inquire about loaner devices while your device is being serviced.

When you call for technical support

We will need to know the serial number of your unit when you call us. The serial number label is located on the bottom of the warming unit.

Section 2: Introduction

Product Description

The Ranger blood/fluid warming system consists of a Model 245 warming unit and a sterile disposable fluid warming set.

The 3M™ Ranger™ blood/fluid warming system includes a warming unit and a disposable warming set. The Ranger blood/fluid warming system is designed to warm blood, blood products, and liquids and deliver these at flow rates from KVO to 500 mL/min. At these flow rates, the device maintains fluid output temperatures ranging from 33°C to 41°C (Note: Output temperatures are dependent on input fluid temperatures and flow rate). It takes less than 2 minutes to warm up to the 41°C set point temperature. The alert points on the Model 245 are 43° and 44°C.

Disposable blood/fluid warming sets are available in: pediatric flow, standard flow, and high flow applications. Warming sets are sterile, not made with natural rubber latex, single-use-only items and are designed to be used with the warming unit.

The Ranger blood/fluid warming unit is designed to be mounted to an I.V. pole. A handle located on the top of the unit makes transport easy. When mounted to the I.V. pole, the unit fits easily above the 3M™ Bair Hugger™ warming unit.

This manual includes operating instructions and unit specifications for the Ranger blood/fluid warming system. For information about using Ranger blood/fluid warming sets with the Ranger blood/fluid warming unit, Model 245, refer to the "Instructions for Use" included with each warming set. The Ranger blood/fluid warming system should only be used in healthcare facilities by trained medical professionals.

Indications for use

The Ranger blood/fluid warmer is intended to warm blood, blood products, and liquids.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency trauma settings, or other areas where blood/fluids are infused.

Explanation of Signal Word Consequences

WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE: Indicates a situation which, if not avoided, could result in property damage only.

WARNING:

1. To reduce the risks associated with hazardous voltage and fire and thermal energy hazards:
 - Do not substitute other devices (i.e. Model 247) for the Ranger blood/fluid warming unit or Ranger blood/fluid warming sets.
 - Do not continue use of the unit if the over-temperature alert continues to sound and the temperature does not return to the set point temperature. Immediately stop fluid flow and discard the warming set. Have the blood/fluid warming unit tested by a biomedical technician or call 3M.
2. To reduce the risks associated with hazardous voltage and fire:
 - Do not modify or service this device, and do not open the warming unit case as there are no user-serviceable parts in the unit.
 - Connect power cord to receptacles marked "Hospital Only," "Hospital Grade," or a reliably grounded outlet.
 - Use only the power cord specified for this product and certified for the country of use.
 - Do not allow the power cord to get wet.
 - Do not use the Ranger blood/fluid warming system when it appears the unit power cord, or warming set is damaged. Use only 3M specified replacement parts.
 - Keep power cord visible and accessible at all times. The plug on the power cord serves as the disconnect device. The wall socket outlet shall be as close as practical and shall be easily accessible.
 - Do not use with a multiple socket outlet or extension cord.
 - Do not use this equipment adjacent to, stacked with or in combination with other equipment without verifying that the total leakage current from the combined equipment does not exceed the safety limits for Type BF equipment and to assure normal operation in the configuration in which it will be used.

- To reduce the risks associated with air embolism and incorrect routing of fluids:
 - Never infuse fluids if air bubbles are present in the fluid line.
 - Ensure all luer connections are tightened.
- Do not mount the Ranger blood fluid warming system higher than patient level during infusion, as air embolism may result.
- To reduce the risks associated with potential blood loss, do not use in combination with an extracorporeal circuit.

CAUTION:

- To reduce the risks associated with cross-contamination:
 - The cleaning tool provides only superficial cleaning, it does not disinfect or sterilize the interior of the unit.
- To reduce the risks associated with impact and facility medical device damage:
 - Clamp the Ranger blood/fluid warming unit to an I.V. pole with a minimum 14" (35.6 cm) radius wheelbase and at a height no higher than 44" (112 cm).
- To reduce the risks associated with environmental contamination:
 - Follow applicable regulations when disposing of this device or any of its electronic components.
- Do not use for direct cardiac application. When using the Ranger blood/fluid warming system with a central venous catheter (CVC), ensure that the catheter tip does not have direct contact with the heart and ensure that all electrical devices connected to or near the patient have an appropriate leakage current rating for the application. If a CVC tip is found to be in direct contact with a patient's heart, the blood/fluid warming unit, Model 245 should be immediately disconnected from the CVC until the CVC is safely repositioned. Failure to follow these precautions may cause cardiac disturbance and/or patient injury.
- Follow the AABB Guidelines for the use of blood warming devices which caution against warming when administering platelets, cryoprecipitate, or granulocyte suspensions.
- The 3M Ranger blood/fluid warming system has been tested to be resistant to both electromagnetic fields (EMI) and electrostatic discharge (ESD). To reduce the risk associated with EMI due to portable and mobile RF communications equipment:
 - Install and put into service the 3M Ranger blood/fluid warming system according to the EMC information provided in the Guidance and Manufacturer's Declaration.
 - Should interference occur move away from the portable or mobile RF communications equipment.

NOTICE:

- To avoid device damage:
 - Do not clean the blood/fluid warming unit with solvents. Damage to the case, label, and internal components may result.
 - Do not immerse the blood/fluid warming unit in cleaning or sterilizing solutions. The unit is not liquid proof.
 - Do not insert metallic instruments in the blood/fluid warming unit.
 - Do not use abrasive materials or solutions to clean the heater plates.
 - Do not allow spills to dry inside the unit, as this may make it more difficult to clean the unit.
- The Ranger blood/fluid warming unit meets medical electronic interference requirements. If radio frequency interference with other equipment should occur, connect the unit to a different power source.


Overview and Operation

The 3M™ Ranger™ blood/fluid warming system consists of a Model 245 warming unit and a sterile disposable fluid warming set.

The warming unit is a compact, lightweight, liquid-resistant device with a clamp located on the side for attachment to an I.V. pole (see Figure 1). A carrying handle on the top of the unit makes it easy to transport.

Located on the front panel you will find:

- Alphanumeric display that indicates the heater temperature during normal operation. In an over-temperature condition, the display alternately flashes a temperature of 43°C or higher and the word "HI." An audible alert also sounds. In an under-temperature condition, the display alternately flashes a temperature of 33°C or lower and the word "LO."
- Alert indicator light that comes on when either an over- or under-temperature condition occurs.


The Ranger blood/fluid warming system consists of a Model 245 warming unit and sterile warming set.

Figure 1

Ranger blood/fluid warming set example

- Inlet line
- Blue inlet pinch clamp
- Fluid warming cassette
- Bubble trap
- Roller Clamp
- Patient line
- Injection port
- White outlet pinch clamp
- Patient connection

Refer to the instructions provided with each warming set for information on use.

Model 245 Product safety features

The following chart describes the safety alert features of the Ranger blood/fluid warming unit.

Alert Type	What To Look For	Description/Cause	Action
Over-temperature alert - 43°C	Alert indicator light illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 43°C or higher and the word "HI."	Heater temperature rose to 43°C because of transient conditions.*	Observe alphanumeric display. If temperature does not drop to 41°C (may take a few minutes), discontinue use of unit. Call 3M.
Under-temperature alert - 33°C	Alert indicator light illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 33°C or lower and the word "LO."	Heater temperature has fallen to 33°C.	Alerts stop when temperature rises above 33°C. Continue use of unit. If the temperature does not rise above 33°C, unplug the unit and call 3M.

* Transient conditions may cause an over-temperature alert condition. These conditions include:

- There was an extreme change in flow rates (e.g., from 500 mL/min to stop flow).
- The unit was turned on and reached the set point temperature of 41°C before the warming cassette was inserted in the device.
- Fluids were prewarmed to above 42°C before being infused.

Alert Type	What To Look For	Description/Cause	Action
Independent back-up safety alert - 44°C	Alert indicator light illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 43°C or higher and the word "HI." Alphanumeric display is dark, alert sounds (backup safety alert still works even if display is dark).	Heater temperature rose to 44°C. Safety backup system is activated at 44°C and unit shuts off power to the heating plates.	TURN OFF THE UNIT AND UNPLUG IT. Do not use the blood/fluid warming unit. Discard the disposable set. Call 3M

Section 3: Instructions for Use

Preparation and setup of the Ranger blood/fluid warming unit

1. Attach the Ranger blood/fluid warming unit to the I.V. pole. Tighten the pole clamp securely (see Figure 3).

CAUTION:

To reduce the risks associated with impact and facility medical device damage:

- Clamp the Ranger blood/fluid warming unit to an I.V. pole with a minimum 14" (35.6 cm) radius wheelbase and at a height no higher than 44" (112 cm).
2. Slide the cassette into the slot in the warming unit. The cassette can only fit into the device one way (see Figure 4, #1).
 3. Prime the warming set. For more information about priming the set, refer to instructions provided with the warming sets (see Figure 4, #2).
 4. Place the bubble trap in the holder.
 5. Plug the power cord into an appropriate outlet. Turn the unit ON (see Figure 4, #3). In a few seconds the alphanumeric display will illuminate. It takes less than two minutes to warm up to the 41°C set point temperature.
 6. Begin infusion. When infusion is completed, remove the warming set and discard according to institutional protocol.


Figure 3


Figure 4

Removing the warming set from the Ranger blood/fluid warming unit

1. Close the inlet clamp proximal to the cassette and open all clamps distal to the cassette.
2. Disconnect the warming set from the fluid source, if applicable.
3. Allow fluid to flow to the patient (this may take 2-3 seconds). Close a distal clamp.
4. Remove the cassette from the warming unit and discard according to institutional protocol.
5. Reconnect the patient I.V. line to the fluid source to continue the infusion without warming.

Transferring the warming set from one Ranger warming unit to another

1. Follow steps 1-3 above, then remove the warming set from the first warming unit.
2. During transport, keep the clamps closed and do not infuse fluids while the cassette is outside of the warming unit.
3. Slide the cassette into the second warming unit.
4. Ensure that air is removed from the tubing.
5. Open the clamps and continue the infusion.

Section 4: Troubleshooting

Condition	Cause	Solution
Nothing illuminates on the warming unit panel.	Unit is not turned on, plugged in, or power cord is not plugged into an appropriate outlet.	Turn unit on. Make sure the power cord is plugged into the power entry module of the warming unit. Make sure the warming unit is plugged into a properly grounded outlet.
	Unit failure.	Check panel fuses. Call 3M.
Alert indicator illuminates and alert sounds, alpha-numeric display alternately flashes a temperature of 43°C or higher and the word "HI."	Temporary over-temperature condition because: An extreme change in flow rates occurred (e.g., from 500 mL/min to stop flow). Unit was turned on and reached set point temperature before cassette was inserted. Fluids were prewarmed to above 42°C before being run through the warming unit.	Open flow to reduce temperature. Alerts will stop when the display reads 41°C. The unit is ready to use. Alerts will stop when the display reads 41°C. The unit is ready to use. Turn off unit and unplug it. Discontinue infusion of fluids. Do not warm fluids before infusing them through the Ranger warming unit.
Alert sounds, alphanumeric display and alert indicator light go dark.	Primary controller failure. Unit will no longer operate.	Power to heating plates will shut off if temperature rises to 44°C (warming unit Rev N and newer) or 46°C (warming units Rev A to M). Turn unit off and unplug it. Discontinue use of unit. Discard disposable set. Alert will continue to sound if you do not unplug unit. Call 3M.
Unit alerts soon after plugging it in (unit does not have to be turned on for this condition to occur). Heater temperature rises to 44°C (warming units Rev N and newer) or 46°C (warming units Rev A to M) and unit shuts down soon after plugging it in (unit does not have to be turned on for this condition to occur).	Test screw on bottom of unit is loose or missing.	Make sure test screw is completely tightened. If it is missing, turn unit off and unplug it. Call 3M.
Alert sounds but unit has been turned off.	Independent backup safety system has been activated.	Unplug unit. Call 3M.
Cannot remove cassette from unit.	Cassette is too full, fluids are still being infused, or clamp is open proximal to the cassette. Warming unit is below patient level, creating excessive back pressure.	Make sure fluid is drained from cassette before sliding out the cassette, that fluids are no longer infusing, and that clamp is closed proximal to the warming cassette. Raise unit above patient level.
Alert indicator light illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 33°C or below and the word "LO."	Under temperature condition caused by very high flow using very cold fluid, or defective heater/relay.	Alert should stop when temperature rises above 33°C. If alert continues, turn unit off, unplug unit and discontinue use. Call 3M.
Alphanumeric display reads "Er 4" or "Open."	Open wire on temperature sensor.	Do not use unit. Call 3M.
Alphanumeric display reads "Er 5" or "Open."	Electrical interference.	Remove the unit. Refer to biomedical technician or call 3M

Section 5: General Maintenance and Storage

Cleaning Instructions

To clean the exterior of the warming unit:

1. Disconnect the warming unit from the power source before cleaning.
2. Cleaning should be performed in accordance with hospital practices for cleaning OR equipment. After every use; wipe the warming unit and any other surfaces that may have been touched. Use a damp, soft cloth and a hospital approved mild detergent, germicidal disposable wipes, disinfecting towelettes, or antimicrobial spray. The following active ingredients are acceptable for use in cleaning the warming unit:
 - Oxidizers (e.g. 10% Bleach)
 - Quaternary Ammonium Compounds (e.g. 3M™ Quat Disinfectant Cleaner)
 - Phenolics (e.g. 3M™ Phenolic Disinfectant Cleaner)
 - Alcohols (e.g. 70% Isopropyl Alcohol)
3. Let air dry.

CAUTION

To reduce the risks associated with cross-contamination:

- The cleaning tool provides only superficial cleaning, it does not disinfect or sterilize the interior of the unit.

To clean the heating plates:

The Ranger hardware cleaning tool, is intended to clean both heating plates of the warming unit. It is not necessary to disassemble the warming unit to use the tool.

1. Unplug the irrigation fluid warming unit.
2. Unfold the cleaning tool. Wet the foam pads with a nonabrasive solution. The following active ingredients are acceptable for use with the cleaning tool:
 - Oxidizers (e.g. 10% Bleach)
 - Quaternary Ammonium Compounds (e.g. 3M™ Quat Disinfectant Cleaner)
 - Phenolics (e.g. 3M™ Phenolic Disinfectant Cleaner)
 - Alcohols (e.g. 70% Isopropyl Alcohol)
3. Insert the tool from the back of the unit and pull the tool all the way out from the front (see Figure 5).
4. Rinse the tool with water and repeat 3 times. Discard the tool according to institutional protocol.
5. Wipe off the unit to remove excess fluid.


Figure 5

To clean resistant, dried-on fluids:

1. Spray a nonabrasive solution inside the slot of the warming unit and let sit for 15-20 minutes.
2. Clean the unit by using the cleaning tool.

NOTICE

1. To avoid device damage:
 - Do not immerse the Ranger unit or accessories in any liquid or subject them to any sterilization process.
 - Do not clean the warming unit with solvents. Damage to the case, label, and internal components may result.
 - Do not insert metallic instruments in the warming unit.
 - Do not use abrasive materials or solutions to clean the heater plates.
 - Do not allow spills to dry inside the unit, as this may make it more difficult to clean the unit.

NOTE: You may use a nonmetal instrument, such as a cotton swab, to clean the upper channels. If you are unable to adequately clean the unit, call 3M technical service.

Storage

Store all components in a cool, dry place when not in use.

Servicing

There are no user-serviceable parts in the Ranger blood/fluid warming unit. All service must be performed by 3M or an authorized service technician. In the USA, call 3M at 1-800-228-3957 for service information. Outside of the USA contact your local 3M representative.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

The following symbols may appear on the product's labeling or exterior packaging.

"OFF" (power)		To indicate disconnection from the mains, at least for main switches, or their positions, and all those cases where safety is involved. Source: IEC 60417-5008
"ON" (power)		To indicate connection to the mains, at least for mains switches, or their positions, and all those cases where safety is involved. Source: IEC 60417-5007
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Caution		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4
CE Mark 2797		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3

Type BF applied part		To identify a type BF applied part complying with IEC 60601-1. Source: IEC 60417-5333
Equipotentiality		To identify the terminals which, when connected together, bring the various parts of an equipment or of a system to the same potential, not necessarily being the earth (ground) potential. Source: IEC 60417-5021
Refer to instruction manual/ booklet		To signify that the instruction manual/booklet must be read. Source: ISO 7010-M002
Fuse		Indicates a replaceable fuse
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Keep dry		Indicates a medical device that needs to be protected from moisture. Source: ISO 15223, 5.3.4
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Protective earth (ground)		To identify any terminal which is intended for connection to an external conductor for protection against electric shock in case of a fault, or the terminal of a protective earth (ground) electrode. Source: IEC 60417, 5019
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle. Source: Directive 2012/19/EC on waste electrical and electronic equipment (WEEE)
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of physician. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Serial number		Indicates the manufacturer's serial number so that a specific medical device can be identified. Source: ISO 15223, 5.1.7
UL Classified		Indicates product was evaluated and Listed by UL for the USA and Canada.
Unique device identifier		Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10

For more information see HCBGregulatory.3M.com

Section 6: Specifications

Guidance and manufacturer's declaration – electromagnetic emissions

The model 245 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 245 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The model 245 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The model 245 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The model 245 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 245 should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	2± kV power lines	2± kV power lines	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	1± kV line to line 2± kV line to gnd	1± kV line to line 2± kV line to gnd	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 6 cycles 70% U_T (30% dip in U_T) for 30 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 6 cycles 70% U_T (30% dip in U_T) for 30 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model 245 requires continued operation during power mains interruptions, it is recommended that the model 245 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The model 245 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 245 should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the model 245, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model 245 is used exceeds the applicable RF compliance level above, the model 245 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the model 245.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the model 245

The model 245 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model 245 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model 245 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.37	0.37	0.74
1	1.17	1.17	2.30
10	3.69	3.69	7.37
100	11.67	11.67	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Physical Characteristics

Warming Unit

4.5 in. (11 cm) high x 7.5 in. (19 cm) wide x 10 in. (25 cm) long; wt.: 7 lb. 7 oz. (3.4 kg)

Classifications

- Protection against electric shock: Class I Medical Electrical Equipment with Type BF applied part.
- Protection against ingress of water: IPX0 (Ordinary Equipment).
- Mode of operation: Continuous operation.

MEDICAL –GENERAL MEDICAL EQUIPMENT AS TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES 60601-1 (2005) + AMD (2012), CAN/CSA-C22.2 No. 60601-1 (2008) + (2014), and IEC 60601-1-6:2010 (Third Edition) + A1:2013; Control No. 4HZ8

Electrical Characteristics	Temperature and Accuracy Characteristics
Device Rating 100-120 VAC, 50/60 Hz 220-240 VAC, 50/60 Hz	Set Point Temperature 41° ± 1.5°C
Maximum Heating Power 900 W	Over-temperature Alert 43° +3/-2°C
Fuses 2 x T10A-H (250V) for 100-120 VAC 2 x T6.3A-H (250V) for 220-240 VAC	Under-temperature Alert 33° ± 2°C
Fuse Type Time delay, high breaking	Over-temperature Cutoff 46° +3/-2°C (warming units Rev N and newer) 46° ± 2°C (warming units Rev A to M)

Leakage Current

Meets leakage current requirements in accordance with IEC 60601-1.

Environmental Conditions

Operating temperature range: 15° to 40°C (59° to 104°F)

Storage and Transport temperature range: -20° to 45°C (-4° to 113°F)

Operating humidity: 10 to 85% RH, noncondensing

Atmospheric pressure range: 50 kPa to 106 kPa

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3M, Ranger and Bair Hugger are trademarks of
Unauthorized use prohibited. All rights reserved.
3M Deutschland GmbH