REF 63000

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Indications for Use

The Bair Hugger family of temperature management systems are indicated for hypothermic patients or normothermic patients for whom induced hypothermia or localized temperature therapy is clinically indicated. In addition, the temperature management systems can be used to provide patient thermal comfort when conditions exist that may cause patients to become too warm or too cold. The temperature management systems can be used with adult and pediatric patients.

The Bair Hugger temperature management system should only be used by trained medical professionals.

- STERILE

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency departments and other departments in the hospital setting where patient temperature management is required.

Explanation of Signal Word Sequences

WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION: indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

Contraindications

- Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

WARNING: To reduce the risk of thermal injury:

- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- Do not place the non-perforated side of the warming blanket on the patient. Always place the perforated side (with the small holes) directly on top of the patient in contact with the patient's skin.
- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 series, 700 series, or 675 warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- Do not place patient securement device (i.e. safety strap or tape) over the warming blanket.
- Do not place the warming blanket directly over a dispersive electrode pad.

WARNING: To reduce the risk of patient injury or death due to altered drug delivery:

- Do not use a warming blanket over transdermal medication patches.

WARNING: To reduce the risk of injury due to interference with ventilation:

- Do not allow the warming blanket or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.

WARNING: To reduce the potential for injury due to patient falls:

- Do not use a warming blanket to transfer or move the patient.

WARNING: To maintain proper aseptic technique:

- Any surgical incision in the lower extremities must be closed before applying the blanket to the lower part of the body.

CAUTION: To reduce the risk of cross-contamination:

- This warming blanket is sterile and is intended for single patient use ONLY. Placing a sheet between the warming blanket and the patient does not prevent contamination of the product.

CAUTION: To reduce the risk of fire:

- This product is classified as Class I Normal Flammability as defined by the Consumer product Safety Commission's flammable fabric regulation, 16 CFR 1610. Follow standard safety protocols when using high intensity heat sources.

CAUTION: To reduce the risk of thermal injury:

- Do not use if primary packaging has been previously opened or is damaged.

CAUTION: To reduce the risk of thermal injury, hyperthermia or hypothermia:


- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

Directions for Use


1. Leave the warming blanket folded up until ready to begin warming therapy.

WARNING: Any surgical incision in the lower extremities must be closed before applying the sterile warming blanket to the lower part of the body.

2. Remove the backing from the adhesive tape strip and adhere the warming blanket to the patient's waist. (Figure A).


3. Unfold the warming blanket and position the clear window to allow access to the groin. Extend the warming blanket so that the end is outside the sterile field (Figure B).

4. Attach the sterile hose of the warming blanket to the Bair Hugger warming unit hose (Figure C). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to insure secure attachment.

WARNING: Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.

5. Tuck the edges of the warming blanket under the patient.
6. Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator Manual for your specific Warming Unit Model.)

CAUTION: Patient Monitoring Recommendations:

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.

- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.
7. Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the hose from the warming blanket and discard the blanket per hospital policy.

Storage/Shelf Life/Disposal

For best results, store at room temperature. Avoid excessive heat and humidity.

For shelf life information, refer to the expiration date on the package.

Dispose of the blanket per facility protocol. If the sterile packaging is damaged or unintentionally opened, discard the product and do not use.

How Supplied

Product is sterile unless package is damaged or opened.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Sterilized using ethylene oxide		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3

Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Caution		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Single sterile barrier system		Indicates a single sterile barrier system. Source: ISO 15223, 5.2.11
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

For more information see, HCBGRegulatory.3M.com

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3M, Bair Hugger and the Bair Hugger logo are trademarks of 3M.
©2020, 3M. Unauthorized use prohibited. All rights reserved.
3M Deutschland GmbH