**3M™ Tegaderm™ Antimicrobial Transparent Dressing**  
**3M™ Tegaderm™ Antimicrobial I.V. Advanced Securement Dressing**

en

**Product Description**

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##### 3M™ Tegaderm™ Antimicrobial I.V. Advanced Securement Dressing (with borders)

Catalog number	Dressing Size	Target Amount of CHG per dressing (mg)
9132	7 cm x 8.5 cm	5.65 mg

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**Glossary Table**

	Manufacturer	Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.
	Authorized Representative in European Community	Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
	Date of manufacture	Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
	Use-by date	Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
	Do not resterilize	Indicates a medical device that is not to be resterilized. Source: ISO 15223, 5.2.6
	Do not use if package is damaged or open	Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8
	Upper limit of temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed. Source: ISO 15223, 5.3.6
	Do not reuse	Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
	Consult instructions for use	Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3

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<b>Rx Only</b>	Rx Only	Indicates that U.S. Federal Law restricts this device to sale by or on the order of physician. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).
<b>Antimicrobial/2%</b>	Antimicrobial Agent- 2% CHG	Indicates product contains 2% Chlorhexidine gluconate as antimicrobial agent.

For more information, see [HCBGRegulatory.3M.com](http://HCBGRegulatory.3M.com)

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**CE 2797**  
 Made in the USA with  
 Globally Sourced Materials

**3M Company**  
 2510 Conway Ave.  
 St. Paul, MN 55144 USA  
 1-800-228-3957 (USA Only)  
[3M.com/Tegaderm](http://3M.com/Tegaderm)

Issue Date: 2020-03  
 34-8724-3503-6

**EC REP 3M Deutschland GmbH**  
**Health Care Business**  
 Carl-Schurz-Str. 1  
 41453 Neuss, Germany

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<b>dZ</b>		Printed Colors – Front: 
<b>Requester:</b> Susan Barker		Printed Colors – Back: 
<b>Creator:</b> deZimnia_27040		Match Colors:
<b>File Name:</b> 34872435036.indd		
<b>Structure #:</b> SS-87661		
<b>Date:</b> 10/16/20		
Scale:  1 Inch		This artwork has been created as requested by 3M. 3M is responsible for the artwork AS APPROVED and assumes full responsibility for its correctness.

**PROOF REQUIRED**

Please refer to 3M Corporate General Spec. 205 for details

**3M™ Tegaderm™ Antimicrobial Transparent Dressing**  
**3M™ Tegaderm™ Antimicrobial I.V. Advanced Securement Dressing**

(en)

**Product Description**

3M™ Tegaderm™ Antimicrobial Transparent Dressing and 3M™ Tegaderm™ Antimicrobial I.V. Advanced Securement Dressing consist of a polyurethane film coated with a transparent chlorhexidine gluconate (CHG) acrylic adhesive. CHG, a broad spectrum antimicrobial/antifungal agent known to inhibit microbial growth has been formulated into the acrylic adhesive.

The transparent film is breathable, allowing oxygen and moisture vapor exchange, yet is impermeable to external contaminants, including fluids (waterproof), bacteria, viruses\*, yeast, and mold. The dressing must remain intact to protect the IV site from external contaminants.

3M™ Tegaderm™ Antimicrobial I.V. Advanced Securement Dressing is bordered, notched and reinforced with soft cloth tape and is designed to provide securement around catheters and other devices.

*In vitro* testing (time kill) demonstrates that 3M™ Tegaderm™ Antimicrobial Transparent Dressing and 3M™ Tegaderm™ Antimicrobial I.V. Advanced Securement Dressing have an antimicrobial effect against a variety of gram-positive bacteria, gram-negative bacteria, yeast and mold in the dressing.

\**In vitro* testing shows that the film of the dressing provides a barrier to viruses 27 nm in diameter or larger while the dressing remains intact without leakage. These results have not been studied with regard to prevention of viral infection. No clinical study has been conducted regarding the ability of the dressing to prevent viral infections.

**Indications for Use**

3M™ Tegaderm™ Antimicrobial Transparent Dressing and 3M™ Tegaderm™ Antimicrobial I.V. Advanced Securement Dressing are intended to be used to cover and protect catheter sites and to secure devices to the skin. Common applications include covering and securing IV catheters, other intravascular catheters and percutaneous devices.

**Warnings**

- DO NOT USE 3M™ TEGADERM™ ANTIMICROBIAL TRANSPARENT DRESSING AND 3M™ TEGADERM™ ANTIMICROBIAL I.V. ADVANCED SECUREMENT DRESSING ON PREMATURE INFANTS OR INFANTS YOUNGER THAN 2 MONTHS OF AGE. USE OF THIS PRODUCT ON PREMATURE INFANTS MAY RESULT IN HYPERSENSITIVITY REACTIONS OR NECROSIS OF THE SKIN.

- FOR EXTERNAL USE ONLY. DO NOT ALLOW CONTACT OF THIS DEVICE WITH EARS, EYES, MOUTH OR MUCOUS MEMBRANES.
- INTENDED FOR SINGLE PATIENT USE ONLY.
- DO NOT REUSE. AS WITH ALL ADHESIVE-BASED PRODUCTS, ADHESIVE EFFECTIVENESS AND FUNCTIONALITY CAN DECLINE AFTER THE FIRST USE AND THE PRODUCT WILL NOT PERFORM AS SPECIFIED. REUSE MAY LEAD TO INFECTION OR OTHER ILLNESS/INJURY.

**Adverse Reactions**

- ADVERSE REACTIONS SUCH AS IRRITATIONS, SENSITIZATION AND GENERALIZED ALLERGIC REACTIONS HAVE BEEN REPORTED WITH THE USE OF CHLORHEXIDINE GLUCONATE.
- IF ALLERGIC REACTIONS OCCUR, DISCONTINUE USE IMMEDIATELY AND, IF SEVERE, CONTACT A PHYSICIAN.

**Precautions**

The dressing should not be placed over infected wounds. This device is not intended to treat, prevent, or reduce catheter-related bloodstream infections (CRBSIs) or other percutaneous device-related infections. This device has not been studied in a randomized clinical study to determine its effectiveness in preventing such infections.

Active bleeding at insertion sites should be stabilized before applying the dressing.


Do not stretch the dressing before applying it to the skin. Dressings applied under tension can cause trauma to the skin.

To ensure good adhesion and to help prevent skin irritation, remove detergent residues and allow all skin preparations and protectants to dry completely before applying the dressing to the skin.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner. For proper use, clinicians should be trained in the use of this device.

**Instructions for Use:**

Failure to follow the manufacturer's instructions may result in complications including skin irritation and/or maceration.

**Site Preparation:**

1. Prepare the site in accordance to facility protocol.
2. Stabilize any active bleeding before applying the dressing.
3. Ensure the skin is clean, dry and free from detergent, lotion and soap residues.
4. Ensure all skin preps and protectants dry completely before applying the dressing. This will promote optimal adhesion to the skin and help minimize any skin irritation.
5. Clipping of hair at site may improve dressing adhesion. Shaving is not recommended.

**Dressing Selection**

Choose a dressing large enough to provide at least one-inch margin of adherence on dry, healthy skin around the catheter site. A summary of safety and clinical performance can be found at <https://ec.europa.eu/tools/eudamed>. This product is to be used by health care providers.

**Dressing Application:**

1. Prep the catheter site according to your facility's protocol. Let all prep solutions dry completely. Remove the liner from the dressing, exposing the adhesive surface of the dressing.
2. Center the transparent portion of the dressing over the insertion site. Do not stretch the dressing during application. Press the transparent portion of the dressing into place.
3. While slowly peeling off the paper frame, smooth down the dressing edges with your fingertips. Apply firm pressure to the entire dressing to enhance dressing adhesion.
4. After the dressing has been applied, apply additional tape to further secure I.V. tubing or to stabilize the catheter.
5. Document dressing change information on label according to facility protocol. Remove label from frame and place on dressing

**Removal:**

1. Remove any tape strips applied to the top of the dressing.
2. Using a low and slow removal technique, start removing the dressing from where the catheter or tubing exits the dressing toward the catheter insertion site. Avoid skin trauma by peeling the dressing back, rather than pulling it up from the skin.
3. If needed, a medical adhesive solvent can be used to help remove the dressing. Continue the low and slow removal method until the dressing is completely removed.
4. Dispose dressing per facility protocol.

**Site Care:**

1. The site should be observed at least daily for signs of infection or other complications. If infection is suspected, remove the dressing, inspect the site directly, and determine the appropriate medical intervention. Infection may be signaled by fever, pain, redness, swelling, or unusual discharge or odor.
2. Change the dressing as necessary, in accordance with facility protocols. At a minimum, dressings should be changed at least every 7 days. Change the dressing more frequently if the site has high levels of exudate or if the dressing integrity has become compromised.
3. The dressing should be changed sooner than 7 days if:
  - the dressing becomes loose, soiled or compromised in any way
  - the site is obscured or no longer visible
  - there is visible drainage underneath the dressing

**Storage/ Shelf Life/Disposal**

For best results, store in a cool, dry place. For shelf life, refer to the expiration date on the package. If the sterile packaging is damaged or unintentionally opened, discard the product, do not use.

**How supplied**

**3M™ Tegaderm™ Antimicrobial Transparent Dressing**

Catalog number	Dressing Size	Target Amount of CHG per dressing (mg)
9124	6 cm x 7 cm	4.04 mg

**3M™ Tegaderm™ Antimicrobial I.V. Advanced Securement Dressing (with borders)**

Catalog number	Dressing Size	Target Amount of CHG per dressing (mg)
9132	7 cm x 8.5 cm	5.65 mg


Please report serious incidents occurring in relation to the device to 3M and the local competent authority. For questions or comments, contact the 3M Health Care Customer Help Line at 1-800-226-3957.

**Glossary Table**

	Manufacturer	Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.
	Authorized Representative in European Community	Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
	Date of manufacture	Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
	Use-by date	Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
	Do not re sterilize	Indicates a medical device that is not to be re sterilized. Source: ISO 15223, 5.2.6
	Do not use if package is damaged or open	Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8
	Upper limit of temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed. Source: ISO 15223, 5.3.6
	Do not reuse	Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
	Consult instructions for use	Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3

	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
	Natural rubber latex is not present	Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
	Medical Device	Indicates the item is a medical device
	Contains medicinal substance	Indicates a medical device that contains or incorporates a medicinal substance
	Unique device identifier	Indicates bar code to scan product information into patient electronic health record
	Importer	Indicates the entity importing the medical device into the EU
	Single sterile barrier system	To identify a single layer that makes up the sterile barrier system (SBS) and differentiate it from layers of protective packaging designed to prevent damage to the sterile barrier system and its contents
	CE Mark 2797	Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
	Green Dot	Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

<b>Rx Only</b>	Rx Only	Indicates that U.S. Federal Law restricts this device to sale by or on the order of physician. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).
	Antimicrobial Agent- 2% CHG	Indicates product contains 2% Chlorhexidine gluconate as antimicrobial agent.

For more information, see [HCBRegulatory.3M.com](http://HCBRegulatory.3M.com)

For further information, please contact your local 3M representative or contact us at [3M.com](http://3M.com) and select your country.

**CE 2797**

Made in the USA with  
Globally Sourced Materials

**3M Company**  
2510 Conway Ave.  
St. Paul, MN 55144 USA  
1-800-228-3957 (USA Only)  
[3M.com/Tegaderm](http://3M.com/Tegaderm)

Issue Date: 2020-03  
34-8724-3503-6

**EC REP 3M Deutschland GmbH**  
**Health Care Business**  
Carl-Schurz-Str. 1  
41453 Neuss, Germany

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