ENGLISH	Cautions: Full Legal Restricts this device to sale or use on the order of a dental professional.
General Information 3M™ ESPE™ Ketac™ Cem Plus is a radiopaque, fluoride-releasing, resin-modified glass ionomer luting cement. It is self-curing with an option for light curing of excess cement. The dental cement consists of a base and catalyst paste packaged in the 3M™ ESPE™ Clicker™ Dispensing System. The paste/paste formulation offers greater convenience over the traditional powder/liquid cement systems. The cement is available in a white shade. The Clicker dispensing system contains 11g of material.	Warranty 3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of product for use. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.
Indications Final cementation of: • Metal and porcelain fused to metal (PFM) crowns and bridges • Metal inlays and onlays • Crowns and bridges made with all-alumina or all-zirconia cores such as 3M™ ESPE™ Lava™ or Procera™ AllCeram • Prefabricated and cast post • Orthodontic bands and appliances • Porcelain fused to metal (PFM), metal, all-alumina or all-zirconia core restorations on implant abutments	Limitation of Liability Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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Instructions for Use

Preparation

1. Prepare teeth using accepted clinical guidelines.

2. Cover areas in close proximity to the pulp by applying small amounts of hard-setting calcium hydroxide material or resin-modified glass ionomer liner (e.g. Tack Light Cure Glass Ionomer Liner/Base) prior to taking an impression for the final restoration.

3. Remove the temporary restoration and thoroughly clean the preparation of any temporary cement residue using an oil-free pumice paste.

4. Rinse the prepared tooth or cavity with a water spray, and dry with air, cotton, or a paper tip.

5. Leave tooth surface moist. Do not over dry.

6. The preparation cavity should be dry enough that the surface has a slightly glossy appearance. As is the case with any permanent cement, over drying can lead to poor operative sensitivity.

7. Try in final restoration and check for fit. Adjust if necessary.

8. Sandblast surfaces to be cemented with 30 or 50 micron aluminum oxide at a pressure of 2 Bar (30 psi) to create a matte surface appearance.

9. Thoroughly clean the bonding surface of the restoration.

10. Keep areas isolated from blood and saliva contamination during cementation process.

This product contains calcium persulfate, an allergic reaction trigger. Allergic reactions can occur in some people.

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Directions	Arbeits- und Aushärtezeit		
1. Application			
1. Remove cap from the Clicker dispenser by holding down the cap lever and sliding the cap off of the dispenser.	Arbeitszeit ab Mischbeginn:	2,5 Minuten	
1.2. For first time use from a new cartridge, dispense a small amount of material to ensure even dispensing. Discard this material.	Überbissgrünig Zement:		
1.2.1. To ensure even dispensing for each consequent use, dispense a small amount of material and discard. If the Clicker dispenser does not dispense evenly, check dispenser opening for blockage; remove any paste plugs.	Kurze Lichthärtung („Tack Cure“)	5 Sekunden pro Oberfläche	Für die Lichthärtung ein konventionelles Polymerisationsgerät verwenden
1.3. Fully depress clicker lever to dispense "1 Click" of cement onto mix pad. Allow paste to fully extrude, then release lever. Repeat dispensing process for additional material. The paste is automatically dispensed in equal volumes. (The actual weight ratio used is 1:3:1)	Selbsthärtende Gel-Phase	2 Minuten nach Einbringung in den Mund	
2. Storing The Used Clicker Dispenser	Aushärtezeit bei Selbsthärtung:	5 Minuten nach Einbringung in den Mund	
2.1. Wipe the dispenser tips clean, independently of one another with gauze or tissue to prevent cross contamination of the two pastes.			
2.2. Replace cap by sliding onto unit securely latched and an audible click is heard.			
2.2.1. Check to ensure that the cap is securely attached to the dispenser by pulling on the cap.			
2.2.2. If cap moves when pulled on, repeat steps 2.2 and 2.2.1 until the cap does not move when pulled on. This ensures that the cap is properly attached to the dispenser.			
3. Mixing The Cement			
3.1. Note: Dispense and mix the cement immediately to avoid water evaporation and drying out of the cement pastes.			
3.2. Using a plastic or metal cement spatula, mix the pastes together for 20 seconds until a uniform color is achieved. Avoid the incorporation of air bubbles.			
3.3. The working time of the mixed cement is 2.5 minutes at 23°C/73°F. Higher temperatures and longer mix times may shorten the working time. Lower temperatures will lengthen the working time.			
3.4. Apply a thin layer of cement to the inside surface of the restoration. The cement may also be applied directly to the tooth surface for inlay/onlay restoration.			
3.5. Seat the restoration with light pressure. Maintain light pressure on the restoration to hold proper positioning during setting process.			
4. Clean Up and Finishing			
4.1. Excess cement is best removed after brief light exposure (5 seconds per surface with a conventional polymerization device) or during self-cured waxy stage (starting 2 minutes after seating in the mouth). Use a suitable instrument (e.g. scaler, explorer) for this process.			
4.2. Finish and polish restoration when material has completely set after 5 minutes from placement in the mouth.			

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3. **Mixing The Cement**
 - 3.1. **Note:** Dispense and mix the cement immediately to avoid water evaporation and drying out of the cement pastes.
 - 3.2. Using a plastic or metal cement spatula, mix the pastes together for 20 seconds until a uniform color is achieved. Avoid the incorporation of air bubbles.
 - 3.3. The working time of the mixed cement is 2.5 minutes at 23°C/73°F. Higher temperatures and longer mix times may shorten the working time. Lower temperatures will lengthen the working time.
 - 3.4. Apply a thin layer of cement to the inside surface of the restoration. The cement may also be applied directly to the tooth surface for inlay/onlay restorations.
 - 3.5. Seat the restoration with light pressure. Maintain light pressure on the restoration to hold proper positioning during setting process.
4. **Clean Up and Finishing**
 - 4.1. Excess cement is best removed after brief light exposure (5 seconds per surface with a conventional polymerization device) or during self-cured waxy stage (starting 2 minutes after seating in the mouth). Use a suitable instrument (e.g. scaler, explorer) for this process.
 - 4.2. Finish and polish restoration when material has completely set after 5 minutes from placement in the mouth.

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4. Clean Up and Finishing

- 4.1. Excess cement is best removed after brief light exposure (5 seconds per surface with a conventional polymerization device) or during self-curing waxy stage (starting 2 minutes after seating in the mouth). Use a suitable instrument (e.g. scaler, explorer) for this process.
- 4.2. Finish and polish restoration when material has completely set after 5 minutes from placement in the mouth.

Storage and Use

This product is designed to be used at room temperature. If stored in a refrigerator, allow product to reach room temperature prior to use. Shelf life at room temperature is 24 months. Ambient temperatures routinely greater than 28°C/80°F may reduce shelf life. See outer package for expiration date.

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References: CDC Guidelines for Disinfection and Sterilization in Healthcare Facilities, 2008, Rutala WJ, 2008, Rutala
--

Working Time From Start Of Mixing:	2.5 Minutes	For Tack Light Curing, Use A Conventional Polymerization Device
Excess Cement:	5 Seconds Per Surface	
Tack Light Cure	2 Minutes After Placement In Mouth	
Self-Curing Gel Phase	5 Minutes After Placement In Mouth	
Self-Curing Setting Time		
Directions		
1. Application		
1. Remove cap from the Clicker dispenser by holding down the cap lever and sliding the cap off of the dispenser.		
1.2. For first time use from a new cartridge, dispense a small amount of material to ensure even dispensing. Discard this material.		
1.2.1. To ensure even dispensing for each consequent use, dispense a small amount of material and discard. If the Clicker dispenser does not dispense evenly, check dispenser opening for blockage; remove any paste plugs.		
1.3. Fully depress clicker lever to dispense "1 Click" of cement onto mix pad. Allow paste to fully extrude, then release lever. Repeat dispensing process until the desired amount of material is automatically dispensed in equal volumes. (The actual weight ratio used is 1:3:1)		
2. Storing The Used Clicker Dispenser		
2.1. Wipe the dispenser tips clean, independently of one another with gauze or tissue to prevent cross contamination of the two pastes.		
2.2. Replace cap by sliding onto unit securely latched and an audible click is heard.		
2.2.1. Check to ensure that the cap is securely attached to the dispenser by pulling on the cap.		
2.2.2. If cap moves when pulled on, repeat steps 2.2 and 2.2.1 until the cap does not move when pulled on. This ensures that the cap is properly attached to the dispenser.		
3. Mixing The Cement		
3.1. Note: Dispense and mix the cement immediately to avoid water evaporation and drying out of the cement pastes.		
3.2. Using a plastic or metal cement spatula, mix the pastes together for 20 seconds until a uniform color is achieved. Avoid the incorporation of air bubbles.		
3.3. The working time of the mixed cement is 2.5 minutes at 23°C/73°F. Higher temperatures and longer mix times may shorten the working time. Lower temperatures will lengthen the working time.		
3.4. Apply a thin layer of cement to the inside surface of the restoration. The cement may also be applied directly to the tooth surface for inlay/onlay restoration.		
3.5. Seat the restoration with light pressure. Maintain light pressure on the restoration to hold proper positioning during setting process.		
4. Clean Up and Finishing		
4.1. Excess cement is best removed after brief light exposure (5 seconds per surface with a conventional polymerization device) or during self-cured waxy stage (starting 2 minutes after seating in the mouth). Use a suitable instrument (e.g. scaler, explorer) for this process.		
4.2. Finish and polish restoration when material has completely set after 5 minutes from placement in the mouth.		

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