3M ESPE

RelyXTM Ultimate ClickerTM

- (en) Adhesive Resin Cement
- (ru) Цемент адгезивный композитный двойного отверждения RelyX Ultimate. Дозирующая система Clicker с цементом.
- (bg) Адхезивен композитен цимент
- (hr) Cement za adhezivno cementiranje
- (hu) Kompozit rögzítőcement
- (pl) Adhezyjny cement żywiczny
- (ro) Cement răšinic adeziv
- (sk) Kompozitový upevňovací cement
- (sl) Adhezivni kompozitni cement
- (cs) Kompozitní upevňovací cement
- (tr) Adeziv Rezin Siman
- (et) Adhesiivne vaiktsement
- (lv) Kompožīta nostiprināšanas cements
- (lt) Dervinis cementas fiksavimui
- (uk) Адгезивний композитний цемент
- (es) Cemento de fijación de composite
- (pt) Cimento de fixação compósito
- (ko) 접착형 레진시멘트
- (zh) 樹脂黏合劑
- (th) สารแอดไฮซีฟ เรชิน ชีเมนต
- (fr) Ciment de scellement composite adhésif

Instructions for Use

Инструкция по применению
Указания за употреба
Upute o upotrebi
Felhasználási tudnivalók
Instrukcja użycia
Instrucțiuni de utilizare
Informácie pre použitie
Navodila za uporabo
Informace pro uživatele
Kullanma Talimatları

Kasutusjuhend

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Adhesive Resin Cement

CONTENTS

| | |
|--|---|
| Product Description | 2 |
| Indications | 2 |
| Contraindications | 2 |
| Precautionary Measures | 2 |
| For Patients and Dental Personnel | 2 |
| For Patients | 2 |
| For Dental Personnel | 2 |
| Precautions during Processing | 2 |
| Pulp Protection | 2 |
| Tooth Preparation | 2 |
| Preparation of the Cavity/Tooth Stump | 2 |
| Preparation of Root Canals | 3 |
| Preparation for Veneers | 3 |
| Preparation for Maryland and Inlay/Onlay Bridges | 3 |
| Preparation for Chairside CAD/CAM Restorations | 3 |
| Trial Fitting of the Restoration | 3 |
| Pre-treatment of Restorations | 3 |
| Pre-treatment of Metal Restorations | 3 |
| Pre-treatment of Etchable Glass Ceramic Restorations | 3 |
| Pre-treatment of Zirconia and Aluminum Oxide Ceramic Restorations | 3 |
| Pre-treatment of Composite and Resin Nano Ceramic (e.g. Lava Ultimate) Restorations | 4 |
| Pre-treatment of Glass Fiber-reinforced Posts | 4 |
| Pre-treatment of Maryland and Inlay/Onlay Bridges | 4 |
| Etching of the Tooth Structure (Optional, depending on material type and indication) | 4 |
| Dosing and Applying Adhesive | 4 |
| Application of the Adhesive | 5 |
| Dosing Cement | 5 |
| Times | 5 |
| Mixing | 5 |
| Application to Tooth or Restoration | 5 |
| Application in the Root Canal | 5 |
| Removal of Excess | 5 |
| Option 1 for best aesthetic results | 6 |
| Option 2 for fast removal of excess | 6 |
| Polymerization and Shaping | 6 |
| After Shaping | 6 |
| Hygiene and Disinfection | 6 |
| Disposal | 6 |
| Notes | 6 |
| Storage and Stability | 6 |
| Customer Information | 6 |
| Warranty | 6 |
| Limitation of Liability | 6 |

Product Description

RelyX™ Ultimate Clicker™ is a dual-curing resin cement for manual mixing supplied in a Clicker dispenser. It is used for the adhesive cementation of indirect restorations and available in various shades. RelyX Ultimate Clicker is used in combination with Single Bond Universal Adhesive. The adhesive can be used either as a "Total Etch" or "Self-Etch" procedure. The additional etching of the tooth structure with Scotchbond™ Universal Etchant increases the adhesive strength of the adhesive to the tooth structure even further. RelyX Ultimate Clicker contains bi-functional (meth)acrylates. The proportion of inorganic fillers is about 43% by volume; the grain size (D 90%) is about 13 µm. The mixing ratio, based on volume, is 1 part base paste : 1 part catalyst.

 These Instructions for Use must be kept for reference for the duration of product use. Please see the pertinent information for use for details on all of the products mentioned below. The product may only be used when the product labeling is clearly readable.

Indications

Scotchbond Universal Etchant:

- Etching of enamel and dentin using the acid etching technique

RelyX Ultimate Clicker in combination with Single Bond Universal as an adhesive for:

- Final cementation of all-ceramic, composite, or metal inlays, onlays, crowns and bridges; 2-3-unit Maryland bridges and 3-unit inlay/onlay bridges
- Final cementation of posts and screws
- Final cementation of all-ceramic, or composite veneers
- Final cementation of all-ceramic, composite, or metal restorations on implant abutments

Indications such as Maryland and inlay/onlay bridges (resin bonded bridges) demand especially high adhesive bond strength. Regardless of the manufacturer of the cement and restoration, these indications may be exposed to a higher risk of decementation. For Resin Nano Ceramic restorations (e.g. Lava™ Ultimate) it is important to maximize bond strength. For the detailed bonding procedure, please refer to the instructions for use of the restoration material.

Contraindications

Maryland and inlay/onlay bridges are contraindicated for patients with bruxism or periodontitis.

Precautionary Measures

For Patients and Dental Personnel

- Single Bond Universal Adhesive: contact with eyes may cause severe eye damage. Wear eye protection to prevent injury. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- RelyX Ultimate Clicker base paste: contains sodium persulfate that may trigger an allergic reaction in certain individuals when inhaled. This product may not be appropriate for use in those individuals with known sensitivity to sulfites since a cross-reaction may occur with sodium persulfate.

- Scotchbond Universal Etchant: contact with eyes or skin may cause burns. Wear suitable protective clothing, gloves, and eye/face protection. In case of contact, flush area with plenty of water and seek medical advice.

For Patients

- RelyX Ultimate Clicker and Single Bond Universal Adhesive contain substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of these products in patients with known acrylate or perester allergies.
- If prolonged contact with oral tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary, and discontinue future use.

For Dental Personnel

- RelyX Ultimate Clicker and Single Bond Universal Adhesive contain substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured paste. If skin contact occurs, wash immediately with soap and water.
- The use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If RelyX Ultimate Clicker or Single Bond Universal Adhesive contacts the glove, remove and discard the gloves. Wash your hands immediately with soap and water and re-glove.
- If allergic reaction occurs, seek medical attention as needed.
- Single Bond Universal Adhesive is highly flammable. Avoid static electricity build-up. Keep away from heat, sparks, open flame, pilot lights, and other ignition sources.

3M MSDSs can be obtained from www.mmm.com, or contact your local subsidiary.

Precautions during Processing

- Seat temporary restorations using a eugenol-free product (e.g., with RelyX™ Temp NE). The use of temporary cements containing eugenol can inhibit the polymerization process of RelyX Ultimate Clicker during the final cementing.
- Do not use hydrogen peroxide (H_2O_2) or sodium bicarbonate ($NaHCO_3$) as the residues can adversely affect the adhesive strength and setting reaction of RelyX Ultimate Clicker.
- For the best possible bond, do not treat the surface of Zirconia, aluminum oxide, or metal restorations with phosphoric acid.

Pulp Protection

To avoid pulpal irritation, cover areas in close proximity to the pulp by applying small amounts of hard-setting calcium hydroxide material. To avoid bite increases during the cementing of the later restoration, the pulp protection should be applied prior to taking an impression for the final restoration.

Tooth Preparation

Preparation of the Cavity/Tooth Stump

- Remove loose preparation debris by spraying with water.

- ▶ Clean the stump or the cavity thoroughly, rinse with a water spray, and lightly air dry in only 2–3 bursts of water-free and oil-free air, or use cotton pellets or similar materials to dry it off. Do not overdry.
- The cavity should be just dry enough that the surface has a slightly glossy appearance. As is the case with any permanent cement, overdrying can lead to post-operative sensitivity.

Do not use substances such as desensitizers, disinfectants, astringents, dentin sealants, rinsing solutions containing EDTA, scanning powder, etc., after the final cleaning. Their residues may have a detrimental effect on the bonding strength and setting reaction of the adhesive.

Preparation of Root Canals

- ▶ Treat the root canal endodontically as usual (root canal filling with guttapercha and removal of the root canal filling, leaving 4 mm of guttapercha apically).
- ▶ Clean the root canal with a 2.5%–5.25% solution of sodium hypochlorite (NaOCl).
- ▶ Rinse immediately with water, dry with paper tips, and continue with the application of the adhesive; see "Dosing and Applying Adhesive".
- The use of a rubber dam during the cementation of posts is recommended.

Preparation for Veneers

- ▶ Prepare and clean the tooth as is usual for veneers.
- ▶ Subsequently, it must be etched, as described under "Etching of the Tooth Structure."

Preparation for Maryland and Inlay/Onlay Bridges

Abutment teeth must have an adequate enamel surface for bonding. They should be healthy or only slightly restored and the periodontal conditions should be good. It is the sole responsibility of the dentist to ensure proper selection of indication and technique. The guidelines of the relevant national professional associations must be observed for such indications.

- ▶ Prepare retentive elements such as cingular rests and/or approximal grooves.
- ▶ Subsequently, it must be etched, as described under "Etching of the Tooth Structure."
- ▶ Use a rubber dam and keep the cavity free of any contamination during cementation.

Preparation for Chairside CAD/CAM Restorations

Please follow the instructions for use of the restoration material. For Resin Nano Ceramic restorations (e.g. Lava Ultimate), please refer to the Lava Ultimate Instructions for Use.

Trial Fitting of the Restoration

- ▶ Remove the temporary restoration and thoroughly clean the tooth structure of any residue from the temporary cement.
- ▶ Make a trial insertion of the final restoration and check its fit, color, and contact points. Use RelyX™ Try-In Paste or similar material for this purpose.

- If a low viscosity silicone is used for the trial fit, the silicone residue must then be removed completely, e.g., with alcohol.
- For glass ceramic restorations, do not check the occlusion until the restoration has been cemented. Breakage could occur if done prior to being cemented in place.

Pre-treatment of Restorations

Avoid any contamination whatsoever of the surfaces to be treated during pre-treatment and until the final cementation.

Pre-treatment of Metal Restorations

Please follow the instructions for use of the restoration material. In the absence of deviating instructions, we recommend the following procedure:

- ▶ Sandblast surfaces to be cemented. Use 30 or 50 μm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance.
- ▶ Clean the blasted surface with alcohol and dry it with water-free and oil-free air.
- ▶ Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl 5%) solution is recommended), rinsed with water, and dried.
- ▶ Use the disposable applicator to apply Single Bond Universal Adhesive to the entire surface to be luted and allow it to react for 20 sec.
- During treatment, rewet the disposable applicator as required.
- ▶ Subsequently direct a gentle stream of air over the liquid for about 5 sec until it no longer moves and the solvent has evaporated completely.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Pre-treatment of Etchable Glass Ceramic Restorations

Please follow the instructions for use of the restoration material. In the absence of deviating instructions, we recommend the following procedure:

- ▶ Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.
- ▶ Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.
- ▶ Use the disposable applicator to apply Single Bond Universal Adhesive to the entire surface to be luted and allow it to react for 20 sec.
- During treatment, rewet the disposable applicator as required.
- ▶ Subsequently direct a gentle stream of air over the liquid for about 5 sec until it no longer moves and the solvent has evaporated completely.

Alternatively, instead of Single Bond Universal Adhesive any silane can be applied in accordance with its Instructions for Use after etching, e.g., RelyX™ Ceramic Primer: allow to react for 5 sec, then air dry completely so that the solvent evaporates completely.

Pre-treatment of Zirconia and Aluminum Oxide Ceramic Restorations

Please follow the instructions for use of the restoration material. In the absence of deviating instructions, we recommend the following procedure:

- ▶ Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance.
- ▶ Clean the blasted surface with alcohol and dry it with water-free and oil-free air.
- ▶ Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl 5%) solution is recommended), rinsed with water, and dried.
- ▶ Use the disposable applicator to apply Single Bond Universal Adhesive to the entire surface to be luted and allow it to react for 20 sec.
 - During treatment, rewet the disposable applicator as required.
- ▶ Subsequently direct a gentle stream of air over the liquid for about 5 sec until it no longer moves and the solvent has evaporated completely.

As an alternative, the restoration surface to be luted can be coated (silicated) with the micro-blasting device CoJet™ Prep and the blast-coating agent CoJet™ Sand from a distance of 2–10 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prep and CoJet Sand. Then blow away any residues of the blasting agent with water-free and oil-free air. Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted, rub in for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec, or use a suitable silane in accordance with the instructions for use.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Pre-treatment of Composite and Resin Nano Ceramic (e.g. Lava Ultimate) Restorations

Please follow the instructions for use of the restoration material. In the absence of deviating instructions, we recommend the following procedure:

- ▶ Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance. Do not etch with HF or phosphoric acid.
- ▶ Clean the blasted surface with alcohol and dry it with water-free and oil-free air.
- ▶ Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl 5%) solution is recommended), rinsed with water, and dried.
- ▶ Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted, rub in for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec. As an alternative, a suitable silane can be used in accordance with the instructions for use.

Pre-treatment of Glass Fiber-reinforced Posts

Please follow the instructions for use of the post to be used. In the absence of deviating instructions, we recommend the following procedure:

- ▶ Clean the post with alcohol and dry it with water-free and oil-free air.
- ▶ Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted, rub in for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec. As an alternative, a suitable silane can be used in accordance with the instructions for use.

Pre-treatment of Maryland and Inlay/Onlay Bridges

In the absence of deviating instructions, we recommend the following procedure for metal, Zirconia, and aluminum oxide surfaces:

- ▶ Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet Prep and the blast-coating agent CoJet Sand from a distance of 2–10 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prep and CoJet Sand.
- ▶ Blow away any residues of the blasting agent with water-free and oil-free air.
- ▶ Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted, rub in for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec. As an alternative, a suitable silane can be used in accordance with the instructions for use.

For Maryland and inlay/onlay bridges made of etchable glass ceramics, please refer to "Pre-treatment of Etchable Glass Ceramic Restorations".

Etching of the Tooth Structure (Optional, depending on material type and indication)

A clinically adequate adhesive bond is achieved by the application of RelyX Ultimate Clicker in combination with Single Bond Universal Adhesive. However, selective etching of the enamel or the etching of the enamel and dentin ("Total Etch Procedure") can produce higher adhesive values.

For veneers, Maryland and inlay/onlay bridges in general as well as for all Lava Ultimate restorations, higher adhesive strengths are necessary; therefore, etching is mandatory with these indications.

- ▶ The application of a rubber dam is recommended.
- ▶ Remove the Scotchbond Universal Etchant syringe from the foil pouch, unscrew the cap, and attach an Etchant Dispensing Tip.
- ▶ Apply the etching gel to the precise area of the enamel and/or dentin and allow to react for 15 sec. Rinse thoroughly for 15 sec with water, using a suction at the same time.
- ▶ Dry with water-free and oil-free air or with cotton pellets; do not overdry.
- ▶ Avoid any contamination of the etched surfaces. If contamination occurs, repeat the etching, rinsing, and drying.
- ▶ Continue immediately with the application of the adhesive.
- ▶ After etching, remove the Etchant Dispensing Tip from the syringe and replace the cap.

Dosing and Applying Adhesive

The dosing bottle cap has a thumb depression so that it can be opened to one side and closed again with one hand. The cap will remain open if the opening angle is >90°.

- ▶ To open, take the bottle in your hand, place your thumb in the thumb depression, and push the cap to the side until it remains in the open position.
- ▶ Hold the bottle upside down in a vertical position and pour the required quantity of Single Bond Universal Adhesive in a mixing well. Protect the adhesive in the mixing well from light.

- ▶ After use, carefully close the bottle again by applying pressure to the lever from above.
 - The cap tangibly snaps into place when pressed closed and a small, even gap between the cap and bottle can be seen all the way around.

Application of the Adhesive

Avoid any contamination with blood, saliva, or sulcus fluid during application. A rubber dam is recommended to keep the treated area sufficiently dry.

- ▶ Use the disposable applicator to apply the adhesive to the entire tooth structure and **rub it in** for 20 sec. Avoid contact between the adhesive and the oral mucosa.
- ▶ If necessary, rewet the disposable applicator during treatment.
- ▶ Blow a gentle stream of air over the liquid for about 5 sec until it no longer moves and the solvent has evaporated completely. Remove any excess in the root canal with a paper tip and blow a gentle stream of air over it for 5 sec as well.
- ▶ When cementing indirect restorations or posts, Single Bond Universal Adhesive is not light-cured, but it can optionally be polymerized for 10 sec using any commonly available light-curing device with a lightwave range of 400–500 nm.

Note: For maximum adhesion of Resin Nano Ceramic restorations (e.g. Lava Ultimate), light cure the adhesive for 10 seconds. For detailed bonding procedure, please refer to the instructions for use of the restoration material.

Dosing Cement

When dosing the first time:

- ▶ Remove the Clicker dispenser from the foil pouch and note the removal date on the Clicker label. Discard the pouch.
- ▶ Remove the cap from the Clicker dispenser.
- ▶ Dispense a small quantity of paste on a mixing pad by pressing down the Clicker lever. This levels the filling quantity of base paste and catalyst paste in the Clicker. Discard the paste which has been dispensed.

When dosing, always:

- ▶ Remove the cap.
- ▶ Dispense the paste on a mixing pad:
 - Press the Clicker lever and hold it down.
 - When the paste stops flowing, wipe off the Clicker openings on the mixing pad. Doing so dispenses the correct volume of the base paste and catalyst paste.
 - Dispense only as much paste on the mixing pad (number of clicks) as required for cementing.
- ▶ Cleaning the Clicker openings:
 - Wipe the Clicker openings with a clean cloth to prevent blockage from cross-contamination of the pastes.
- ▶ Replace the cap by pressing until it locks into place noticeably.

Times

The processing and setting times of RelyX Ultimate Clicker depend on the ambient and oral temperature. The times shown are based on conditions relevant for practice. As is

the case with every composite cement, the setting of RelyX Ultimate Clicker slows down significantly at room temperature. RelyX Ultimate Clicker is a dual-curing cement and therefore also sensitive to natural or artificial light. The working time is significantly reduced during application under operating lights.

| | RelyX Ultimate Clicker min:sec | Single Bond Universal Adhesive min:sec |
|---|---|---|
| Working time from start of mixing | 02:30 | — |
| Light-curing in a lightwave range of 400–500 nm | | (optional) 00:10 |
| - single surface, from occlusal | 00:20 | — |
| - any other surface additional | 00:20 | — |
| - RelyX Fiber Post posts, from occlusal | 00:40 | — |
| Self-curing | | |
| - Setting time after start of mixing | 06:00 | — |

Mixing

- ▶ Mix base paste and catalyst paste into a homogenous paste within 20 seconds using a spatula. Avoid incorporating air bubbles.
- ▶ Then cover the cement with a light protection cover to protect it from intensive light.

Application to Tooth or Restoration

- ▶ Protect the working area from contamination with water, blood, saliva, and sulcal fluid during the application and setting phase.
- ▶ Wet the entire area of the cavity walls, the cavity floor, and, if necessary, the surface of the restoration to be cemented evenly with RelyX Ultimate Clicker.
- ▶ Seat the restoration and stabilize long enough for the cement to fully set.

For veneers: owing to the relatively short processing time of RelyX Ultimate Clicker (dual curing begins after 2 min 30 sec), we recommend cementing a maximum of 2 veneers at the same time.

Application in the Root Canal

Do not use Lentulo spirals to insert the cement in the root canal as this can excessively accelerate setting.

- ▶ Spread cement to the post.
 - As an alternative, the cement can be inserted into the root canal with the aid of an applicator tip, e.g., Centrix® NeedleTubes (Centrix is a registered trademark of the U.S.-based company Centrix); this will minimize the entrapment of air bubbles.
- ▶ Immediately place the post in the pretreated root canal; apply moderate pressure to hold it in position. We recommend rotating the post slightly during insertion to avoid the inclusion of air bubbles.

Removal of Excess

Notes on oxygen inhibition: As is the case with all composite materials, RelyX Ultimate Clicker is subject to oxygen inhibi-

tion, i.e., the upper layer (about 50 µm) which is in contact with atmospheric oxygen during the polymerization does not harden. Leaving sufficient excess permits the removal of the uncured layer during shaping/polishing without leaving behind a deficit.

Option 1 for best aesthetic results

- ▶ Remove the excess before polymerization, e.g., using a sponge pellet, while using a suitable instrument to hold the restoration in position.
- ▶ Cover the margins of the restoration with a glycerin gel to prevent oxygen inhibition.
- ▶ Immediately light-cure or wait for the self-curing to finish (see Polymerization and Shaping).

Option 2 for fast removal of excess

- ▶ Remove cement excess either after brief light-curing (about 1–2 seconds using a commonly available polymerization device), or use a suitable instrument (e.g., scaler) to remove the cement excess during the self-curing (from 3 minutes after start of mixing in the "gel phase"). Larger volumes of excess are easier to remove.
 - If the excess is removed during the self-curing, a suitable instrument must be used to hold the restoration in position.

Polymerization and Shaping

- ▶ We recommend light-curing the cement through the restoration when doing ceramic and composite work. Select the exposure times appropriate for the number of surfaces (please refer to "Times"). Polymerize the cement through the post when using translucent posts. The exposure time depends on the translucency of the post being used; for RelyX Fiber Post, it is 40 sec.
- ▶ Shape any remaining uneven areas on the margins of restorations and polish the marginal area with diamond devices, aluminum oxide-coated discs (e.g., Sof-Lex™), and diamond polishing paste.
- ▶ Then check the occlusion.

After Shaping

- ▶ Carefully check the sulcus of the treated teeth and the surrounding areas; remove any cement residues still remaining. This is especially important when cementing to implant abutments.

Hygiene and Disinfection

- ▶ Clean any contaminated RelyX Ultimate Clicker dispensers, Single Bond Universal Adhesive bottles, and Scotchbond Universal Etchant syringes with commonly used agents and disinfect them in accordance with the commonly used hygienic regulations.

Disposal

- Scotchbond Universal Etchant syringes must be completely emptied before being disposed of for recycling. Small amounts of Etching Gel can be disposed of with normal waste water when well diluted with water or when neutralized.

- Dispose of empty RelyX Ultimate Clicker dispensers and Single Bond Universal Adhesive bottles with the normal household waste.

Notes

- RelyX Ultimate Clicker and Single Bond Universal Adhesive can prematurely polymerize when exposed to natural or artificial lighting. Avoid intensive light exposure during application.

Storage and Stability

Store RelyX Ultimate Clicker in the foil pouch at 15°–25° C / 59°–77° F. After removal from the foil pouch, use RelyX Ultimate Clicker within 6 months and before the lapse of the expiration date. Constantly high humidity influences the setting and must be avoided.

Store Single Bond Universal Adhesive and Scotchbond Universal Etchant at 0°–25° C / 36°–77° F.

Do not use after the expiration date.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

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- อีกทางเลือกหนึ่งคือสามารถใช้ซีเมนต์ลงในคลองรากฟันโดยใช้ปลายสำหรับจีดวัสดุ เช่น Centrix® Needle Tubes (Centrix is a registered trademark of the U.S.-based company Centrix) ซึ่งจะช่วยให้มีฟองอากาศในเนื้อวัสดุอย่างมาก
3M, ESPE, Clicker, CoJet, RelyX, Scotchbond and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH.
© 2019, 3M. All rights reserved.