ENGLISH

General Information

3M™ ESPE™ RelyX™ Luting Plus is a radiopaque, fluoride-releasing, resin-modified glass ionomer luting cement. It is self-curing with an option for tack light curing of excess cement. The dental cement consists of a base and catalyst paste packaged in the 3M™ ESPE™ Clicker™ Dispensing System. The paste/paste formulation offers greater convenience over the traditional powder/liquid cement systems. The cement is available in a white shade. The Clicker dispensing system contains 11g of material.

Indications

Final cementation of:

- Metal and porcelain fused to metal (PFM) crowns and bridges
- Metal inlays and onlays
- Crowns and bridges made with all-alumina or all-zirconia cores such as 3M™ ESPE™ Lava™ or Procera® AllCeram
- Prefabricated and cast post
- Orthodontic bands and appliances
- Porcelain fused to metal (PFM), metal, all-alumina or all-zirconia core restorations on implant abutments

Precautionary Information for Patients:


This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

This product contains potassium persulfate, which may produce an allergic reaction. Potassium persulfate may trigger an allergic respiratory reaction in certain individuals. This product may not be appropriate for use in those individuals with known sensitivity to sulfites, since a cross-reaction may occur with potassium persulfate.

Precautionary Information for Dental Personnel:

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M SDSs can be obtained by visiting www.3M.com or contact your local subsidiary.

3M ESPE
RelyX™ Luting Plus
Resin Modified Glass Ionomer Cement
Ciment en verre ionomère modifié à la résine
Cemento de ionómero de vidrio modificado con resina

Instructions for Use

Preparation

1. Prepare teeth using accepted clinical guidelines.
2. Cover areas in close proximity to the pulp by applying small amounts of hard-setting calcium hydroxide material or resin-modified glass ionomer liner (e.g. 3M™ ESPE™ Vitrebond™ Light Cure Glass Ionomer Liner/Base) prior to taking an impression for the final restoration.
3. Remove the temporary restoration and thoroughly clean the preparation of any temporary cement residue using an oil-free pumice paste.
4. Rinse the prepared tooth or cavity with a water spray, and dry with air, cotton, or a paper tip.
 - 4.1. Leave tooth surface moist. Do not over dry.
 - 4.2. The preparation or cavity should be just dry enough that the surface has a slightly glossy appearance. As is the case with any permanent cement, over drying can lead to post operative sensitivity.

5. Try in final restoration and check for fit. Adjust if necessary.
6. Sandblast surfaces to be cemented with 30 or 50 micron aluminum oxide at a pressure of 2 Bar (30 psi) to create a matte surface appearance.
7. Thoroughly clean the bonding surface of the restoration.
8. Keep area isolated from blood and saliva contamination during cementation process.

Working and Setting Times

Working Time From Start Of Mixing:	2.5 Minutes	
Excess Cement:		
Tack Light Cure	5 Seconds Per Surface	For Tack Light Curing, Use A Conventional Polymerization Device
Self-Curing Gel Phase	2 Minutes After Placement In Mouth	
Self-Curing Setting Time:	5 Minutes After Placement In Mouth	

Directions

1. Application

- 1.1. Remove cap from the Clicker dispenser by holding down the cap lever and sliding the cap off of the dispenser.
- 1.2. For first time use from a new cartridge, dispense a small amount of material to ensure even dispensing. Discard this material.
 - 1.2.1. To ensure even dispensing for each consequent use, dispense a small amount of material and discard. If the Clicker dispenser does not dispense evenly, check dispenser opening for blockage; remove any paste plugs.
- 1.3. Fully depress clicker lever to dispense “1 Click” of cement onto mix pad. Allow paste to fully extrude, then release lever. Repeat dispensing process for additional material. The paste is automatically dispensed in equal volumes. (The actual weight ratio dispensed is 1.3:1)

2. Storing The Used Clicker Dispenser

- 2.1. Wipe the dispenser tips clean, independently of one another with gauze or tissue to prevent cross contamination of the two pastes.
- 2.2. Replace cap by sliding onto dispenser until securely latched and an audible click is heard.**
 - 2.2.1. Check to ensure that the cap is securely attached to the dispenser by pulling on the cap.
 - 2.2.2. If cap removes when pulled on, repeat steps 2.2 and 2.2.1 until the cap does not remove when pulled on. This ensures that the cap is properly attached to the dispenser.

3. Mixing The Cement

- 3.1. **Note:** Dispense and mix the cement immediately to avoid water evaporation and drying out of the cement pastes.
- 3.2. Using a plastic or metal cement spatula, mix the pastes together for 20 seconds until a uniform color is achieved. Avoid the incorporation of air bubbles.
- 3.3. The working time of the mixed cement is 2.5 minutes at 23°C/73°F. Higher temperatures and longer mix times may shorten the working time. Lower temperatures will lengthen the working time.
- 3.4. Apply a thin layer of cement to the inside surface of the restoration. The cement may also be applied directly to the tooth surface for inlay/onlay restorations.
- 3.5. Seat the restoration with light pressure. Maintain light pressure on the restoration to hold proper positioning during setting process.

4. Clean Up and Finishing

- 4.1. Excess cement is best removed after brief light exposure (5 seconds per surface with a conventional polymerization device) or during self-curing waxy stage (starting 2 minutes after seating in the mouth). Use a suitable instrument (e.g. scaler, explorer) for this process.
- 4.2. Finish restoration and check occlusion when material has completely set after 5 minutes from placement in the mouth.

Storage and Use

This product is designed to be used at room temperature. If stored in cooler allow product to reach room temperature prior to use. Shelf life at room temperature is 24 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.

Do not allow the pastes to dry out. Store with cap securely attached to the dispenser.

Cleaning and Disinfection

Disinfect the product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control (CDC) and endorsed by the American Dental Association. Centers for Disease Control and Prevention. Summary of Infection Prevention Practices in Dental Settings: Basic Expectations for Safe Care. Atlanta, GA: Centers for Disease Control and Prevention, US Dept of Health and Human Services; October 2016.

References:

CDC Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008. Rutala WA, Weber DJ, Healthcare Infection Control Practices Advisory Committee (HICPAC).

Disposal – See the Safety Data Sheet (available at www.3M.com or through your local subsidiary) for disposal information.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Caution: U.S. Federal Law restricts this device to sale or use on the order of a dental professional.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and

3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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Procera  is a registered trademark of Nobel Biocare AB.
3M, ESPE, RelyX, Clicker, Lava and Vitrebond are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada.
  3M 2018. All rights reserved.