3M™ Steri-Drape™ loban™ 2 Specialty Drapes

Indication for Use
3M™ Steri-Drape™ loban™ 2 Specialty Drapes are used to create a sterile field for the surgical procedure. They are intended for external use only.

Contraindications

Do not use on patients with known sensitivity to iodine.

Warning

• Do not defibrillate through plastic drapes as arching or ineffective defibrillation may result.

• Do not re-sterilize.

• Do not ethylene oxide sterilize.

Directions for Use

1. Disinfect the surgical or operative site according to hospital policy or physician-prescribed routine.
2. Make sure skin is dry.
3. Apply drape to thoroughly dry, intact skin; smooth with a sterile towel.
4. To remove an adherent area, gently peel back at 180-degree angle from the skin. This means folding the adherent back on itself and then pulling it off the skin (see Figure A and B). During the removal process, hold onto the adhered area with the thumb and hand as close as possible to the junction of the skin and the drape. This minimizes the stretching. Skin of elderly patients is typically fragile and thus requires even greater care when applying and removing adhesive products such as drapes.

Storage Conditions

For best results, store at room temperature. Avoid excessive heat and humidity. For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

Product Details
The 3M™ Steri-Drape™ loban™ 2 Specialty Drape products meet the requirements of EN 13795-2:2019 Surgical clothing and drapes – Requirements and test methods – Part 1: Surgical drapes and gowns.

High Performance Critical Area Catalog Numbers: 6665, 6677, 6677P, 6678, 6681, 6682, 6687, 6687P, 6697, 6697P, 6697CA

Standard Performance Critical Area Catalog Numbers: 6617, 6617P001, 6619, 6619P001, 6657, 6658, 6659

*Exclusions: EN 13795:2019 does not cover requirements for incise drapes or films. For this reason, the incise portion of the drape has not been tested. EN 13795:2019 covers only the critical areas of the drape. The critical area is defined as the fabric or plastic portion of the drape and does not include the product seams or other components that may fall within the critical area such as Tube and Cord Organizers and Fluid Pouches.

Explanation of Symbols

REF	Catalogue Number	LOT	Batch code
			Manufacturer
			Date of manufacture
			STERILE R Sterilized using irradiation
			Do not reuse
			Use by date

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are trademarks of 3M.
© 2009 3M. All rights reserved.