**Sterile**

**Product Description**

3M™ Steri-Drape™ and 3M™ Steri-Drape™ 2 Incise Drapes are single use incise drapes that include an adhesive incise area. They are constructed of a film material that is capable of inhibiting strike through to reduce the risk of exposure to blood and body fluids. The sterile drape is applied to the prepped surgical site and the surgeon incises through the adhesive film and skin. The drape works by immobilizing bacteria on the skin and by creating a sterile surface at the beginning of surgery.

**Indications for Use**

3M™ Steri-Drape™ and Steri-Drape™ 2 Incise Drapes are indicated for use as an incise drape to cover patient's skin. It is intended for external use only. This product is intended to be used by trained medical professionals in a clinical/ surgical setting.

**Cautions:**

**Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury**

**CAUTION:** To reduce the risks associated with arcing, injury, or ineffective defibrillation, do not defibrillate through plastic incise drapes.

**CAUTION:** To reduce the risk of skin irritation or skin stripping:

- Do not stretch the drape or apply the drape under tension.
- Take greater care when applying and removing drapes on those with fragile skin.
- Follow recommended removal instructions

**CAUTION:** To reduce the risk of infection due to product or package degradation do not re-sterilize.

**CAUTION:** To avoid the risk of infection related to cross contamination and to ensure product integrity, do not reuse. Product is single use only.

**CAUTION:** To reduce the risk of electric shock if used in conjunction with potentially flammable products, follow standard safety protocols regarding ignition sources and product instructions regarding dry and dissipation times, if applicable.

**Directions for Use**

**SITE PREPARATION**

1. Disinfect the operative site and allow prep to fully dry.

**APPLICATION**

2. Hold incise drape a few inches above the skin. Grasp edges of non-adhesive film handle. Firmly peel off paper liner (see figure a).
3. When "stop" symbol appears, discontinue peeling (see figure b).
4. Using a sterile towel, press down firmly on the film, first contacting the skin along the intended incision line to ensure good skin contact (see figure c). Firmly smooth film into place working away from the intended incision line to achieve wrinkle free adhesion. Apply to site WITHOUT TENSION. When applying over a knee, bend the knee 20 degrees during application.
5. Peel the remainder of the drape off the paper liner and complete drape application.

**REMOVAL**

6. To remove, grasp the drape as close as possible to the junction of the drape to the skin. Place one hand on the skin for counter-traction. This prevents severe stretching.
7. Slowly and gently peel back drape at a 180-degree angle from the skin. This means folding the film back on itself and maintaining that angle while pulling the film away from the skin (see figure d & e).

**Storage/ Shelf Life/Disposal**

For best results store at room temperature. Avoid excessive heat and humidity. For shelf life information, refer to the expiration date on the package.

Dispose of the drape and its packaging per facility policy.

If the sterile packaging is damaged or unintentionally opened, discard the product and do not use.

**How supplied**

Product is sterile unless package is damaged or open.

3M™ Steri-Drape™ Incise Drape catalog numbers: 1035, 1040, 1050, 1051  
3M™ Steri-Drape™ 2 Incise Drape catalog numbers: 2040, 2045, 2050, 2051

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

**Symbol Glossary**

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Sterilized using irradiation		Indicates a medical device that has been sterilized using irradiation. Source: ISO 15223, 5.2.4
Do not resterilize		Indicates a medical device that is not to be resterilized. Source: ISO 15223, 5.2.6

**REF** 1040, 1050, 1051, 2040, 2045, 2050, 2051  
Made in USA with globally sourced materials

**3M Company**  
2510 Conway Ave.  
St. Paul, MN 55144 USA  
1-800-228-3957 (USA only)  
3M.com

**REF** 1035  
Made in Mexico  
**3M Company**  
2510 Conway Ave.  
St. Paul, MN 55144 USA  
1-800-228-3957 (USA only)  
3M.com

**EC REP** **3M Deutschland GmbH**  
Health Care Business  
Carl-Schurz-Str. 1  
41453 Neuss, Germany

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Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Caution		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Single sterile barrier system		Indicates a single sterile barrier system. Source: ISO 15223, 5.2.11

For more information see, [HCBGRegulatory.3M.com](http://HCBGRegulatory.3M.com)

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3

<p>9300 Winnetka Ave. N Brooklyn Park, MN 55445 (763) 488-5700</p>		<b>CUST ID#</b>	<b>PRINT METHOD</b>	<b>SGS#</b>	<b>ROUND</b>
				5272930-2	9
<b>Production Art</b> <b>Description:</b> 3M™ Steri-Drape™ Inzise Drape  <b>File Name:</b> 34872607337.indd <b>GA:</b> <b>Software:</b> InDesign CC 2019 <b>SGS Contact:</b> <b>Email:</b> <b>Printer:</b> <b>Printer Location:</b> <b>Supplier:</b>	<b>Item Spec#:</b> 34-8726-0733-7  <b>Supersedes#:</b>  <b>Structure#:</b>  <b>Cat/Product#:</b>  <b>Reference:</b>  <b>Requester:</b> Frank Boeshart  <b>Die # / Doc. Size:</b> 3.5" x 12.25" (36 pages saddle stitched)	<b>DATE</b>	<b>CHANGE LOG</b>	<b>INT</b>	
		08.27.20	New Build	SA	
		11.04.20	Update as per annotations	SA	
		11.13.20	Update as per Annotation	VP	
		12.10.20	Back from OMNI. Updated charts as they needed to be the most current 12-2020 Excel sheet	kat	
		12.16.20	Update as per symbol Added	SA	
		12.22.20	Updated document Size	SV	
		01.04.21	Final Release & English only.	Jer	
		01.12.21	Add page numbers and re-release	kat	
<b>INKS</b>  <b>PROCESS</b> <b>BLACK</b>	<b>MATCH</b> 	<b>SPECIAL</b> 	<b>LIMITATION OF LIABILITY</b> – The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. Should the supplied materials be found to be incorrect please contact an SGS representative immediately. The liability of SGS shall be limited to correcting its own product. SGS will not be responsible for any loss or damage including, but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost. <b>UPC NOTE</b> – Because of the differences in scanning equipment, SGS cannot guarantee the scanability of UPC barcodes. The Uniform Code Council recommends a UPC symbol be printed black on white at 100% with no truncation. All UPC codes are test scanned and verified before leaving our facility.		

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3M and Steri-Drape are trademarks of 3M. Used under license in Canada. © 2020, 3M. Unauthorized use prohibited. All rights reserved. 3M et Steri-Drape sont des marques de commerce de 3M. Utilisées sous licence au Canada. © 2020, 3M. Toute utilisation non autorisée est interdite. Tous droits réservés. Issue Date: 2020-05\_2 34-8726-0733-7