
en Indications for Use

The Ranger blood/fluid warmer is intended to warm blood, blood products, and liquids.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency trauma settings, or other areas where blood/fluids are infused.

Explanation of Signal Word Consequences
WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.
CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

WARNING: To reduce the risks associated with air embolism:

- Never infuse fluids if air bubbles are present in the fluid line.
- Bubble trap does not eliminate air automatically. See instructions on "Removing Air from the Disposable Set".
- Ensure all luer connections are tightened.
- Do not mount the Ranger blood fluid warming system higher than patient level, as air embolism may result.
- Always keep bubble trap in a vertical position during infusion

WARNING: To reduce the risks of biohazard exposure, cross-contamination or infection:

- Do not reuse. This product is designed for single patient use only.
- Do not use if caps are not in place.
- Ensure that all luer connections are securely tightened prior to priming set.
- Do not exceed 300 mm/Hg pressure; ensure inline hand pressure pumps or other pumps do not exceed 300 mm/Hg pressure.

WARNING: To reduce the risk of leaking and loss of therapy:

- Do not reuse. This product is designed for single patient use only.
- Ensure that all luer connections are securely tightened prior to priming set.

WARNING: To reduce the risks associated with potential blood loss:

- Do not use in combination with an extracorporeal circuit.

WARNING: To reduce the risk associated with patient exposure to leakage current that exceeds Type BF equipment limits:

- Do not use in combination with other medical equipment without verifying that the total leakage current from the combined equipment does not exceed the safety limits for Type BF equipment

CAUTION: To reduce the risk of leaking and loss of therapy:

- Do not exceed 300 mmHg pressure; ensure inline hand pressure pumps or other pumps do not exceed 300 mm/Hg pressure.
- Do not remove warming cassette from warming unit while in use

CAUTION: To reduce the risks associated with warming platelets, cryoprecipitate, or granulocyte suspensions:

- Follow the AABB Guidelines for the use of blood warming devices which caution against warming when administering platelets, cryoprecipitate, or granulocyte suspensions.

CAUTION: To reduce the risks associated with potential material toxicity:

- Do not use in combination with an extracorporeal circuit as the Ranger blood/fluid warming system is not designed for use with a circulating blood system; therefore, biocompatibility testing has not been completed to assess genotoxicity.

Instructions for Use

Do not use if package has been previously opened or damaged. Refer to the 3M™ Ranger™ blood/fluid warming unit model 245 Operator's Manual for set up and use of the warming system.




Figure 1: Illustration of Disposable Set

- A) Inlet line
- B) Blue inlet pinch clamp
- C) Fluid warming cassette
- D) Bubble trap
- E) Roller Clamp
- F) Patient line
- G) Injection port
- H) White outlet pinch clamp
- I) Patient connection

Figure 2

Peel To Open

Peel To Open


Figure 3

Instructions

Set-up and Priming of the warming set

1. Adjust warming unit height to patient level or lower.
2. Open warming set pouch where indicated by the arrows (Fig. 2)
3. Tighten all luer connections
4. Before priming, slide the fluid warming cassette (Fig. 1, C) into the slot in the warming unit (Fig. 3, #1).
5. Connect the short inlet line (Fig. 1, A), on the Ranger warming set to the fluid source.
6. Prime all tubing to purge air from the line (Fig. 3, #2). Invert the bubble trap (Fig. 1, D) while priming and fill.
7. Turn the bubble trap right side up and prime the patient line (Fig. 1, F).
8. Place the bubble trap, right side up, into the holder on the warming unit.
9. Close all clamps.
10. Turn the warming unit on (Fig. 3, #3).
11. Monitor fluid lines to ensure that they are air-free.

WARNING: Bubble trap does not eliminate air automatically. See instructions on "Removing Air from the Disposable Set".

The fluid warming set is now ready for use.

Removing Air from the Disposable Set

1. Close the clamp (Fig. 1, H) between the injection port (Fig 1, G) and the patient connection (Fig. 1, I).
2. Invert the bubble trap.
3. Insert a syringe into the injection port (Fig. 1, G) and aspirate air until the bubble trap and I.V. line are free of air.
4. Place the bubble trap right side up in the holder.
5. Open the clamp and continue the infusion.

Removing the Disposable Set

1. Close the inlet clamp proximal to the fluid warming cassette. (Fig. 1, B)
2. Open the outlet clamp distal to the cassette (Fig. 1, E, H).
3. Allow fluid to flow into the patient (this may take 2-3 seconds).
4. Remove the fluid warming cassette from the warming unit and discard according to institutional protocol.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Sterilized using ethylene oxide		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
Do not use if package is damaged or open		Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8
Sterile fluid path		Indicates the presence of a sterile fluid path within the medical device in cases when other parts of the medical device, including the exterior, might not be supplied sterile. Source: ISO 15223, 5.2.9
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2

Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device.
Unique device identifier		Indicates bar code to scan product information into patient electronic health record
Importer		Indicates the entity importing the medical device into the EU.
Single sterile barrier system		Indicates a single sterile barrier system. Source: ISO 15223, 5.2.11
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)

For more information see, HCBGRegulatory.3M.com

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