- ⓔn Universal Electrosurgical Pads
- ⓕr Plaques électrochirurgicales universelles
- ⓓe Universalneutralelektroden
- ⓔt Piastra universală per elettrochirurgia
- ⓔs Almohadillas electroquirúrgicas universales
- ⓓl Universele elektrochirurgische pads
- ⓔv Universella elektrokirurgiska dynor
- ⓓa Elektrokirurgiske universalelektroder
- ⓓo Elektrokirurgiske universalplater
- ⓕi Universaalit sähkökirurgiset diatermialevyt
- ⓔt Placas Universais para Eletrocirurgia, Séries
- ⓔl Γενικά Ηλεκτροχειρουργικά Επιθέματα
- ⓔl Uniwersalne elektrody bierne
- ⓕu Univerzális elektrosebészeti érintkezőtápaszok
- ⓔs Univerzální elektrochirurgické podložky
- ⓔk Univerzálne elektrochirurgické podložky
- ⓔl Univerzalne elektrokirurške blazinice
- ⓔt Universaalsed elektrokirurgilised padielektroodid
- ⓔv Universālie elektroķirurgiskie paliktņi
- ⓔt universalūs elektrochirurginiai elektrodai
- ⓔo Plăcuțe electrochirurgicale universale
- ⓔu Универсальные электрохирургические электроды
- ⓕr Univerzalni elektrokirurški jastučići
- ⓔg Универсални електрохирургични подложки
- ⓔr Univerzalni elektrohirurški jastučići
- ⓔr Üiversal Elektro Cerrahi Pedleri
- ⓔa 汎用電気手術用パッド
- ⓔh 通用型电刀负极板
- ⓔr الضمادات الجراحية الكهربائية العام
- ⓔq PIASTRAT elektrokirurgjikalë universale
- ⓔk Универзални електрохируршки подлоги

3M™ Universal Electrosurgical Pads 9160, 9160F, 9165E

en

General Use

Read and save this document. Make sure everyone who will use this product knows and understands all information contained within this document and AORN recommended practices for electrosurgery. READ WARNING

Intended Use

3M™ Universal Electrosurgical Pads (9160, 9160F, 9165E) are designed to work with most electrosurgical units (ESU's) for virtually every surgical application where electrosurgery is utilized to provide a safe return path for electrosurgical current. Split style Universal Electrosurgical Pads are for use with ESUs that have a CQMS (i.e. REM, ARM, NESSY, etc.). The 3M™ Universal Electrosurgical Pads are designed to be used on any patient where full skin contact and a suitable placement site can be obtained. There are no patient weight restrictions for use of this product. Use of this product for unintended applications could lead to an unsafe condition.

They are intended to be used by healthcare professionals in hospitals and surgical centers

Product Description

3M™ Universal Electrosurgical Pads (i.e. grounding pads, neutral electrodes) consist of a conductive adhesive area surrounded by a nonconductive border adhesive. The pads are supplied pre-corded or non-corded, have fluid resistant non-woven backing, and have a conductor area of 15in² (97cm²). The pads are single use, disposable and non-sterile.

WARNING

Improper use of Universal Electrosurgical Pads can cause electrosurgical burns or pressure necroses. For patient safety, follow all of the instructions below. Failure to follow any of these instructions increases the risk of electrosurgical burns or pressure necroses.

Instructions for Safe Use

1. Use Appropriate Pads, Equipment, and Accessories

- Does the electrosurgical generator have a Contact Quality Monitoring System (e.g. REM™, ARM™, NESSY™)?
 - If **YES** - use **split-style** Universal Pads (9160/9160F/9165E).
- Use 3M™ cables and adapters as required with 3M™ Universal Pads.
- Check expiration date on package. 3M™ Universal Pads are safe to use for 14 days after package is opened.
- To avoid risk of alternate site burns, do not allow the patient to contact earth grounded metal or items which are capacitively coupled to earth ground.
- To avoid increased risk of burns and infection related to cross contamination, do not reuse pads.

2. To Reduce the Risk of Burns, Do Not Overload the Universal Pad with Too Much Current

- Do not activate the electrosurgical device or active accessory for more than 60 seconds in any 2-minute period, as this will overload the Universal Pad with current and may result in a patient burn.
- Any combination of **high power, long activation time, and a conductive irrigant** (e.g., saline) may overload the Universal Pad with current and may result in a patient burn. To reduce this risk:
 - Use non-conductive solutions unless specific medical reasons indicate otherwise. Conductive fluids (e.g., blood or saline) in direct contact with or near an active electrode may carry electrical current and/or heat away from target tissues, which may lead to unintended burns to the patient.
 - Use the lowest possible power setting.
 - Use short activation times. If long activation is necessary, allow time between activations to allow the tissue under patient plate to cool.
 - Use two Universal Pads with a Y-adapter.
 - If you do not receive the desired surgical effect, stop and verify the correct distention/irrigation solution and good Universal Pad contact before proceeding with electrosurgery or increasing the power setting.

3. Select an Appropriate Site

To reduce the risk of burns and pressure necroses:

- Select a smooth, well-vascularized, muscular area close to surgical site that allows full Universal Pad-to-skin contact.
- Site must be clean, dry, and free of hair. Remove hair at application site.
- Locate Universal Pad closer to the surgical site than to the ECG electrodes.
- Remove metal jewelry from patient.
- Avoid placement over bony prominences, metal prostheses, or scar tissue.
- Avoid placement such that current flows through a metal prosthesis or conductive implant. For patients with implanted devices, contact device manufacturer for precautions to avoid interference.
- Avoid placement over surgical prep solutions containing iodine (Betadine, Povidone-iodine etc.)
- Do not apply Universal Pad where fluids may pool.
- Do not apply Universal Pad over injection site.
- Select a suitable site remote from any warming device.
- Do not place the Universal Pad under the patient. Weight bearing sites have restricted blood flow and may reduce the performance of the Universal Pad.
- Do not place the Universal Pad across skin gaps such as the gap between the buttocks or gap between the arms and body.

4. Pad Application


To reduce the risk of burns and pressure necroses:

- Inspect Universal Pad, cord, and cable. Do not use if cut, modified, or damaged.
- Remove clear liner from Universal Pad before applying to patient.
- Apply one end of Universal Pad and smoothly press to other end. Avoid air entrapment.

- Avoid stretching or folding either Universal Pad or patient's skin.
- Smooth down the Universal pad edges after application to ensure complete Universal Pad adherence.
- Do not use electrode gel or barrier cream under Universal Pad.
- Do not wrap Universal Pad completely around a limb. Do not overlap the Universal Pad edges.
- Do not place the Universal Pad over compromised skin.
- Do not reposition Universal Pad after initial application. If patient is repositioned, confirm full pad-to-skin contact and integrity of all connections.
- Do not place compression stocking or device over Universal Pad.
- Do not coil or wrap cord or cable around patient's limb or metal object.
- Do not allow cord or cable to lie on or under patient.
- Do not place cable clamp under patient.
- Use ECG cables with RF suppressors/chokes to prevent electrosurgical current from flowing through the ECG electrodes. Any ECG electrodes should be placed as far as possible from the surgical electrodes. Needle monitoring electrodes are not recommended.

Use of two Universal Pads with Y-adapter:

- Patients with dry skin, adipose tissue, and/or poor vascularization may generate a high impedance alarm and may require two Universal Pads.
- Do not plug cords into Y-adapter until after each Universal Pad has been applied.
- Preferred placement of each pad is bilaterally (i.e., left and right side) equally distant from surgical site.
- For information about Y-adapters in Europe, please contact your local 3M representative. In the US and other countries outside of Europe, use the 3M™ 1157-C Y-adapter.

5. Pad removal

- Do not remove by pulling on cable or cord.
- Start at corner. Peel back slowly at 180-degree angle to prevent skin trauma.

Notice to manufacturers and users of radio frequency (RF) medical devices and active accessories:

All 3M™ Universal Electrosurgical Pads conform to Section 201.15.101.5 of the ANSI/AAMI/IEC 60601-2-2:2017 Standard, which specifies that a grounding pad must be capable of carrying a current of 700 milliamperes (mA) for a continuous period of 60 seconds. Manufacturers of RF medical devices or accessories should not recommend 3M™ Universal Electrosurgical Pads for use with RF medical devices or accessories that can deliver a current load to the Universal Pad that exceeds these limits. 3M™ certifies that the 3M™ 9100 Series Split Style Universal Electrosurgical Pads, when used in accordance with these Instructions for Use, meet the requirements of the ANSI/AAMI/IEC 60601-2-2:2017 standard for compatibility when used with high frequency (HF) electrosurgical generators with CQM systems that operate with both a maximum impedance limit (not to exceed 150 ohms) and a differential (dynamic) impedance limit (not to exceed 40%).

For questions regarding compatibility of 3M™ Universal Electrosurgical Pads with specific generators, in the USA, please contact 3M at 1-800-228-3957. Outside the USA, please contact your 3M representative.

Shelf Life:

For shelf life refer to the expiration date printed on each package.

Disposal:

Dispose of contents/container in accordance with the local/regional/national/international regulations.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

Symbols Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Refer to instruction manual/booklet		To signify that the instruction manual/booklet must be read. ISO 7010-M002 • This symbol is blue on the packaging
Does not contain natural rubber latex		Indicates there is no presence of dry natural rubber nor natural rubber latex as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
CE Mark 2797		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
UKCA Mark 0086		Indicates conformity to all applicable regulations and/or directives in the United Kingdom (UK), for products placed on the market in Great Britain (GB) with UK Approved Body involvement. (0086 to be replaced with UK Approved Body ID number)
Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch

For more information, see HCBGRegulatory.3M.com

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REM is a trademark of Valleylab.
ARM is a trademark of Conmed/Aspen Labs.
NESSY is a trademark of ERBE.
3M is a trademark of 3M. Used under license in Canada.
© 2022, 3M. Unauthorized use prohibited. All rights reserved.
3M Deutschland GmbH