MD

	Importer		Medical Device
	Importateur		Dispositif médical
	Importeur in die EU		Medizinprodukt
	Importatore		Dispositivo medico
	Importador		Producto sanitario
	Importeur		Medisch hulpmiddel
	Importör		Medicinteknisk produkt
	Importer		Medicinsk udstyr
	Importer		Medisinsk utstyr
	Maahantuaja		Lääkinnällinen laite
	Importador		Dispositivo médico
	Εισαγωγέας		Ιατροτεχνολογικό προϊόν
	Importer		Wyrób medyczny
	Importőr		Orvostechnikai eszköz
	Dovozce		Zdravotnický prostředek
	Importér		Zdravotnícka pomôcka
	Uvoznik		Medicinski pripomoček
	Maaletooja		Meditsiiniseade
	Importētājs		Medicīniska ierīce
	Importuotojas		Medicinos priemonė
	Importator		Dispozitiv medical
	Імпортёр		Медичний пристрій
	Uvoznik		Medicinski proizvod
	вносител		медицинско изделие
	Uvoznik		Medicinsko sredstvo
	İthalatçı		Tıbbi cihaz
	Импорттаушы		Медициналық құрылғы
	ผู้นำเข้า		เครื่องมือแพทย์
	Importuesi		Pajisje Mjekësore
	Увозник		Медицинско средство

CH	REP
----	-----

	Swiss Authorised Representative
	Représentant autorisé en Suisse
	Bevollmächtigter in der Schweiz
	Rappresentante svizzero autorizzato
	Representante autorizado de Suiza
	Zwitserse gemachtigde vertegenwoordiger
	Schweizisk auktoriserad representant
	Schweizisk autoriseret repræsentant
	Autorisert representant i Sveits
	Valtuutettu edustaja Sveitsissä
	Representante autorizado suíço
	Εξουσιοδοτημένος αντιπρόσωπος Ελλάδας
	Autoryzowany przedstawiciel w Szwajcarii
	Svájci meghatalmazott képviselő
	Švýcarský zplnomocněný zástupce
	Autorizovaný zástupca pre Švajčiarsko
	Pooblaščení zastopnik v Švici
	Šveitsi volitatud esindaja
	Šveices pilnvarots pārstāvis
	Igaliotas atstovas Šveicarijoje
	Reprezentant autorizat elvețian
	Уповноважений представник у Швейцарії
	Ovlašteni predstavnik za Švicarsku
	Упълномощен представител в Швейцария
	Ovlašćeni predstavnik u Švajcarskoj
	İsviçre Yetkili Temsilcisi
	Швейцарияның уәкілетті өкілі
	ตัวแทนที่ได้รับอนุญาตจากสวิส
	Përfaqësuesi zviceran i autorizuar
	Овластен претставник на Швајцарија

3M™ Tegaderm™ Chlorhexidine Gluconate gel pad device + I.V. Port Dressing

en

Product Description

3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad Device is used to protect vascular access sites. The Chlorhexidine Gluconate Gel Pad Device is breathable and transparent, allowing continuous site observation.

The transparent film provides an effective barrier against external contamination including fluids (waterproof), bacteria, viruses* and yeast.

The antimicrobial CHG Gel Device consists of an integrated gel pad containing 2% w/w Chlorhexidine Gluconate (CHG), a well-known antiseptic agent with broad spectrum antimicrobial and antifungal activity. The gel absorbs fluid.

**In vitro* testing shows that the transparent film of the gel pad device provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage. The barrier to viruses is due to the physical properties of the dressing, rather than the ancillary properties of CHG.

In vitro testing (time kill and zone of inhibition) demonstrates that the CHG gel pad of the device has an antimicrobial effect against a variety of gram-positive and gram-negative bacteria and yeast.

Intended Use / Intended Purpose:

3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad Device can be used to protect catheter sites. Common applications include protecting intravascular catheters and percutaneous devices. Tegaderm™ CHG is intended to reduce skin colonization and catheter colonization and to suppress regrowth of microorganisms commonly related to bloodstream infections. Tegaderm™ CHG is intended to reduce catheter-related bloodstream infections (CRBSI) in patients with central venous or arterial catheters.

Warnings

DO NOT PUNCTURE CHG GEL WITH NONCORING NEEDLE.

DO NOT USE TEGADERM™ CHG DRESSINGS ON PREMATURE INFANTS OR INFANTS YOUNGER THAN 2 MONTHS OF AGE. USE OF THIS PRODUCT ON PREMATURE INFANTS MAY RESULT IN HYPERSENSITIVITY REACTIONS OR NECROSIS OF THE SKIN.

THE SAFETY AND EFFECTIVENESS OF TEGADERM™ CHG DRESSINGS HAS NOT BEEN EVALUATED IN CHILDREN UNDER 18 YEARS OF AGE. FOR EXTERNAL USE ONLY. DO NOT ALLOW THIS PRODUCT TO CONTACT EARS, EYES, MOUTH OR MUCOUS MEMBRANES.

DO NOT USE THIS PRODUCT ON PATIENTS WITH KNOWN HYPERSENSITIVITY TO CHLORHEXIDINE GLUCONATE. THE USE OF CHLORHEXIDINE GLUCONATE CONTAINING PRODUCTS HAS BEEN REPORTED TO CAUSE IRRITATIONS, SENSITIZATION, AND GENERALIZED ALLERGIC REACTIONS. IF ALLERGIC REACTIONS OCCUR, DISCONTINUE USE IMMEDIATELY, AND IF SEVERE, CONTACT A PHYSICIAN.

Hypersensitivity reactions associated with topical use of chlorhexidine gluconate have been reported in several countries. The most serious reactions (including anaphylaxis) have occurred in patients treated with lubricants containing chlorhexidine gluconate, which were used during urinary tract procedures. Caution should be used when using chlorhexidine gluconate containing preparations, and the patient should be observed for the possibility of hypersensitivity reactions.

Precautions

3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad Device should not be placed over infected wounds.

It is not intended to be used as a treatment of percutaneous device-related infections.

In the case of clinical wound infection, systemic antibacterials should be used if indicated.

Any active bleeding at the insertion site should be stabilized before applying the device.

Do not stretch the device during application. Mechanical skin trauma may result if the device is applied with tension.

The skin should be clean, dry, and free of detergent residue to prevent skin irritation and to ensure good adhesion.

Do not reuse. Reuse may result in compromising product integrity and lead to device failure.

Use care and cover the CHG gel pad device to protect from water when showering.

The CHG gel pad device has not been tested for use during radiation therapy.


Clinical Trial Results: A randomized, controlled clinical trial consisting of 1879 subjects with 4163 central venous and arterial catheter insertion sites was conducted at 11 hospitals (1). Results showed that the use of Tegaderm™ CHG resulted in a statistically significant 60% reduction in the incidence of catheter-related bloodstream infections (P=0.02). Study results also demonstrate a statistically significant reduction in skin colonization (P<0.001) and catheter colonization (P<0.0001) in the chlorhexidine vs. non-chlorhexidine group.

Variable	Non-chlorhexidine vs. chlorhexidine dressings (941 patients/2055 catheters vs. 938 patients/2108 catheters)
Catheter-related bloodstream infection	
Incidence densities (n per 1000 catheter-days)	1.3 vs. 0.5
Hazard ratio	0.402 [0.186 to 0.868], <i>P</i> =0.02
Catheter colonization	
Incidence densities (n per 1000 catheter-days)	10.9 vs. 4.3
Hazard ratio	0.412 [0.306 to 0.556], <i>P</i> <0.0001

(1) Timsit JF et al. Randomized Controlled Trial of Chlorhexidine Dressing and Highly Adhesive Dressing for Preventing Catheter-Related Infections in Critically Ill Adults *Am. J. Respir. Crit. Care Med.* 2012; 186:1272-1278

A summary of safety and clinical performance (SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI. The EUDAMED public website is <https://ec.europa.eu/tools/eudamed>. The product basic UDI-DI is 06082238401010000000039AM.

Instructions for Use: Refer to packaging figures and those included in package insert. Failure to follow the manufacturer's instructions for use may result in complications including skin irritation and/or maceration.

Site Preparation: Prepare port site according to institution protocol. Clipping of hair at the site may improve CHG gel pad device adhesion. Shaving is not recommended. The skin should be clean, dry and free of detergent residue.

Any active bleeding at the insertion site should be stabilized before application of CHG gel pad device.

Post Needle Insertion Application of CHG Gel Pad Device (Figures A1 to A5)

- Open package and remove the sterile CHG gel pad device using aseptic technique.
- Prepare port site and insert noncoring needle per facility protocol, allowing for a space under the noncoring needle for the CHG gel pad device of about ½ cm. Allow all preps and protectants to dry completely before applying CHG gel pad device to prevent skin irritation and to ensure good adhesion.
- Under green tab labeled **1**, grasp the soft cloth tabs. Peel and remove top liner and discard (Figure A1).
- Do not unfold the liner. Pull the soft cloth tabs apart to widen the slit of the gel (Figure A2).
- Using the soft cloth tabs center the device under the noncoring needle (Figure A3).
- Lift and grasp the soft cloth tabs and remove liner labeled **2** in direction of arrow (Figure A4).
- If needed, using the soft cloth tabs, bring slit edges closer together. Do not stretch device during application. Mechanical skin trauma may result if the device is applied with tension (Figure A5).
- Use firm pressure to smooth down device border and enhance adhesion.
- See application instructions for 3M™ I.V. Port Dressing.

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
3M and Tegaderm are trademarks of 3M.
Used under license in Canada.
© 2022, 3M. Unauthorized use prohibited.
All rights reserved.

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Issue Date: 2022-09, rev. 4
34-8728-3429-5

EC REP 3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss, Germany

Simultaneous Application of the CHG Gel Pad Device and Noncoring Needle (Figure B1 to B5)

1. Open package and remove the sterile CHG gel pad device using aseptic technique.
2. Prepare port site and pre-fill noncoring needle system according to facility protocol. Allow all preps and protectants to dry completely before application to prevent skin irritation and to ensure good adhesion.
3. Under green tab labeled 1, grasp the soft cloth tabs. Peel and remove top liner and discard (Figure B1).
4. Do not unfold device liner. Leave protective sheath in place over needle. Orient the CHG gel pad device facing upward and slide needle all the way into base of slit (Figure B2).
5. Wrap the CHG gel pad device around noncoring needle. Remove protective sheath and insert needle into port site per facility protocol (Figure B3).
6. Before removing liner, ensure needle is positioned all the way into base of slit. The CHG gel pad device may be adjusted under noncoring needle by grasping soft cloth tabs. Press gently on soft cloth tabs and remove liner labeled 2 in direction of arrow (Figure B4).
7. If needed, using soft cloth tabs, bring slit edges closer together. Do not stretch device during application. Mechanical skin trauma may result if device is applied with tension (Figure B5).
8. Use firm pressure to smooth down device border and enhance adhesion.
9. See application instructions for 3M™ I.V. Port Dressing.

Site Care:

1. The site should be observed daily for signs of infection or other complications. If infection is suspected, remove the CHG gel pad device, inspect the site directly, and determine appropriate medical intervention. Infection may be signaled by fever, pain, redness, swelling, or unusual odor or discharge.
2. Inspect the CHG gel pad device daily and change as necessary, in accordance with facility protocol; at a minimum of every 7 days, per Centers for Disease Control and Prevention (CDC) recommendations. CHG gel pad device changes may be needed more frequently with highly exudative sites.
3. Use care and cover the CHG gel pad device to protect from water when showering. The 3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad Device should be changed as necessary:
 - If the device becomes loose, soiled or compromised in any way
 - If the site is obscured or no longer visible
 - If there is visible drainage outside the CHG gel pad
 - If the CHG gel pad device appears to be saturated or overly swollen*

*Note: To test if the CHG gel pad is fully saturated, lightly press down on a corner of the gel pad with your finger. If the CHG gel pad remains displaced once your finger is removed, the CHG device should be changed.

Note: Tegaderm™CHG gel pad device is not designed to absorb large quantities of blood or fluid.

Removal

CHG Gel Pad Device: See removal instructions for 3M™ I.V. Port Dressing. Remove non-coring needle per facility protocol. Remove CHG gel pad device using soft cloth tabs. A medical adhesive solvent can be used to help remove the device border. If needed, use sterile alcohol swabs or wipes, or sterile solutions (i.e., sterile water or normal saline) to facilitate removal of the gel pad. Continue the low and slow removal method until the CHG gel pad device is completely removed.

Care should be taken not to dislodge catheters or other medical devices when the CHG gel pad device is removed. Support the skin and catheter while removing the CHG gel pad device.

Note: Dispose of product per facility protocol.

Storage/ Shelf Life/Disposal

For best results, store in a cool, dry place. For shelf life, refer to the expiration date on the package.

Sterility of the product is guaranteed unless individual package is damaged or open.

If the sterile packaging is damaged or unintentionally opened, discard the product, do not use.

If you have any questions or comments, contact the 3M Health Care Customer Help Line at 01509 611611 or go to www.3M.com.

How supplied / Ordering Information

Catalog #	CHG Gel Pad Device Size	Average Amount of CHG per CHG Gel Pad Device (mg based on gel pad size)
1665R 	CHG Gel Pad Device Dimensions 6.2 cm x 4.9 cm (2 2/16 in x 1 15/16 in)	30
1664R 	CHG Gel Pad Device Dimensions 6.2 cm x 4.9 cm (2 2/16 in x 1 15/16 in)	30

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) of local regulatory authority.

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Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Sterilized using irradiation		Indicates a medical device that has been sterilized using irradiation. Source: ISO 15223, 5.2.4
Do not resterilize		Indicates a medical device that is not to be resterilized. Source: ISO 15223, 5.2.6
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2

			
en Manufacturer fr Fabricant de Hersteller it Produttore es Fabricante nl Fabrikant sv Tillverkare da Producent no Produsent fi Valmistaja pt Fabricante el Κατασκευαστής pl Wytwórca hu Gyártó cs Výrobce sk Výrobca sl Proizvajalec et Tootja lv Ražotājs lt Gamintojas ro Producător uk Виробник hr Proizvođač bg Производител sr Proizvođač tr Üretici kk Өндіруші th ผู้ผลิต sq Prodhuesi mk Производител	en Date of manufacture fr Date de fabrication de Herstellungsdatum it Data di produzione es Fecha de fabricación nl Productiedatum sv Tillverkningsdatum da Fremstillingsdato no Produksjonsdato fi Valmistuspäivä pt Data de fabrico el Ημερομηνία κατασκευής pl Data produkcji hu Gyártási időpont cs Datum výroby sk Dátum výroby sl Datum izdelave et Tootmise kuupäev lv Ražošanas datums lt Gamybos data ro Data de fabricație uk Дата виготовлення hr Datum proizvodnje bg Дата на производство sr Datum proizvodnje tr Üretim Tarihi kk Өндірілген күні th วันผลิต sq Data e prodhimit mk Датум на производство	en Use-by date fr A utiliser avant de Verwendbar bis it Utilizzabile fino al es Fecha de caducidad nl Te gebruiken tot sv Bäst före datum da Anvendes inden no Utløpsdato fi Viimeinen käyttöpäivä pt Data de validade el Ημερομηνία ανάλωσης pl Użyć do daty hu Lejárati idő cs Použitelné do sk Použitelný do sl Uporabno do et Kõlblikkusaeg lv Izmantojams līdz lt Naudoti iki ro Data limită de utilizare uk Термін при-датності hr Rok valjanosti bg Да се използва преди sr Rok upotrebe tr Son kullanma tarihi kk Дейін пайдаланыңыз th ใช้ภายในวันที่ sq Afati i përdorimit mk Употребливо до	en Batch code fr Batch code de Fertigungslosnummer, Charge it Numero di lotto es Código de lote nl Lotnummer sv Partikod da Batchkode no Batchkode fi Eräkoodi pt Código do lote el Κωδικός παρτίδας pl Kod partii hu Tételszám cs Číslo šarže sk Číslo šarže sl Številka serije et Partii number lv Sērijas numurs lt Partijos kodas ro Cod de lot uk Номер партії hr Serijski broj bg Код на партида sr Šifra partije tr Seri kodu kk Партия коды th รหัสชุด sq Numri i ngarkesës mk Сериски код

Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
Does not contain natural rubber latex		Indicates there is no presence of dry natural rubber nor natural rubber latex as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Single sterile barrier system		Indicates a single sterile barrier system. Source: ISO 15223, 5.2.11
Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

For more information, see HCBGregulatory.3M.com

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